



January 22, 2018

Techlab, Inc.  
Donna Link  
Director Regulatory and Compliance  
2001 Kraft Drive  
Blacksburg, Virginia 24060

Re: K173217

Trade/Device Name: Campylobacter Quik Chek  
Regulation Number: 21 CFR 866.3110  
Regulation Name: Campylobacter fetus serological reagents  
Regulatory Class: Class I  
Product Code: LQP  
Dated: September 29, 2017  
Received: October 25, 2017

Dear Donna Link:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

 Ribhi Shawar -S

For

Uwe Scherf, M.Sc., Ph.D.

Director

Division of Microbiology Devices

Office of In Vitro Diagnostics

and Radiological Health

Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K173217

Device Name

CAMPYLOBACTER QUIK CHEK

Indications for Use (Describe)

The CAMPYLOBACTER QUIK CHEK™ test is a rapid membrane enzyme-linked immunosorbent assay for the qualitative detection of a Campylobacter-specific antigen in human fecal specimens. The CAMPYLOBACTER QUIK CHEK™ test is designed to detect *C. jejuni* and *C. coli* from patients with signs and symptoms of gastroenteritis. The test is intended for use with preserved fecal specimens in transport media and unpreserved fecal specimens. Test results should be considered in conjunction with clinical findings and patient history.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## **CAMPYLOBACTER QUIK CHEK™ 510(k) SUMMARY**

This summary of 510(k) safety and effectiveness is being submitted in accordance with the requirements of 21 CFR 807.92.

### **Applicant/Contact Information:**

Date Prepared: January 22, 2018  
Name: TECHLAB, Inc.  
Address: 2001 Kraft Drive  
Corporate Research Center  
Blacksburg, VA 24060 USA

Contact Person: Donna T. Link  
Phone Number: 540-953-1664  
Email: dlink@techlab.com

### **1.1 Manufacturing Facility Address**

TECHLAB, Inc.  
20 Corporate Drive  
Radford, VA 24141 USA

### **1.2 Product and Trade Name of the Device**

*CAMPYLOBACTER QUIK CHEK™*

### **1.3 Common Name or Classification Name**

*Campylobacter* spp. detection test

### **1.4 Classification and Regulation**

Class I  
21 CFR 866.3110; *Campylobacter fetus* serological reagents

### **1.5 Product Code**

LQP – *Campylobacter* spp.

### **1.6 Panel**

83 Microbiology

### **1.7 Reason for Premarket Notification**

The development of a new rapid membrane enzyme immunoassay for the qualitative detection of *Campylobacter* spp. in a single use cassette.

## **Intended Use**

The *CAMPYLOBACTER QUIK CHEK*™ test is a rapid membrane enzyme-linked immunosorbent assay for the qualitative detection of a *Campylobacter*-specific antigen in human fecal specimens. The *CAMPYLOBACTER QUIK CHEK*™ test is designed to detect *C. jejuni* and *C. coli* from patients with signs and symptoms of gastroenteritis. The test is intended for use with preserved fecal specimens in transport media and unpreserved fecal specimens. Test results should be considered in conjunction with clinical findings and patient history.

## **Explanation**

Worldwide, *Campylobacter* species are the most common cause of bacterial gastroenteritis, with 400-500 million cases of diarrhea each year. Infants in developing countries are at even greater risk, as are travelers to those countries. *Campylobacter*-associated gastroenteritis is estimated to affect nearly 1 million people a year in the USA. In approximately 1 of 1000 cases, *Campylobacter jejuni* is closely linked to the subsequent development of Guillian-Barre Syndrome, an acute auto-immune paralysis. *C. jejuni* infection has also been associated with reactive arthritis in both children and adults. When individuals with severe symptoms of gastroenteritis seek medical help, the clinician is faced with multiple possible causes that can present with similar clinical features (e.g., diarrhea, nausea, vomiting, fever, abdominal pain) but that require very different, often conflicting, types of treatment.

For *Campylobacter*, the current standard for identification is bacterial culture followed by microscopic examination of the organisms. Although this traditional method is straightforward, it has two major limitations. First, pathogenic species of *Campylobacter* are microaerophilic or strictly anaerobic, so that exposure of culture or feces to environmental oxygen leads to death or inactivation of the bacteria. Thus, during transport or storage of specimens under aerobic conditions, the number of viable organisms can decrease, leading to potentially inaccurate culture results. Second, *Campylobacter* species are slow-growing, requiring from 48-72 hours before reaching a point where the culture can safely be reported as negative. Such delays can leave the clinician in a quandary and the patient with non-specific, ineffective, or even inappropriate treatment.

The *CAMPYLOBACTER QUIK CHEK*™ test allows detection of *Campylobacter jejuni* and *Campylobacter coli*, the species most commonly associated with human disease, in less than 30 minutes. Furthermore, the *CAMPYLOBACTER QUIK CHEK*™ test does not rely on bacterial viability, and can be performed on the bench-top with samples that have been exposed to air.

## **Device Description**

The *CAMPYLOBACTER QUIK CHEK*™ test uses antibodies that recognize a *Campylobacter*-specific antigen in human fecal samples. The device contains a *Reaction Window* with two vertical lines of immobilized antibodies. The test line ("T") contains antibodies against a *Campylobacter*-specific antigen. The control line ("C"), contains anti-IgG antibodies. The *Conjugate* consists of antibodies to a *Campylobacter*-specific antigen coupled to horseradish peroxidase. To perform the test, a fecal specimen is added to a tube containing a mixture of *Diluent* and *Conjugate*. The diluted sample-conjugate mixture is added to the *Sample Well* and the device is allowed to incubate at room temperature for 15 minutes. During the incubation, the *Campylobacter*-specific antigens in the sample bind to the antibody-peroxidase conjugate. The antigen-antibody complexes migrate through a filter pad to a membrane where they are captured by the immobilized anti-*Campylobacter* antibodies in the line. The *Reaction Window* is subsequently washed with *Wash Buffer*, followed by the addition of *Substrate*. After a 10-minute incubation, the "T" reaction is examined visually for the appearance of a vertical blue line. A blue line indicates a positive test. A positive "C" reaction, indicated by a vertical blue line, monitors/confirms that the sample and reagents were added correctly, the reagents were active at the time of performing the assay, and that the sample migrated properly through the

*Membrane Device*. It also confirms the reactivity of the other reagents associated with the assay and that the results are valid.

### **Materials Provided**

- **Membrane Devices** – 25, each pouch contains 1 device
- **Conjugate (2.5 mL)** – Antibody to a *Campylobacter*-specific antigen coupled to horseradish peroxidase in a buffered protein solution
- **Diluent (22 mL)** – Buffered protein solution with graduated dropper assembly
- **Positive Control (2 mL)** – *Campylobacter*-specific antigen in a buffered protein solution
- **Wash Buffer (12 mL)** – Buffered solution with graduated dropper assembly
- **Substrate (3.5 mL)** – Solution containing tetramethylbenzidine
- **Disposable plastic transfer pipettes** – graduated at 25 µL, 100 µL, 200 µL, 300 µL, 400 µL and 500 µL

The predicate device (ImmunoCard STAT!® CAMPY) and the *CAMPYLOBACTER QUIK CHEK*™ test both detect *Campylobacter* spp. (*C. jejuni* and *C. coli*) in fecal specimens and are substantially equivalent in principle. The following table shows a comparison of both devices.

| Similarities        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Item                | Device<br>K173217                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | Predicate<br>K090700                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Indications for use | The <i>CAMPYLOBACTER QUIK CHEK</i> ™ test is a rapid membrane enzyme-linked immunosorbent assay for the qualitative detection of a <i>Campylobacter</i> -specific antigen in human fecal specimens. The <i>CAMPYLOBACTER QUIK CHEK</i> test is designed to detect <i>C. jejuni</i> and <i>C. coli</i> from patients with symptoms of gastroenteritis. The test is intended for use with unpreserved human fecal specimens and fecal specimens that are in transport media. Test results should be considered in conjunction with clinical findings and patient history. | ImmunoCard STAT! CAMPY is an immunochromatographic rapid test for the qualitative detection of specific <i>Campylobacter</i> antigens in human stool. ImmunoCard STAT! CAMPY detects <i>C. jejuni</i> and <i>C. coli</i> in human stool, where stool may be either unpreserved or preserved in Cary-Blair-based transport media. Test results are to be used in conjunction with information available from the patient clinical evaluation and other diagnostic procedures. |
| Measured analyte    | Detection of <i>Campylobacter</i> -specific antigens ( <i>C. jejuni</i> and <i>C. coli</i> )                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | Detection of <i>Campylobacter</i> -specific antigens ( <i>C. jejuni</i> and <i>C. coli</i> )                                                                                                                                                                                                                                                                                                                                                                                 |
| Type of Test        | Qualitative                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | Same                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Controls            | Positive and negative control included in the kit Internal "C" Control line                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | Same                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Format              | Single Use Membrane Cassette                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | Same                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Target Population   | Persons suspected of having <i>Campylobacter</i> infection                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              | Same                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| Storage             | Refrigerated (2°C – 8°C)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | Same                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |

There are no differences between the subject device and the predicate(s) with respect to indications and intended use.

| Differences     |                                                       |                                               |
|-----------------|-------------------------------------------------------|-----------------------------------------------|
| Item            | Device<br>K173217                                     | Predicate<br>K090700                          |
| Specimen Type   | Fecal specimens in Cary-Blair and C&S Transport Media | Fecal specimens in Cary-Blair Transport Media |
| Time to Result  | <30 minutes                                           | < 25 minutes                                  |
| Antibody Format | Monoclonal/Polyclonal                                 | Monoclonal/Monoclonal                         |

### **Summary of Performance Data**

#### **Prospective Study**

The performance of the *CAMPYLOBACTER QUIK CHEK*™ test was evaluated at 4 independent sites. Prospective incoming fecal specimens were collected and tested by culture and the *CAMPYLOBACTER QUIK CHEK*™ test. The following table shows a summary of the clinical performance of the *CAMPYLOBACTER QUIK CHEK*™ test for all 4 sites combined. The results of the study show that the *CAMPYLOBACTER QUIK CHEK*™ test exhibited a sensitivity of 97.1%, and a specificity of 99.1% with culture.

#### **Age and Gender Distribution**

Age information was available for 1552 patients. The ages ranged from less than 1 year to 100 years. Of the 1552 patients, 15.7% were ≤ 18 years. The gender identification was 38.7% females and 61.3% males. No difference in test performance was observed based on patient age or gender.

#### ***CAMPYLOBACTER QUIK CHEK*™ test versus Culture**

| N = 1552                                  | Culture Positive | Culture Negative |
|-------------------------------------------|------------------|------------------|
| <i>CAMPYLOBACTER QUIK CHEK</i> ™ Positive | 34               | 13*              |
| <i>CAMPYLOBACTER QUIK CHEK</i> ™ Negative | 1**              | 1504             |

|             |       | 95% Confidence Limits |
|-------------|-------|-----------------------|
| Sensitivity | 97.1% | 85.5% - 99.9%         |
| Specificity | 99.1% | 98.5% - 99.5%         |

The 14 discrepant specimens were further characterized by additional testing at TECHLAB. This testing included an FDA-cleared commercial Microassay well EIA, an FDA-cleared commercial molecular test, in-house PCR (detecting the 16s rRNA gene of *Campylobacter* spp., and species-specific identification), and bidirectional sequencing.

\* Nine of the 13 specimens that were culture negative and *CAMPYLOBACTER QUIK CHEK*™ test positive were confirmed to be positive for *C. jejuni* with all test methods.

Two of the 13 specimens that were culture negative and *CAMPYLOBACTER QUIK CHEK*™ test positive were confirmed to be positive with the commercial EIA, in-house PCR, and bidirectional sequencing.

One of the 13 specimens that was culture negative and *CAMPYLOBACTER QUIK CHEK*™ test positive was confirmed to be positive with an FDA-cleared commercial molecular test, in-house PCR and bidirectional sequencing.

One specimen that was culture negative and *CAMPYLOBACTER QUIK CHEK*™ test positive was confirmed to be positive for *C. upsaliensis* (an important pathogen) by species-specific PCR and sequencing.

\*\*The one specimen that was culture positive and *CAMPYLOBACTER QUIK CHEK*™ test negative was confirmed to be negative for *C. jejuni* or *C. coli* with all test methods.

### **Retrospective Study**

Supplemental testing was performed on 30 retrospective positive specimens. The patient ages ranged from less than 11 months to 74 years. All retrospective specimens were *Campylobacter* spp. culture positive and were further characterized as *Campylobacter* spp. positive by an FDA-cleared commercial Microassay well EIA, an FDA-cleared commercial molecular test, in-house PCR (detecting the 16s rRNA gene of *Campylobacter* spp., and species-specific identification), and bidirectional sequencing. These specimens were then tested in the *CAMPYLOBACTER QUIK CHEK*™ test. All 30 specimens tested positive for *Campylobacter* spp. by all methods, yielding 100% correlation with all test methods.

### **Reproducibility**

The reproducibility of the *CAMPYLOBACTER QUIK CHEK*™ test was determined using 8 fecal specimens that were coded to prevent their identification during testing. Testing was performed at 2 independent laboratories and on-site at TECHLAB, Inc. The samples included 2 negative samples, 2 high negative samples, 2 low positive samples and 2 moderate positive samples. The samples were tested in triplicate twice a day over a 5-day period by multiple technicians at each site using 2 different kit lots. A positive and negative control was run with each panel of the masked samples. The results from each laboratory were submitted to TECHLAB, Inc. and compared with in-house results. The results were consistent among the different locations, and exhibited a correlation of 100%. The samples produced the expected results 100% of the time.

### **Analytical Sensitivity**

The analytical sensitivity of the test was determined by using *C. jejuni* and *C. coli* whole organism culture preparations in a sample matrix. The concentration of *C. jejuni* and *C. coli* organisms in fecal matrix at which specimens are positive by the *CAMPYLOBACTER QUIK CHEK*™ test 95% of the time is the assay Limit of Detection (LoD).

The Limit of Detection (LoD) for the *CAMPYLOBACTER QUIK CHEK*™ test with raw fecal samples was established at  $8.39 \times 10^4$  CFU/mL (1271 CFU/test) for *C. jejuni*. For specimens in Protocol™ Cary Blair media, the LoD was established at  $1.78 \times 10^5$  CFU/mL (2781 CFU/test) for *C. jejuni*. For specimens in Protocol™ C&S media, the LoD was established at  $7.25 \times 10^4$  CFU/mL (1133 CFU/test) for *C. jejuni*.

The Limit of Detection (LoD) for the *CAMPYLOBACTER QUIK CHEK*™ test with raw fecal samples was established at  $7.70 \times 10^5$  CFU/mL (11667 CFU/test) for *C. coli*. For specimens in Protocol™ Cary Blair media, the LoD was established at  $2.22 \times 10^6$  CFU/mL (34688 CFU/test) for *C. coli*. For specimens in Protocol™ C&S media, the LoD was established at  $1.56 \times 10^6$  CFU/mL (24375 CFU/test) for *C. coli*.

**Analytical Specificity (Cross Reactivity)**

The CAMPYLOBACTER QUIK CHEK™ test was evaluated for cross-reactivity with common intestinal organisms and viruses listed below. None of the organisms or viruses were shown to interfere with the performance of the CAMPYLOBACTER QUIK CHEK™ test.

|                                             |                                                 |                                      |
|---------------------------------------------|-------------------------------------------------|--------------------------------------|
| <i>Acinetobacter baumannii</i>              | <i>Aeromonas hydrophila</i>                     | <i>Bacillus cereus</i>               |
| <i>Bacillus subtilis</i>                    | <i>Bacteroides fragilis</i>                     | <i>Campylobacter concisus</i>        |
| <i>Campylobacter fetus</i>                  | <i>Candida albicans</i>                         | <i>Citrobacter freundii</i>          |
| <i>Clostridium bifermentans</i>             | <i>Clostridium difficile</i>                    | <i>Clostridium perfringens</i>       |
| <i>Edwardsiella tarda</i>                   | <i>Enterobacter cloacae</i>                     | <i>Enterococcus faecalis</i>         |
| <i>Escherichia coli</i>                     | <i>Escherichia coli</i> EIEC                    | <i>Escherichia coli</i> EPEC         |
| <i>Escherichia coli</i> ETEC                | <i>Escherichia coli</i> O157:H7 (non-toxigenic) |                                      |
| <i>Escherichia coli</i> O157:H7 (toxigenic) | <i>Escherichia fergusonii</i>                   | <i>Escherichia hermanii</i>          |
| <i>Helicobacter pylori</i>                  | <i>Klebsiella pneumoniae</i>                    | <i>Lactobacillus acidophilus</i>     |
| <i>Lactococcus lactis</i>                   | <i>Listeria monocytogenes</i>                   | <i>Peptostreptococcus anaerobius</i> |
| <i>Plesiomonas shigelloides</i>             | <i>Porphyromonas asaccharolytica</i>            | <i>Prevotella melaninogenica</i>     |
| <i>Proteus vulgaris</i>                     | <i>Pseudomonas aeruginosa</i>                   | <i>Pseudomonas fluorescens</i>       |
| <i>Salmonella enterica typhimurium</i>      | <i>Serratia marcescens</i>                      | <i>Shigella dysenteriae</i>          |
| <i>Shigella flexneri</i>                    | <i>Shigella sonnei</i>                          | <i>Staphylococcus aureus</i>         |
| <i>Staphylococcus aureus</i> (Cowan's)      | <i>Streptococcus agalactiae</i>                 | <i>Staphylococcus epidermidis</i>    |
| <i>Vibrio parahaemolyticus</i>              | <i>Yersinia enterocolitica</i>                  |                                      |
| Adenovirus Type 1, 2, 3, 5, 40, 41          | Human Coronavirus                               | Coxsackievirus B2, B3, B4, B5        |
| Echovirus 9, 11, 18, 22, 33                 | Enterovirus 68, 69, 70, 71                      | Norovirus                            |
| Human Rotavirus                             |                                                 |                                      |

*Campylobacter* species that were shown to be reactive with the CAMPYLOBACTER QUIK CHEK™ test.

*C. helveticus* (strain 54661) was found to be positive at  $3.08 \times 10^6$  CFU/mL (4 x LoD of *C. coli*)

*C. lari* (strain 23947) was found to be positive at  $3.08 \times 10^6$  CFU/mL (4 x LoD of *C. coli*)

*C. upsaliensis* (strain 14913) was found to be positive at  $1.54 \times 10^6$  CFU/mL (2 x LoD of *C. coli*)

**Inclusivity Study**

The specificity of the CAMPYLOBACTER QUIK CHEK™ test was evaluated using several strains of *Campylobacter jejuni* and *Campylobacter coli*. All strains listed generated positive results when tested.

*C. coli* strains: 11283, 10956, 17755, 36994, 53138

*C. jejuni* sub-species *jejuni* strains: 11284, 6951, 12081, 29411, 38106

*C. jejuni* sub-species *doylei* strain: 24567

**Interfering Substances (U.S. Formulation)**

The following substances had no effect on positive or negative CAMPYLOBACTER QUIK CHEK™ test results analyzed at the concentrations indicated:

Barium sulfate (5% w/v), Benzalkonium Chloride (1% w/v), Ciprofloxacin (0.25% w/v), Ethanol (1% w/v), Hog gastric mucin (3.5% w/v), Human blood (40% v/v), Hydrocortisone (1% w/v), Imodium® (5% v/v), Kaopectate® (5% v/v), Leukocytes (0.05% w/v), Maalox® Advanced (5% v/v), Mesalazine (10% w/v), Metronidazole (0.25% w/v), Mineral Oil (10% w/v), Mylanta® (4.2 mg/mL), Naproxen Sodium (5% w/v), Nonoxynol-9 (40% w/v), Nystatin (1% w/v), Palmitic Acid/Fecal Fat (40% w/v), Pepto-Bismol® (5% v/v), Phenylephrine (1% w/v), Polyethylene glycol 3350 (10% w/v), Prilosec OTC® (5 µg/mL), Sennosides (1% w/v), Simethicone (10% w/v), Steric Acid/Fecal Fat (40% w/v), Tagamet® (5 µg/mL), TUMS® (50 µg/mL), Human Urine (5% v/v), and Vancomycin (0.25% w/v).

### **Prozone**

To ensure that a high concentration of *Campylobacter* antigen does not interfere with a positive reaction in the CAMPYLOBACTER QUIK CHEK™ test, high positive samples were prepared by spiking a negative fecal pool at a concentration possibly observed in clinical specimens. A total of 5 different dilutions of *C. jejuni* and *C. coli* whole organism culture preparation, up to and including the clinically observed high concentration, were prepared and tested in triplicate. The results demonstrated that there was no overall prozone effect, that elevated levels of antigen did not affect the detection of the antigen.

### **CONCLUSION**

The conclusions drawn from the nonclinical and clinical tests demonstrate that the CAMPYLOBACTER QUIK CHEK™ test is as safe and as effective and performs as well or better than standard culture, and is equivalent to the predicate device in performance. The information submitted in this premarket notification is complete and supports a substantial equivalence decision.