



September 13, 2018

Q'Apel Medical LLC
% Michele Lucey
President
Lakeshore Medical Device Consulting LLC
128 Blye Hill Landing
Newbury, New Hampshire 03255

Re: K181000

Trade/Device Name: SelectFlex 072 Neurovascular Access System
Regulation Number: 21 CFR 870.1250
Regulation Name: Percutaneous Catheter
Regulatory Class: Class II
Product Code: DQY
Dated: August 10, 2018
Received: August 15, 2018

Dear Michele Lucey:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/CombinationProducts/GuidanceRegulatoryInformation/ucm597488.htm>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,


Xiaolin Zheng -S

for Carlos L. Peña, PhD, MS
Director
Division of Neurological
and Physical Medicine Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K181000

Device Name

SelectFlex 072 Neurovascular Access System

Indications for Use (Describe)

The SelectFlex 072 Neurovascular Access System is indicated for the introduction of interventional devices into the peripheral and neurovasculature.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

K181000
510(k) SUMMARY
As required by 21 CFR 807.92

Submitter Information:

Submitter's Name: Q'Apel Medical LLC
Address: 330 Wilshire Boulevard
Santa Monica CA 90401
Telephone: 310-395-3950
Fax: 310-395-3950
Contact Person: Michele Lucey
Telephone: 603-748-1374

Date Prepared: September 12, 2018

Device Trade Name: SelectFlex 072 Neurovascular Access System

Classification: Class II

Product Code(s): DQY

Regulation Number(s): 870.1250

Predicate Devices: Penumbra Benchmark Intracranial Access System, (K142321)

Reference Devices: Depuy Synthes Envoy Guiding Catheter STR (K965247, K982770)
Depuy Synthes Envoy XB (Extra Backup Support) Guiding Catheter (K140307)
Cook Flexor® = Shuttle Guiding Sheath (K142819)

Indication for Use:

The SelectFlex 072 Neurovascular Access System is indicated for the introduction of interventional devices into the peripheral and neurovasculature.

Device Description

The SelectFlex 072 Neurovascular Access System is a sterile, single-use intravascular catheter used to facilitate access target vasculature during interventional procedures. The system is composed of a the SelectFlex 072 Catheter, a 3cc Inflation Syringe, and a 7Fr Peel Away Introducer. The 072 SelectFlex Catheter has a usable length of 105cm. The SelectFlex 072 Catheter has variable stiffness along its length and has a dual mode stiffness mechanism on the distal portion of the catheter that is activated by the user allowing the device to transition between tracking and support modes. The distal end of the SelectFlex 072 Catheter has a hydrophilic coating.

Performance Testing Summary:

The SelectFlex 072 Neurovascular access system has successfully completed the following relevant performance testing demonstrating that the device is suitable for its intended use.

Test Description	Test Method	Results
Packaging Integrity (sterile barrier)	Tested per ISO 11607-1 and -2	PASS All samples met the pre-determined acceptance criteria
Visual Surface Requirements	Visual inspection of catheter surfaces	PASS All samples met the pre-determined acceptance criteria
Dimensional Verification	Device dimensions were measured to confirm conformance to the product specification	PASS All samples met the pre-determined acceptance criteria
Liquid Leakage Under Pressure	Tested per ISO 10555-1:2013 Annex C	PASS All samples met the pre-determined acceptance criteria
Hub Aspiration Air Leakage	Tested per ISO 10555-1:2013 for Hub Aspiration Air Leakage	PASS All samples met the pre-determined acceptance criteria
Simulated Use/Usability	Device preparation, delivery, access, was evaluated in a challenging neurovascular model. Simulated use testing includes a usability assessment with multiple physicians	PASS All samples met the pre-determined acceptance criteria
Flex Fatigue	Tested per ISO 10555-1: 2013 for Flexural Fatigue	PASS All samples met the pre-determined acceptance criteria
Tip Deflection	The amount of tip deflection under load was evaluated and compared to the predicate device	PASS All samples met the pre-determined acceptance criteria
Inflation Fatigue	Tested per ISO 10555-1: 2013 for Inflation Fatigue – 20 inflation cycles	PASS All samples met the pre-determined acceptance criteria
Burst Volume	Tested per ISO 10555-1: 2013 for Inflation Fatigue – tested to 2x the inflation volume	PASS All samples met the pre-determined acceptance criteria
Torque Test	Tested per ISO 10555-1: 2013 for Torque Testing	PASS All samples met the pre-determined acceptance criteria
Flow Rate	Tested per ISO 10555-1: 2013 for Flow rate compared to the predicate device	PASS All samples met the pre-determined acceptance criteria
ISO 594-2: Conical Luer Fittings	Tested per ISO 594-2 for Conical Fittings	PASS All samples met the pre-determined acceptance criteria

Test Description	Test Method	Results
Corrosion Resistance	Tested per ISO 10555-1 Annex A for corrosion resistance	PASS All samples met the pre-determined acceptance criteria
Particulate Count	Effluent tested per AAMI TIR42, USP 788 using multiple insertion and withdrawal cycles	PASS All samples met the pre-determined acceptance criteria
Coating Integrity, Lubricity, Durability	Tested in consideration of FDA CTQ for Hydrophilic Coated vascular catheters	PASS All samples met the pre-determined acceptance criteria
Peak Tensile Testing	Tested per ISO 10555-1 for tensile strength including all bonds/joints	PASS All samples met the pre-determined acceptance criteria
Device Removal in Support and Tracking Modes	Removal force in both modes was compared to a reference device	PASS All samples met the pre-determined acceptance criteria
Radiopacity	Radiopacity was evaluated during simulated use testing confirming visualization under fluoroscopy	PASS All samples met the pre-determined acceptance criteria
Bacterial Endotoxin	Testing was performed in accordance with USP 40, NF 35, 2017. <85> Bacterial Endotoxins Test	PASS All samples met the pre-determined acceptance criteria
Shelf Life	Accelerated aging studies were completed to confirm the stability of the product and package	PASS All samples met the pre-determined acceptance criteria

No animal or clinical studies were required to demonstrate substantial equivalence.

Biocompatibility Testing Summary

Categorized as Externally Communicating Device, Circulating Blood, Limited Contact (≤ 24 hours), per ISO 10993-1, the following testing was conducted:

Test Name	Test Method	Results
Cytotoxicity	Tested in accordance with ISO 10993-5, Biological Evaluation of Medical Devices – Part 5: Tests for <i>in vitro</i> toxicity, Neutral Red Uptake Method	Pass Noncytotoxic according to the predetermined acceptance criteria
Intracutaneous Irritation	Tested in accordance with ISO 10993-10, Biological Evaluation of Medical Devices – Part 10: Tests for Irritation and Skin Sensitization	Pass Test requirements for intracutaneous reactivity were met according to the predetermined acceptance criteria
Sensitization	Tested in accordance with ISO 10993-10, Biological Evaluation of Medical Devices – Part 10 Tests for Irritation and Skin Sensitization, Kligman Maximization Test	Pass did not elicit a sensitization response according to the predetermined acceptance criteria
Systemic Toxicity	Tested in accordance with ISO 10993-11, Biological Evaluation of Medical Devices – Part 11: Tests for Systemic Toxicity	Pass Test requirements for systemic toxicity were met according to the predetermined acceptance criteria

Material Mediated Pyrogenicity	Tested in accordance with ISO 10993-11, Biological Evaluation of Medical Devices – Part 11: Tests for Systemic Toxicity and USP 40 <151> Pyrogen Test	Pass Nonpyrogenic, met the predetermined acceptance criteria
Hemolysis	Tested in accordance with ASTM F756-17, Standard Practice for Assessment of Hemolytic Properties of Materials and ISO 10993-4, Biological Evaluation of Medical Devices – Part 4: Selection of Tests for Interactions with Blood, Tests for Hemolytic Properties, Direct and Indirect Methods	Pass Non-hemolytic, met the predetermined acceptance criteria
<i>In Vitro</i> Hemocompatibility	Tested in accordance with ISO 10993-4, Biological Evaluation of Medical Devices – Part 4: Selection of Tests for Interactions with Blood, Hemocompatibility, Direct Contact Method	Pass Not expected to result in adverse effects <i>in vivo</i> , met the predetermined acceptance criteria
Complement Activation	Tested in accordance with ISO 10993-4, Biological Evaluation of Medical Devices – Part 4: Selection of Tests for Interactions with Blood, SC5b-9 Complement Activation	Pass Does not activate the complement system, met the predetermined acceptance criteria
Un-activated Partial Thromboplastin Time	Tested in accordance with ISO 10994-4, Biological Evaluation of Medical Devices – Part 4: Selection of Tests for Interactions with Blood and ASTM F2382-04, Standard Test Method for Assessment of Intravascular Medical Device Materials on Partial Thromboplastin Time (UPTT)	Pass Does not have an effect on coagulation of human plasma, met the predetermined acceptance criteria
Thrombogenicity	Tested in accordance with ISO 10994-4, Biological Evaluation of Medical Devices – Part 4: Selection of Tests for Interactions with Blood	Pass Demonstrates similar thromboresistance characteristics as the control device, met the predetermined acceptance criteria

Predicate Device Comparison

The following table provides a comparison of the key characteristics of the SelectFlex 072 Neurovascular Access Catheter to the predicate device.

Substantial Equivalence Comparison Chart			
Feature	Subject Device	Predicate device	Comparison
	SelectFlex 072 Neurovascular Access System	Penumbra Benchmark™ Intracranial Access System	
Material	Commonly used medical grade plastics, stainless steel, nitinol	Commonly used medical grade plastics and stainless steel	Similar, difference does not raise new questions regarding safety and efficacy
Outer Diameter	.095-in	0.081-0.083-in	Similar, Reference devices include a diameter range above and below the subject device.
Inner Dimension	.072-in	0.070-in	Similar, Reference devices include a diameter range above and below the subject device.
Fr Designation	7 Fr	6 Fr	Similar, Reference devices includes a diameter range above and below the subject device.
Usable Length	105cm	95cm, 105cm	Same
Tip Shape	Straight	Straight & Multipurpose	Same
Injection Port	Yes	Yes	Same

Substantial Equivalence Comparison Chart			
Feature	Subject Device	Predicate device	Comparison
	SelectFlex 072 Neurovascular Access System	Penumbra Benchmark™ Intracranial Access System	
Radiopaque	Distal Tip has a radiopaque marker band, stainless steel reinforcement and nitinol scaffold in the catheter shaft renders the shaft visible on fluoroscopy	Distal Tip has a radiopaque marker band, stainless steel reinforcement in the catheter shaft renders the shaft visible on fluoroscopy	Same
Coating	Hydrophilic Coating – Distal Portion	Hydrophilic Coating – Distal Portion	Same
Reinforced Shaft	Stainless steel reinforced shaft	Stainless steel reinforced shaft	Same
Variable Stiffness Mechanism	Variable durometer catheter shaft construction to deliver a flexible distal tip and transition to a stiffer proximal section; ten transition segments between the distal tip and proximal shaft. & User selected tracking and support modes	Variable durometer catheter shaft construction to deliver a flexible distal tip and transition to a stiffer proximal section; ten transition segments between the distal tip and proximal shaft.	Similar, the range in catheter stiffness for the subject device is encompassed in the combined range reported for the predicate devices and reference devices with the same intended use. Difference in mechanism does not raise new questions regarding safety and effectiveness.
Accessories Supplied	7 Fr Introducer Sheath 3cc syringe	5F Penumbra Select Catheter – offering includes option for Benchmark and Select Introducer Device as a kit.	Different, there are no special accessories required for use with the predicate device. Difference does not raise new questions regarding safety and effectiveness.
Guidewire Compatibility	0.035 – 0.038-in	0.035 - 0.038-in	Same
How Supplied	Sterile, single use	Sterile, single use	Same
Sterilization Method	ETO	ETO	Same
Sterility Assurance Level	10 ⁻⁶	10 ⁻⁶	Same

Substantial Equivalence:

This summary demonstrates that the SelectFlex 072 Neurovascular Access System and the predicate device have the same intended use, similar technological characteristics, materials, and principals of operation. It can therefore be concluded that the SelectFlex 072 Neurovascular Access System is substantially equivalent to the predicate device.