



October 5, 2018

Tornier Inc.
Renee Stoffel
Senior Regulatory Affairs Specialist
10801 Nesbitt Avenue South
Bloomington, Minnesota 55437

Re: K181587

Trade/Device Name: ORTHOLOC™ SPS Shoulder Plating System
Regulation Number: 21 CFR 888.3030
Regulation Name: Single/Multiple Component Metallic Bone Fixation Appliances And Accessories
Regulatory Class: Class II
Product Code: HRS, HWC
Dated: August 30, 2018
Received: August 31, 2018

Dear Ms. Stoffel:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for

devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/CombinationProducts/GuidanceRegulatoryInformation/ucm597488.htm>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Laurence D. Coyne -S

For Mark N. Melkerson
Director
Division of Orthopedic Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K181587

Device Name

ORTHOLOC™ SPS Shoulder Plating System

Indications for Use (Describe)

The ORTHOLOC™ SPS is intended for fractures and fracture dislocations, osteotomies, and non-unions of the proximal humerus, particularly in osteopenic bone.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



Date Prepared: October 4, 2018

Administrative Information

Name: Tornier, Inc.
 Address: 10801 Nesbitt Avenue South
 Bloomington, MN 55437
 United States of America

 Contact Person: Renee Stoffel
 Title: Sr. Regulatory Affairs Specialist
 Phone: 952-683-7471
 Fax: 952-426-7601

Device Information

Name of Device: ORTHOLOC™ SPS Shoulder Plating System
 Common Name (s): Humeral Plating System
 Classification Name: Single/multiple component metallic bone fixation appliances and accessories
 Regulatory Class: II
 Regulation Number: 888.3030
 Product Code: HRS, HWC

Predicate Device Information

Primary Predicate: K011815, SYNTHES LCP PROXIMAL HUMERUS PLATES

 Additional Predicates: K000684, SMALL FRAGMENT DYNAMIC COMPRESSION LOCKING (DCL) SYSTEM
 K011335, SYNTHES ONE-THIRD TUBULAR DCL PLATE
 K041860, SYNTHES (USA) LCP PROXIMAL HUMERUS PLATES, LONG
 K082625, SYNTHES (USA) 3.5MM LCP PERIARTICULAR PROXIMAL HUMERUS PLATES
 K100146 and K102352, EVOLVE EPS ORTHOLOC
 K111432, EVOLVE (R) TRIAD (TM) PLATING SYSTEM AND EVOLVE (R) TRIAD (TM) BONE SCREWS

Device Description

ORTHOLOC™ SPS is a plating system designed to provide anatomical reduction and stable primary fixation of the proximal humerus using implanted plates and screws. The plates are available in three designs: Standard, Posterior, and Greater Tuberosity (GT) and are manufactured from stainless steel per ASTM F138. The system also includes locking cortical and cancellous screws manufactured from stainless steel per ASTM F2229 and non-locking cortical screws manufactured from stainless steel per ASTM F138. The system is provided non-sterile and requires steam sterilization prior to use.

Tornier, Inc. ORTHOLOC™ SPS

**Indications for Use**

The ORTHOLOC™ SPS is intended for fractures and fracture dislocations, osteotomies, and non-unions of the proximal humerus, particularly in osteopenic bone.

Comparison of Technological Characteristics with the Predicate Devices

The ORTHOLOC™ SPS is substantially equivalent in intended use, materials, design, and sterilization method to the predicate devices. The design differences do not raise new issues of safety or effectiveness and are supported by performance testing.

Non-clinical Performance Testing

Non-clinical bench testing (mechanical testing) was performed to demonstrate substantial equivalence to the predicate devices.

- Plate Bending and Fatigue to ASTM F382-17
- Screw Torsion, Driving Torque, and Axial Pullout to ASTM F543-17
- Engineering Analysis of the Greater Tuberosity Plate
- Cross Sectional Analysis comparing the subject and predicate plates

Clinical Testing

No clinical studies were performed.

Conclusions

The ORTHOLOC™ SPS does not raise new questions of safety or effectiveness. Differences in technological characteristics have been addressed with performance testing. The results of performance testing for the ORTHOLOC™ SPS system support substantial equivalence to the predicate devices.