



March 18, 2019

Beckman Coulter
Radha Goolabsingh
Director, Regulatory Affairs
11800 SW 147th Ave
Miami, Florida 33196-2500

Re: K181599

Trade/Device Name: Unicel DxH 800 Cellular Analysis System with Early Sepsis Indicator
Application

Regulation Number: 21 CFR 866.3215

Regulation Name: Device to detect and measure non-microbial analyte(s) in human clinical specimens
to aid in assessment of patients with suspected sepsis.

Regulatory Class: Class II

Product Code: QFS, GKZ

Dated: June 15, 2018

Received: June 18, 2018

Dear Radha Goolabsingh:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/CombinationProducts/GuidanceRegulatoryInformation/ucm597488.htm>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Steven R. Gitterman -S for

Uwe Scherf, Ph.D.
Director
Division of Microbiology Devices
Office of In Vitro Diagnostics
and Radiological Health
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K181599

Device Name

UniCel DxH 800 Coulter Cellular Analysis System with Early Sepsis Indicator Application

Indications for Use (Describe)

The Unicel DxH 800 Coulter Cellular Analysis System with Early Sepsis Indicator Application is the quantitative measurement of Monocyte Distribution Width (MDW). The Early Sepsis Indicator is intended for use with adult patients presenting to the emergency department, on whom a white cell differential test has been ordered.

MDW is measured from a (K2EDTA) whole-blood venous sample within 2 hours of collection. MDW values greater than 20.0 together with other laboratory findings and clinical information, aids in identifying patients with sepsis or at increased risk of developing sepsis within the first 12 hours of hospital admission.

MDW values greater than 20.0 should be interpreted in association with other clinical information and diagnostic testing, as a proportion of patients without sepsis may have an elevated MDW value at baseline.

MDW values less than or equal to 20.0 cannot rule out sepsis or the development of sepsis within 12 hours of hospital admission. The Early Sepsis Indicator should not be used as the sole basis to determine the absence of sepsis.

The predictive value of the Early Sepsis Indicator for identifying sepsis in patients with hematological abnormalities has not been established.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Section 5.0: – 510(k) Summary

510(k) Summary for Unicel DxH 800 Coulter Cellular Analysis System with Early Sepsis Indicator Application**510(k) Owner / Submitter Information**

Name: Beckman Coulter Inc.
Address: 11800 SW 147th Ave., Miami, FL 33196
Cellular: 786-423-4575
Desk Phone #: (305) 380-2584
Fax: (786) 639-2584
Contact Person: Radha Goolabsingh
Email Address: rgoolabsingh@beckman.com

Date Submitted: June 15, 2018

Device Name and Classification

Trade Name: Unicel DxH 800 Coulter Cellular Analysis System with Early Sepsis Indicator Application

Common Name: Early Sepsis Indicator

Classification: Class II (Special Controls)

Classification Name: Device to detect and measure non-microbial analyte(s) in human clinical specimens to aid in assessment of patients with suspected sepsis (21 CFR 866.3215)

Product Code: QFS

Panel: 83- Division of Microbiology Devices

Unicel DxH 800 Cellular Analysis System - Last 510(k) Clearance: K140911, 09/05/2014

Predicate Device Information

Predicate Product	510(k) Number	Date Cleared	Classification	21 CFR	Product Code
VIDAS BRAHMS PCT (PCT)	K162827	02/23/2017	Class II (Special Controls)	866.3215	PMT

Device Description**Systems Overview**

The Early Sepsis Indicator (ESI) requires the use of the UniCel DxH 800 Coulter Cellular Analysis System (DxH 800) and its reagents, controls and calibrators last cleared under 510(k) K140911.

The UniCel DxH 800 Coulter Cellular Analysis System contains a quantitative, automated hematology analyzer (DxH 800) designed for *in vitro* diagnostic use in screening patient populations by clinical laboratories. The system provides a Complete Blood Count (CBC), Leukocyte 5 Part Differential (Diff), Reticulocyte (Retic), Nucleated Red Blood Cell (NRBC) on whole blood, as well as, Total Nucleated Count (TNC), and Red Cell Count (RBC) on Body Fluids (cerebrospinal, serous and synovial). This submission adds a new parameter, Monocyte Distribution Width (MDW) to those mentioned above. This parameter has been shown to aid in the early detection of Sepsis in emergency room patients.

The system consists of two primary components, the workstation and the DxH 800 analyzer as shown in Figure 1. DxH 800 System Configuration. The primary function of the DxH 800 analyzer is to process samples and provide results to the workstation. The primary functions of the workstation are: user interface, system control, results processing and storage and external communications. The analyzer runs embedded code and the workstation runs Microsoft Windows 7 Operating System (OS).

Optionally, the workstation can be connected to:

- A printer for creating reports;
- A Laboratory Information System (LIS) for receiving test orders and releasing results; and
- Pro-Service Remote Management System (RMS) which provides secure access to the DxH800 workstation for BCI Service Personnel to perform troubleshooting, system monitoring and for assisting customers.

The version of the system proposed in this submission is DxH800 v3.8.0 This system is an update to DxH800 v3.0.0 (see K140911) to add a new feature.

Neither the system architecture nor the top-level software architecture have been modified from DxH800 v3.0 to add the Early Sepsis Indicator feature.

The DxH 800 v.3.8.0 user interface and labeling will be offered in English in the U.S. and in French, German, Italian, Spanish, Japanese, Chinese and English outside the U.S.




Figure 1. DxH 800 System Configuration

Additionally, the UniCel DxH 800 Coulter Cellular Analysis System with Early Sepsis Indicator Application uses the currently cleared controls and calibrator provided in Table 1 below.

Table 1: Summary of Device Components, Product Description, Proposed Classification and Proposed Predicate:

Product	510(k) Clearance Status	Description	Classification	Proposed Predicate
DxH 800 Hematology Analyzer	Existing Product (cleared under K140911) Modified- Addition of new parameter (MDW) requiring Software algorithm modifications	The UniCel DxH 800 COULTER Cellular Analysis System (DxH 800) is intended for In Vitro Diagnostic Use in clinical laboratories. The DxH 800 System provides automated complete blood count, leukocyte differential, nucleated red blood cell (NRBC) enumeration and reticulocyte analysis as well as an automated method for enumeration of the Total Nucleated Cells (TNC) and Red Blood Cells (RBC) in body fluids.	Class II, GKZ, 864.5220, Hematology, Automated Differential Cell Counter Early Sepsis Indicator Application: 866.3215, Product code (QFS)	VIDAS B.R.A.H.M.S. PCT (PCT) (K16287) (Paper Predicate)
COULTER DxH Diff Pack	Existing Product (Included in K140911) No change	The COULTER DxH Diff Pack is an erythrolytic reagent and a leukocyte preservative used to perform a five-part differential analysis using VCSn technology on the DxH 800/DxH 600.	Class I Exempt, GGK, 864.8540, Hematology, Red cell lysing reagent	N/A
COULTER DxH Diluent	Existing Product (Included in K140911) No change	COULTER DxH Diluent provides the ability to analyze portions of the diluted blood samples for different blood cell types and measure hemoglobin on the DxH 800/600. Acts as a rinsing agent on all DxH Systems.	Class I Exempt, GIF, 864.8200, Hematology, Blood cell diluent	N/A
COULTER DxH Cleaner	Existing Product (Included in K140911) No change	DxH Cleaner is a cleaning agent for use on the DxH System components that come in contact with blood samples.	Class I Exempt, JCB, 864.4010, Hematology, General Purpose Reagent	N/A
COULTER 6C Plus Cell Control	Existing Product (6C Cell Control) (cleared under K081822)	COULTER 6C Plus Cell Control is an integrated control that enables monitoring of system performance for CBC, Diff, NRBC, MDW and LDW parameters. The COULTER 6C Plus Cell Control	Class II, Exempt, JPK, 864.8625, Hematology, Quality Control Mixture	COULTER 6C Cell Control (K081822)

Product	510(k) Clearance Status	Description	Classification	Proposed Predicate
	Modified - to provide two new parameters (LDW and MDW) for control monitoring of the calculation of cellular dispersion	proposed for use with the Sepsis Application is the current COULTER 6C Cell Control (K081822), assayed to include values for the MDW and LDW parameters. LDW is not intended for the reporting of patient results.		
COULTER Latron CP-X Control	<u>Existing Product</u> (cleared under K885028) Modified - with tighter expected ranges for control monitoring of the differential volume measurement	COULTER LATRON CP-X is a suspension of stable polystyrene particles of uniform size with a diameter CV of $\leq 3.0\%$. Latron CP-X is used by customers as part of the daily quality control procedure to monitor the stability of the electrical processing and the fluidic flow rate systems used to measure the volume (size), conductivity and light scattering characteristics of cells as they pass through the flow cell. It is also used by Beckman Coulter manufacturing and service engineers to calibrate the differential volume measurements	Class II, GKL, 864.5220, Hematology, Automated Differential Cell Counter	COULTER Latron CP-X Control (K885028)
*Coulter S-CAL Calibrator	<u>Existing Product</u> (cleared under K840794) No change	The current COULTER S-CAL Calibrator used with the DxH 800 is intended for customers to determine the calibration factors for directly measured CBC parameters; it is not required differential parameters	Class II, Exempt, KRY, 864.8175, Calibrator for Platelet Counting	N/A

*Note: The current COULTER S-CAL Calibrator used with the DxH 800 is intended for customers to determine the calibration factors for directly measured CBC parameters; it is not required for differential parameters. Since the MDW and LDW parameters are derived from the monocyte and lymphocyte differential, the use of a calibrator by the customer is not required for the Sepsis Application.

Intended Use:
UniCel DxH 800 Coulter Cellular Analysis System with Early Sepsis Indicator Application:

The Unicel DxH 800 Coulter Cellular Analysis System with Early Sepsis Indicator Application is the quantitative measurement of Monocyte Distribution Width (MDW). The Early Sepsis Indicator is intended for use with adult patients presenting to the emergency department, on whom a white cell differential test has been ordered.

MDW is measured from a (K2EDTA) whole-blood venous sample within 2 hours of collection. MDW values greater than 20.0 together with other laboratory findings and clinical information, aids in identifying patients with sepsis or at increased risk of developing sepsis within the first 12 hours of hospital admission.

MDW values greater than 20.0 should be interpreted in association with other clinical information and diagnostic testing, as a proportion of patients without sepsis may have an elevated MDW value at baseline.

MDW values less than or equal to 20.0 cannot rule out sepsis or the development of sepsis within 12 hours of hospital admission. The Early Sepsis Indicator should not be used as the sole basis to determine the absence of sepsis.

The predictive value of the Early Sepsis Indicator for identifying sepsis in patients with hematological abnormalities has not been established.

Indication for Use:

See Intended Use above.

Overall Comparison between Early Sepsis Indicator and the Predicate:

Table 2 provides an overall comparison of the differences and similarities against the Predicate Device (VIDAS BRAHMS PCT (PCT) K162827).

Table 2: Overall Comparison between ESI and the Predicate

Item	DxH 800 Coulter Cellular Analysis System with Early Sepsis Indicator Application	Predicate Device (Paper Predicate) VIDAS BRAHMS PCT (PCT) K162827
Intended Use and Indications for Use	The Unicel DxH 800 Coulter Cellular Analysis System with Early Sepsis Indicator Application is the quantitative measurement of Monocyte Distribution Width	VIDAS B·R·A·H·M·S PCT (PCT) is an automated test for use on the instruments of the VIDAS family for the determination of human procalcitonin in human serum or plasma (lithium heparinate) using the ELFA (Enzyme-Linked Fluorescent Assay) technique.

Item	DxH 800 Coulter Cellular Analysis System with Early Sepsis Indicator Application	Predicate Device (Paper Predicate) VIDAS BRAHMS PCT (PCT) K162827
	(MDW). The Early Sepsis Indicator is intended for use with adult patients presenting to the emergency department, on whom a white cell differential test has been ordered. MDW is measured from a (K2EDTA) whole-blood venous sample within 2 hours of collection. MDW values greater than 20.0 together with other laboratory findings and clinical information, aids in identifying patients with sepsis or at increased risk of developing sepsis within the first 12 hours of hospital admission. MDW values greater than 20.0 should be interpreted in association with other clinical information and diagnostic testing, as a proportion of patients without sepsis may have an elevated MDW value at baseline. MDW values less than or equal to 20.0 cannot rule out sepsis or the development of sepsis within 12 hours of hospital admission. The Early Sepsis Indicator should not be used as the sole basis to determine the absence of sepsis. The predictive value of the Early Sepsis Indicator for identifying sepsis in patients with hematological	Used in conjunction with other laboratory findings and clinical assessments, VIDAS BRAHMS PCT is intended for use as follows: • to aid in the risk assessment of critically ill patients on their first day of ICU admission for progression to severe sepsis and septic shock, • to aid in assessing the cumulative 28-day risk of all-cause mortality for patients diagnosed with severe sepsis or septic shock in the ICU or when obtained in the emergency department or other medical wards prior to ICU admission, using a change in PCT level over time • to aid in decision making on antibiotic therapy for patients with suspected or confirmed lower respiratory tract infections (LRTI) defined as community-acquired pneumonia (CAP), acute bronchitis, and acute exacerbation of chronic obstructive pulmonary disease (AECOPD) – in an inpatient setting or an emergency department, • to aid in decision making on antibiotic discontinuation for patients with suspected or confirmed sepsis.

Item	DxH 800 Coulter Cellular Analysis System with Early Sepsis Indicator Application	Predicate Device (Paper Predicate) VIDAS BRAHMS PCT (PCT) K162827
	abnormalities has not been established.	
Specimen	Venous Whole Blood	Human Serum or Plasma (lithium heparinate)
Analyte	MDW Parameter	Procalcitonin (PCT)
Automated	Automated Hematology Analyzer	Automated Assay
Assay Technique	Whole blood analysis and cellular differentiation on the DxH800	ELFA (Enzyme-Linked Fluorescent Assay) technique
Assay Principle	Coulter Principle: Volume, Conductivity, Light Scatter Analysis (VCSn) Technology using : Aperture impedance (DC) Conductivity (RF) Laser Light Scatter (Multiple angles) Laser Light Absorbance	Immunoassay based on sandwich principle
Detection method	VCSn	Fluorescence (ELFA) of 4-methy-umbelliferyl measured at 450 nm
Assay/Parameter Duration	≥90 specimens per hour (40 sec/cycle) for CBC and Diff with NRBC cycle (approximately 40 seconds per sample)	Approximately 20 minutes
Combination Devices	N/A	Instruments of the VIDAS family: VIDAS, miniVIDAS or VIDAS 3
Antibodies	N/A	Conjugate: Alkaline phosphatase-labeled mouse monoclonal anti-human procalcitonin Solid phase: Mouse monoclonal anti-procalcitonin immunoglobulins coated on interior of the SPR Sample
Sample Volume	Automatic, cap-piercing 165 µL	200 µL

Table 3: Summary of the MDW Parameter Performance Testing

Study	Testing Approach	FDA Guidance Documents	Standards/ References	Testing Results
Performance - Bench Testing				
Effect of Incomplete Lysing on MDW	To evaluate the effect of incomplete lysing on the MDW parameter results	None	None	The addition of the incomplete lysis flagging condition causes less than 1.2% of the clinical specimens to have the MDW inhibited. At this rate of flagging the clinical utility of the MDW is not diminished.
Specimen Stability	To characterize the specimen age profile for the MDW parameter on the DxH 800 system using venous whole blood stored at room temperature.	None	CLSI EP25-A: Evaluation of Stability of In Vitro Diagnostic Reagents; Approved Guideline – September 2009	MDW trends as sample ages. The trend exhibits a quadratic pattern and tends to plateau after 4 hours. Consistent with published literature on monocyte morphological changes, this study shows MDW values tend to increase over the first 3 hours from the draw time and then plateau from 4 to 8 hours post draw at room temperature

Study	Testing Approach	FDA Guidance Documents	Standards/ References	Testing Results
Comparability – Specimen Tube Diameter	To evaluate the effect of mixing using the smaller diameter 8 mm and 11 mm tubes as compared to the standard 13 mm tube on the MDW parameter results.	None	CLSI H26-A2: Validation, Verification and Quality Assurance of Automated Hematology Analyzers; Approved Standard – Second Edition - June 2010 EP09-A3: Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline - Third Edition – August 2013	The MDW reported from 8 mm and 11 mm diameter tubes is equivalent to the MDW reported from the 13 mm diameter tube.
Comparability – Sampling modes	To assess the impact on MDW parameter results due to sampling via the Closed Vial Cassette sampling mode and the Open Vial Single-Tube sampling mode on the DxH 800.	None	CLSI H26-A2: Validation, Verification and Quality Assurance of Automated Hematology Analyzers; Approved Standard – Second Edition - June 2010 CLSI EP09-A3: Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline - Third Edition - August 2013	Testing of the MDW parameter in the Open Vial Single-Tube sampling mode was found to be equivalent to the Closed Vial Cassette sampling mode.

Interfering Substances	To assess the effect of interfering substances (lipemia, bilirubin and hemoglobin) on the measurement of the MDW parameter on the DxH 800 system.	None	CLSI H26-A2: Validation, Verification and Quality Assurance of Automated Hematology Analyzers - June 2010 CLSI C56-A: Hemolysis, Icterus, and Lipemia/Turbidity Indices as Indicators of Interference in Clinical Laboratory Analysis; Approved Guideline - July 2012 CLSI EP07-A2: Interference testing in Clinical Chemistry; Approved Guideline – Second Edition - November 2005	Lipids (Triglycerides) at a maximum concentration of 1500 mg/dL do not interfere with MDW measurement. Conjugated bilirubin at the maximum concentration of 40 mg/dL does not interfere with the MDW parameter. Unconjugated bilirubin at the maximum concentration of 20 mg/dL does not interfere with the MDW parameter. Free HGB at a maximum concentration of 1.87 g/dL does not interfere with MDW measurement. Hemolysis at a maximum of 2.02 g/dL does not interfere with MDW measurement.
-------------------------------	---	------	--	--

Study	Testing Approach	FDA Guidance Documents	Standards/ References	Testing Results
Performance - Clinical Testing				
System Precision - Repeatability and Reproducibility	To assess the system precision of the MDW parameter using COULTER 6C Plus Cell Controls on the DxH 800 system.	Class II Special Controls Guidance Document: Premarket Notifications for Automated Differential Cell Counters for Immature or Abnormal Blood Cells; Final Guidance for Industry and FDA December 4, 2001 – Section 13: Specimens.	CLSI EP05-A3 Evaluation of precision of Quantitative Measurement Procedures ; Approved Guideline – Third Edition	Repeatability and Reproducibility study met precision specification of $\leq 10\%$ CV and $\leq 15\%$ CV, respectively by site and all sites combined.
System Precision- Repeatability Profile	To demonstrate the system precision performance of MDW on the DxH 800	Class II Special Controls Guidance Document: Premarket Notifications for Automated Differential Cell Counters for Immature or Abnormal Blood Cells; Final Guidance for Industry and FDA December 4, 2001 – Section 13: Specimens.	CLSI H26-A2 Validation, Verification, and Quality Assurance of Automated Hematology Analyzers, Approved Standard – 2 nd Edition; June 2010 CLSI EP05-A3 Evaluation of precision of Quantitative Measurement Procedures ; Approved Guideline – Third Edition	System repeatability performance acceptance criteria were $\leq 10\%$ for all parameters.

Study	Testing Approach	FDA Guidance Documents	Standards/ References	Testing Results
Adult Reference Interval	To assess the normal reference interval of MDW for apparently healthy adult males and females. A total of 146 specimens were collected consisting of 70 females and 76 males to establish the normal reference interval of MDW.	None	CLSI EP28-A3c, Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline, (10/2010)	The adult normal reference interval assessed for MDW on the DxH 800 system across all sites was found to be between the lower limit of 13.98 and the upper limit of 21.28. Data was included in the product labeling, however the upper limit of normal will not be used for MDW decision-making.

Study	Testing Approach	FDA Guidance Documents	Standards/ References	Testing Results
Clinical Accuracy	To assess the diagnostic ability of MDW. A multi-center prospective cohort study was conducted to clinically validate the performance of Monocyte Volume Distribution Width (MDW) parameter on the DxH 800 hematology analyzer, as an aid in the early detection of sepsis in adult patients presenting to the emergency department (ED).	Class II Special Controls Guidance Document: Premarket Notifications for Automated Differential Cell Counters for Immature or Abnormal Blood Cells; Final Guidance for Industry and FDA December 4, 2001 Statistical Guidance on Reporting Results from Studies Evaluating Diagnostic Tests; Final Guidance for Industry and FDA Staff; March 2007	CLSI EP12-A2: User Protocol for Evaluation of Qualitative Test Performance; Approve Guideline. Second Edition guideline. CLSI EP24-A2: Assessment of the Diagnostic Accuracy of Laboratory Tests Using Receiver Operating Characteristic Curves; Approved Guideline – Second Edition. 2011.	The study results support the proposed intended use and indications for use. The study validated the two (2) predefined cut-offs and demonstrated that the lower limit of sensitivity and specificity were within acceptance criteria for both cut-offs. However, the study demonstrated that an MDW cut-off of 20.0 units provided an optimum diagnostic ability by balancing the ability to detect positive patients (sensitivity) and negative patients (specificity) and this is the cut-off selected.

Substantial Equivalence Conclusion:

The conclusions drawn from the testing discussed in **Section 20**, which demonstrate clinical accuracy (sensitivity and specificity) of the MDW parameter and together with multiple literature citations (as provided in this submission), supports the use of MDW as an aid in the early detection of adult patients with or developing sepsis.

Furthermore, it can be concluded that the submitted information in this premarket notification is complete, offers a reasonable assurance of safety and effectiveness and supports a substantial equivalence determination for the DxH 800 Coulter Cellular Analysis System with Early Sepsis Indicator Application.