



October 4, 2018

Pixxgen Corporation
% Daniel Kamm
Principal Engineer
Kamm & Associates
8870 Ravello Ct
Naples, Florida 34114

Re: K182533

Trade/Device Name: PIXX 1717; PIXX 1417; PIXX 1212; Digital Diagnostic X-ray Receptor Panels
Regulation Number: 21 CFR 892.1680
Regulation Name: Stationary x-ray system
Regulatory Class: Class II
Product Code: MQB
Dated: September 12, 2018
Received: September 14, 2018

Dear Daniel Kamm:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal

statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/CombinationProducts/GuidanceRegulatoryInformation/ucm597488.htm>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

A handwritten signature in blue ink, appearing to read "Rob 2. Mills", is written over a large, light blue "FDA" watermark.

Robert A. Ochs, Ph.D.

Director

Division of Radiological Health

Office of In Vitro Diagnostics

and Radiological Health

Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K182533

Device Name

PIXX 1717; PIXX 1417; PIXX 1212; Digital Diagnostic X-Ray Receptor Panels

Indications for Use (Describe)

Indicated for use in general radiographic images of human anatomy. It is intended to replace radiographic film/screen systems in all general-purpose diagnostic procedures, excluding fluoroscopic, angiographic, and mammographic applications

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary
510(k) Number K182533
PIXXGEN Corporation
5F, SMART BAY, 123, Beolmal-ro, Dongan-gu,
Anyang-si, Gyeonggi-do, 14056, Korea
Tel +82.70.4846.8888
Date Prepared: September 28, 2018
Contact: Young Kim, President

1. Identification of the Device:

Proprietary-Trade Names: PIXX 1717; PIXX 1417; PIXX 1212; Digital Diagnostic X-Ray Receptor Panels

Device: Solid State X-Ray Imager (Flat Panel/Digital Imager)

Regulation Name: Stationary x-ray system.

Product Code: MQB

Regulation Number: 892.1680

Device Class: 2

2. Equivalent legally marketed device: PIXX 1717; PIXX 1417; PIXX 1212; Digital Diagnostic X-Ray Receptor Panels.

Manufacturer: PIXXGEN Corporation

510(k) Number: K180976

Device: Solid State X-Ray Imager (Flat Panel/Digital Imager)

Regulation Name: Stationary x-ray system.

Product Code: MQB

Regulation Number: 892.1680

Device Class: 2

3. Indications for Use (intended use) Indicated for use in general radiographic images of human anatomy. It is intended to replace radiographic film/screen systems in all general-purpose diagnostic procedures, excluding fluoroscopic, angiographic, and mammographic applications. (Rx Only)

4. Description of the Device:

The PIXX series are digital radiography systems, featuring an integrated flat panel digital detector (FPD). They are designed to perform digital radiographic examinations as a replacement for conventional film. This integrated platform provides the benefits of PACS with the advantages of digital radiography for a filmless environment and improves cost effectiveness. The major functions and principle of operation of the updated panels are the same as our previous panel retaining the Wi-Fi wireless features. The modification involves adding battery power capability to the devices in our previous submission K180976.

5. Safety and Effectiveness, comparison to predicate device. The results of bench and test laboratory results indicates that the new devices are as safe and effective as the predicate devices.

6. Substantial Equivalence Chart

	K180976, PIXX 1717, PIXX 1417, PIXX 1212	PIXX 1717	PIXX 1417	PIXX 1212
Intended Use	Indicated for use in general radiographic images of human anatomy. It is intended to replace radiographic film/screen systems in all general-purpose diagnostic procedures, excluding fluoroscopic, angiographic, and mammographic applications	UNCHANGED	UNCHANGED	UNCHANGED
Configuration	This submission is for the Digital Panel and Software only, no generator or stand provided.	UNCHANGED	UNCHANGED	UNCHANGED
Pixel Pitch	140um	SAME	SAME	SAME
Limiting Resolution	Over 3 lp/mm	SAME	SAME	SAME
DQE(CSI)	At 2 lp/mm 26.5%	SAME	SAME	SAME
MTF(CSI)	At 2 lp/mm 44%	SAME	SAME	SAME
DQE(GOS)	At 2 lp/mm 21	SAME	SAME	SAME
MTF(GOS)	At 2 lp/mm 35%	SAME	SAME	SAME
DQE and MTF Comments	Adding battery power capability did not change the panel performance.			
A/D Conversion	16 bits	SAME	SAME	SAME
Active Area	17 x 17 inch 14 x 17 inch 12 x 12 inch	17 x 17 inch	14 x 17 inch	12 x 12 inch
Dimensions(m m)/ Weights(Kg)	460(W)x461(L)x15(H)/3.0Kg 385(W) x 460(L) x 15(H)/ 2.8Kg 308.5(W) x 319.5(L) x 15(H)/ 1.9Kg	SAME	SAME	SAME
Pixels	3,072x3,072 2,560 x 3,072 2,048 x 2,048	3,072 x 3,072	2,560 x 3,072	2,048 x 2,048
Software	Outputs a DICOM image	SAME	SAME	SAME
DICOM	Yes	Yes	Yes	Yes
Scintillator	CsI orGOS	UNCHANGED	UNCHANGED	UNCHANGED
Interface	Wired: Gigabit Ethernet (1000BASE-T) Wireless:IEEE802.11ac, backward compatible	SAME	SAME	SAME
Power source	AC Line only	AC Line and/or Rechargeable Lithium Battery; 5 hours/300 images	AC Line and/or Rechargeable Lithium Battery; 5 hours/300 images	AC Line and/or Rechargeable Lithium Battery 8 hours/480 images

	K180976, PIXX 1717, PIXX 1417, PIXX 1212	PIXX 1717	PIXX 1417	PIXX 1212
Standards	Electrical Safety per IEC 60601-1:2012 and EMC per IEC 60601-1-22007+AC:2010 as well as IEEE802.11ac. Meets FCC requirements.	SAME plus IEC 62133 Battery safety	SAME plus IEC 62133 Battery safety	SAME plus IEC 62133 Battery safety

- 7. Summary of Bench Testing Conducted:** IEC Standards were employed for: Electrical Safety and Electromagnetic Compatibility, and Battery Safety Tests. Standards met:
Electrical safety per: IEC/UL 60601-1, Medical Electrical Equipment – Part 1: General Requirements for Safety. (General)
Electromagnetic Compatibility per IEC 60601-1-2, Collateral Standard: Electromagnetic compatibility Requirements and tests.
IEC 62133 Edition 2.0 2012-12 Secondary cells and batteries containing alkaline or other non-acid electrolytes – Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications
Risk Analysis was conducted in accordance with ISO 14971:2012 and EN 62304. The software remains the same as in K180976. MTF and DQE measurements were not needed because the panel characteristics are identical to our predicate cleared earlier this year.
Battery life testing was conducted to confirm that the PIXX 1212 runs for 8 hours/480 images and for the PIXX 1717/1417 for 5 hours/300 images.
Cybersecurity precautionary labeling was added per the FDA guidance: Content of Premarket Submissions for Management of Cybersecurity in Medical Devices, Guidance for Industry and Food and Drug Administration Staff, Document Issued on: October 2, 2014.
- 8. Summary of Clinical Testing:** Clinical images were not required because the panel technology did not change from K180976.
- 9. Conclusion:** After analyzing bench, clinical image, and external laboratory testing to applicable standards, it is the conclusion of Pixxgen that the PIXX 1717; PIXX 1417; PIXX 1212; Digital Diagnostic X-Ray Receptor Panels are as safe and effective as the predicate device, have few technological differences, and has the identical indications for use, thus rendering them substantially equivalent to the predicate device.