



February 26, 2019

Medos International SARL
% Elizabeth Messana
Regulatory Affairs Specialist II
DePuy Synthes Mitek, a Johnson and Johnson Company
325 Paramount Drive
Raynham, Massachusetts 02767

Re: K182941

Trade/Device Name: HEALIX ADVANCE™ SP PEEK Anchor
Regulation Number: 21 CFR 888.3040
Regulation Name: Smooth or threaded metallic bone fixation fastener
Regulatory Class: Class II
Product Code: MBI
Dated: January 16, 2019
Received: January 17, 2019

Dear Ms. Messana:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's

requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/CombinationProducts/GuidanceRegulatoryInformation/ucm597488.htm>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

 Laurence D. Coyne -S

For Mark N. Melkerson
Director
Division of Orthopedic Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

DEPARTMENT OF HEALTH AND HUMAN SERVICES
Food and Drug AdministrationForm Approved: OMB No. 0910-0120
Expiration Date: 06/30/2020
See PRA Statement below.**Indications for Use**

510(k) Number (if known)

K182941

Device Name

HEALIX ADVANCE™ SP PEEK Anchor

Indications for Use (Describe)

The HEALIX ADVANCE SP PEEK Anchor is indicated for use in the following procedures for reattachment of soft tissue to bone:

- Shoulder: Rotator cuff repair, biceps tenodesis, deltoid repair
- Knee: Posterior oblique repair, medial collateral ligament (MCL) reconstruction, lateral collateral ligament reconstruction, iliotibial (IT) band tenodesis
- Foot/Ankle: Achilles tendon repair

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)☐ Over-The-Counter Use (21 CFR 801 Subpart C)**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(K) SUMMARY
HEALIX ADVANCE™ SP PEEK Anchor
Date Prepared: 1/9/2019

Submitter's Name and Address DePuy Synthes Mitek Sports Medicine
a Johnson & Johnson company
 325 Paramount Drive
 Raynham, MA 02767

On behalf of:
 Medos International SARL
 Chemin-Blanc 38, Le Locle Neuchatel
 CH 2400, Switzerland

Contact Person Elizabeth Messana Telephone: 508-828-3150
 Regulatory Affairs Specialist II Email: emessan1@its.jnj.com

DePuy Synthes Mitek Sports Medicine
a Johnson & Johnson company
 325 Paramount Drive
 Raynham, MA 02767

Name of Medical Device Proprietary Name: HEALIX ADVANCE™ SP PEEK Anchor
 Classification Name: Smooth or threaded metallic bone fixation fastener
 Product Code: MBI
 Common Name: Suture Anchor

Substantial Equivalence The HEALIX ADVANCE SP PEEK Anchor is substantially equivalent to:

- K101679, K082810, K071177, K061863, K051219 Arthrex PushLock Suture Anchors

Reference Devices:

- K161638, SPEEDTRAP™ Graft Preparation System
- K171114, HEALIX ADVANCE™ KNOTLESS PEEK Anchor
- K120824, K090530, Stryker ReelX STT™ Suture Anchor System

Device Classification	HEALIX ADVANCE SP PEEK Anchor is classified as: Fastener, Fixation, Nondegradable, Soft Tissue, classified as Class II, product code MBI, regulated under 21 CFR 888.3040.
Device Panel	Orthopedic Devices
Device Description	The HEALIX ADVANCE SP PEEK Anchor is a two-piece “self-punching” design consisting of a cannulated, threaded anchor and dilator which are permanently implanted into the body. The anchor-dilator are preloaded on a disposable inserter shaft with handle, held in place by a non-absorbable stay suture comprised of Ultra-High Molecular Weight Polyethylene (UHMWPE) and Green poly (ethylene terephthalate) (PET). During implantation the stay suture may be discarded or incorporated into the repair at the surgeon’s discretion. The proposed HEALIX ADVANCE SP PEEK Anchor is offered in two sizes, 4.9 mm and 5.5 mm and is molded from Polyetheretherketone (PEEK) material. The HEALIX ADVANCE SP PEEK Anchor is provided sterile via Ethylene Oxide (EO) sterilization and is for single use only.
Technological Characteristics	The suture anchor design, anchor materials, principal of operation and indications for use are similar when compared to the predicate Arthrex PushLock Suture Anchors (K101679, K082810, K071177, K061863, K051219). The suture anchor design, anchor materials, principal of operation and indications for use are similar when compared to the reference device-ReelX STT™ Suture Anchor by Stryker Corp. (K120824, K090530). The stay-suture design and material are identical to the winding suture present in the reference device- SPEEDTRAP™ Graft Preparation System (K161638). The PEEK material is identical to that present in the reference device-HEALIX ADVANCE™ KNOTLESS PEEK Anchor (K171114).

Indications for Use	The HEALIX ADVANCE SP PEEK Anchor is indicated for use in the following procedures for reattachment of soft tissue to bone: • Shoulder: Rotator cuff repair, biceps tenodesis, deltoid repair • Knee: Posterior oblique repair, medial collateral ligament (MCL) reconstruction, lateral collateral ligament reconstruction, iliotibial (IT) band tenodesis • Foot/Ankle: Achilles tendon repair
Non-clinical Testing	Verification activities were performed on the proposed device and / or its predicates. Performance testing included evaluation of fixation strength, insertion torque, torque to failure and cyclic migration at time zero. Ethylene Oxide Sterilization was validated according to ANSI/AAMI/ISO 11135: 2014 to a SAL of 1×10^{-6}. EO residuals were tested per AAMI/ANSI/ISO 10993-7:2008 The proposed device has been determined to be non-pyrogenic per the requirements set forth in ANSI/AAMI ST-72:2011, United States Pharmacopeia (USP), and European Pharmacopeia (EP) using the bacterial endotoxin testing (BET) method.
Safety and Performance	Results of performance testing have demonstrated that the proposed devices are suitable for their intended use. Based on similarities in the indications for use, technological characteristics, and performance in comparison to the predicate devices, the proposed HEALIX ADVANCE SP PEEK Anchor has shown to be substantially equivalent to the predicate devices under the Federal Food, Drug and Cosmetic Act.
