



Siemens Healthcare Diagnostics Inc.  
Julie Warren  
Regulatory Affairs Professional  
500 GBC Dr. P.O. Box 6101, Mailstop 514  
Newark, DE 19714

Re: K221605  
Trade/Device Name: Emit® II Plus Buprenorphine Assay  
Regulation Number: 21 CFR 862.3650  
Regulation Name: Opiate Test System  
Regulatory Class: Class II  
Product Code: DJG  
Dated: March 16, 2023  
Received: March 17, 2023

Dear Julie Warren:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR

803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Joseph A. Kotarek -S  
Digitally signed by  
Joseph A. Kotarek -S  
Date: 2023.07.25  
10:40:32 -04'00'

Joseph Kotarek, Ph.D.  
Branch Chief  
Division of Chemistry  
and Toxicology Devices  
OHT7: Office of In Vitro Diagnostics  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)  
K221605

Device Name  
Emit® II Plus Buprenorphine Assay

### Indications for Use (Describe)

The Emit® II Plus Buprenorphine Assay is a homogeneous enzyme immunoassay with a 5 ng/mL cutoff. The assay is intended for use in laboratories for the qualitative and/or semiquantitative analyses of buprenorphine in human urine. Emit® II Plus assays are designed for use with a number of chemistry analyzers.

The semiquantitative mode is for the purpose of enabling laboratories to determine an appropriate dilution of the specimen for confirmation by a confirmatory method such as Liquid Chromatography/Mass Spectrometry (LC/MS) or permitting laboratories to establish quality control procedures.

The Emit® II Plus Buprenorphine Assay provides only a preliminary analytical test result. A more specific alternative chemical method(s) must be used to obtain a confirmed analytical result. Gas Chromatography/Mass Spectrometry (GC/MS) or LC/MS are the preferred confirmatory methods. Other chemical confirmation methods are available. Clinical consideration and professional judgment should be applied to any drug-of-abuse test result, particularly when preliminary positive results are used.

For Professional Use.

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.

For in vitro diagnostic use.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*“An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number.”*

This 510(k) Summary of Safety and Effectiveness is being submitted in accordance with the requirements of Safe Medical Device Act of 1990 and 21 CFR 807.92.

The assigned 510(k) number is: K221605

**1. Submitter**

Siemens Healthcare Diagnostics Inc.  
 500 GBC Dr.  
 P.O. Box 6101  
 Newark, DE 19714

Primary Contact Person: Julie Warren  
 Date of Preparation: July 7, 2023

**2. Device Information**

Trade Name	Emit® II Plus Buprenorphine Assay
Common Name	Buprenorphine Assay
Device	Enzyme Immunoassay, Opiates
Review Panel	Toxicology
Product Code	DJG
Submission Type	510(k)
Regulation Number	21 CFR § 862.3650
Device Class	2

**3. Purpose of Submission**

The purpose of this premarket notification is for the Emit® II Plus Buprenorphine Assay on the DxC 500 AU Clinical Chemistry Analyzer.

**4. Predicate Device**

Predicate Device Name: Emit® II Plus Buprenorphine Assay  
 510(k) Number: K150606

## 5. Device Description

### **Emit® II Plus Buprenorphine Assay**

The Emit® II Plus Buprenorphine Assay is a homogeneous enzyme immunoassay technique used for the analysis of specific compounds in human urine. The assay is based on competition between drug in the specimen and drug labeled with the recombinant glucose-6-phosphate dehydrogenase (rG6PDH) for antibody binding sites. Enzyme activity decreases upon binding to the antibody, so the drug concentration in the specimen can be measured in terms of enzyme activity. Active enzyme converts nicotinamide adenine dinucleotide (NAD) to NADH in the presence of glucose-6-phosphate (G6P), resulting in an absorbance change that is measured spectrophotometrically. Endogenous serum G6PDH does not interfere because the coenzyme NAD functions only with the bacterial (*Leuconostoc mesenteroides*) enzyme employed in the assay.

The Emit® II Plus Buprenorphine Assay reagents are provided liquid, ready to use and may be used directly from the refrigerator. The product is sold in three (3) kit sizes: 28 mL, 115 mL, and 1000 mL. Reagents 1 and 2 are provided as a matched set. They should not be interchanged with components of kits with different lot numbers.

#### **Antibody/Substrate Reagent 1**

Mouse monoclonal antibodies to buprenorphine (0.53 µg/mL)\*, NAD (6.9 mM), G6P (10.9 mM), bovine serum albumin, preservatives, and stabilizers.

\*The antibody titer and enzyme conjugate activity may vary from lot to lot.

#### **Enzyme Reagent 2**

Norbuprenorphine labeled with bacterial rG6PDH (0.50 µg/mL), HEPES buffer, bovine serum albumin, preservatives, and stabilizers, where the antibody titer and enzyme conjugate activity may vary from lot to lot.

**6. Intended Use/Indications for Use**

**Emit® II Plus Buprenorphine Assay**

The Emit® II Plus Buprenorphine Assay is a homogeneous enzyme immunoassay with a 5 ng/mL cutoff. The assay is intended for use in laboratories for the qualitative and/or semiquantitative analyses of buprenorphine in human urine. Emit® II Plus assays are designed for use with a number of chemistry analyzers.

The semiquantitative mode is for the purpose of enabling laboratories to determine an appropriate dilution of the specimen for confirmation by a confirmatory method such as Liquid Chromatography/Mass Spectrometry (LC/MS) or permitting laboratories to establish quality control procedures.

The Emit® II Plus Buprenorphine Assay provides only a preliminary analytical test result. A more specific alternative chemical method(s) must be used to obtain a confirmed analytical result. Gas Chromatography/Mass Spectrometry (GC/MS) or LC/MS are the preferred confirmatory methods. Other chemical confirmation methods are available. Clinical consideration and professional judgment should be applied to any drug-of-abuse test result, particularly when preliminary positive results are used.

For Professional Use.

**Caution:** Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.

For *in vitro* diagnostic use.

**7. Comparison of Technological Characteristics with the Predicate Device**

The similarities and differences between the proposed device Emit II Plus Buprenorphine Assay on Beckman Coulter DxC AU 500 Clinical Chemistry Analyzer and the predicate device, Emit II Plus Buprenorphine Assay, were compared side by side.

Feature	Predicate Device: Emit II Plus Buprenorphine Assay	Proposed Device: Emit II Plus Buprenorphine Assay
Intended Use	A homogeneous enzyme immunoassay with a 5 ng/mL cutoff. The assay is intended for use in laboratories for the qualitative and/or	Same

Feature	Predicate Device: Emit II Plus Buprenorphine Assay	Proposed Device: Emit II Plus Buprenorphine Assay
	semiquantitative analyses of buprenorphine in human urine. Preliminary analytical test result.	
Analyte	Buprenorphine	Same
Cutoff Level	5 ng/mL	Same
Measurement	Qualitative, semi-quantitative	Same
Methodology	EMIT® technology: Homogeneous enzyme immunoassay	Same
Detection	Absorbance change measured Spectrophotometrically at 340nm	Same
Instrument	Analyzers must be capable of maintaining a constant reaction temperature, pipette specimens/reagents, mix thoroughly, measure enzyme rates precisely and time the reaction accurately.	Beckman Coulter DxC 500 AU Analyzer
Sample Type	Human Urine	Same
Reagent form and storage	Liquid – ready to use, on-board storage	Same
Calibrator	0, 2.5, 5, 15, and 25 ng/mL	Same
Control Levels	Positive: 7 ng/mL Negative: 3 ng/mL	Same

## 8. Performance Data

The following performance studies were conducted to demonstrate substantial equivalence of the Emit II Plus Buprenorphine Assay on the DXC 500 AU Clinical Chemistry Analyzer to the Emit II Plus Buprenorphine Assay previously cleared under premarket notification 510(k), K150606.

### 8.1 Method Comparison

The method comparison study was evaluated in accordance with CLSI EP09c 3rd Edition Measurement Procedure Comparison and Bias Estimation Using Patient Samples Evaluation of Precision Performance of Quantitative Measurement Methods and CLSI EP 12-A2 User Protocol for Evaluation of Qualitative Test Performance to determine accuracy and verification of the 5 ng/mL cutoff.

One hundred twenty (120) native urine individual patient samples were collected from external suppliers. The samples were received frozen, and thawed before testing. The samples were distributed across the assay range, 21% of which were within -50% of the cutoff and 19% of which were within +50% of the cutoff.

Each specimen was tested on the DxC 500 AU and on LC/MS/MS. The results were analyzed using a 2x2 box plot to assess qualitative agreement relative to the 5 ng/mL cutoff value.

**Table 1. Method Comparison - Qualitative Analysis**

		LC/MS/MS	
		Neg	Pos
Emit II Plus Buprenorphine on DxC 500 AU	Neg	48	0
	Pos	9	63

Result: 92.5% Agreement

**Table 2. Method Comparison – Discordant Results**

Sample ID	DxC 500 AU (ng/mL)	Buprenorphine (ng/mL)	Norbuprenorphine (ng/mL)	Buprenorphine+ Norbuprenorphine (ng/mL)	DxC 500 AU Pos/Neg	LC/MS/MS Pos/Neg
BUP35	6.1	0.0	4.2	4.2	Pos	Neg
BUP45	6.0	0.0	4.6	4.6	Pos	Neg
BUP51	6.1	0.0	4.3	4.3	Pos	Neg
BUP60	5.6	0.0	3.8	3.8	Pos	Neg
BUP62	6.1	0.0	3.9	3.9	Pos	Neg
BUP87	5.5	0.0	4.3	4.3	Pos	Neg
BUP95	5.6	0.0	4.5	4.5	Pos	Neg
BUP101	6.3	0.0	4.8	4.8	Pos	Neg
BUP109	5.0	0.0	4.3	4.3	Pos	Neg

One hundred twenty (120) samples were tested. Relative to the 5 ng/mL cutoff, forty-eight (48) samples were negative by Emit II Plus Buprenorphine Assay and LC/MS/MS, and sixty-three (63) were positive. Nine (9) discordant samples were observed as positive by Emit II Plus Buprenorphine and negative by LC/MS/MS.

## 8.2 Recovery

The recovery was evaluated in accordance with the following testing protocol. Nine (9) urine samples were prepared by spiking buprenorphine into negative urine pools. The following levels were tested: 2.0, 3.0, 4.0, 5.0, 8.0, 12.0, 18.0, 22.0, and 25.0 ng/mL.

Experimental design was (1) reagent lot, one (1) calibrator lot, one (1) instrument, one (1) calibration with two (2) replicate measurements each. Samples were assayed in five (5) replicates. Samples were additionally analyzed by LC/MS/MS to determine the reference value. The nominal value was determined gravimetrically based on the stock concentration.

**Table 3. Recovery of Emit II Plus Buprenorphine - Semi-Quantitative Analysis**

Sample ID	Nominal Value (ng/mL)	DxC 500 AU Value (ng/mL)	% Recovery vs Nominal	LC/MS/MS (ng/mL)	% Recovery vs LC/MS/MS
BUP1	2.0	1.9	95	2.2	86
BUP2	3.0	2.9	97	3.2	91
BUP3	4.0	4.1	103	4.4	93
BUP4	5.0	4.9	98	5.1	96
BUP5	8.0	7.6	95	8.8	86
BUP6	12.0	11.6	97	13.1	89
BUP7	18.0	18.1	101	18.1	100
BUP8	22.0	22.9	104	22.9	100
BUP9	25.0	24.7	99	25.6	96

## 8.3 Precision

### Repeatability and Within-Lab Precision

The precision was evaluated in accordance with CLSI EP05-A3 *Evaluation of Precision Performance of Quantitative Measurement Methods*.

Eleven (11) urine samples were prepared by spiking buprenorphine into negative urine pools, aliquoted and frozen prior to the start of testing. The following levels were tested: 0, 1.25, 2.5, 3.0, 3.75, 5.0, 6.25, 7.0, 7.5, 8.75, and 10.0 ng/mL.

Experimental design was a 20x2x2 design using (1) reagent lot, one (1) calibrator lot, one (1) instrument, one (1) calibration, twenty (20) test days and two (2) runs per test

day with two (2) replicate measurements each. Each day, two runs were performed, with a minimum of two (2) hours in between.

**Table 4. Precision of Emit II Plus Buprenorphine Qualitative Analysis**

Sample Concentration (ng/mL)	% of the Cutoff	# of Results	Results (Pos/Neg)	Mean (mA/min)
0.0	-100%	80	80 Negative	280
1.25	-75%	80	80 Negative	297
2.5	-50%	80	80 Negative	316
3.0	-40%	80	80 Negative	323
3.75	-25%	80	80 Negative	333
5.0	0%	80	21 Negative 59 Positive	363
6.25	+25%	80	80 Positive	380
7.0	+40%	80	80 Positive	392
7.5	+50%	80	80 Positive	402
8.75	+75%	80	80 Positive	427
10.0	+100%	80	80 Positive	444

**Table 5. Precision of Emit II Plus Buprenorphine Semi-Quantitative Analysis**

Sample Concentration (ng/mL)	% of the Cutoff	# of Results	Results (Pos/Neg)	Mean (ng/mL)	Repeatability		Within-Lab Precision	
					SD	%CV	SD	%CV
0.0	-100%	80	80 Negative	0.1	0.06	-	0.11	-
1.25	-75%	80	80 Negative	1.2	0.08	6.6	0.08	6.6
2.5	-50%	80	80 Negative	2.3	0.09	3.7	0.11	4.6
3.0	-40%	80	80 Negative	2.7	0.07	2.8	0.15	5.7
3.75	-25%	80	80 Negative	3.3	0.08	2.5	0.15	4.4
5.0	0%	80	24 Negative 56 Positive	5.1	0.08	1.6	0.18	3.5
6.25	+25%	80	80 Positive	6.1	0.08	1.3	0.21	3.4
7.0	+40%	80	80 Positive	6.8	0.11	1.7	0.24	3.5
7.5	+50%	80	80 Positive	7.3	0.08	1.1	0.21	2.9
8.75	+75%	80	80 Positive	8.8	0.10	1.1	0.24	2.8
10.0	+100%	80	80 Positive	9.7	0.10	1.1	0.37	3.8

### **Reproducibility Precision**

The reproducibility was evaluated in accordance with the following protocol.

Three (3) urine samples were tested: Emit II Plus Specialty Drug Control Negative (3 ng/mL) Emit II Plus Specialty Drug Cal/Ctrl (Cutoff Calibrator), and the Emit II Plus Specialty Drug Control Positive (7 ng/mL).

Experimental design was a 3x5x5 design; 3 (Sites/Instruments) × 5 (days) × 5 (replicates per day) for a total seventy-five (75) replicates per reagent lot. To capture lot-to-lot variability, three (3) Emit II Plus Buprenorphine reagent lots were processed in parallel in this study design. Testing was completed by three (3) operators. Each testing day, one (1) run was performed on three (3) instruments.

**Table 6. Total Reproducibility Precision of Emit II Plus Buprenorphine Assay**

Sample	N	Mean (ng/mL)	SD	%CV
Negative Control	225	3.1	0.17	5.3
Cutoff Calibrator	225	5.0	0.18	3.7
Positive Control	225	7.1	0.21	3.0

**Table 7. Reproducibility Precision of Emit II Plus Buprenorphine Assay**

			Repeatability		Between Day		Between Lot		Between Instrument/Site	
Sample	N	Mean (ng/mL)	SD	%CV	SD	%CV	SD	%CV	SD	%CV
Negative Control	225	3.1	0.11	3.6	0.06	1.8	0.03	0.8	0.11	3.4
Cutoff Calibrator	225	5.0	0.11	2.2	0.07	1.4	0.00	0.0	0.13	2.6
Positive Control	225	7.1	0.12	1.7	0.08	1.2	0.05	0.7	0.15	2.1

#### 8.4 Limit of Detection

The limits of detection, Limit of Blank (LoB) and Limit of Detection (LoD) were evaluated in accordance with CLSI EP17-A2: *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures*.

Blank samples were drug-free negative urine, and the low samples were drug-free urine spiked with buprenorphine.

Experimental design was (2) reagent lots, one (1) calibrator lot, one (1) instrument, one (1) calibration, three (3) testing days, four (4) blank samples with minimum of sixty (60) replicate measurements (across all blank samples, days, instrument) and four (4) samples, each containing a low level of measurand (LoD) with minimum of sixty (60) replicate measurements (across all low-level samples, days, instrument).

**Table 8. Limit of Detection Summary**

Sample	Limit of Detection Lot 1 (ng/mL)	Limit of Detection Lot 2 (ng/mL)	Limit of Detection Final (ng/mL)	Claimed Limit of Detection (ng/mL)
LoD Sample	0.7	0.8	0.8	≤ 1.0

### 8.5 Specificity and Cross-reactivity

The specificity was evaluated in accordance with CLSI EP07-ED3 *Interference Testing Clinical Chemistry*.

The experimental design was one (1) reagent lot, one (1) calibrator lot, one (1) instrument, one (1) sample pool for each compound tested. The samples were assayed in five (5) replicates.

**Table 9. Specificity - Buprenorphine Metabolites**

Compound	Concentration Tested (ng/mL)	Mean (ng/mL)	Cutoff (ng/mL)	% Cross-reactivity
Buprenorphine	5 ng/mL	4.9	5.0	98%
Norbuprenorphine	5 ng/mL	5.3	5.0	106%
Buprenorphine Glucuronide	1,000 ng/mL	1.4	5.0	0.10%
Norbuprenorphine Glucuronide	1,000 ng/mL	1.4	5.0	0.10%

**Table 10. Cross-reactivity - Structurally Related Compounds**

Description	Concentration @ Testing	Mean Rate (mA/min)	Cutoff Rate (mA/min)	Response vs. Cutoff (Pos/Neg)	Mean (ng/mL)	%Cross-reactivity to 5 ng/mL Cutoff
6-Acetylcodeine	100,000 ng/mL	131	167	Neg	-0.6	<0.1%
6-Acetylmorphine	100,000 ng/mL	137	167	Neg	0.4	<0.1%
Codeine	100,000 ng/mL	135	167	Neg	0.1	<0.1%
Dextromethorphan	100,000 ng/mL	137	167	Neg	0.3	<0.1%
Dihydrocodeine	100,000 ng/mL	136	167	Neg	0.3	<0.1%
Ethyl Morphine	100,000 ng/mL	136	167	Neg	0.3	<0.1%
Heroin	100,000 ng/mL	133	167	Neg	0.0	<0.1%
Hydrocodone	100,000 ng/mL	136	167	Neg	0.2	<0.1%

Description	Concentration @ Testing	Mean Rate (mA/min)	Cutoff Rate (mA/min)	Response vs. Cutoff (Pos/Neg)	Mean (ng/mL)	%Cross-reactivity to 5 ng/mL Cutoff
Hydromorphone	100,000 ng/mL	137	167	Neg	0.4	<0.1%
Levorphanol	100,000 ng/mL	141	167	Neg	1.0	<0.1%
Morphine	100,000 ng/mL	135	167	Neg	0.1	<0.1%
Morphine 3-glucuronide	100,000 ng/mL	137	167	Neg	0.4	<0.1%
Morphine 6-glucuronide	100,000 ng/mL	136	167	Neg	0.3	<0.1%
Nalorphine	100,000 ng/mL	136	167	Neg	0.3	<0.1%
Naloxone	100,000 ng/mL	138	167	Neg	0.6	<0.1%
Naltrexone	100,000 ng/mL	137	167	Neg	0.4	<0.1%
Norcodeine	100,000 ng/mL	135	167	Neg	0.1	<0.1%
Normorphine	100,000 ng/mL	136	167	Neg	0.2	<0.1%
Noroxycodone	100,000 ng/mL	135	167	Neg	0.2	<0.1%
Noroxymorphone	100,000 ng/mL	136	167	Neg	0.2	<0.1%
Oxycodone	100,000 ng/mL	136	167	Neg	0.3	<0.1%
Oxymorphone	100,000 ng/mL	137	167	Neg	0.4	<0.1%

**Table 11. Cross-reactivity - Structurally Un-Related Compounds**

Compound	Conc. Tested (µg/mL)	-40% Control Mean Semi-Quant (ng/mL)	-40% Control Response vs. Cutoff (Pos/Neg)	+40% Control Mean Semi-Quant (ng/mL)	+40% Control Response vs. Cutoff (Pos/Neg)
10,11-Dihydrocarbamazepine	85.0	2.8	Neg	6.9	Pos
Acetaminophen	1000.0	2.5	Neg	7.0	Pos
Acetylsalicylic Acid	1500.0	2.5	Neg	7.0	Pos
Amitriptyline	100.0	3.1	Neg	7.5	Pos
Amoxicillin	500.0	2.5	Neg	6.9	Pos
AZT (Zidovudine)	2000.0	2.5	Neg	6.8	Pos
Brompheniramine	75.0	2.7	Neg	6.9	Pos
Caffeine	1000.0	2.7	Neg	7.0	Pos
Captopril	500.0	2.6	Neg	6.9	Pos

Compound	Conc. Tested (µg/mL)	-40% Control Mean Semi-Quant (ng/mL)	-40% Control Response vs. Cutoff (Pos/Neg)	+40% Control Mean Semi-Quant (ng/mL)	+40% Control Response vs. Cutoff (Pos/Neg)
Chlordiazepoxide	100.0	2.6	Neg	7.1	Pos
Chlorpromazine	10.0	2.9	Neg	7.2	Pos
Cimetidine	1000.0	2.9	Neg	7.4	Pos
Clomipramine	2.5	2.8	Neg	7.3	Pos
Clonidine	1000.0	2.9	Neg	7.2	Pos
Cyclobenzaprine	125.0	3.1	Neg	7.6	Pos
d-Amphetamine	700.0	2.9	Neg	7.3	Pos
Desipramine	800.0	3.3	Neg	8.0	Pos
Diazepam	100.0	3.1	Neg	7.4	Pos
Digoxin	0.01	2.9	Neg	6.9	Pos
Diphenhydramine	1000.0	4.6	Neg	8.8	Pos
d-Methamphetamine	500.0	2.9	Neg	7.2	Pos
Doxepin	100.0	3.0	Neg	7.6	Pos
EDDP (2-Ethylidene-1.5dimethyl-3.3-duogenylidine)	1000.0	3.6	Neg	8.3	Pos
EMDP	100.0	2.7	Neg	7.0	Pos
Enalapril	500.0	2.7	Neg	7.3	Pos
Fluoxetine	500.0	3.0	Neg	7.8	Pos
Glutethimide	500.0	2.8	Neg	7.4	Pos
Haloperidol	100.0	2.7	Neg	7.4	Pos
Hydroxyzine (dihydrochloride)	500.0	3.3	Neg	8.0	Pos
Ibuprofen	1000.0	2.9	Neg	7.4	Pos
Imipramine	200.0	3.2	Neg	7.6	Pos
Ketamine	100.0	2.6	Neg	7.0	Pos
Ketorolac Tromethamine	400.0	2.3	Neg	6.7	Pos
LAAM (L-a-Acetylmethadol)	25.0	2.7	Neg	7.4	Pos
L-Cotinine	100.0	2.6	Neg	7.0	Pos

Compound	Conc. Tested (µg/mL)	-40% Control Mean Semi-Quant (ng/mL)	-40% Control Response vs. Cutoff (Pos/Neg)	+40% Control Mean Semi-Quant (ng/mL)	+40% Control Response vs. Cutoff (Pos/Neg)
Levofloxacin	100.0	2.8	Neg	7.1	Pos
Levothyroxine (L-Thyroxine)	50.0	2.6	Neg	7.1	Pos
Lidocaine	1000.0	2.7	Neg	7.1	Pos
Lormetazepam	1.0	2.5	Neg	7.0	Pos
LSD	10.0	2.5	Neg	7.1	Pos
MDMA (Ecstasy)	1000.0	2.6	Neg	7.3	Pos
Meperidine	800.0	2.9	Neg	7.5	Pos
Methadone	500.0	3.3	Neg	7.9	Pos
Methaqualone	600.0	3.0	Neg	7.6	Pos
NAPA (N-Acetylprocainamide)	400.0	2.6	Neg	6.8	Pos
Naproxen	1000.0	2.4	Neg	6.6	Pos
Nicotinic Acid	500.0	2.5	Neg	7.2	Pos
Nifedipine	500.0	2.5	Neg	6.4	Pos
Nordiazepam	100.0	2.7	Neg	7.2	Pos
Nortriptyline	250.0	2.6	Neg	7.4	Pos
Oxazepam	300.0	2.4	Neg	6.6	Pos
Perphenazine	150.0	2.9	Neg	7.6	Pos
Phencyclidine	900.0	4.6	Neg	8.9	Pos
Phenobarbital	500.0	2.5	Neg	7.1	Pos
Phenelzine (sulfate)	100.0	2.6	Neg	6.0	Pos
Phenytoin	1000.0	2.4	Neg	6.0	Pos
Procainamide (HCl)	1000.0	2.6	Neg	7.0	Pos
Procyclidine (HCl)	800.0	3.2	Neg	8.0	Pos
Promethazine (HCl)	100.0	3.7	Neg	8.1	Pos
Propoxyphene	1000.0	3.5	Neg	8.0	Pos
Protriptyline (HCl)	200.0	3.3	Neg	7.9	Pos
Pseudoephedrine	1000.0	2.7	Neg	7.0	Pos

Compound	Conc. Tested (µg/mL)	-40% Control Mean Semi-Quant (ng/mL)	-40% Control Response vs. Cutoff (Pos/Neg)	+40% Control Mean Semi-Quant (ng/mL)	+40% Control Response vs. Cutoff (Pos/Neg)
Quinacrine (HCl)	900.0	4.8	Neg	9.0	Pos
Ranitidine	1000.0	2.6	Neg	7.1	Pos
Ritalin	1000.0	2.7	Neg	7.5	Pos
Salicylic Acid	500.0	2.5	Neg	7.1	Pos
Scopolamine	500.0	2.5	Neg	6.9	Pos
Secobarbital	1000.0	2.7	Neg	7.5	Pos
Tapentadol	100.0	2.4	Neg	7.0	Pos
THC (11-nor- $\Delta^9$ -THC-9-COOH)	100.0	2.4	Neg	6.7	Pos
Thioridazine	100.0	4.1	Neg	8.6	Pos
Tramadol	1000.0	3.0	Neg	7.8	Pos
Trazodone	5.0	2.4	Neg	6.8	Pos
Trimethoprim	1000.0	2.5	Neg	7.1	Pos
Tripolidine (zymine)	50.0	2.6	Neg	7.2	Pos
Tyramine	100.0	2.5	Neg	6.9	Pos
Verapamil	500.0	2.8	Neg	7.6	Pos
Zolpidem	100.0	2.7	Neg	7.1	Pos

### 8.6 Interference

Interference was evaluated in accordance with CLSI EP07-ED3 *Interference Testing Clinical Chemistry*.

Urine testing samples were prepared by spiking exogenous materials (metabolites and structurally related compounds) into urine samples at  $\pm$  40% positive and negative control levels (3 ng/mL and 7 ng/mL respectively).

The experimental design was one (1) reagent lot, one (1) calibrator lot, one (1) instrument, one (1) sample pool for each compound tested. The samples were assayed in five (5) replicates.

**Table 12. Interference – Endogenous and Exogenous Substances**

<b>Interferent Sample</b>	<b>Interference concentration</b>	<b>-40% Control Response vs. Qual Cutoff (Pos/Neg)</b>	<b>-40% Control Sample Mean (ng/mL)</b>	<b>+40% Control Response vs. Qual Cutoff (Pos/Neg)</b>	<b>+40% Control Sample Mean (ng/mL)</b>
Acetone	1.0 g/dL	Neg	2.6	Pos	6.3
Ascorbic Acid	1.5 g/dL	Neg	2.9	Pos	6.3
Bilirubin (Conjugated)	2.0 mg/dL	Neg	2.2	Pos	6.2
Bilirubin (Unconjugated)	2.0 mg/dL	Neg	2.4	Pos	5.9
Creatinine	0.5 g/dL	Neg	2.7	Pos	6.3
Ethanol	1.0 g/dL	Neg	2.5	Pos	6.4
IgG	0.5 g/dL	Neg	2.5	Pos	7.1
Glucose	2.0 g/dL	Neg	2.5	Pos	6.1
Hemoglobin	115 mg/dL	Neg	3.0	Pos	6.7
Human Serum Albumin	0.5 g/dL	Neg	2.9	Pos	6.7
Oxalic Acid	0.1 g/dL	Neg	2.5	Pos	6.5
Riboflavin	7.5 mg/dL	Neg	2.5	Pos	6.4
Sodium Chloride	6.0 g/dL	Neg	2.5	Pos	6.6
Urea	6.0 g/dL	Neg	2.4	Pos	6.2
Sodium Azide	1% w/v	Neg	2.5	Pos	6.4
Sodium Fluoride	1% w/v	Neg	2.5	Pos	6.1
Galactose	1.0 g/dL	Neg	2.5	Pos	5.8

**Table 13. Interference – pH and Specific Gravity**

<b>Interferent Sample</b>	<b>-40% Control Response vs. Qual Cutoff (Pos/Neg)</b>	<b>-40% Control Mean (ng/mL)</b>	<b>+40% Control Response vs. Qual Cutoff (Pos/Neg)</b>	<b>+40% Control Mean (ng/mL)</b>
pH 3	Neg	2.9	Pos	7.2
pH 4	Neg	2.9	Pos	7.5
pH 5	Neg	3.0	Pos	7.3
pH 6	Neg	3.0	Pos	7.3
pH 7	Neg	2.6	Pos	7.0
pH 8	Neg	2.4	Pos	6.9
pH 9	Neg	2.6	Pos	6.2
pH 10	Neg	2.6	Pos	6.0
pH 11	Neg	2.5	Pos	6.2
sG 1.002	Neg	2.7	Pos	6.4
sG 1.005	Neg	2.9	Pos	6.8
sG 1.010	Neg	2.9	Pos	7.0
sG 1.015	Neg	2.9	Pos	7.2
sG 1.020	Neg	2.9	Pos	7.2
sG1.025	Neg	3.0	Pos	7.2
sG 1.030	Neg	2.6	Pos	7.0
sG 1.035	Neg	2.5	Pos	6.9

**9. Conclusion**

Based on a comparison of the technological characteristic similarities and differences, the Emit® II Plus Buprenorphine Assay (K150606) and the Emit® II Plus Buprenorphine Assay have equivalent intended use and technology attributes. Performance studies were completed to demonstrate the performance of the proposed device to the predicate device were substantially equivalent. Through the accuracy comparison of both methods to the reference LC/MS/MS as well as performance studies of precision, recovery, and specificity, the substantial equivalence is supported.