



January 4, 2024

Koya Medical, Inc.
% Alex Chang
Regulatory Affairs Consultant
BioDesign Regulatory Services, LLC
16185 Los Gatos Blvd.
Los Gatos, California 95032

Re: K223228
Trade/Device Name: Dayspring
Regulation Number: 21 CFR 870.5800
Regulation Name: Compressible Limb Sleeve
Regulatory Class: Class II
Product Code: JOW
Dated: May 1, 2023
Received: May 2, 2023

Dear Alex Chang:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Eric E. Richardson -S 

for Nicole Gillette

Assistant Director

DHT2B: Division of Circulatory Support,
Structural and Vascular Devices

OHT2: Office of Cardiovascular Devices

Office of Product Evaluation and Quality

Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K223228

Device Name

Dayspring

Indications for Use (Describe)

The Koya Dayspring system is a prescription only wearable compression system that is intended for use in a clinic or home setting by medical professionals and patients who are under medical supervision to increase lymphatic flow in the treatment of many conditions such as:

- Lymphedema
- Primary lymphedema
- Post mastectomy edema
- Edema following trauma and sports injuries
- Post immobilization edema
- Venous insufficiency
- Reducing wound healing time
- Treatment and assistance in healing stasis dermatitis, venous stasis ulcers, or arterial and diabetic leg ulcers
- Lipedema
- Phlebolymphe

The Dayspring system is developed on a wearable compression technology platform, which is designed to provide mobility for patients.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

K223228

510(k) Summary

This 510(k) summary of safety and effectiveness information is submitted in accordance with the requirements of SMDA 1990 and 21 CFR 807.92.

I. SUBMITTER: Koya Medical, Inc.
2461 Peralta St. Oakland CA 94607 USA
Establishment Registration: 3017424826

CONTACT: Jarren Baldwin
VP Operations & Technology
P: 415 851 0337
E-mail: jarren@koyamedical.com

DATE PREPARED: May 1 2023

II. DEVICE:

TRADE NAME: DAYSPRING
CLASSIFICATION NAME: SLEEVE, LIMB, COMPRESSIBLE
DEVICE CLASSIFICATION: CLASS II
PRODUCT CODE: JOW

III. PREDICATE DEVICES:

Primary Predicate
Manufacture: Koya Medical Inc
Trade Name: Dayspring
510(k): K210885

Secondary Predicate
Manufacture: Tactile System Technology
Trade Name: Flexitouch Plus System (PD32-G3)
510(k): K203178

IV. DEVICE DESCRIPTION:

The Dayspring system consists of two main components: a controller and garment. The garment is powered by an active smart compression technology that is calibrated, instant-acting, and silent. This technology uses a Nickel Titanium (Ni-Ti) shape-memory alloy programmed by the controller. A liner or lite weight clothing is worn under the garment to prevent direct patient contact with the garment. The garment is wrapped around the patient's affected area so that the device fits snugly. The device has independently controlled sections in each limb. The controller can be programmed to provide graduated sequential compression therapy to the affected area over a range of 0-100 mmHg. The device is powered by a rechargeable Lithium-ion battery pack. The device was developed to provide patients with untethered access and a functional range of motion and mobility.

V: INDICATIONS FOR USE:

The Dayspring system is a prescription only wearable compression system that is intended for use in a clinic or home setting by medical professionals and patients who are under medical supervision to increase lymphatic flow in the treatment of many conditions such as:

- Lymphedema
- Primary lymphedema
- Post mastectomy edema
- Edema following trauma and sports injuries
- Post immobilization edema
- Venous insufficiency
- Reducing wound healing time
- Treatment and assistance in healing stasis dermatitis, venous stasis ulcers, or arterial and diabetic leg ulcers
- Lipedema
- Phlebolympheidema

The Dayspring system is developed on a wearable compression technology platform, which is designed to provide mobility for patients.

VI. COMPARISON OF TECHNOLOGICAL CHARACTERISTICS WITH THE PREDICATE DEVICE:

Feature	Subject Device	Primary Predicate (K210885)	Additional Predicate (K203178)
Indications for use	The Koya Dayspring system is a prescription only wearable compression system that is intended for use in a clinic or home setting by medical professionals and patients who are under medical supervision to increase lymphatic flow in the treatment of many conditions such as: • Lymphedema • Primary lymphedema • Post mastectomy edema • Edema following trauma and sports injuries • Post immobilization edema • Venous insufficiency • Reducing wound healing time • Treatment and assistance in healing stasis dermatitis, venous stasis ulcers, or arterial and diabetic leg ulcers • Lipedema • Phlebolymphedema The Dayspring system is developed on a wearable compression technology platform, which is designed to provide mobility for patients.	The Koya Dayspring system is a prescription only wearable compression system that is intended for use in a clinic or home setting by medical professionals and patients who are under medical supervision to increase lymphatic flow in the treatment of many conditions such as: • Lymphedema • Primary lymphedema • Post mastectomy edema • Edema following trauma and sports injuries • Post immobilization edema • Venous insufficiency • Reducing wound healing time • Treatment and assistance in healing stasis dermatitis, venous stasis ulcers, or arterial and diabetic leg ulcers • Lipedema • Phlebolymphedema The Dayspring system is developed on a wearable compression technology platform, which is designed to provide mobility for patients.	The Flexitouch Plus System and garments for legs, arms, trunk, and chest are intended for use by medical professionals and patients who are under medical supervision to increase lymphatic flow in the treatment of many conditions such as: • Lymphedema • Primary lymphedema • Post mastectomy edema • Edema following trauma and sports issues • Post immobilization edema • Venous insufficiency • Reducing wound healing time • Treatment and assistance in healing stasis dermatitis, venous stasis ulcers, arterial ulcers, and diabetic leg ulcers • Lipedema • Phlebolymphedema The Flexitouch Plus System and garments for the head and neck are intended for use by medical professionals and patients who are

			under supervision for the treatment of head and neck lymphedema.
Electrical Requirements	Rechargeable Li-ion Battery Pack, with 0.800 A input from 90-264 VAC 50/60 Hz to AC Adapter, with output voltage of 25.0V DC and 3.0A	Rechargeable Li-ion Battery Pack, with 0.800 A input from 90-264 VAC 50/60 Hz to AC Adapter, with output voltage of 25.0V DC and 3.0A	100-240 VAC 50/60 Hz to AC Adapter with output voltage of 12.0V DC and 3.0A
Output	Sequential calibrated gradient Pressure	Sequential calibrated gradient Pressure	Sequential calibrated gradient Pressure
Mechanism of Action	Exertion of sequential pressure to affected area	Exertion of sequential pressure to affected area	Exertion of sequential pressure to affected area
Principles of Operation	Lithium-ion battery powered integrated shape memory alloy channels creating compressive pressure	Lithium-ion battery powered integrated shape memory alloy channels creating compressive pressure	Electrically powered integrated pneumatic air channels creating compressive pressure
Device Total Pressure Range	0-100 mmHg	0-100 mmHg	0-100 mmHg
Controller Unit			
Controller unit size and weight	3.4" x 5.2" x 1.5" 0.80 lbs	3.4" x 5.2" x 1.5" 0.80 lbs	8"x10"x8" 6.2 lbs
Controller Enclosure Material	All plastic construction	All plastic construction	All plastic construction
User Interface	Pushbuttons. Also available is Bluetooth Low Energy (BLE) Module for communication with mobile application on mobile device	Pushbuttons. Also available is Bluetooth Low Energy (BLE) Module for communication with mobile application on mobile device	Pushbuttons Mobile application or BLE not available
Software/Hardware	Analog and digital electronic with microprocessor	Analog and digital electronic with microprocessor	Analog and digital electronic with microprocessor

Appliances	• Upper Extremity (Arm) • Lower Extremity (Leg + Abdomen)	• Upper Extremity (Arm) • Lower Extremity (Leg)	• Upper Extremity (Arm + Chest) • Lower Extremity (Leg + Abdomen) • Head and Neck
Appliance Material	Nylon fabric with velcro straps	Nylon fabric with velcro straps	Nylon fabric with velcro straps
Stockinette/Liner	Class I biocompatible liner provided with the unit	Class I biocompatible liner provided with the unit	Class I biocompatible liner provided with the unit

VII. PERFORMANCE DATA:

The following performance data were provided in support of substantial equivalence determination.

Biocompatibility Testing

The subject device contains the addition of an abdominal garment that consists of a static and active compression wrap that is placed over lite-weight clothing. The arm and leg garments evaluated in previous 510(k) applications are applied over patient-contacting circular-knit liners that have been evaluated for biocompatibility per ISO 10993-1, ISO 10993-5, and ISO 10993-10. Previous applications submitted for the Dayspring system demonstrated that all patient contacting material is compliant to all applicable biocompatibility safety standards.

Sterilization & Shelf-life Testing

The subject device is non-sterile, and components are unlikely to deteriorate with age. Accelerated 1-year accelerated shelf-life testing was performed and shown to support shelf stability.

Electrical Safety and Electromagnetic Compatibility (EMC)

The subject device was evaluated based on the following applicable performance and safety standards: IEC 60601-1:2012, IEC 60601-1-11:2015 and IEC 60601-1-2:2014. Results demonstrated that the subject device was compliant to all applicable performance and safety standards.

Software Verification and Validation Testing

The subject device includes embedded firmware in the controller which has the ability to connect with custom software installed on a mobile device. The system software exhibits a moderate level of concern.

Software lifecycle planning and documentation as well verification and validation testing were performed in accordance with IEC 62304:2015 and as recommended by the following FDA Guidance documents for Industry and FDA Staff:

- Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices
- Guidance for Industry, FDA Reviewers and Compliance on Off-The-Shelf Software Use in Medical Devices
- Guidance for Industry and Food and Drug Administration Staff Content of Premarket Submissions for Management of Cybersecurity in Medical Devices
- Guidance for Industry, FDA Reviewers and Compliance on Postmarket Management of Cybersecurity in Medical Devices

Mechanical Testing

The subject device was evaluated based on the following benchtop performance tests.

Performance Study	Overview	Status
Packaging Transportation Simulation	ASTM D4169	Completed / Pass
Pressure Verification Testing	Pressure capable of delivering 0-100 mmHg compression pressures.	Completed / Pass

Animal Study

Animal performance testing was not required to demonstrate safety and effectiveness of the device.

Clinical Study

Clinical testing was not required to demonstrate the safety and effectiveness of the Dayspring system. Instead, substantial equivalence is based upon benchtop performance testing.

VIII. CONCLUSION:

The data included in this submission demonstrates that the modified Dayspring is substantially equivalent to the primary predicate device, Dayspring (K210885).