



Spineway
% J.D. Webb
Consultant
The OrthoMedix Group, Inc.
4313 West 3800 South
West Haven, Utah 84401

October 10, 2023

Re: K231658

Trade/Device Name: VEOS Spinal Fixation System
Regulation Number: 21 CFR 888.3070
Regulation Name: Thoracolumbosacral Pedicle Screw System
Regulatory Class: Class II
Product Code: NKB, KWP
Dated: October 5, 2023
Received: October 6, 2023

Dear J. D. Webb:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Eileen

Cadel -S

Digitally signed
by Eileen Cadel -S
Date: 2023.10.10
16:21:03 -04'00'

for

Colin O'Neill, M.B.E.

Assistant Director

DHT6B: Division of Spinal Devices

OHT6: Office of Orthopedic Devices

Office of Product Evaluation and Quality

Center for Devices and Radiological Health

Enclosure

Indications for Use

Submission Number (if known)

K231658

Device Name

VEOS Spinal Fixation System

Indications for Use (Describe)

The VEOS Spinal Fixation System is intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of acute and chronic instabilities or deformities of the thoracic, lumbar and sacral spine (from T1 to S1).

The VEOS Spinal Fixation System is intended for noncervical pedicle fixation for the following indications : degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies); spondylolisthesis; trauma (i.e., fracture or dislocation); spinal stenosis; curvatures (i.e., scoliosis, kyphosis, and/or lordosis); tumor; pseudoarthrosis; and failed previous fusion in skeletally mature patients.

When used in a posterior percutaneous approach with MIS instrumentation, the VEOS Spinal Fixation System is intended for noncervical pedicle fixation for the following indications: degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies); spondylolisthesis; trauma (i.e., fracture or dislocation); spinal stenosis; curvatures (i.e., scoliosis, kyphosis, and/or lordosis); tumor; pseudoarthrosis; and failed previous fusion in skeletally mature patients.

When used for posterior non-cervical pedicle screw fixation in pediatric patients, The VEOS Spinal Fixation System is indicated as an adjunct to fusion to treat adolescent idiopathic scoliosis. Pediatric pedicle screw fixation is limited to a posterior approach.

This system is intended to be used with autograft and/or allograft.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

Contact Details

[21 CFR 807.92\(a\)\(1\)](#)

Applicant Name	Spineway
Applicant Address	7 allée Moulin Berger Ecully 69130 France
Applicant Contact Telephone	+33614845362
Applicant Contact	Mr. Maurad Ben Mokhtar
Applicant Contact Email	rd.mbm@spineway.com
Correspondent Name	The OrthoMedix Group, Inc.
Correspondent Address	4313 West 3800 South West Haven UT 84401 United States
Correspondent Contact Telephone	512-590-5810
Correspondent Contact	Mr. J.D. Webb
Correspondent Contact Email	jdwebb@orthomedix.net

Device Name

[21 CFR 807.92\(a\)\(2\)](#)

Device Trade Name	VEOS Spinal Fixation System
Common Name	Thoracolumbosacral pedicle screw system
Classification Name	Thoracolumbosacral Pedicle Screw System
Regulation Number	888.3070
Product Code	NKB

Legally Marketed Predicate Devices

[21 CFR 807.92\(a\)\(3\)](#)

Predicate #	Predicate Trade Name (Primary Predicate is listed first)	Product Code
K161387	Mont Blanc and Mont Blanc MIS Spinal Systems	NKB
K180179	Firebird Spine System, including the JANUS® Fenestrated Screws	NKB
K954696	Smith Nephew - Rogozinski	KWP
K121728	Tiger Spine System	NKB

Device Description Summary

[21 CFR 807.92\(a\)\(4\)](#)

The Spineway VEOS Spinal Fixation System is composed of implant devices made from a titanium alloy Ti6Al4V-ELI per ISO 5832-3 and Cobalt-Chrome per ISO 5832-12. All implant components are provided sterile or non-sterile. The implants are to be implanted from the posterior approach. The screws are available as monoaxial, polyaxial, and polyaxial reduction screws in multiple diameters and length.

Titanium rods are available in two diameters and multiple lengths. CoCr rods are available in one diameter and two length. CoCr rods are only available non-sterile. Transverse connectors are available in various sizes to attach to the two parallel rods. Associated instrumentation to complete the procedure is provided.

Screws for MIS applications are available as polyaxial cannulated screws in multiple diameters and lengths. Rods for MIS applications are available as pre-bent rods and straight rods.

Associated instrumentation to complete the procedure is provided.

Intended Use/Indications for Use

[21 CFR 807.92\(a\)\(5\)](#)

The VEOS Spinal Fixation System is intended to provide immobilization and stabilization of spinal segments in skeletally mature patients as an adjunct to fusion in the treatment of acute and chronic instabilities or deformities of the thoracic, lumbar and sacral spine (from T1 to S1).

The VEOS Spinal Fixation System is intended for noncervical pedicle fixation for the following indications : degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies); spondylolisthesis; trauma (i.e., fracture or dislocation); spinal stenosis; curvatures (i.e., scoliosis, kyphosis, and/or lordosis); tumor; pseudoarthrosis; and failed previous fusion in skeletally mature patients.

When used in a posterior percutaneous approach with MIS instrumentation, the VEOS Spinal Fixation System is intended for noncervical pedicle fixation for the following indications: degenerative disc disease (defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies); spondylolisthesis; trauma (i.e., fracture or dislocation); spinal stenosis; curvatures (i.e., scoliosis, kyphosis, and/or lordosis); tumor; pseudoarthrosis; and failed previous fusion in skeletally mature patients.

When used for posterior non-cervical pedicle screw fixation in pediatric patients, The VEOS Spinal Fixation System is indicated as an adjunct to fusion to treat adolescent idiopathic scoliosis. Pediatric pedicle screw fixation is limited to a posterior approach.

This system is intended to be used with autograft and/or allograft.

Indications for Use Comparison

[21 CFR 807.92\(a\)\(5\)](#)

The VEOS Spinal Fixation System and the primary predicate have substantially equivalent indications for use.

Technological Comparison

[21 CFR 807.92\(a\)\(6\)](#)

Intended Use:

The VEOS Spinal Fixation System and the predicates are intended to help provide immobilization and stabilization of spinal segments as an adjunct to fusion of the thoracic, lumbar, and/or sacral spine.

Material:

The VEOS Spinal Fixation System and the predicate pedicle screws are manufactured from Ti6Al4V ELI. Rods are made from Ti6Al4V ELI or CoCr alloy.

Design/geometry:

The VEOS Spinal Fixation System and the predicates are top loading systems, polyaxial and monoaxial screws. They all include standard and reduction pedicle screws, straight, lordotic, and pre-bent rods. They also include transverse connectors. The VEOS and Janus screws have cannulations and fenestrations.

Dimensions:

VEOS sizes fit within range of predicate screws.

There are no differences between the VEOS Spinal Fixation System, the Mont Blanc / MIS Mont Blanc, and Firebird® Spinal Fixation System including JANUS Fenestrated Screws that would affect the substantial equivalence of the devices.

Non-Clinical and/or Clinical Tests Summary & Conclusions

[21 CFR 807.92\(b\)](#)

The structural competence of the VEOS Spinal Fixation System was evaluated by performing a series of tests complying with -ASTM F1717 "Standard Test Methods for Spinal Implant Constructs in a Vertebrectomy Model"

-ASTM F1798 " Standard Test Method for Evaluating the Static and Fatigue Properties of Interconnection Mechanisms and Subassemblies Used in Spinal Arthrodesis Implants."

Clinical testing was not performed.

Spineway considers the VEOS Spinal Fixation System to be substantially equivalent to the predicate devices listed above. This conclusion is based upon the devices' similarities in principles of operation, technology, materials, and indications for use.