



November 14, 2023

Saremco Dental AG
% Nevine Erian
Principal Regulatory Consultant
BQC Consulting, LLC
24341 Barbados Dr.
Dana Point, California 92629

Re: K232050

Trade/Device Name: saremco print BRIDGETEC
Regulation Number: 21 CFR 872.3690
Regulation Name: Tooth Shade Resin Material
Regulatory Class: Class II
Product Code: EBF, EBG, PZY
Dated: August 16, 2023
Received: August 17, 2023

Dear Nevine Erian:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device"

(<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



Bobak
Shirmohammadi -S

For Michael E. Adjodha, M. ChE., CQIA
Assistant Director
DHT1B: Division of Dental and
ENT Devices

OHT1: Office of Ophthalmic, Anesthesia,
Respiratory, ENT and Dental Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K232050

Device Name

saremco print BRIDGETEC

Indications for Use (Describe)

saremco print BRIDGETEC is a light-curing 3D-printed material intended as an indirect restorative for both anterior and posterior restorations, including occlusal surfaces. The saremco print BRIDGETEC material is used for fabricating permanent restorations such as inlays, onlays, veneers, full crown and 3-part bridge restorations. saremco print BRIDGETEC can also be used for the fabrication of artificial teeth and temporary crowns & bridges.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D) ☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

Submitter **Saremco Dental AG**
 Gewerbestrasse 4
 CH-9445 Rebstein, Switzerland
 Establishment Reg. No. N/A

Contact Dr. Christoph Evers
 Head of R&D
Phone +41 71 775 80 90
Fax +41 71 775 80 99
E-Mail christoph.evers@saremco.ch

Official Correspondent Nevine Erian
 Regulatory Consultant
BQC Consulting, LLC
 24341 Barbados Dr.
 Dana Point, CA 92629
Phone (949) 370-7155
Fax (949) 240-8829
E-mail nevine@bqcconsulting.com

Date Prepared August 15, 2023

Trade/Device Name	saremco print BRIDGETEC
Common Name	Crown, Bridge and Denture Teeth Material
Classification Name & Regulation Number	Tooth Shade Resin Material – 21 CFR 872.3690
Primary Product Code	EBF
Other Product Codes	EBG & PZY

Predicate Devices

CROWNTEC (Saremco Dental AG) – K213343 – **Primary Predicate**

TERA HARZ (Graphy Inc.) – K202846 – **Reference Device**

Device Description

saremco print BRIDGETEC is a resin-based material for the production of permanent and temporary restorations, veneers and artificial teeth on 3D printers. The material is an alternative to traditional heat-cured and auto polymerization resins. It requires a computer-aided and manufacturing (CAD/CAM) system including scanner, design software, additive printer and a post-curing unit.

Statement of Intended Use

saremco print BRIDGETEC is a methacrylate-based material for 3D printing of dental restorations.

Statement of Indication for Use

saremco print BRIDGETEC is a light-curing 3D-printed material intended as an indirect restorative for both anterior and posterior restorations, including occlusal surfaces. The saremco print BRIDGETEC material is used for fabricating permanent restorations such as inlays, onlays, veneers, full crown and 3-part bridge restorations. saremco print BRIDGETEC can also be used for the fabrication of artificial teeth and temporary crowns & bridges.

Material Composition

saremco print BRIDGETEC is composed of methacrylate-based resins, dental glass-filler, photo initiators and pigments.

Technological Characteristics

saremco print BRIDGETEC is a light-curing resin for 3D printing.

Non-Clinical Performance Testing

saremco print BRIDGETEC was tested and met the applicable requirements of the following ISO standard:

- ISO 4049:2019 – Dentistry – Polymer-based restorative materials
- ISO 10477:2020 – Dentistry – Polymer-based crown and veneering materials
- ISO 22112:2017 – Dentistry – Artificial teeth for dental prostheses

Bench test results allowed us to conclude that saremco print BRIDGETEC meets its intended use.

Biocompatibility

saremco print BRIDGETEC meets the biocompatibility requirements of the following standards:

- ISO 10993-1:2018 – Biological Evaluation of Medical Devices – Part 1: Evaluation and Testing Within a Risk Management Process
- ISO 7405:2018 Dentistry – Evaluation of Biocompatibility of Medical Devices Used in Dentistry

Clinical Performance Data

The performance of methacrylate-based polymer resins in the clinical environment has been well established. No human clinical testing was performed to support the substantial equivalence of saremco print BRIDGETEC.

Substantial Equivalence

The technical characteristics of saremco print BRIDGETEC is substantially equivalent to the predicate devices.

Material

saremco print BRIDGETEC is a resin-based material as the predicate devices.

Physical Properties

saremco print BRIDGETEC has similar physical properties as the predicate devices.

Comparison of saremco print BRIDGETEC to Predicate Devices

<i>Attribute</i>	saremco print BRIDGETEC	saremco print CROWNTEC	TERA HARZ
Indications			
Intended for both anterior and posterior restorations, including occlusal surfaces.	Yes	Yes	Yes
Used for fabricating permanent restorations such as inlays, onlays, veneers, full crown and 3-part bridge restorations.	Yes	Yes	Yes
Used for permanent 3-part bridge restorations	Yes	No	Yes
Used for the fabrication of artificial teeth	Yes	Yes	No
Used for the fabrication of temporary crowns & bridges.	Yes	Yes	Yes
Physical Properties			
Product State	Liquid	Liquid	Liquid
Material Type	Methacrylate based resin	Methacrylate based resin	Methacrylate based resin
Flexural Strength	Avg. ≥ 135 MPa	Avg. ≥ 135 MPa	Avg. 148.73 MPa
Performance Testing	ISO 4049 ISO 10477 ISO 22112	ISO 4049 ISO 10477 ISO 22112	ISO 4049 ISO 10477

Attribute	saremco print BRIDGETEC	saremco print CROWNTEC	TERA HARZ
Material Attributes			
Chemical Characterization	Methacrylate-based resins, dental glass-filler, photo initiators and pigments.	Methacrylate-based resins, dental glass-filler, photo initiators and pigments.	Polyurethane Resin, Methacrylate, Dimethacrylate, Phosphine oxide, Butylated hydroxytoluene and Pigments
Shelf Life	2 years	3 years	1 year
Storage Conditions	Do not expose to light	Do not expose to light	Do not expose to light
Storage Temperature	4-28°C (39-82°F)	4-28°C (39-82°F)	15-25°C (60-77°F)
Technical Attributes			
Product Code(s)	EBF, EBG & PZY	EBF, EBG & PZY	EBF & EBG
Fabrication Method	Additive 3D printing	Additive 3D printing	Additive 3D printing
UV Laser Wavelength	385 nm	385 or 405 nm	385 or 405 nm
Layer Thickness when Printing	50 µm	50 µm	50 or 100 µm
Polymerization Method	UV Curing Light	UV Curing Light	UV Curing Light
Sterile	No	No	No
Single Use	No	No	No

The differences in the physical properties and chemical compositions between saremco print BRIDGETEC and the predicate devices do not impact safety and effectiveness, as the finished clinical product is a custom-fitted biocompatible restoration regardless of the material variation.

Conclusion

Information provided in this application demonstrates that saremco print BRIDGETEC is substantially equivalent to the predicate devices. saremco print BRIDGETEC shares the same indications for use, similar material composition, similar physical properties and technological characteristics as the predicate devices.