



March 22, 2024

Medivators (A Subsidiary of STERIS Corporation)

Nick Wang

Senior Manager, Regulatory Affairs

3150 Pollok Drive

Conroe, Texas 77303

Re: K232067

Trade/Device Name: Defendo Fujifilm 500/600 Single Use Air/Water Valve, Defendo Fujifilm 500/600 Single Use Suction Valve, Defendo Fujifilm 500/600 Single Use Biopsy Valve, Endogator Fujifilm 500/600 Single Use Connector

Regulation Number: 21 CFR 876.1500

Regulation Name: Endoscope And Accessories

Regulatory Class: Class II

Product Code: ODC

Dated: February 22, 2024

Received: February 22, 2024

Dear Nick Wang:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Shanil P. Haugen -S

Shanil P. Haugen, Ph.D.

Assistant Director

DHT3A: Division of Renal, Gastrointestinal,
Obesity and Transplant Devices

OHT3: Office of Gastrorenal, ObGyn,
General Hospital, and Urology Devices

Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K232067

Device Name

Defendo Fujifilm 500/600 Single Use Air/Water Valve, Defendo Fujifilm 500/600 Single Use Suction Valve, Defendo Fujifilm 500/600 Single Use Biopsy Valve, Endogator Fujifilm 500/600 Single Use Connector

Indications for Use (Describe)

The DEFENDO™ FUJIFILM™ 500/600 Single Use Air/Water Valve is intended to be used to control the air/water function on an endoscope during a GI endoscopic procedure.

The DEFENDO™ FUJIFILM™ 500/600 Single Use Suction Valve is intended to be used to control the suction function on an endoscope during a GI endoscopic procedure.

The DEFENDO™ FUJIFILM™ 500/600 Single Use Biopsy Valve is intended for covering the endoscope biopsy port during an endoscopy procedure. The single use Biopsy Valve provides access for endoscopic device passage and exchange, help maintain sufflation, and minimize leakage of biomaterial from the biopsy port throughout the endoscopic procedure.

The ENDOGATOR™ FUJIFILM™ 500/600 Single Use Connector is intended to be used in conjunction with Endogator Irrigation Tubing to provide irrigation via sterile water during GI endoscopic procedures when used with an irrigation pump (or cautery unit).

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

Defendo Fujifilm 500/600 Single Use Valve Set

Manufacturer: Medivators
3150 Pollok Drive
Conroe, TX 77303

Submitter: Nick Wang, Ph.D., RAC
Senior Manager, Regulatory Affairs
Phone: 440-392-7482

Summary Date: 14 March 2024

**1. Device Name**

Trade Name: Defendo Fujifilm 500/600 Single Use Valves Set

Device Classification: Class II

Common/usual Name: Defendo Fujifilm 500/600 Single Use Air/Water Valve, Defendo Fujifilm 500/600 Single Use Suction Valve, Defendo Fujifilm 500/600 Single Use Biopsy Valve, Endogator Fujifilm 500/600 Single Use Connector

Classification Name: Endoscope and Accessories

Classification Number: 21 CFR 876.1500

Product Code: ODC

2. Predicate Device(s)

Primary predicate is listed first:

K102409 - Defendo Disposable Air/Water Valve Model 100304

K102581 - Defendo Disposable Suction Valve Model 100305

K090851 - Defendo Biopsy Valves

K220395 - Endogator Endoscopy Irrigation Tubing

3. Description of Device

The Defendo Fujifilm 500/600 Single Use Valves Set, the subject of this 510(k), are a family of ready-to-use, sterile, single use valves designed to be used on the Fujifilm 500/600 series GI endoscopes as accessories. The specific valves include the Defendo Fujifilm 500/600 Single Use Air/Water Valve, the Defendo Fujifilm 500/600 Single Use Suction Valve, the Defendo Fujifilm 500/600 Single Use Biopsy Valve, and the Endogator Fujifilm 500/600 Single Use Connector.

- The Air/Water valve allows the end-user to control and maintain air or CO₂ insufflation down the endoscope's accessory channel and control the water used to wash the lens of the endoscope during an endoscopic procedure.
- The Suction valve allows the user to control suction through the endoscope's accessory channel and serves as a fluid-removal conduit to aspirate fluids from the patient. This allows the clinician to remove fluids/soil from the GI tract during an endoscopic procedure.

STERIS Traditional 510(k) PREMARKET NOTIFICATION

Defendo Fujifilm 500/600 Single Use Valves Set

- The Biopsy valve allows instruments/devices to be passed through it while maintaining insufflation, allowing suction, and guiding instruments through the biopsy channel during endoscopic procedures.
- The Connector allows attachment between the irrigation tube set via its luer connection and the compatible endoscope at its dedicated auxiliary water connection. This is sometimes also referred to as the auxiliary water jet valve.

These valves are supplied in various kits. The kits are combinations of two or more of the different valves listed above.

4. Intended Use

The DEFENDO™ FUJIFILM™ 500/600 Single Use Air/Water Valve is intended to be used to control the air/water function on an endoscope during a GI endoscopic procedure.

The DEFENDO™ FUJIFILM™ 500/600 Single Use Suction Valve is intended to be used to control the suction function on an endoscope during a GI endoscopic procedure.

The DEFENDO™ FUJIFILM™ 500/600 Single Use Biopsy Valve is intended for covering the endoscope biopsy port during an endoscopy procedure. The single use Biopsy Valve provides access for endoscopic device passage and exchange, help maintain sufflation, and minimize leakage of biomaterial from the biopsy port throughout the endoscopic procedure.

The ENDOGATOR™ FUJIFILM™ 500/600 Single Use Connector is intended to be used in conjunction with Endogator Irrigation Tubing to provide irrigation via sterile water during GI endoscopic procedures when used with an irrigation pump (or cautery unit).

5. Technological Characteristic Comparison Table

Table 1. Proposed/Predicate Device Technological Characteristics Comparison Table – Air / Water Valve

Feature	Defendo Fujifilm 500/600 Single Use Air/Water Valve (Subject Device)	Defendo Disposable Air/Water Valve (Predicate K102409)	Comparison
Intended use	The DEFENDO™ FUJIFILM™ 500/600 Single Use Air/Water Valve is intended to be used to control the	The Defendo Disposable Air/Water Valve is intended to be used to control the air/water function	Same

STERIS Traditional 510(k) PREMARKET NOTIFICATION **Defendo Fujifilm 500/600 Single Use Valves Set**

Feature	Defendo Fujifilm 500/600 Single Use Air/Water Valve (Subject Device)	Defendo Disposable Air/Water Valve (Predicate K102409)	Comparison
	air/water function on an endoscope during a GI endoscopic procedure.	on an endoscope during a GI endoscopic procedure.	
Sterile/Non-sterile	Sterile	Sterile	Same
Sterilization Method	EtO	EtO	Same
Sterilization Assurance Level	10 ⁻⁶	10 ⁻⁶	Same
Usage	Single-Use	Single-Use	Same
Main Device Dimensions (lengths/widths)	Air/Water Valve: Length: ~32.0 mm Diameter: ~16.5 mm	Air/Water Valve: Length: ~46 mm Diameter: ~17 mm	Similar
Target Population	Patients undergoing an endoscopic procedure	Patients undergoing an endoscopic procedure	Same
Energy Used/Delivered	None	None	Same
Method of Application	Manual actuation	Manual actuation	Same
Compatible Endoscopes	Fujifilm 500/600 series GI endoscope	Olympus endoscope	Reason for this 510(k)
Packaging	Supplied as a kit in PETG Tray with Tyvek Lid	Supplied individually in Tyvek Pouch	Package systems all maintain sterile barrier

Table 2. Proposed/Predicate Device Technological Characteristics Comparison Table – Suction Valve

Feature	Defendo Fujifilm 500/600 Single Use Suction Valve (Subject Device)	Defendo Disposable Suction Valve (Predicate K102581)	Comparison
Intended Use	The DEFENDO™ FUJIFILM™ 500/600 Single Use Suction Valve is intended to be used to control	The Defendo Disposable Suction Valve is intended to be used to control the suction function on	Same

STERIS Traditional 510(k) PREMARKET NOTIFICATION **Defendo Fujifilm 500/600 Single Use Valves Set**

Feature	Defendo Fujifilm 500/600 Single Use Suction Valve (Subject Device)	Defendo Disposable Suction Valve (Predicate K102581)	Comparison
	the suction function on an endoscope during a GI endoscopic procedure.	an endoscope during a GI endoscopic procedure.	
Sterile/Non-sterile	Sterile	Sterile	Same
Sterilization Method	EtO	EtO	Same
Sterilization Assurance Level	10 ⁻⁶	10 ⁻⁶	Same
Usage	Single Use	Single Use	Same
Main Device Dimensions (lengths/widths)	Suction Valve: Length: ~29.1 mm Diameter: ~17.1 mm	Suction Valve: Length: ~27 mm Diameter: ~18 mm	Similar, the slight differences are due to the differences in the endoscope port.
Target Population	Patients undergoing an endoscopic procedure	Patients undergoing an endoscopic procedure	Same
Energy Used/Delivered	None	None	Same
Method of Application	Manual actuation	Manual actuation	Same
Compatible Endoscopes	Fujifilm 500/600 series GI endoscope	Olympus endoscope	Reason for this 510(k)
Packaging	Supplied as a kit in PETG Tray with Tyvek Lid	Supplied individually in Tyvek Pouch	Package systems all maintain sterile barrier

STERIS Traditional 510(k) PREMARKET NOTIFICATION **Defendo Fujifilm 500/600 Single Use Valves Set**

Table 5-3. Proposed/Predicate Device Technological Characteristics
Comparison Table – Biopsy Valve

Feature	Defendo Fujifilm 500/600 Single Use Biopsy Valve (Subject Device)	Defendo Biopsy Valve (for Olympus and Fujinon Endoscopes) (Predicate K090851)	Comparison
Intended Use	The DEFENDO™ FUJIFILM™ 500/600 Single Use Biopsy Valve is intended for covering the endoscope biopsy port during an endoscopy procedure. The single use Biopsy Valve provides access for endoscopic device passage and exchange, help maintain sufflation, and minimize leakage of biomaterial from the biopsy port throughout the endoscopic procedure.	The Defendo Biopsy Valve is intended for covering the endoscope biopsy port during an endoscopy procedure. The valve provides access for endoscopic device passage and exchange, help maintain sufflation, and minimize leakage of biomaterial from the biopsy port throughout the endoscopic procedure.	Same
Sterile/Non-sterile	Sterile	Sterile	Same
Sterilization Method	EtO	EtO	Same
Sterilization Assurance Level	10 ⁻⁶	10 ⁻⁶	Same
Usage	Single Use	Single Use	Same
Main Device Dimensions (lengths/widths)	Length: ~17.8 mm Diameter: ~15.8 mm	Length: ~17 mm Diameter: ~17 mm	Similar, minor differences did not affect performance (confirmed through testing)
Target Population	Patients undergoing an endoscopic procedure	Patients undergoing an endoscopic procedure	Same
Energy Used/Delivered	None	None	Same

STERIS Traditional 510(k) PREMARKET NOTIFICATION
Defendo Fujifilm 500/600 Single Use Valves Set

Feature	Defendo Fujifilm 500/600 Single Use Biopsy Valve (Subject Device)	Defendo Biopsy Valve (for Olympus and Fujinon Endoscopes) (Predicate K090851)	Comparison
Compatible Endoscopes	Olympus and Fujifilm endoscope	Olympus and Fujinon endoscopes	Same
Packaging	Continue to be supplied individually in Tyvek Pouch (predicate) and supplied as a set in PETG Tray with Tyvek Lid (new).	Supplied individually in Tyvek Pouch	Package systems all maintain sterile barrier

Table 5-4. Proposed/Predicate Device Technological Characteristics
Comparison Table – Connector

Feature	Endogator Fujifilm 500/600 Single Use Connector (Subject Device)	Endogator Connector (K092429)	Comparison
Intended Use	The ENDOGATOR™ FUJIFILM™ 500/600 Single Use Connector is intended to be used in conjunction with Endogator Irrigation Tubing to provide irrigation via sterile water during GI endoscopic procedures when used with an irrigation pump (or cautery unit).	The Endogator system is intended to provide irrigation via sterile water supply during GI endoscopic procedure when sued in conjunction with an irrigation pump (or cautery unit).	Similar (Language modified for the subject device to reflect the connector provided does not include the Endogator tubing)
Sterile/Non-sterile	Sterile	Sterile	Same
Sterilization Method	EtO	EtO	Same
Sterilization Assurance Level	10 ⁻⁶	10 ⁻⁶	Same
Usage	Single Use	Single Use	Same

STERIS Traditional 510(k) PREMARKET NOTIFICATION

Defendo Fujifilm 500/600 Single Use Valves Set

Feature	Endogator Fujifilm 500/600 Single Use Connector (Subject Device)	Endogator Connector (K092429)	Comparison
Main Device Dimensions (lengths/widths)	Length: ~22 mm Diameter: ~12.7 mm	Length: ~22 mm Diameter: ~12.7 mm	Same
Target Population	Patients undergoing an endoscopic procedure	Patients undergoing an endoscopic procedure	Same
Energy Used/Delivered	None	None	Same
Compatible Endoscopes	Fujifilm 500/600 series GI endoscope	None specified (commonly used with Fujifilm endoscopes)	Same
Packaging	Continue to be supplied individually in Tyvek Pouch (predicate) and supplied as a set in PETG Tray with Tyvek Lid (new).	Supplied individually in Tyvek Pouch	Package systems all maintain sterile barrier

4. Non-Clinical Performance Testing

Non-clinical performance testing was performed on the subject device to support the substantial equivalence determination. The following tests were performed:

Air / Water Valve

- Gas (Air and CO₂) flow rate testing
- Water flow rate testing
- Backflow prevention testing
- Mechanical testing
 - Depression Force
 - Valve Application Force
 - Valve Removal Force
 - Force to Disassemble via Axial Pull
- Button Cycling

Suction Valve

- Suction Bypass
- Suction Rate (water and soil)

STERIS Traditional 510(k) PREMARKET NOTIFICATION

Defendo Fujifilm 500/600 Single Use Valves Set

- Mechanical testing
 - Depression Force
 - Valve Application Force
 - Valve Removal Force
 - Force to Disassemble via Axial Pull
- Valve Operation / Procedure Duration Test

Biopsy Valve

- Leak Testing
- Insufflation Testing

Auxiliary Waterjet Valve (Connector)

- Flow Durability Testing
- Torque Test

In addition to performance testing, ISO 11607/ASTM 4169 testing was completed to ensure the packaging of the subject device kit maintained sterile barrier after shipping/transit.

5. Biocompatibility

The biocompatibility of the subject devices was assessed in accordance with the FDA guideline “Use of International Standard ISO 10993-1, “Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management process”. Biocompatibility testing conducted on the patient-contacting portions of the subject device in accordance with the ISO 10993 standard. The test result shows that the subject devices are biocompatible.

6. Conclusion

Based on the intended use, technological characteristics, non-clinical performance testing, and biocompatibility assessment, the subject device has shown to be substantially equivalent to the predicate and having met the acceptance criteria based on its indications for use.