



November 20, 2024

Shenzhen Ruiyi Business Technology Co., Ltd.
% Libray Zhang
Official Correspondent
Shanghai Spica Management Consulting Co., Ltd.
609 Room, No.133 Shengang Avenue, Pudong New District
Shanghai, 201306
China

Re: K232965

Trade/Device Name: Leg Massager (Models: RP-ALM070H, RP-ALM071H)
Regulation Number: 21 CFR 890.5650
Regulation Name: Powered Inflatable Tube Massager
Regulatory Class: Class II
Product Code: IRP, IRT
Dated: September 30, 2024
Received: September 30, 2024

Dear Libray Zhang:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

All medical devices, including Class I and unclassified devices and combination product device constituent parts are required to be in compliance with the final Unique Device Identification System rule ("UDI Rule"). The UDI Rule requires, among other things, that a device bear a unique device identifier (UDI) on its label and package (21 CFR 801.20(a)) unless an exception or alternative applies (21 CFR 801.20(b)) and that the dates on the device label be formatted in accordance with 21 CFR 801.18. The UDI Rule (21 CFR 830.300(a) and 830.320(b)) also requires that certain information be submitted to the Global Unique Device Identification Database (GUDID) (21 CFR Part 830 Subpart E). For additional information on these requirements, please see the UDI System webpage at <https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/unique-device-identification-system-udi-system>.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory->

[assistance/contact-us-division-industry-and-consumer-education-dice](#)) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Heather L. Dean -S

Heather Dean, PhD
Assistant Director, Acute Injury Devices Team
DHT5B: Division of Neuromodulation and
Physical Medicine Devices
OHT5: Office of Neurological and
Physical Medicine Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K232965

Device Name

Leg Massager (Models:RP-ALM070H, RP-ALM071H)

Indications for Use (Describe)

Leg Massager (Models:RP-ALM070H, RP-ALM071H) is intended for home to temporarily relieve minor muscle aches and/or pains, and to temporarily increase circulation to the treated areas in people who are in good health.

Type of Use (Select one or both, as applicable)

☐ Prescription Use (Part 21 CFR 801 Subpart D)

☒ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

K232965

Type of Submission

Traditional

Date Prepared

Jun 27, 2023

Submission Sponsor

Manufacturer Name

Shenzhen Ruiyi Business Technology Co., Ltd.

Address

Qianhai Complex A201, Qianwan Road 1, Qianhai

Shenzhen-Hong Kong Cooperation Zone,

Shenzhen, China

Tel

86-18566621145

Email

derric@renpho.com

Contact Person

derric bai

Application Correspondent

Company Name

Shanghai Spica Management Consulting Co., Ltd.

Address

609 Room, No.133 Shengang Avenue, Pudong New

District, Shanghai, China

Tel

86-13020102321

Email

Libray@spicagloble.com

Contact Person

Libray Zhang

Device Identification

Trade Name

Leg Massager (Models:RP-ALM070H, RP-ALM071H)

Regulation Number

21 CFR 890.5650

Classification Name

Massager, Powered Inflatable Tube

Device Classification

Class II

Panel

Physical Medicine

Product Code

IRP

Previous Submissions

None

Indications for Use

Leg Massager (Models:RP-ALM070H, RP-ALM071H) is intended for home to temporarily relieve minor muscle aches and/or pains, and to temporarily increase circulation to the treated areas in people who are in good health.

Device Description

510(k) Summary

Leg Massager (Models:RP-ALM070H, RP-ALM071H) is consist of air pressure sensor, air pump, sleeves etc working together as one unit. The air pump is connected to the dedicated sleeves via a series of hoses. The compression massage direction is from limb end to body center by inflating the air chambers sequentially and then deflating as one cycle, the pressure can be adjusted to avoid any discomfort to the patient. The sleeve works under the action of sensor and microprocessor.

Leg Massager (Models:RP-ALM070H, RP-ALM071H), in medical market, it has a sequential squeezing from distal to proximal, thus help to improve circulation to the treated areas.

Predicate Device Information

Sponsor	Shenzhen Dongjilian Electronics Co., Ltd.
Trade/Device Name	Air Compression Therapy Device
510(K) number	K193354
Regulation Number	21 CFR 890.5650

Sponsor	Medella Health Limited
Trade/Device Name	FLOWpresso
510(K) number	K223729
Regulation Number	21 CFR 890.5650

510(k) Summary

Table 6A: Summary of Comparison

	Subject Device	Predicate Device	Differences Discussion
Device name	Leg Massager (Models:RP-ALM070H, RP-ALM071H)	Air Compression Therapy Device	N/A
510(k) number	K232965	K193354	N/A
Manufacturer	Shenzhen Ruiyi Business Technology Co., Ltd.	Shenzhen Dongjilian Electronics Co., Ltd.	N/A
Product regulation	21 CFR 890.5650	21 CFR 890.5650	Same
Classification name	Massager, Powered Inflatable Tube	Massager, Powered Inflatable Tube	Same
Regulation class	2	2	Same
Product code	IRP	IRP	Same
Indications for use	Leg Massager (Models:RP-ALM070H, RP-ALM071H) is intended for home to temporarily relieve minor muscle aches and/or pains, and to temporarily increase circulation to the treated areas in people who are in good health.	The Air Compression Therapy Device is indicated for the temporary relief of minor muscle aches and pains and for temporary increase in circulation to the treated areas in people who are in good health. The Air Compression Therapy Device simulates kneading and stroking of tissues by using an inflatable garment.	Same


510(k) Summary

Rx or OTC	OTC	OTC	Same
Pressure range	0~240mmHg	0~240mmHg	Same
Inflation Time	3-30s	3-30s	Same
Deflation Time	1-5s	1-5s	Same
Treatment time	20 minutes	20 minutes	Same
Standard	ES 60601-1; IEC60601-1-2; ISO 10993-5; ISO 10993-10; IEC 60601-1-11	ES 60601-1; IEC60601-1-2; ISO 10993-5; ISO 10993-10; IEC 60601-1-11	Same
Mode of compression	Sequential	Sequential/ Peristaltic	Same
Power source	100~240V 50/60Hz	100~240V 50/60Hz	Same
Power consumption	24W	12W	Similar
Size and appearance of sleeves (leg part)	Thighs: One size: 11*16*24.9inch+8*33.5inch	Leg: One size: 73*26cm	


510(k) Summary

Photo	RP-ALM070H:  RP-ALM071H: 		Similar
Housing materials	Molded ABS enclosure	Molded ABS enclosure	Same
Number of chambers	3	3	Same
Work mode	RP-ALM070H: Six models (3 combine massage modes and 3 Separate Massage Modes) C1: Massage full legs. C2: Massage feet and calves. (It can be used individually, and don't need	Mode 1: Starting with the foot chamber and progressing up the thigh chamber, each section compresses and the pressure gradually rises to the	Although the subject device provides 6 kinds of work mode, the Mode F1, F2, F3, M1, M2, M3 are the similar with predicate device (K193354), while the other work modes of subject device just have difference about inflatable order of the different chambers. The treatment pressure range are the

510(k) Summary

	to connect the air hose of thighs wraps.) C3: Massage feet, calves, thighs by turn. It will be turn off When you press the Button C again. T: Massages thighs C: Massages calves F : Massages feet. It will be turn off when you press the Button S again. RP-ALM071H: Five models(3 combine massage modes and 2 Separate Massage Modes): Combine 1: Massage feet and calves Combine 2: Massage from calves to feet Combine 3: Massage from feet to calves. Press it again to turn off combine message function. F: Massage feet. C: Massage calves. Press it again to turn off separate massage function.	pre-determined air pressure level, then decompresses and the air pressure drops. Once the thigh section decompresses, the cycle begins again. Mode 1 follows this pressure sequence:  Mode 2: Starting with the foot chamber and progressing up the thigh, each section compresses and the pressure gradually rises to the pre-determined air pressure level, holds the air until the entire garment is compressed. All three sections then decompress simultaneously and the air pressure drops, then cycle begins again. Mode 2 follows this pressure	same under different work modes, so the difference of pressure range would not raise adversely impact on safety and effectiveness.
--	---	--	---

510(k) Summary

		sequence:  Mode 3: include two stage, stage 1: it work according to the method of mode 1, after the stage 1 is completed, it go to stage 2(working according to the method of mode 2) without interruption time until finish the stage 2, then enter next cycle without interruption. Mode1 $\longleftrightarrow$ Mode2 The pressure sequence of mode 3 combines mode 1 and mode 2	
Safety feature	Button on display allows user to stop or pause therapy session at any time	Button on display allows user to stop or pause therapy session at any time	Same
Technology	Compressor and valve system which sequentially inflates inflatable chambers	Compressor and valve system which sequentially inflates inflatable chambers	Same

510(k) Summary

Operating environment	Temperature: 5°C- 40°C, Humidity: 15%- 90%	Temperature: 5°C - 40°C, Humidity:5%- 90% non-condensing	Same
Transportation and Storage environment	Temperature: - 25°C~70°C Humidity:15%-90% non condensing Atmospheric Pressure:75kPa-106kPa	Temperature: - 20°C~55°C; Humidity:5%-90% non condensing Atmospheric Pressure:75kPa-106kPa	Similar

Table 6B: Summary of Comparison

	Subject Device	Predicate Device	Differences Discussion
Device name	Leg Massager (Models:RP-ALM070H, RP-ALM071H)	FLOWpresso	N/A
510(k) number	K232965	K223729	N/A
Manufacturer	Shenzhen Ruiyi Business Technology Co., Ltd.	Medella Health Ltd	N/A
Product regulation	21 CFR 890.5650	21 CFR 890.5650	Same
Classification name	Massager, Powered Inflatable Tube	Massager, Powered Inflatable Tube	Same
Regulation class	2	2	Same
Product code	IRP, IRT	IRP, IRT	Same
Indications for Use	Leg Massager (Models:RP-ALM070H, RP-ALM071H) is intended for home to temporarily relieve minor muscle aches	The Flowpresso Device combines heat and compression therapy. Flowpresso is intended to treat post	Similar, although the predicate device is intended to treat post traumatic and post-surgical medical and/or surgical

510(k) Summary

	and/or pains, and to temporarily increase circulation to the treated areas in people who are in good health.	traumatic and post-surgical medical and/or surgical conditions for which localized thermal therapy are indicated. Flowpresso is intended to be used by, or on the order of, licensed health care professionals based in rehabilitation facilities, outpatient clinics and athletic training settings.	conditions for which localized thermal therapy are indicated, the subject device is ,the working principles and expected outcomes of the two are consistent. All of them relieve pain by combining heat therapy with air compression on the affected area (subject device only covers the thighs, calves, and feet, predictive device also includes the back/hip), allowing patients to relax.
Number of patients that can be treated at one time	One	One	Same
Treatment Time	20 minutes	01-40 minutes	Different, Although the “Treatment Time” of the subject device is different than the predicate devices, the treatment time of the subject device is included in the predicate device.
Heat Therapy	Not exceeding 45°C	Default: 86°F – 104°F	Similar, the maximum temperature of the subject device is slightly higher than that of the predicate device. The heating program of the subject device can be selected to be turned on or not, and all have passed IEC testing. Small differences do not affect safety and effectiveness.

510(k) Summary

Compression Pressure Levels	0~240mmHg	5-130mm/Hg	Different, pressure range is dependent on patient comfort levels and are adjusted in real-time.
Static or Intermittent Pressure	Both	Both	Same
Power Down	Available	Available	Same
Heating Mechanism	Thermotherapy	Thermotherapy	Same
User Interface	Touch Screen	Touch Screen	Same
Line voltage	100~240V	110-220V	Similar
Line Frequency	50/60Hz	50/60 Hz	Same
Electrical Safety Standards	ES 60601-1; IEC60601-1-2; IEC 60601-1-11	ANSI/AAMI ES60601-1:2005 + A1:2012 IEC 60601-1:2005 + A1:2012 EN60601-1:2006 +A1:2013	Similar
Operating Temperature	5°C- 40°C	60°F – 80°F (16°C – 27°C)	Different, Although the “Operating Temperature” of the subject device is different than the predicate devices, they comply with ANSI/AAMI ES60601-1, so the difference does not affect the safety and effectiveness.
Storage Temperature	- 25°C~70°C	33°F – 122°F (1°C - 50°C)	Different, Although the “Storage Temperature” of the subject device is different than the predicate devices, they comply with ANSI/AAMI ES60601-1, so

510(k) Summary

			the difference does not affect the safety and effectiveness.
Operating Humidity	15%- 90%	Below 60% Noncondensing	Different, Although the “Operating Humidity” of the subject device is different than the predicate devices, they comply with ANSI/AAMI ES60601-1, so the difference does not affect the safety and effectiveness.
Storage Humidity	15%-90%	Below 60% Noncondensing	Different, Although the “Storage Humidity” of the subject device is different than the predicate devices, they comply with ANSI/AAMI ES60601-1, so the difference does not affect the safety and effectiveness.
Operating Atmospheric Pressure and Altitude	75kPa-106kPa	700 hPa – 1060 hPa (corresponds to a max. elevation of 9,842 ft. 6 in (3000m))	Similar
Types of Garments	Sleeves for thighs, calves, and feet	Various anatomical thermal/compression garments for: Back/Hip, arm, leg, feet	Different, the type of garments of the predicate device exceeds that of the subject device, which does not include the back/ip part. So the risk of subject device is smaller than that of predicate device.
Patient Contacting Material	All encompassed with a Nylon with a Polyurethane laminate material.	Needle cotton, non-conductive and non-woven lining. All encompassed with a PVC material.	Different, Although the material of the subject device is different than the predicate devices, they comply with ISO10993-5 and ISO10993-10, so the difference does not

510(k) Summary

			affect the safety and effectiveness.
Multi-Patient Use and Single Patient Use Wraps	Multi-Patient Use	Multi-Patient Use	Same
Biocompatibility	Cytotoxicity testing per ISO 10993-5 Sensitization testing per ISO 10993-10 Irritation testing per ISO 10993-10	Cytotoxicity testing per ISO 10993-5 Sensitization testing per ISO 10993-10 Irritation testing per ISO 10993-10	Same
Sterile/Non Sterile	Non-sterile only	Non-sterile only	Same
Cleaning Disinfection Validation of Labeling	Yes- for Multi-Patient use reusable wraps	Yes- for Multi-Patient use reusable wraps	Same
Human Factors testing to confirm intended users have found instructions for cleaning and disinfection easy to use	Yes- for Multi-Patient use reusable wraps	Yes- for Multi-Patient use reusable wraps	Same
Expected Life of	Based on frequency of use and	Based on frequency of use and	Same

510(k) Summary

Garments	continued functional performance	continued functional performance	
Validation of repeated cleaning and disinfection for reusable garments	Yes- for Multi-Patient use reusable wraps	Yes- for Multi-Patient use reusable wraps	Same

510(k) Summary

Summary of the technological characteristics of the device

The device meets all the applicable technical requirements of :

IEC 60601-1-11: 2015 - Medical electrical equipment - Part 1-11: General Requirements for Basic Safety and Essential Performance - Collateral Standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment

IEC 60601-1-2: 2014 - Medical electrical equipment - Part 1-2: General Requirements for Basic Safety and Essential Performance - Collateral Standard: Electromagnetic Compatibility

ANSI AAMI ES60601-1:2005/(R)2012 and A1:2012, C1:2009/(R)2012 and A2:2010/(R)2012 - Medical electrical equipment - Part 1: General requirements for basic safety and essential performance

IEC TS 60601-4-2: Medical electrical equipment - Part 4-2: Guidance and interpretation - Electromagnetic immunity: performance of medical electrical equipment and medical electrical systems

ISO 10993-5: 2009 - Biological Evaluation of Medical Device - Part 5: Tests for in vitro Cytotoxicity

ISO 10993-10: 2010 - Biological Evaluation of Medical Devices - Part 10: Tests for irritation and skin sensitization

Conclusion

Based on the indications for use, technological characteristics, and non-clinical performance data, “Leg Massager (Models:RP-ALM070H, RP-ALM071H)” is as safe, as effective, and performs as well as the legally marketed predicate devices. Therefore, the subject device is substantially equivalent to the predicate device.