



May 17, 2024

Smedtrum Medical Technology Co., Ltd.
Crimson Wu
Senior Regulatory Engineer
1F., No. 8, Ln. 97, Wugong Rd., Xinzhuang Dist.
New Taipei City, 248016
Taiwan

Re: K240482

Trade/Device Name: Intense Pulsed Light System (ST-690)
Regulation Number: 21 CFR 878.4810
Regulation Name: Laser Surgical Instrument For Use In General And Plastic Surgery And In Dermatology
Regulatory Class: Class II
Product Code: ONF
Dated: February 16, 2024
Received: February 20, 2024

Dear Crimson Wu:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>).

Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Tanisha L. Hithe -S
Digitally signed by
Tanisha L. Hithe -S
Date: 2024.05.17
15:10:25 -04'00'

Tanisha Hithe
Assistant Director
DHT4A: Division of General Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K240482

Device Name

Intense Pulsed Light System (ST-690)

Indications for Use (Describe)

The Intense Pulsed Light System is intended for medical use in the treatment of the following conditions:

1. Moderate inflammatory acne vulgaris;
2. Benign pigmented epidermal lesions including dyschromia, hyperpigmentation, melasma, ephelides (freckles);
3. Benign cutaneous vascular lesions including port wine stains, hemangiomas, facial truncal and leg telangiectasias, erythema of rosacea, leg veins, spider angiomas and venous malformations;
4. Permanent hair reduction-long-term stable reduction in number of hairs re-growing after a treatment regimen.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."



Smedtrum Medical Technology Co., Ltd.
1F., No. 8, Ln. 97, Wugong Rd., Xinzhuang Dist.,
New Taipei City 248016, Taiwan
TEL : +886 (02) 2298 9578 FAX : +886 (02) 2298 9426

510(k) Summary: K240482

Intense Pulsed Light System

I. SUBMITTER

Smedtrum Medical Technology Co., Ltd.
1F., No. 8, Ln. 97, Wugong Rd.,
Xinzhuang Dist., New Taipei City 248016,
Taiwan (R.O.C.)

Contact Person

Crimson Wu

Position: Senior Regulatory Engineer

Tel: +886-2-2298-9578, Ext. 301

Fax: +886-2-2298-9426

E-mail: crimsonwu@smedtrum.com

Date of preparation: May 17th, 2024

II. DEVICE

Trade Name:	Intense Pulsed Light System
Common or Usual Name:	Intense Pulsed Light System
Product code:	ONF
Classification Name:	Powered Light Based Non-Light Surgical Instrument With Thermal Effect 21 C.F.R. § 878.4810, Device Class II

III. PREDICATE DEVICE

Manufacturer	Smedtrum Medical Technology Co., Ltd.
Trade Name:	Intense Pulsed Light System
Common or Usual Name:	Laser surgical instrument for use in general and plastic surgery and in dermatology
Classification Name:	Powered Light Based Non-Light Surgical Instrument With Thermal Effect 21 C.F.R. § 878.4810, Device Class II
Premarket Notification:	K231394 Aug 9 th , 2023

IV. DEVICE DESCRIPTION

Intense Pulsed light (IPL) System is a type of intensive, broadband, coherent light source which has a wavelength spectrum of 420 nm -1200 nm. There are six optical filters that can be using in this system listed in table. With these special properties, the IPL System has a wide application in non-ablative therapies based on theory of human



skin tissue's selective absorption.

Spectra	Application areas	Expected Working Life
420 -1200nm	Acne	100,000 shots
510 -1200nm	Acne, vascular, pigment	80,000 shots
560 -1200nm	Vascular, pigmentation	70,000 shots
610 -1200nm	Hair removal	60,000 shots
640 -1200nm	Hair removal	50,000 shots
610 – 980nm (SHR)	Hair removal	30,000 shots

V. INDICATIONS FOR USE

The Intense Pulsed Light System is intended for medical use in the treatment of the following conditions:

1. Moderate inflammatory acne vulgaris;
2. Benign pigmented epidermal lesions including dyschromia, hyperpigmentation, melasma, ephelides (freckles);
3. Benign cutaneous vascular lesions including port wine stains, hemangiomas, facial truncal and leg telangiectasias, erythema of rosacea, leg veins, spider angiomas and venous malformations;
4. Permanent hair reduction-long-term stable reduction in number of hairs re-growing after a treatment regimen.

COMPARISON OF TECHNOLOGICAL CHARACTERISTICS WITH THE PREDICATE DEVICE

Feature	Proposed device	Predicate device (K231394)
<i>Device Name</i>	Intense Pulsed Light System	Intense Pulsed Light System
<i>Product Code</i>	ONF	ONF
<i>Regulation No.</i>	21 CFR 878.4810	21 CFR 878.4810
<i>Device Class</i>	Class II	Class II

Feature	Proposed device	Predicate device (K231394)
<i>Indication for Use</i>	The Intense Pulsed Light System is intended for medical use in the treatment of the following conditions: 1. Moderate inflammatory acne vulgaris; 2. Benign pigmented epidermal lesions including dyschromia, hyperpigmentation, melasma, ephelides (freckles); 3. Benign cutaneous vascular lesions including port wine stains, hemangiomas, facial truncal and leg telangiectasias, erythema of rosacea, leg veins, spider angiomas and venous malformations; 4. Permanent hair reduction-long-term stable reduction in number of hairs re-growing after a treatment regimen.	The Intense Pulsed Light System is intended for medical use in the treatment of the following conditions: 1. Moderate inflammatory acne vulgaris; 2. Benign pigmented epidermal lesions including dyschromia, hyperpigmentation, melasma, ephelides (freckles); 3. Benign cutaneous vascular lesions including port wine stains, hemangiomas, facial truncal and leg telangiectasias, erythema of rosacea, leg veins, spider angiomas and venous malformations; 4. Permanent hair reduction-long-term stable reduction in number of hairs re-growing after a treatment regimen.
<i>Light Source</i>	Intense Pulsed light (Xenon Flash Lamp)	Intense Pulsed light (Xenon Flash Lamp)
<i>Wavelength</i>	420 – 1200 nm	420 – 1200 nm
<i>Light Delivery System</i>	Handpiece	Handpiece
<i>Filter</i>	420 -1200nm: Acne; 510 -1200nm: Acne, vascular, pigment; 560 -1200nm: Vascular, pigment; 610-1200nm: Hair removal; 640-1200nm: Hair removal; 610-980nm: Hair removal;	420 -1200nm: Acne; 510 -1200nm: Acne, vascular, pigment; 560 -1200nm: Vascular, pigment; 610-1200nm: Hair removal; 640-1200nm: Hair removal; 610-980nm: Hair removal;
<i>Fluence</i>	420 -1200nm: 4.1-34.8 J/cm ² ; 510 -1200nm: 3.8-31.6 J/cm ² ; 560 -1200nm: 3.9-30.2 J/cm ² ; 610-1200nm: 3.9-27.8 J/cm ² ; 640-1200nm: 3.3-24.9 J/cm ² ; 610-980nm: 3.2-23.3 J/cm ² ;	420 -1200nm: 4.1-34.8 J/cm ² ; 510 -1200nm: 3.8-31.6 J/cm ² ; 560 -1200nm: 3.9-30.2 J/cm ² ; 610-1200nm: 3.9-27.8 J/cm ² ; 640-1200nm: 3.3-24.9 J/cm ² ; 610-980nm: 3.2-23.3 J/cm ² ;
<i>Pulsed Energy density</i>	4.1~35 J/cm ² (10~50 Level)	4.1~35 J/cm ² (10~50 Level)
<i>Pulsed width</i>	5~20 ms	5~20 ms

<i>Feature</i>	Proposed device	Predicate device (K231394)
<i>Pulsed duration</i>	5~50 ms	5~50 ms
<i>Spot Size</i>	12 mm × 35 mm & 15 mm × 50 mm	12 mm × 35 mm & 15 mm × 50 mm
<i>Cooling</i>	water + air + TEC	water + air + TEC
<i>Appearance</i>	Height 451 mm	Height 1190 mm
	Width 468 mm	Width 590 mm
	Depth 630 mm	Depth 550 mm
	Weight Approx. 45 kg	Weight Approx. 68 kg

VI. PERFORMANCE DATA

The Intense Pulsed Light System has been determined through engineering testing to verify optical energy output and electrical safety.

Electrical safety and electromagnetic compatibility

The test results demonstrated that the proposed device complies with the following standards:

IEC 60601-1: 2005+corr.1:2006+Corr.2:2007+A1:2012 Medical electrical equipment
- Part 1: General requirements for basic safety and essential performance

IEC 60601-1-2:2020 Medical electrical equipment -Part 1-2: General requirements for basic safety and essential performance-Collateral Standard: Electromagnetic disturbances - Requirements and tests

IEC 62471 First edition 2006-07 Photobiological safety of lamps and lamp systems

IEC 60601-2-57 Edition 1.0 2011-01 Medical Electrical Equipment - Part 2-57: Particular requirements for the basic safety and essential performance of non-laser light source equipment intended for therapeutic, diagnostic, monitoring and cosmetic/aesthetic use

Bench testing

Intense Pulsed Light system is an Light source equipment. The spatial variation of the LS equipment output over the treatment area shall not deviate from the average irradiance or radiant exposure by more than $\pm 20\%$

The product fulfills the requirements of IEC 60601-2-57.



Smedtrum Medical Technology Co., Ltd.
1F., No. 8, Ln. 97, Wugong Rd., Xinzhuang Dist.,
New Taipei City 248016 , Taiwan
TEL : +886 (02) 2298 9578 FAX : +886 (02) 2298 9426

VII. CONCLUSION

The Intense Pulsed Light System (ST-690) has the same intended use, same indications for use, the same technological characteristics, the same energy used, and the same operating principles as its predicates. The non-clinical data and performance testing reports in this submission demonstrate that Intense Pulsed Light System meets the expected performance requirements. ST-690 is designed in Tabletop using without caster, but predicate device is designed in Standing using with moving by caster. ST-690 can only setup one handpiece on the device. There is no big difference in intended use or indications for use for the Intense Pulsed Light System. Any difference between the subject and predicate device do not raise new issues of safety or effectiveness. Based on above analysis, the Intense Pulsed Light System is as safe, as effective, and performs as well as the cited predicate device.