



June 21, 2024

Guangzhou Decheng Biotechnology Co., Ltd.

% Joe Shia

Regulatory Consultant

LSI International Inc.

504 East Diamond Ave., Suite H

Gaithersburg, Maryland 20877

Re: K240643

Trade/Device Name: MissLan[®] Early Detection Pregnancy Test Strip; MissLan[®] Early Detection Pregnancy Test Cassette; MissLan[®] Early Detection Pregnancy Test Midstream

Regulation Number: 21 CFR 862.1155

Regulation Name: Human Chorionic Gonadotropin (HCG) Test System

Regulatory Class: Class II

Product Code: LCX

Dated: May 20, 2024

Received: May 20, 2024

Dear Joe Shia:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the

Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Paula V. Caposino -S

Paula Caposino, Ph.D.

Deputy Division Director

Division of Chemistry

and Toxicology Devices

OHT7: Office of In Vitro Diagnostics

Office of Product Evaluation and Quality

Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known)

K240643

Device Name

MissLan® Early Detection Pregnancy Test Strip

MissLan® Early Detection Pregnancy Test Cassette

MissLan® Early Detection Pregnancy Test Midstream

Indications for Use (Describe)

MissLan® Early Detection Pregnancy Test Strip is used for qualitative detection of Human Chorionic Gonadotropin (HCG) in human urine. This test is intended for use as an aid in early detection of pregnancy, in some cases as early as six (6) days before the day of the missed period, (i.e. as early as five (5) days before the day of the expected period). It is intended for use by people who would like to find out whether they are pregnant in a home environment. Only for use outside the body. For over-the-counter use.

Important note regarding positive results:

Because this test detects low levels of hCG, it is possible that this test may give positive results even if you are not pregnant. All results should be confirmed by your healthcare provider, especially when making decisions about future medical care.

This device is intended for home-use only.

MissLan® Early Detection Pregnancy Test Cassette is used for qualitative detection of Human Chorionic Gonadotropin (HCG) in human urine. This test is intended for use as an aid in early detection of pregnancy, in some cases as early as six (6) days before the day of the missed period, (i.e. as early as five (5) days before the day of the expected period). It is intended for use by people who would like to find out whether they are pregnant in a home environment. Only for use outside the body. For over-the-counter use.

Important note regarding positive results:

Because this test detects low levels of hCG, it is possible that this test may give positive results even if you are not pregnant. All results should be confirmed by your healthcare provider, especially when making decisions about future medical care.

This device is intended for home-use only.

MissLan® Early Detection Pregnancy Test Midstream is used for qualitative detection of Human Chorionic Gonadotropin (HCG) in human urine. This test is intended for use as an aid in early detection of pregnancy, in some cases as early as six (6) days before the day of the missed period, (i.e. as early as five (5) days before the day of the expected period). It is intended for use by people who would like to find out whether they are pregnant in a home environment. Only for use outside the body. For over-the-counter use.

Important note regarding positive results:

Because this test detects low levels of hCG, it is possible that this test may give positive results even if you are not pregnant. All results should be confirmed by your healthcare provider, especially when making decisions about future medical care.

This device is intended for home-use only.

Type of Use (Select one or both, as applicable)

☐ Prescription Use (Part 21 CFR 801 Subpart D)

☒ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) SUMMARY

K240643

- 1. Date:** June 21, 2024
- 2. Submitter:** Guangzhou Decheng Biotechnology Co., Ltd.
Floor 3/4/5/7, Building A1, No.12, Nanyun 1st
Road, Science City, Huangpu District,
Guangzhou, Guangdong, 510663, P.R. China
- 3. Contact person:** Joe Shia
LSI International Inc.
504 East Diamond Ave., Suite H
Gaithersburg, MD 20877
Telephone: 240-505-7880
Email: shiajl@yahoo.com
- 4. Device Name:** MissLan® Early Detection Pregnancy Test Strip
MissLan® Early Detection Pregnancy Test Cassette
MissLan® Early Detection Pregnancy Test Midstream
- Classification:** Class II
Product Code: LCX
CFR: 862.1155
- 5. Predicate Devices:** Wondfo One Step HCG Urine Pregnancy Test Midstream,
Wondfo One Step HCG Urine Pregnancy Test Strip,
Wondfo One Step HCG Urine Pregnancy Test Cassette
(K150022)

6. Intended Use

MissLan® Early Detection Pregnancy Test Strip is used for qualitative detection of Human Chorionic Gonadotropin (HCG) in human urine. This test is intended for use as an aid in early detection of pregnancy, in some cases as early as six (6) days before the day of the missed period, (i.e. as early as five (5) days before the day of the expected period). It is intended for use by people who would like to find out whether they are pregnant in a home environment. Only for use outside the body. For over-the-counter use.

Important note regarding positive results:

Because this test detects low levels of hCG, it is possible that this test may give positive results even if you are not pregnant. All results should be confirmed by your healthcare provider, especially when making decisions about future medical care.

This device is intended for home-use only.

MissLan® Early Detection Pregnancy Test Cassette is used for qualitative detection of Human Chorionic Gonadotropin (HCG) in human urine. This test is intended for use as an aid in early detection of pregnancy, in some cases as early as six (6) days before the day of the missed period, (i.e. as early as five (5) days before the day of the expected period). It is intended for use by people who would like to find out whether they are pregnant in a home environment. Only for use outside the body. For over-the-counter use.

Important note regarding positive results:

Because this test detects low levels of hCG, it is possible that this test may give positive results even if you are not pregnant. All results should be confirmed by your healthcare provider, especially when making decisions about future medical care.

This device is intended for home-use only.

MissLan® Early Detection Pregnancy Test Midstream is used for qualitative detection of Human Chorionic Gonadotropin (HCG) in human urine. This test is intended for use as an aid in early detection of pregnancy, in some cases as early as six (6) days before the day of the missed period, (i.e. as early as five (5) days before the day of the expected period). It is intended for use by people who would like to find out whether they are pregnant in a home environment. Only for use outside the body. For over-the-counter use.

Important note regarding positive results:

Because this test detects low levels of hCG, it is possible that this test may give positive results even if you are not pregnant. All results should be confirmed by your healthcare provider, especially when making decisions about future medical care.

This device is intended for home-use only.

7. Device Description

MissLan® Early Detection Pregnancy Test will be sold in Strip, Cassette and Midstream format. The Strip format is a single test strip. The Cassette format consists of a single test strip assembled in a plastic housing. The Midstream format consists of a single test strip assembled in a plastic housing with an absorbent tip, and is designed to be tested in dip or midstream mode.

MissLan® Early Detection Pregnancy Test Strip, MissLan® Early Detection Pregnancy Test Cassette and MissLan® Early Detection Pregnancy Test Midstream each contains a pouch with the device and instructions, and in addition, cassette format is packaged with pipette dropper and urine collection cup.

The devices utilize a combination of antibodies to detect hCG in urine as well as to serve as a run control. Each device contains mouse monoclonal anti- β -hCG antibody colloidal gold conjugate pre-dried on the sample pad. Mouse monoclonal

anti- α -hCG antibody (on the Test Line) and goat anti mouse IgG polyclonal antibody (on the Control Line) are coated and immobilized on a nitrocellulose membrane. The result is displayed to the user in the test window as two lines for a 'Pregnant' positive result and one line for a 'Not Pregnant' negative result.

8. Substantial Equivalence Information

Similarities		
Item	Candidate device	Predicate device
Intended use	Early detection of pregnancy	Same
Specimen	Urine	Same
Assay technical	Immunochromatographic assay	Same
Sensitivity	10 mIU/mL	Same
Results	Qualitative	Same
Device format	Strip, Cassette, Midstream	Same
Differences		
Item	Device	Predicate
Target user	For over-the-counter use	Prescription use and OTC use
Read time	3-10 minutes	5 minutes

9. Test Principle

MissLan® Early Detection Pregnancy Test Strip, MissLan® Early Detection Pregnancy Test Cassette and MissLan® Early Detection Pregnancy Test Midstream use lateral flow immunoassay for in vitro qualitative detection of Human Chorionic Gonadotropin (HCG) in human urine. If hCG is present in the sample, it will reach the Test Zone (T) of the membrane and form a colored line. When the test is performed properly, a colored line will always appear in the Control Zone (C). The test result is shown in the result window and read visually in 3 minutes of urine addition. Two distinct colored lines, one in the Test Zone (T) and another in the Control Zone (C) indicate a positive test result (pregnant). Absence of a colored line in the Test Zone (T) and only a colored line in the Control Zone (C) indicates a negative test result (not pregnant). Absence of a colored line in the Control Zone (C) even in the presence of a colored line in the Test Zone (T) indicates an invalid test result.

10. Performance Characteristics

A. Analytical performance

a. Precision/Reproducibility/Sensitivity

Negative female urine was spiked with hCG standard (Traceable to the 5th WHO) to hCG concentrations of 0, 3, 3.5, 5, 6.5, 7.5, 10, 15 and 25 mIU/mL. Each sample

was tested in 10 replicates per day for 5 days for each device lot. Total of three device lots for each format were tested and each operator tested one lot separately. Tests were performed by three different operators for each sample concentration. The midstream format was performed using both dip and in-stream methods. The results are summarized in the table below:

Strip format

hCG Concentration (mIU/mL)	Operator 1		Operator 2		Operator 3		Total result		% Negative	% Positive
	Lot I		Lot II		Lot III					
	-	+	-	+	-	+	-	+		
0	50	0	50	0	50	0	150	0	100%	0%
3	50	0	50	0	50	0	150	0	100%	0%
3.5	44	6	45	5	45	5	134	16	89.33%	10.67%
5	36	14	37	13	36	14	109	41	72.67%	27.33%
6.5	27	23	23	27	26	24	76	74	50.67%	49.33%
7.5	10	40	13	37	13	37	36	114	24.00%	76.00%
10	0	50	0	50	0	50	0	150	0%	100%
15	0	50	0	50	0	50	0	150	0%	100%
25	0	50	0	50	0	50	0	150	0%	100%

Cassette format

hCG Concentration (mIU/mL)	Operator 1		Operator 2		Operator 3		Total result		% Negative	% Positive
	Lot I		Lot II		Lot III					
	-	+	-	+	-	+	-	+		
0	50	0	50	0	50	0	150	0	100%	0%
3	50	0	50	0	50	0	150	0	100%	0%
3.5	45	5	45	5	45	5	135	15	90.00%	10.00%
5	37	13	35	15	35	15	107	43	71.33%	28.67%
6.5	23	27	26	24	24	26	73	77	48.67%	51.33%
7.5	11	39	13	37	10	40	34	116	22.67%	77.33%
10	0	50	0	50	0	50	0	150	0%	100%
15	0	50	0	50	0	50	0	150	0%	100%
25	0	50	0	50	0	50	0	150	0%	100%

Midstream format (in-stream method)

hCG Concentration (mIU/mL)	Operator 1		Operator 2		Operator 3		Total result		% Negative	% Positive
	Lot I		Lot II		Lot III					
	-	+	-	+	-	+	-	+		
0	50	0	50	0	50	0	150	0	100%	0%
3	50	0	50	0	50	0	150	0	100%	0%
3.5	45	5	45	5	45	5	135	15	90.00%	10.00%
5	37	13	36	14	37	13	110	40	73.33%	26.67%
6.5	24	26	26	24	24	26	74	76	49.33%	50.67%
7.5	11	39	12	38	14	36	37	113	24.67%	75.33%
10	0	50	0	50	0	50	0	150	0%	100%
15	0	50	0	50	0	50	0	150	0%	100%
25	0	50	0	50	0	50	0	150	0%	100%

Midstream format (dip method)

hCG Concentration (mIU/mL)	Operator 1		Operator 2		Operator 3		Total result		% Negative	% Positive
	Lot I		Lot II		Lot III					
	-	+	-	+	-	+	-	+		
0	50	0	50	0	50	0	150	0	100%	0%
3	50	0	50	0	50	0	150	0	100%	0%
3.5	45	5	44	6	45	5	134	16	89.33%	10.67%
5	37	13	36	14	35	15	108	42	72.00%	28.00%
6.5	23	27	24	26	26	24	73	77	48.67%	51.33%
7.5	12	38	14	36	9	41	35	115	23.33%	76.67%
10	0	50	0	50	0	50	0	150	0%	100%
15	0	50	0	50	0	50	0	150	0%	100%
25	0	50	0	50	0	50	0	150	0%	100%

MissLan® Early Detection Pregnancy Test exhibited reproducibility of results.

Based on the above results, the sensitivity of MissLan® Early Detection Pregnancy Test is demonstrated to be 10 mIU/mL.

b. Linearity/assay reportable range:

Linearity is not applicable since this is a qualitative test.

The test device was evaluated for high dose or hook effect.

Hook effect test:

Negative urine samples were spiked with varying hCG concentrations (6,250 mIU/mL, 12,500 mIU/mL, 25,000 mIU/mL, 50,000 mIU/mL, 100,000 mIU/mL, 200,000 mIU/mL and 500,000 mIU/mL). The results demonstrated that no hook effect was observed at hCG concentration up to 500,000 mIU/mL.

c. Traceability, Stability, Expected values (controls, calibrators, or methods):

Traceability:

MissLan® Early Detection Pregnancy Test is calibrated against reference material traceable to WHO International Standard 5th edition, NIBSC code 07/364.

Stability:

The shelf-life of the MissLan® Early Detection Pregnancy Test at 2~30°C is 36 months based on real time stability data.

d. Analytical specificity

Interfering substance:

To evaluate potential interfering substances of the MissLan® Early Detection Pregnancy Test, urine samples containing 0 mIU/mL, 3 mIU/mL and 10 mIU/mL hCG were spiked with interfering substances to the concentrations listed below. No interference effect was observed at the tested concentrations.

<i>Substance</i>	<i>Concentration</i>
Albumin	2000 mg/dL

<i>Substance</i>	<i>Concentration</i>
Bilirubin	40 mg/dL
Glucose	2000 mg/dL
Hemoglobin	1000 mg/dL
Uric acid	23.5 mg/dL
Ketone	20 mg/dL
β -hydroxybutyrate	2000 mg/dL
Pregnanediol	1.5 mg/dL
Acetaminophen	20 mg/dL
Acetylsalicylic acid	80 mg/dL
Amoxicillin	20 mg/dL
Ampicillin	20 mg/dL
Ascorbic acid	20 mg/dL
Atropine	20 mg/dL
Benzoyllecgonine	10 mg/dL
Caffeine	20 mg/dL
Cannabinol	10 mg/dL
Codeine	6 μ g/dL
EDTA	80 mg/dL
Ephedrine	20 mg/dL
Ethanol	1%
Folic acid	0.03 mg/dL
Salicylic acid	20 mg/dL
Gentisic acid	20 mg/dL
Ibuprofen	40 mg/dL
Phenothiazine	20 mg/dL
Phenylpropanolamine	20 mg/dL
Tetracycline	20 mg/dL
Thiophene	20 mg/dL
Vitamin B1	80 mg/dL

Cross-reactivity:

To evaluate cross-reactivity, negative and positive urine samples (0 mIU/mL, 3 mIU/mL and 10 mIU/mL hCG) were spiked with potential cross reactants (500 mIU/mL hLH, 1000 mIU/mL hFSH, 1000 μ IU/mL hTSH). No cross-reactivity was observed at tested concentration.

Effects of hCG β -core fragment:

To evaluate the effect of the hCG β -core fragment, negative urine samples (0 mIU/mL and 3 mIU/mL hCG) and positive urine samples (10 mIU/mL and 20,000 mIU/mL hCG) were spiked with hCG β -core fragment (hCG β cf) at concentrations of 50,000 pmol/L, 125,000 pmol/L, 250,000 pmol/L and 500,000 pmol/L. The performance of MissLan[®] Early Detection Pregnancy Test was not affected by hCG β -core fragment concentrations up to 500,000 pmol/L.

Effects of urine pH:

To evaluate the effect of urine pH on the results of MissLan® Early Detection Pregnancy Test, urine samples containing 0 mIU/mL, 3 mIU/mL and 10 mIU/mL hCG were tested at pH values of 4, 5, 6, 7, 8 and 9. The results indicated that urine pH ranges between 4 and 9 does not affect the performance of MissLan® Early Detection Pregnancy Test.

Effects of urine density:

To evaluate the effect of urine density on the results of MissLan® Early Detection Pregnancy Test, urine samples containing 0 mIU/mL, 3 mIU/mL and 10 mIU/mL hCG were tested at density values of 1.000, 1.005, 1.010, 1.015, 1.020, 1.025, 1.030 and 1.035. The results indicated that urine with a relative density range of 1.000 to 1.035 does not affect the performance of MissLan® Early Detection Pregnancy Test.

B. Method comparison study

Method comparison with predicate device:

The performance of the candidate device was compared to the predicate device. Urine samples were collected from 400 women aged 18 to 51 at three clinical sites. Approximately half of the subjects were pregnant in the early stage of less than 5 weeks. Samples were randomly collected at various time throughout the day.

All samples were tested by three professionals at each site using the candidate device and by one professional at each site using the predicate device. 100 samples were tested for strip format, and 100 samples were tested for cassette format; for the midstream format, both dip testing (100 samples) and in-stream testing (100 samples) were evaluated. The samples were blinded and randomized before being tested by professionals. All results are summarized in the table below.

Candidate device		Predicate device		
		Positive	Negative	Total
Strip format	Positive	49	0	49
	Negative	0	51	51
	Total	49	51	100

Candidate device		Predicate device		
		Positive	Negative	Total
Cassette format	Positive	51	0	51
	Negative	0	49	49
	Total	51	49	100

Candidate device		Predicate device		
		Positive	Negative	Total
Midstream format	Positive	48	0	48

(dip method)	Negative	0	52	52
	Total	48	52	100

Candidate device		Predicate device		
		Positive	Negative	Total
Midstream format (in-stream method)	Positive	49	0	49
	Negative	0	51	51
	Total	49	51	100

The conformity between MissLan[®] Early Detection Pregnancy Test (strip / cassette / midstream) and the predicate device is 100%.

C. Clinical studies

1. Clinical Sensitivity:

Not applicable.

2. Clinical Specificity:

Not applicable.

3. Other Clinical Supportive Data (When 1. and 2. Are Not Applicable):

3.1 Lay person study

First study:

A lay user study was performed at three sites with a total of 400 females with diverse educational and occupational backgrounds and ages ranging from 18 to 51 years old. Each subject tested her own urine sample using the device according to the instructions for use and provided a sample for professional testing. 100 subjects tested the strip format, 100 subjects tested the cassette format, 100 subjects tested the midstream format using dip method, and 100 subjects tested the midstream format using in-stream method.

The results are summarized in the table below.

Strip format

Strip format		Professional		
		Positive	Negative	Total
Layperson	Positive	49	0	49
	Negative	0	51	51
	Total	49	51	100

Cassette format

Cassette format		Professional		
		Positive	Negative	Total
Layperson	Positive	51	0	51
	Negative	0	49	49

	Total	51	49	100
--	-------	----	----	-----

Midstream format

Midstream format (dip method)		Professional		
		Positive	Negative	Total
Layperson	Positive	48	0	48
	Negative	0	52	52
	Total	48	52	100

Midstream format (in-stream method)		Professional		
		Positive	Negative	Total
Layperson	Positive	49	0	49
	Negative	0	51	51
	Total	49	51	100

From the above tables, the lay person results showed 100% positive and 100% negative conformity with the professional results.

Second study:

Negative urine sample pools were spiked with 3 mIU/mL, 3.5 mIU/mL, 6.5 mIU/mL and 10 mIU/mL hCG. All aliquots were blind labeled by the person who prepared the samples and didn't take part in the sample testing. A total of 300 laypersons performed the testing, of which 100 using the strip format, 100 using the cassette format and 100 using the midstream format. Each layperson tested 4 blind-labeled samples (3 mIU/mL, 3.5 mIU/mL, 6.5 mIU/mL and 10 mIU/mL), and an aliquot of each of the spiked urine samples was also tested by the professional.

Results of spiked samples are summarized in the table below.

Strip format:

hCG Concentration (mIU/mL)	Layperson result			Professional result			Percent Agreement
	No. of Positive	No. of Negative	% positive	No. of Positive	No. of Negative	% positive	
3	0	100	0%	0	100	0%	100%
3.5	9	91	9%	10	90	10%	95%
6.5	49	51	49%	51	49	51%	96%
10	100	0	100%	100	0	100%	100%

Cassette format:

hCG Concentration (mIU/mL)	Layperson result			Professional result			Percent Agreement
	No. of Positive	No. of Negative	% positive	No. of Positive	No. of Negative	% positive	

3	0	100	0%	0	100	0%	100%
3.5	10	90	10%	11	89	11%	95%
6.5	48	52	48%	52	48	52%	96%
10	100	0	100%	100	0	100%	100%

Midstream format:

hCG Concentration (mIU/mL)	Layperson result			Professional result			Percent Agreement
	No. of Positive	No. of Negative	% positive	No. of Positive	No. of Negative	% positive	
3	0	100	0%	0	100	0%	100%
3.5	10	90	10%	11	89	11%	95%
6.5	49	51	49%	52	48	52%	97%
10	100	0	100%	100	0	100%	100%

The results above show that MissLan® Early Detection Pregnancy Test can be used by the laypersons and get the correct results.

3.2 Early pregnancy detection study

A study was performed to evaluate the early pregnancy detection rate of MissLan® Early Detection Pregnancy Test. A total of 650 early pregnancy urine samples from day -8 to +1 relative to the day of the expected menstrual period were collected from 65 pregnant women. Each sample was tested with each format of the MissLan® Early Detection Pregnancy Test. For the midstream format, both dip and in-stream methods were evaluated.

The specific detection rates per day are summarized in the table below.

Strip format:

Test date	EMP-8	EMP-7	EMP-6	EMP-5	EMP-4	EMP-3	EMP-2	EMP-1	EMP	EMP+1
# of positive results of B-ultrasound	65	65	65	65	65	65	65	65	65	65
# of positive results of HCG	1	8	25	50	63	65	65	65	65	65
Detection rate	1.54%	12.31%	38.46%	76.92%	96.92%	100%	100%	100%	100%	100%

Cassette format:

Test date	EMP-8	EMP-7	EMP-6	EMP-5	EMP-4	EMP-3	EMP-2	EMP-1	EMP	EMP+1
# of positive results of B-ultrasound	65	65	65	65	65	65	65	65	65	65
# of positive results of HCG	1	8	25	50	63	65	65	65	65	65
Detection rate	1.54%	12.31%	38.46%	76.92%	96.92%	100%	100%	100%	100%	100%

Midstream format (in-stream method):

Test date	EMP-8	EMP-7	EMP-6	EMP-5	EMP-4	EMP-3	EMP-2	EMP-1	EMP	EMP+1
# of positive	65	65	65	65	65	65	65	65	65	65

results of B-ultrasound										
# of positive results of HCG	1	8	25	50	63	65	65	65	65	65
Detection rate	1.54%	12.31%	38.46%	76.92%	96.92%	100%	100%	100%	100%	100%

Midstream format (dip method):

Test date	EMP-8	EMP-7	EMP-6	EMP-5	EMP-4	EMP-3	EMP-2	EMP-1	EMP	EMP+1
# of positive results of B-ultrasound	65	65	65	65	65	65	65	65	65	65
# of positive results of HCG	1	8	25	50	63	65	65	65	65	65
Detection rate	1.54%	12.31%	38.46%	76.92%	96.92%	100%	100%	100%	100%	100%

The detection rate of MissLan[®] Early Detection Pregnancy Test (strip/cassette/midstream) in EMP-5 (i.e. five days before the expected menstrual period, 6 days before the missed menstrual period) was 76.92%, in EMP-4 was 96.92%, and in EMP-1 was 100%.

3.3 Specificity Study to Determine False-Positive Results Rate

900 urine samples were collected from non-pregnant females in pre-menopausal (18-40 years old), peri-menopausal (41-55 years old) and post-menopausal (>55 years old) age groups at three sites. 300 participants for each age group. In each age group, 100 participants tested test strip, 100 participants tested test cassette, 50 participants tested test midstream using dip method, and 50 participants tested test midstream using in-stream method. Three batches of each format were used in the study. No false positive results were observed for all age groups.

Strip format:

Group	Lot I	Lot II	Lot III	Total result
Pre-menopausal	0+/34-	0+/33-	0+/33-	0+/100-
Peri-menopausal	0+/33-	0+/34-	0+/33-	0+/100-
Post-menopausal	0+/33-	0+/33-	0+/34-	0+/100-

Cassette format:

Group	Lot I	Lot II	Lot III	Total result
Pre-menopausal	0+/34-	0+/33-	0+/33-	0+/100-
Peri-menopausal	0+/33-	0+/34-	0+/33-	0+/100-
Post-menopausal	0+/33-	0+/33-	0+/34-	0+/100-

Midstream format:

Group	Lot I	Lot II	Lot III	Total result
Pre-menopausal	0+/34-	0+/33-	0+/33-	0+/100-
Peri-menopausal	0+/33-	0+/34-	0+/33-	0+/100-
Post-menopausal	0+/33-	0+/33-	0+/34-	0+/100-

11. Conclusion

Based on the test principle and performance characteristics of the device including precision, cut-off, interference, specificity, method comparison, lay-user study and early pregnancy detection study of the devices, it can be concluded that MissLan[®] Early Detection Pregnancy Test is substantially equivalent to the predicate device.