



June 24, 2024

Won Tech Co., Ltd.
Hyun Sik Yoon
Team Leader
64 Techno8-ro 64, Yuseong-gu
Daejeon, 34028
Korea, South

Re: K241144

Trade/Device Name: Picosecond Nd:YAG Laser (PICOCAREMAJESTY)
Regulation Number: 21 CFR 878.4810
Regulation Name: Laser Surgical Instrument For Use In General And Plastic Surgery And In Dermatology
Regulatory Class: Class II
Product Code: GEX
Dated: April 15, 2024
Received: April 25, 2024

Dear Hyun Sik Yoon:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>).

Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Tanisha L. Hithe-S
Digitally signed by
Tanisha L. Hithe-S
Date: 2024.06.24
15:13:13 -04'00'

Tanisha Hithe
Assistant Director
DHT4A: Division of General Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

Submission Number (if known)

K241144

Device Name

Picosecond Nd:YAG Laser (PICOCAREMAJESTY)

Indications for Use (Describe)

The PICOCAREMAJESTY is indicated for the following at the specified wavelength:

- 1064 nm: Removal of tattoos for all skin types (Fitzpatrick I - VI) to treat the following tattoo colors: black, brown, green, blue and purple.
- 532 nm: Removal of tattoos for Fitzpatrick skin types I - III to treat the following tattoo colors: red, yellow and orange.
- Treatment of benign pigmented lesions for Fitzpatrick skin types I-IV.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

K241144

[As required by 21 CFR 807.92]

1. Date Prepared [21 CFR 807.92(a)(a)]

June 6th, 2024

2. Submitter's Information & Contact Person [21 CFR 807.92(a)(1)]

- Name of Manufacturer: WON TECH Co., Ltd.
- Address: 64 Techno 8-ro, Yuseong-gu, Daejeon, 34028,
Republic of Korea
- Contact Name: Hyun Sik Yoon
- Telephone No.: +82-10-6750-5346
- Fax No.: +82-70-7836-0110
- Email Address: yoonhs21@wtlaser.com

3. Trade Name, Common Name, Classification [21 CFR 807.92(a)(2)]

Common name: Nd:YAG Laser System

Trade name: PICOCAREMAJESTY

Classification Description	21 CFR Section	Product Code
Powered Laser Surgical Instrument	878.4810	GEX

As stated in 21 CFR, parts 878.4810, this generic types of devices has been classified as Class II.

4. Identification of Predicate Device(s) [21 CFR 807.92(a)(3)]

The identified predicate devices within this submission are shown as follow:

#1 Predicate device

- 510(k) Number: K234104
- Applicant: WON TECH Co., Ltd.
- Classification Name: Powered Laser Surgical Instrument
- Trade Name: PICOANDY

#2 Predicate device

- 510(k) Number: K201773
- Applicant: WON TECH Co., Ltd.
- Classification Name: Powered Laser Surgical Instrument
- Trade Name: PICOHI

5. Description of the Device [21 CFR 807.92(a)(4)]

The PICOCAREMAJESTY is the solid state laser capable of delivering energy at wavelengths of 1064nm or 532nm at short durations. The device system consists of a main unit, an articulated arm, a handpiece and a foot switch. The laser output is delivered to the skin through the articulated arm delivery system terminated by the handpiece. The fluence (energy density) and frequency are controlled from the LCD display/Touch Pad located on the front of the main unit. The LCD display is used to obtain feedback from the system, such as the number of pulses delivered or spot size selected.

For treatment, the user can select the appropriate fluence value. (The treatment parameters recommended is provided in the user Manual) The energy is changed automatically in accordance with the selected fluence (down) button and then, the fluence is increased / decreased by 0.1 J/cm².

The selectable fluence values are 0.2 to 10.0 J/cm² at 1064nm and 0.1 to 2.5 J/cm² at 532nm

6. Indications for Use [21 CFR 807.92(a)(5)]

The PICOCAREMAJESTY is indicated for the following at the specified wavelength:

- 1064 nm: Removal of tattoos for all skin types (Fitzpatrick I-VI) to treat the following tattoo colors: black, brown, green, blue and purple.

- 532 nm: Removal of tattoos for Fitzpatrick skin types I-III to treat the following tattoo colors: red, yellow and orange.
- Treatment of benign pigmented lesions for Fitzpatrick skin types I-IV.

7. Determination of Substantial Equivalence [21 CFR 807.92(a)(6) and 21 CFR 807.92(b)]

There are no significant differences between the PICOCAREMAJESTY and the predicate devices that would adversely affect the use of the product. It is substantially equivalent to this device in design, function, and technical characteristics.

	Proposed Device	Predicate Device #1	Predicate Device #2	SE decision
K Number	-	K234104	K201773	-
Manufacturer	WON TECH Co., Ltd.	WON TECH Co., Ltd.	Hironic Co., Ltd	-
Model	PICOCAREMAJESTY	PICOANDY	PICOHI	-
Indications for Use	The PICOCAREMAJESTY is indicated for the following at the specified wavelength: • 1064 nm: Removal of tattoos for all skin types (Fitzpatrick I-VI) to treat the following tattoo colors: black, brown, green, blue and purple. • 532 nm: Removal of tattoos for Fitzpatrick skin types I-III to treat the following tattoo colors: red, yellow and orange. • Treatment of benign pigmented lesions for Fitzpatrick skin types I-IV.	The PICOANDY is indicated for the following at the specified wavelength: • 1064 nm: Removal of tattoos for all skin types (Fitzpatrick I-VI) to treat the following tattoo colors: black, brown, green, blue and purple. • 532 nm: Removal of tattoos for Fitzpatrick skin types I-III to treat the following tattoo colors: red, yellow and orange. • Treatment of benign pigmented lesions for Fitzpatrick skin types I-IV.	The PICOHI is indicated for the following at the specified wavelength: • 1064 nm: Removal of tattoos for all skin types (Fitzpatrick I-VI) to treat the following tattoo colors: black, brown, green, blue and purple. • 532 nm: Removal of tattoos for Fitzpatrick skin types I-III to treat the following tattoo colors: red, yellow and orange. • Treatment of benign pigmented lesions for Fitzpatrick skin types I-IV.	Same
Anatomical site	Skin and subcutaneous tissue	Skin and subcutaneous tissue	Skin and subcutaneous tissue	Same
Wavelength	1064 & 532nm	1064 & 532nm	1064 & 532nm	Same

	Proposed Device	Predicate Device #1	Predicate Device #2	SE decision
Pulse Width	275~300ps	350 ~550 ps	275~300ps	Same with Predicate Device #2
Pulse Energy	1. Zoom 1064nm : 30 ~ 500mJ $\pm 10\%$ 2. Zoom 532nm : 10 ~ 250mJ $\pm 10\%$ 3. Collimation Mode: 38 ~ 500mJ $\pm 10\%$ 4. MLA (1064nm): 30 ~ 500mJ $\pm 10\%$ 5. MLA (532nm): 10 ~ 250mJ $\pm 10\%$ 6. 1064nm DOE : 30 ~ 450mJ $\pm 10\%$ 7. 532nm DOE : 20 ~ 200mJ $\pm 10\%$	1. Zoom 1064nm : 30 ~ 550mJ $\pm 10\%$ 2. Zoom 532nm : 10 ~ 200mJ $\pm 10\%$ 3. Collimation Mode: 38 ~ 500mJ $\pm 10\%$ 4. MLA (1064nm): 30 ~ 550mJ $\pm 10\%$ 5. MLA (532nm): 10 ~ 200mJ $\pm 10\%$ 6. 1064nm DOE : 30 ~ 550mJ $\pm 10\%$ 7. 532nm DOE : 20 ~ 200mJ $\pm 10\%$	1. Zoom 1064nm : 10 ~ 508mJ $\pm 20\%$ 2. Zoom 532nm : 4 ~ 255mJ $\pm 20\%$ 3. Collimation Mode: 7 ~ 471mJ $\pm 20\%$ 4. ZMLA (1064nm): 5 ~ 508mJ $\pm 20\%$ 5. ZMLA (532nm): 2 ~ 255mJ $\pm 20\%$ 6. 1064nm DOE : 10 ~ 500mJ $\pm 20\%$ 7. 532nm DOE : 10 ~ 250mJ $\pm 20\%$	Diffenence However, the pulse energy of PICOAREMAJESTY is included in those of the Predicate Device #2
Spot Size	2 to 10 mm	2 to 10 mm	1.5 to 13 mm	Same
Pulse Repetition	Max.10Hz	Max.10Hz	Max.10Hz	Same
Laser Delivery Type	Articulated Arm with Handpiece	Articulated Arm with Handpiece	Articulated Arm with Handpiece	Same
Handpiece	Zoom(1064nm, 532nm) Collimated(1064nm) MLA (1064nm, 532nm) DOE (1064nm, 532nm)	Zoom(1064nm, 532nm) Collimated(1064nm) MLA (1064nm, 532nm) DOE (1064nm, 532nm)	Zoom(1064nm, 532nm) Collimated(1064nm) ZMLA (1064nm, 532nm) DOE (1064nm, 532nm)	Same
Patient Contact Material	Aluminum powder (Handpiece)	Aluminum powder (Handpiece)	Aluminum powder (Handpiece)	Same

Non-Clinical Test Summary [21 CFR 807.92(b)(1)]

1) Electrical Safety, Electromagnetic Compatibility Testing

Bench tests were conducted to verify that the proposed device met all design specifications. The test results demonstrated that the proposed device complies with the following standards:

Standard (Edition)	Standard Title
--------------------	----------------

AASI AAMI ES60601-1:2005/(R)2012 and A1:2012	Medical electrical equipment – Part 1: General requirements for basic safety and essential performance
ANSI AAMI IEC 60601-1-2:2014	Medical electrical equipment - Part 1-2: General requirements for safety and essential performance - Collateral standard: Electromagnetic compatibility - Requirements and tests
IEC 60601-1-6 Edition 3.1 2013	Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability
IEC 60601-2-22 Edition 3.1 2012-10	Medical electrical equipment - Part 2-22: Particular requirements for basic safety and essential performance of surgical, cosmetic, therapeutic and diagnostic laser equipment
IEC 60825-1:2014 (Third Edition)	Safety of laser products - Part 1: Equipment classification and requirements

-

2) Software Validation

The PICOCAREMAJESTY contains MODERATE level of concern software. Software was designed and developed according to a software development process and was verified and validated.

The software information is provided in accordance with FDA guidance: The content of premarket submissions for software contained in medical devices, on May 11, 2005.

3) Biocompatibility

Part	Material	Patient Contact	Duration of Contact by ISO 10993-1	Bio-compatibility
Handpiece Tip of PICOANDY	Aluminium Powder (Cas No. 7429-90-5)	Intact Skin	Limited (< 24 hours)	Yes

- The material of handpiece for the PICOCAREMAJESTY is same to previous registered material of own product (PICOCARE Family <K181272>, ICOCAREMAJESTY <K212127>, and PICOANDY <K234104>) in formulation, processing, and sterilization, and no other chemicals have been added. Thus, biocompatibility is not necessary for this device, as this material has been deemed biocompatible in others previously cleared submissions.

4) Performance Testing

The performance of the PICOCAREMAJESTY has been defined as follows. - Laser wavelength: 1064nm±10%, 532nm±10% - Laser output power: 1064nm mode: (30 to 500)mJ, 532nm mode: (20 to 250)mJ - Pulse Width: 1064nm mode: 300-400 ps, 532nm mode : 300-400 ps - Pulse repetition rate: 1064nm

64 Techno 8-ro, Yuseong-gu, Daejeon, Republic of Korea

TEL: +82 (42) 934 6800 FAX: +82 (42) 934 9491

mode:1-10Hz, 532nm mode : 1-10Hz - Radiation diameter: 1064nm mode: (2 to 10mm) step: 1mm, 532nm mode: (2 to 10mm) step: 1mm

5) Cleaning Validation

Cleaning validation is not needed, because PICOAREMAJESTY is not in contact with patients

Clinical Test Summary [21 CFR 807.92(b)(2)]

No clinical studies were considered necessary and performed.

Conclusion [21 CFR 807.92(b)(3)]

In according with the Federal Food & drug and cosmetic Act, 21 CFR Part 807, and based on the information provided in this premarket notification WON TECH Co., Ltd. concludes that the PICOAREMAJESTY is substantially equivalent to predicate devices as described herein.