



August 22, 2024

Lijia LLC
% Jennifer Day
Regulatory Correspondent
Prime Path Medtech
1321 Upland Dr
Suite 6792
Houston, Texas 77043

Re: K241369
Trade/Device Name: LIJIA Night Guard
Regulatory Class: Unclassified
Product Code: MQC
Dated: May 14, 2024
Received: May 14, 2024

Dear Jennifer Day:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,



Bobak
Shirmohammadi -S

For Michael E. Adjodha, M. ChE., RAC, CQIA
Assistant Director
DHT1B: Division of Dental and
ENT Devices
OHT1: Office of Ophthalmic, Anesthesia,
Respiratory, ENT, and Dental Devices
Office of Product Evaluation and Quality

Center for Devices and Radiological Health

Enclosure

Indications for Use

Submission Number (if known)

K241369

Device Name

LIJIA Night Guard

Indications for Use (Describe)

The LIJIA Night Guard is indicated for protection of teeth and restorations against grinding and clenching.

Type of Use (Select one or both, as applicable)



Prescription Use (Part 21 CFR 801 Subpart D)



Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(k) Summary

A summary of 510(k) safety and effectiveness information in accordance with the requirements of 21 CFR 807.92.

Submitter: LIJIA, LLC

Company Contact Person: Ben Li
Address: 1520 Brookhollow Dr, Santa Ana, CA 92705
Phone: 310-866-0323
Email: b.jaminli@gmail.com

Submission Correspondent: Jennifer Day, Regulatory Affairs Consultant
Address: 1321 Upland Dr. Suite 6792
Houston, TX 77043
Phone: 314-809-1919
Email: jday@primepathmedtech.com

Date Prepared: 5/9/2024

Proprietary Name: LIJIA Night Guard

Common Name: Mouth Guard

Product Code: MQC

Device Classification: Unclassified

Predicate Device: JS Dental Lab Mouth Guard

Reference Device: FBCC Night Guard

Device Description:

The LIJIA Night Guard is a mouth guard used as a barrier between teeth for nighttime teeth grinding by creating physical separation between upper and lower tooth surfaces preventing tooth damage caused by bruxism (e.g., grinding and clenching).

Indications for Use:

The LIJIA Night Guard is indicated for protection of teeth and restorations against grinding and clenching.

Comparison to Predicate Devices:

Table 1. Predicate and Reference Device Comparison

Specification	Subject Device: <i>LIIJA Night Guard</i>	Predicate Device: <i>JS Dental Lab Mouth Guard (K210011)</i>	Reference Device: <i>FBCC Night Guard (K221369)</i>	Comparison Result
Product Code	MQC	MQC	OBR	Same
Classification Panel	Dental	Dental	Dental	Same
Device Class	Unclassified	Unclassified	Unclassified	Same
OTC or Rx	Rx	Rx	OTC	Same
Anatomical Sites	Worn on maxillary or mandibular teeth	Worn on maxillary or mandibular teeth	Worn on maxillary or mandibular teeth	Same
Sterile	Non-Sterile	Non-sterile	Non-Sterile	Same
Patient Removable?	Yes	Yes	Yes	Same
Indication for Use	The LIIJA Night Guard is indicated for protection of teeth and restorations against grinding and clenching.	The Subject Device (The JS Dental Mouth Guard) is intended for protection against bruxism and teeth grinding. They create a barrier between the upper and lower dentition to protect the patient's overall occlusion.	The FBCC Night Guard is indicated for protection of teeth and restorations against grinding and clenching.	Same as reference

Indications for Use:

The indication for use is very similar to the predicate device, the JS Dental Lab Mouth Guard, and the same as the the reference device, the FBCC Night Guard. All of the devices are indicated for the protection of teeth against grinding and clenching.

Technological Features:

The LIIJA Night Guard uses the same mechanism of action as the predicate product and reference product. It is intended to be worn on the teeth to protect teeth against grinding and clenching and to mitigate the pains associated with temporomandibular disorders.

Non-Clinical Performance Testing

The use of thermoplastic materials for Night Guards has been well documented in scientific literature. However, durability testing was completed on these night guards.

An internal manufacturing validation was performed to test the dimensional accuracy of the manufacturing process for the Night Guards.

Clinical Performance Testing

The technological characteristics, indications for use, material, and manufacturing and processes are the same or similar to the predicate device and therefore, no clinical studies were deemed necessary to demonstrate the safety and effectiveness of the subject device.

Conclusion

Based on the similarities in indications for use, technological characteristics, and non-clinical performance testing the Night Guards are substantially equivalent to the predicate device.