



January 15, 2025

Cepheid
Jennifer Motto
Regulatory Affairs Principal
904 Caribbean Drive
Sunnyvale, California 94089

Re: K242109

Trade/Device Name: Xpert Xpress CoV-2 plus

Regulation Number: 21 CFR 866.3981

Regulation Name: Device To Detect And Identify Nucleic Acid Targets In Respiratory Specimens
From Microbial Agents That Cause The SARS-Cov-2 Respiratory Infection And
Other Microbial Agents When In A Multi-Target Test

Regulatory Class: Class II

Product Code: QQX

Dated: July 18, 2024

Received: July 19, 2024

Dear Jennifer Motto:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801 and Part 809); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

All medical devices, including Class I and unclassified devices and combination product device constituent parts are required to be in compliance with the final Unique Device Identification System rule ("UDI Rule"). The UDI Rule requires, among other things, that a device bear a unique device identifier (UDI) on its label and package (21 CFR 801.20(a)) unless an exception or alternative applies (21 CFR 801.20(b)) and that the dates on the device label be formatted in accordance with 21 CFR 801.18. The UDI Rule (21 CFR 830.300(a) and 830.320(b)) also requires that certain information be submitted to the Global Unique Device Identification Database (GUDID) (21 CFR Part 830 Subpart E). For additional information on these requirements, please see the UDI System webpage at <https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/unique-device-identification-system-udi-system>.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory->

[assistance/contact-us-division-industry-and-consumer-education-dice](#)) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Anna M. Mielech -S

Anna Mielech, PhD.
Deputy Branch Chief (Acting)
Viral Respiratory and HPV Branch
Division of Microbiology Devices
OHT7: Office of In Vitro Diagnostics
Office of Product Evaluation and Quality
Center for Devices and Radiological Health Enclosure

Indications for Use

510(k) Number (if known)

K242109

Device Name

Xpert Xpress CoV-2 plus

Indications for Use (Describe)

The Xpert Xpress CoV-2 plus test, performed on the GeneXpert Xpress System, is a rapid real-time RT-PCR test intended for the qualitative detection of SARS-CoV-2 RNA in nasopharyngeal and anterior nasal swab specimens collected from individuals with signs and symptoms of respiratory tract infection.

The Xpert Xpress CoV-2 plus test is intended for use as an aid in the diagnosis of COVID-19 if used in conjunction with other clinical, epidemiologic, and laboratory findings. Positive results are indicative of the presence of SARS-CoV-2 RNA. Positive results do not rule out bacterial infection or co-infection with other pathogens.

Negative results do not preclude SARS-CoV-2 infection. The results of this test should not be used as the sole basis for diagnosis and patient management decisions.

Type of Use (Select one or both, as applicable)

☒ Prescription Use (Part 21 CFR 801 Subpart D)

☐ Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

**510(k) Summary
for
Xpert Xpress CoV-2 *plus***

TABLE OF CONTENTS

1	510(K) SUMMARY	3
1.1	DEVICE DESCRIPTION.....	4
1.2	DEVICE INTENDED USE	5
1.3	SUBSTANTIAL EQUIVALENCE.....	5
1.4	PERFORMANCE STUDIES	8
1.5	CONCLUSIONS	28

**1 510(k) Summary**

As required by 21 CFR Section 807.92(c).

Submitted by: Cepheid
904 Caribbean Drive
Sunnyvale, CA 90489
Phone number: (408) 541-4191

Contact: Jennifer Motto, M.S., RAC

Date of Preparation: January 14, 2025

Device:

Trade name	Xpert® Xpress CoV-2 <i>plus</i>
Common name	Xpert Xpress CoV-2 <i>plus</i>
Type of Test	Qualitative real-time reverse transcription polymerase chain reaction (RT-PCR) and detection test

Regulation number:	21 CFR 866.3981
Classification name:	Multi-target respiratory specimen nucleic acid test including SARS-CoV-2 and other microbial agents, Respiratory Specimen Nucleic Acid SARS-CoV-2 Test,

Primary Product code	QQX
----------------------	-----

Secondary Product code	OOI
------------------------	-----

Classification Advisory Panel	Microbiology (83)
-------------------------------	-------------------

Prescription Use	Yes
------------------	-----

Predicate Device Assay:	Xpert Xpress CoV-2 <i>plus</i> (K230440)
-------------------------	--

1.1 Device Description

The Xpert Xpress CoV-2 *plus* test is a rapid, automated *in vitro* diagnostic test for the qualitative detection of viral RNA from SARS-CoV-2 in nasopharyngeal swab (NPS) and anterior nasal swab (NS) specimens collected from individuals with signs and symptoms of respiratory tract infection.

The Xpert Xpress CoV-2 *plus* test is performed on GeneXpert Xpress System, which consist of a GeneXpert IV instrument that executes sample preparation, nucleic acid amplification and real-time fluorescent signal detection for the tests, and a GeneXpert Hub with preloaded GeneXpert Xpress software for running the tests and viewing the test results. The GeneXpert Hub accessory integrates the computer, touchscreen monitor and barcode scanner. Each of the GeneXpert modules in the GeneXpert IV instrument can perform independent sample preparation and testing. The GeneXpert Xpress System requires the use of single-use disposable cartridges that hold the RT-PCR reagents and host sample purification, nucleic acid amplification, and detection of the target sequences. Because the cartridges are self-contained, cross-contamination between samples is minimized.

The Xpert Xpress CoV-2 *plus* test includes reagents for the detection of viral RNA from SARS-CoV-2 in NPS and NS specimens. The primers and probes in the Xpert Xpress CoV-2 *plus* test are designed to amplify and detect unique sequences in the genes that encode the following SARS-CoV-2 proteins: nucleocapsid (N2), envelope (E), and RNA-dependent RNA polymerase (RdRP). A Sample Processing Control (SPC) and a Probe Check Control (PCC) are also included in the cartridge serving as internal controls. The SPC is present to control for adequate processing of the sample and to monitor for the presence of potential inhibitor(s) in the RT-PCR reaction. The SPC also ensures that the RT-PCR reaction conditions (temperature and time) are appropriate for the amplification reaction and that the RT-PCR reagents are functional. The PCC verifies reagent rehydration, PCR tube filling, and confirms that all reaction components are present in the cartridge including monitoring for probe integrity and dye stability.

The Xpert Xpress CoV-2 *plus* test is designed for use with NPS or NS specimen collected with nylon flocked swabs and placed into viral transport medium (VTM), Universal Transport Medium (UTM) or eNAT®. Examples of the ancillary specimen collection kits that are compatible for use with the Xpert Xpress CoV-2 *plus* test include:

- **Nasopharyngeal Sample Collection Kit for Viruses**
 - Copan UTM® 3C057N (Flexible Minitip Flocked Swab with UTM® Medium w/o Beads)

- Becton Dickinson Universal Viral Transport Kit P/N 220531 (Flexible Minitip Flocked Swab with UVT Medium)
- Copan eNAT® Molecular Collection and Preservation Medium P/N 6U074S01 (Flexible Minitip Flocked Swab with eNAT® Medium)
- **Nasal Sample Collection Kit for Viruses**
 - Copan UTM® 3C064N (Regular Flocked Swab with UTM® Medium w/o Beads)
 - Copan eNAT® Molecular Collection and Preservation Medium P/N 6U073S01 (Regular Flocked Swab with eNAT® Medium)
- **Alternatively, swabs and transport media can be obtained separately:**
 - Nylon flocked swab (Copan P/N 502CS01, 503CS01)
 - Viral transport medium, 3 mL (Copan P/N 3C047N, BD Universal Transport Medium, Remel M4RT, Remel M5)

These sample collection kits for viruses allow NPS and NS specimens from patients to be collected, preserved and transported to laboratory prior to analysis with the Xpert Xpress CoV-2 *plus* test.

1.2 Device Intended Use

The Xpert® Xpress CoV-2 *plus* test, performed on the GeneXpert® Xpress System, is a rapid real-time RT-PCR test intended for the qualitative detection of SARS-CoV-2 RNA in nasopharyngeal and anterior nasal swab specimens collected from individuals with signs and symptoms of respiratory tract infection.

The Xpert Xpress CoV-2 *plus* test is intended for use as an aid in the diagnosis of COVID-19 if used in conjunction with other clinical, epidemiologic, and laboratory findings. Positive results are indicative of the presence of SARS-CoV-2 RNA. Positive results do not rule out bacterial infection or co-infection with other pathogens.

Negative results do not preclude SARS-CoV-2 infection. The results of this test should not be used as the sole basis for diagnosis and patient management decisions.

1.3 Substantial Equivalence

Table 1 shows the similarities and differences between the subject device and the previously cleared predicate device [K230440].

Table 1. Comparison of Similarities and Differences Between Xpert Xpress CoV-2 *plus* and the Predicate Device

Attribute	Subject Device	Predicate Device – K230440
	Xpert® Xpress CoV-2 <i>plus</i> , Performed on the GeneXpert Xpress System	Xpert® Xpress CoV-2 <i>plus</i> , Performed on the GeneXpert Instrument Systems
Regulation	Same	21 CFR 866.3981 Devices to detect and identify nucleic acid targets in respiratory samples from microbial agents that cause the SARS-CoV-2 respiratory infection and other microbial agents when in a multi-analyte test
Product Code	Same	QQX Respiratory Specimen Nucleic Acid SARS-CoV-2 Test
Device Class	Same	II (Special Controls)
Technology/ Detection	Same	Real-time reverse transcription polymerase chain reaction (RT-qPCR)
Intended Use	The Xpert Xpress CoV-2 plus test, performed on the GeneXpert Xpress System, is a rapid real-time RT-PCR test intended for the qualitative detection of SARS-CoV-2 RNA in nasopharyngeal and anterior nasal swab specimens collected from individuals with signs and symptoms of respiratory tract infection. The Xpert Xpress CoV-2 plus test is intended for use as an aid in the diagnosis of COVID-19 if used in conjunction with other clinical, epidemiologic, and laboratory findings. Positive results are indicative of the presence of SARS-CoV-2 RNA. Positive results do not rule out bacterial infection or co-infection with other pathogens. Negative results do not preclude SARS-CoV-2 infection. The results of this test should not be used as the sole basis for diagnosis and patient management decisions.	The Xpert Xpress CoV-2 plus test, performed on the GeneXpert Dx and GeneXpert Infinity Systems, is a rapid real-time RT-PCR test intended for the qualitative detection of SARS-CoV-2 RNA in nasopharyngeal and anterior nasal swab specimens collected from individuals with signs and symptoms of respiratory tract infection. The Xpert Xpress CoV-2 plus test is intended for use as an aid in the diagnosis of COVID-19 if used in conjunction with other clinical, epidemiologic, and laboratory findings. Positive results are indicative of the presence of SARS-CoV-2 RNA. Positive results do not rule out bacterial infection or co-infection with other pathogens. Negative results do not preclude SARS-CoV-2 infection. The results of this test should not be used as the sole basis for diagnosis and patient management decisions.
Assay Targets	Same	SARS-CoV-2 (E, N2, RdRP)
Specimen Type	Same	Nasopharyngeal swab (NPS) Anterior nasal swab (NS)
Transport Media	Same	Universal Transport Medium (UTM)/Viral Transport Medium (VTM)



Attribute	Subject Device	Predicate Device – K230440
	Xpert® Xpress CoV-2 <i>plus</i> , Performed on the GeneXpert Xpress System	Xpert® Xpress CoV-2 <i>plus</i> , Performed on the GeneXpert Instrument Systems
		eNAT
Test Format	Same	Single Use
Automation	Same	Automated Nucleic Acid Extraction, Detection and Results Interpretation
Assay Results	Same	Qualitative
Internal Control	Same	- Sample Processing Control (SPC) - Probe Check Control (PCC)
Time to Result	Same	~30 minutes
Instrument Systems	Cepheid GeneXpert® Xpress System	Cepheid GeneXpert® Instrument Systems

The following performance data (analytical and clinical) were provided in support of the substantial equivalence determination.

1.4 Performance Studies

1.4.1 Analytical Performance

Analytical Sensitivity/Limit of Detection (LoD) – Clinical NPS Matrix

The analytical sensitivity of the Xpert Xpress CoV-2 *plus* test was first estimated by testing limiting dilutions of a NATtrol inactivated SARS-CoV-2 virus in pooled negative clinical NPS-UTM/VTM matrix, using two reagent lots and following the guidance in Clinical and Laboratory Standards Institute (CLSI) document EP17-A2. LoD was estimated by considering each target gene (E, N2, and RdRP) in addition to the overall positivity rate for the Xpert Xpress CoV-2 *plus* test. The estimated LoD value as determined by Probit regression analysis was based on the weakest target gene (N2) and verified using two lots of Xpert Xpress CoV-2 *plus* reagents in replicates of 20 per lot for two clinical NPS matrices (UTM/VTM and eNAT). The concentration levels with observed hit rates greater than or equal to 95% in the LoD determination study were 403, 200 and 70 copies/mL for the N2 target, RdRP target and E target, respectively. The claimed LoD is 403 copies/mL (**Tables 2 and 3**).

Table 2. Xpert Xpress CoV-2 *plus* Limit of Detection – Verification of NATtrol SARS-CoV-2 LoD in NPS-UTM/VTM Matrix

Virus Strain	Reagent Lot	LoD (copies/mL)	Positive Results out of # Valid Replicates	% Positive	Mean E Ct	Mean N2 Ct	Mean RdRP Ct
NATtrol™ SARS-CoV-2	Reagent Lot 1	403	20/20	100%	34.3	38.3	36.6
	Reagent Lot 2	403	20/20	100%	34.1	37.7	36.5

Table 3. Xpert Xpress CoV-2 *plus* Limit of Detection – Verification of NATtrol SARS-CoV-2 LoD in NPS-eNAT Matrix

Virus Strain	Reagent Lot	LoD (copies/mL)	Positives out of # Valid Replicates	% Positive	Mean E Ct	Mean N2 Ct	Mean RdRP Ct
NATtrol™ SARS-CoV-2	Reagent Lot 1	403	20/20 ^a	100%	33.4	36.8	35.5
	Reagent Lot 2	403	20/20	100%	33.5	36.8	35.5

a. One of 20 replicates tested reported INVALID. The run was successfully repeated to obtain 20 valid replicates.

Analytical Sensitivity/Limit of Detection (LoD) – Clinical NS-UTM/VTM Matrix

The analytical sensitivity of the Xpert Xpress CoV-2 *plus* test was first estimated by testing limiting dilutions of a NATtrol inactivated SARS-CoV-2 virus in pooled negative clinical NS-UTM/VTM matrix, using two reagent lots and following the guidance in Clinical and Laboratory Standards Institute (CLSI) document EP17-A2. LoD was estimated by considering each target gene (E, N2, and RdRP) in addition to the overall positivity rate for the Xpert Xpress CoV-2 *plus* test. The LoD point estimates and 95% upper and lower

confidence intervals (CI) for individual SARS-CoV-2 targets (E, N2 and RdRP) in clinical NS-UTM/VTM matrix were determined using Probit regression analysis and the point estimated LoD values were 97, 462 and 157 copies/mL for the E target, N2 target and RdRP target, respectively. The estimated LoD value of the weakest target gene (N2) was verified using two lots of Xpert Xpress CoV-2 *plus* reagents in replicates of 20 per lot for clinical NS-UTM/VTM matrix. The claimed LoD is 462 copies/mL (**Table 4**).

Table 4. Xpert Xpress CoV-2 *plus* Limit of Detection – Verification of NATtrol SARS-CoV-2 LoD in NS-UTM/VTM Matrix

Virus Strain	Reagent Lot	LoD (copies/mL)	Positive Results out of # of Valid Replicates	% Positive	Mean E Ct	Mean N2 Ct	Mean RdRP Ct
NATtrol™ SARS-CoV-2	Reagent Lot 1	462	20/20	100%	34.4	38.9	36.8
	Reagent Lot 2	462	20/20	100%	34.1	37.6	36.1

Analytical Sensitivity (Limit of Detection) – WHO First International Standard SARS-CoV-2 RNA in clinical NS-UTM/VTM matrix

The analytical sensitivity of the Xpert Xpress CoV-2 *plus* test was estimated by testing limiting dilutions of WHO First International Standard for SARS-CoV-2 RNA in pooled negative clinical NS-UTM/VTM matrix, using one reagent lot. The study was performed in a random and blinded fashion.

The lowest concentration level of the First WHO International Standard for SARS-CoV-2 RNA that generated ≥95% positive results (19/20 positive results) for the weakest target gene (N2) by the Xpert Xpress CoV-2 *plus* test was used to confirm the analytical sensitivity or Limit of Detection (LoD). The verified LoD for the First WHO International Standard for SARS-CoV-2 RNA tested in NS-UTM/VTM matrix with the Xpert Xpress CoV-2 *plus* test was 1000 IU/mL. The positivity rates for the overall SARS-CoV-2 results and the gene targets, including the mean Ct values for E, N2 and RdRP at the verified LoD are presented in **Table 5**.

Table 5. LoD Concentration, Positivity Rate and Mean Ct Values for SARS-CoV-2 for the First WHO International Standard for SARS-CoV-2 RNA in Clinical NS-UTM/VTM Matrix

Virus	LoD Concentration	SARS-CoV-2 Positive Results / Test Replicates	E Target		N2 Target		RdRP Target	
			Positive Results / Test Replicates	Mean Ct	Positive Results / Test Replicates	Mean Ct	Positive Results / Test Replicates	Mean Ct
1 st WHO International Standard for SARS-CoV-2 RNA Virus	1000 IU/mL	20/20	20/20	34.3	19/20	39.5	20/20	36.5

Analytical Reactivity (Inclusivity)

SARS-CoV-2 *in silico* Analyses

The inclusivity of Xpert Xpress CoV-2 *plus* primers was initially evaluated on June 30, 2022 using *in silico* analysis of the assay amplicons in relation to 11,650,640 SARS-CoV-2 sequences available in the GISAID gene database for three targets, E, N2 and RdRP. The 11,650,640 SARS-CoV-2 sequences were separated into the lineages of interest based on the Pango Lineage assigned to each genome by GISAID, and those with ambiguous nucleotides were removed. Thus, the following inclusivity analyses focuses on the combined, non-ambiguous sequences from the variants of interest and variants of concern as of June 30, 2022. These constituted 10,469,612 sequences for the E target, 10,587,381 sequences for the N2 target and 10,333,656 sequences for the RdRP target. **Table 6** summarizes the effective predicted inclusivity for E, N2 and RdRP amplicons for the variants of interests and concern.

Table 6. Predicted Inclusivity for E, N2 and RdRP Amplicons for SARS-CoV-2 Variants of Interests and Concern

SARS-CoV-2 Target Amplicon	Exact Match	1 Mismatch ^a	2 or More Mismatches	Predicted Inclusivity
E	10,420,248 of 10,469,612 total (99.5%)	48,562 (0.5%)	802 (0.01%)	100%
N2	10,386,068 of 10,587,381 total (98.1%)	196,336 (1.9%)	4,977 (0.05%)	99.95%
RdRP	10,247,146 of 10,333,656 total (99.2%)	85,373 (0.8%)	1,137 (0.01%)	100%

a. Single-nucleotide mismatches are predicted to not impact the performance of the test.

The *in silico* inclusivity of the Xpert Xpress CoV-2 *plus* probe oligonucleotides for E, N2 and RdRP were also assessed against the top 20 most frequent matches in the GISAID EpiCoV sequence database as of June 15, 2022, which constituted 10,310,839 for the E target, 10,428,014 for the N2 target and 10,178,602 for the RdRP target. For each of the probe oligonucleotides used in the Xpert Xpress CoV-2 *plus* test, **Table 7** summarizes the number sequences as well as the corresponding percentage of sequences from this dataset with exact match, 1 mismatch/insertion, and 2 or more mismatches/insertions in the alignment.

Table 7. Predicted Inclusivity for E, N2 and RdRP Probes for SARS-CoV-2 Variants of Interests and Concern

SARS-CoV-2 Target Probe	Exact Match	1 Mismatch/Insertion ^a	2 or More Mismatches/Insertions	Predicted Inclusivity
E	10,300,688 of 10,310,839 total (99.9%)	9,853 (0.1%)	22 (0.0002%)	100%
N2	10,351,581 of 10,428,014 total (99.3%)	72,957 (0.7%)	0 (0%)	100%
RdRP	0	10,140,254 of 10,178,602 total (99.6%)	37,492 (0.4%)	99.6%

a. Single-nucleotide mismatches/insertions are predicted to not impact the performance of the test.

In addition, the *in silico* inclusivity of the Xpert Xpress CoV-2 *plus* primer/probe sequences for E, N2 and RdRP was assessed against SAR-CoV-2 sequences submitted to the GISAID EpiCoV database up to 180 days prior to and including July 31, 2023 for SARS-CoV-2 isolates circulating in the United States. There were no instances where a specific primer/probe sequence was observed with greater than 0.5% frequency of mismatch during this 180-day period. The analysis predicted that the Xpert Xpress CoV-2 *plus* test will detect all circulating variants/lineages of SARS-CoV-2 evaluated during this 180-day period.

An updated inclusivity of Xpert Xpress CoV-2 *plus* was evaluated using *in silico* analysis of the assay amplicons and probes against the 2,691,091 SARS-CoV-2 sequences available in the GISAID gene database, for variants that are currently circulating, between August 1st 2022 and September 4th 2023 for three targets: E, N2 and RdRP.

The 2,691,091 SARS-CoV-2 sequences were separated into the lineages of interest based on the Pango Lineage assigned to each genome by GISAID, and those with ambiguous nucleotides were removed. Thus, the following inclusivity analyses focuses on the combined, non-ambiguous sequences from the variants of interest and variants of concern between August 1st 2022 and September 4th, 2023. These constituted 2,519,206 sequences for the E target, 2,551,690 sequences for the N2 target and 2,476,229 sequences for the RdRP target. **Table 8** shows the number of sequences in each amplicon alignment and summarizes the effective predicted inclusivity for E, N2 and RdRP amplicons for the variants of interests and concern.

Table 8. Predicted Inclusivity for E, N2 and RdRP Amplicons for SARS-CoV-2 Variants of Interests and Concern

SARS-CoV-2 Target Amplicon	Exact Match	1 Mismatch ^a	2 or More Mismatches	Predicted Inclusivity
E	2,503,319 of 2,519,206 total (99.4%)	15,694 (0.6%)	193 (0.007%)	100%
N2	2,495,212 of 2,551,690 total (97.8%)	55,765 (2.2%)	713 (0.028%)	99.95%
RdRP	2,448,754 of 2,476,229 total (98.9%)	27,322 (1.1%)	153 (0.006%)	100%

a. Single-nucleotide mismatches are predicted not to impact the performance of the test.

The *in silico* inclusivity of the Xpert Xpress CoV-2 *plus* probe oligonucleotides for E, N2 and RdRP was also assessed against the top 20 most frequent matches in the GISAID EpiCoV sequence database between August 1st 2022 and September 4th, 2023, which constituted 2,519,206 for the E target, 2,551,690 for the N2 target and 2,476,229 for the RdRP target. For each of the probe oligonucleotides used in the Xpert Xpress CoV-2 *plus* test, **Table 9** summarizes the number sequences as well as the corresponding percentage of

sequences from this dataset with exact match, 1 mismatch/insertion, and 2 or more mismatches/insertions in the alignment.

Table 9. Predicted Inclusivity for E, N2 and RdRP Probes for SARS-CoV-2 Variants of Interests and Concern

SARS-CoV-2 Target Probe	Exact Match	1 Mismatch/Insertion ^a	2 or More Mismatches/Insertions	Predicted Inclusivity
E	2,516,153 of 2,519,206 total (99.9%)	2997 (0.1%)	0 (0.0%)	100%
N2	2,528,096 of 2,551,690 total (99.1%)	12728 (0.5%)	0 (0%)	99.6%
RdRP	0	2,465,482 of 2,476,229 total (99.6%)	10418 (0.4%)	99.6%

a. Single-nucleotide mismatches/insertions are predicted to not impact the performance of the test.

Based on the built-in redundancy of the Xpert Xpress CoV-2 *plus* test's SARS-CoV-2 amplification system (i.e., 3 independent targets, only 1 of 3 must be detected to assign a positive result), it is not anticipated that any of the evaluated sequences would be missed by the Xpert Xpress CoV-2 *plus* test.

SARS-CoV-2 Wet-Testing

In addition to the *in silico* analysis of the SARS-CoV-2 primers and probes for inclusivity, the inclusivity of the Xpert Xpress CoV-2 *plus* test was evaluated by bench testing against multiple strains of SARS-CoV-2 at levels near the analytical LoD (~3x LoD). A total of 61 strains comprised of 23 intact SARS-CoV-2 viral particles (18 inactivated viral cultures and 5 BSL-3 live viral cultures), and 38 SARS-CoV-2 *in vitro* RNA transcripts representing variant strains were tested in this study with the Xpert Xpress CoV-2 *plus* test. Three replicates were tested for each strain. All SARS-CoV-2 strains tested positive in all three replicates. Results are shown in **Table 10**.

Table 10. Analytical Reactivity (Inclusivity) of the Xpert Xpress CoV-2 *plus* Test

Count No.	SARS-CoV-2 Strain	Tested Concentration	SARS-CoV-2 Test Results	Number of Positive Results Obtained out of the Total Number of Valid Replicates		
				E	N2	RdRP
1	2019-nCoV/Italy-INMI1 ^a	5 TCID ₅₀ /mL	POS	3/3	3/3	3/3
2	England/204820464/2020 ^a	0.5 TCID ₅₀ /mL	POS ^c	3/3	3/3	3/3
3	Hong Kong/VM20001061/2020 ^a	0.25 TCID ₅₀ /mL	POS	3/3	3/3	3/3
4	South Africa/KRISP-K005325/2020 ^a	0.25 TCID ₅₀ /mL	POS	3/3	3/3	3/3

Count No.	SARS-CoV-2 Strain	Tested Concentration	SARS-CoV-2 Test Results	Number of Positive Results Obtained out of the Total Number of Valid Replicates		
				E	N2	RdRP
5	USA/CA_CDC_5574/2020 ^a	0.25 TCID ₅₀ /mL	POS	3/3	3/3	3/3
6	USA/MD-HP01542/2021 ^a	1.3e3 genome equivalents/mL	POS	3/3	3/3	3/3
7	USA/GA-EHC-2811C/2021 ^a	1.3e3 genome equivalents/mL	POS	3/3	3/3	3/3
8	NY-Wadsworth-103677-01/2020 ^a	0.15 TCID ₅₀ /mL	POS	3/3	3/3	3/3
9	NY-Wadsworth-21006055-01/2021 ^a	0.04 TCID ₅₀ /mL	POS	3/3	3/3	3/3
10	NY-Wadsworth-21025952-01/2021 ^a (Isolate 1)	0.625 TCID ₅₀ /mL	POS	3/3	3/3	3/3
11	NY-Wadsworth-21018781-01/2021 ^a (Isolate 2)	0.625 TCID ₅₀ /mL	POS	3/3	3/3	3/3
12	NY-Wadsworth-21033899-01/2021 ^a	1.25 TCID ₅₀ /mL	POS	3/3	3/3	3/3
13	NY-Wadsworth-33126-01/2020 ^a	0.15 TCID ₅₀ /mL	POS	3/3	3/3	3/3
14	USA/CA-Stanford-15_S02/2021 ^a	0.04 TCID ₅₀ /mL	POS	3/3	3/3	3/3
15	USA/PHC658/2021 ^a	25 TCID ₅₀ /mL	POS	3/3	3/3	3/3
16	Japan/TY7-503/2021 ^a	0.75 TCID ₅₀ /mL	POS	3/3	3/3	3/3
17	hCoV-19/USA/MD-HP30386/2022 ^a	27.5 copies/mL	POS	3/3	3/3	3/3
18	hCoV-19/USA/COR-22-063113/2022 ^a	19 copies/mL	POS	3/3	3/3	3/3
19	Australia/VIC01/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
20	Belgium/ULG/10004/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
21	England/205041766/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
22	England/MILK-9E05B3/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
23	England/SHEF-C05B2/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
24	France/HF2393/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
25	Hong Kong/HKU-211129-001/2021 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
26	Iceland/5/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
27	India/CT-ILSGS00361/2021 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
28	India/MH-NCCS-P1162000182735/2021 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
29	India/MH-SEQ-221 S66 R1 001/2021 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
30	Japan/Hu_DP_Kng_19-020/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
31	Japan/IC-0564/2021 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
32	Portugal/PT9543/2021 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3

Count No.	SARS-CoV-2 Strain	Tested Concentration	SARS-CoV-2 Test Results	Number of Positive Results Obtained out of the Total Number of Valid Replicates		
				E	N2	RdRP
33	South Africa/KRISP-EC-K005299/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
34	Taiwan/NTU02/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
35	USA/CA9/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
36	USA/CA-PC101P/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
37	USA/CA-CZB-12943/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
38	USA/MN2-MDH2/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
39	USA/NY-MSHSPSP-PV24650/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
40	USA/TX1/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
41	USA/WA2/2020 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
42	USA/WA-CDC-UW21061750277/2021 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
43	Wuhan-Hu-1 ^b	1.2e3 copies/mL	POS	3/3	3/3	3/3
44	Germany/BavPat1/2020 ^c	1.2e3 genome copies/mL	POS	3/3	3/3	3/3
45	Singapore/2/2020 ^c	300 genome equivalents/mL	POS	3/3	3/3	3/3
46	USA/AZ1/2020 ^c	300 genome equivalents/mL	POS	3/3	3/3	3/3
47	USA/CA1/2020 ^c	1.2e3 genome equivalents/mL	POS	3/3	3/3	3/3
48	USA/CA2/2020 ^c	1.2e3 genome equivalents/mL	POS	3/3	3/3	3/3
49	USA/CA3/2020 ^c	400 genome equivalents/mL	POS	3/3	3/3	3/3
50	USA/CA4/2020 ^c	1.2e3 genome equivalents/mL	POS	3/3	3/3	3/3
51	USA/IL1/2020 ^c	150 genome equivalents/mL	POS	3/3	3/3	3/3
52	USA/New York-PV08001/2020 ^c	1.2e3 genome copies/mL	POS	3/3	3/3	3/3
53	USA/New York-PV08410/2020 ^c	300 genome equivalents/mL	POS	3/3	3/3	3/3
54	USA/New York-PV08449/2020 ^c	300 genome equivalents/mL	POS	3/3	3/3	3/3
55	USA/New York-PV09158/2020 ^c	300 genome equivalents/mL	POS	3/3	3/3	3/3
56	USA/WI1/2020 ^c	1.2e3 genome equivalents/mL	POS	3/3	3/3	3/3
57	hCoV-19/USA/MD-HP38861/2022 ^d	0.29 TCID ₅₀ /mL	POS	3/3	3/3	3/3

Count No.	SARS-CoV-2 Strain	Tested Concentration	SARS-CoV-2 Test Results	Number of Positive Results Obtained out of the Total Number of Valid Replicates		
				E	N2	RdRP
58	hCoV-19/USA/MD-HP38960/2022 ^d	0.043 TCID ₅₀ /mL	POS	3/3	3/3	3/3
59	hCoV-19/USA/CO-CDPHE-2102544747/2021 ^d	0.028 TCID ₅₀ /mL	POS	3/3	3/3	3/3
60	hCoV-19/USA/MD-HP38288/2022 ^d	0.24 TCID ₅₀ /mL	POS	3/3	3/3	3/3
61	hCoV-19/USA/MD-HP40900/2022 ^d	0.72 TCID ₅₀ /mL	POS	3/3	3/3	3/3

- Inactivated viral culture fluid (intact viral particles)
- Synthetic, *in vitro* RNA transcript, TWIST controls
- Genomic RNA
- BSL-3 live viral culture fluid (intact viral particles)
- One of 3 replicates reported ERROR. The run was successfully repeated to obtain 3 valid replicates.

Analytical Specificity (Exclusivity)

In Silico Analyses

The analytical specificity/cross-reactivity of the Xpert Xpress CoV-2 *plus* included an evaluation of the SARS-CoV-2 test primer and probes with potentially cross-reactive microorganisms by *in silico* analysis. The analysis was conducted by mapping the primers and probes of Xpert Xpress CoV-2 *plus* individually to the microorganism sequences downloaded from the NCBI database. E primers and probes are not specific for SARS-CoV-2 and will detect Human and Bat SARS-coronavirus. No potential unintended cross reactivity with other organisms listed in **Table 11** is expected based on the *in silico* analysis.

Table 11. Microorganisms Analyzed in the *in silico* Analysis for the SARS-CoV-2 Target

Microorganisms from the Same Genetic Family	High Priority Organisms
Human coronavirus 229E	Adenovirus (e.g., C1 Ad. 71)
Human coronavirus OC43	Cytomegalovirus
Human coronavirus HKU1	Enterovirus (e.g., EV68)
Human coronavirus NL63	Epstein-Barr virus
SARS-coronavirus	Human Metapneumovirus (hMPV)
MERS-coronavirus	Influenza A & B
Bat coronavirus	Measles
	Mumps
	Parainfluenza virus 1-4
	Parechovirus
	Respiratory syncytial virus
	Rhinovirus
	<i>Bacillus anthracis</i> (Anthrax)

Microorganisms from the Same Genetic Family	High Priority Organisms
	<i>Bordetella pertussis</i>
	<i>Bordetella parapertussis</i>
	<i>Chlamydia pneumoniae</i>
	<i>Chlamydia psittaci</i>
	<i>Corynebacterium diphtheriae</i>
	<i>Coxiella burnetii</i> (Q-Fever)
	<i>Escherichia coli</i>
	<i>Fusobacterium necrophorum</i>
	<i>Haemophilus influenzae</i>
	<i>Lactobacillus</i> sp.
	<i>Legionella non-pneumophila</i>
	<i>Legionella pneumophila</i>
	<i>Leptospira</i>
	<i>Moraxella catarrhalis</i>
	<i>Mycobacterium tuberculosis</i>
	<i>Mycoplasma genitalium</i>
	<i>Mycoplasma pneumoniae</i>
	<i>Neisseria elongata</i>
	<i>Neisseria meningitidis</i>
	<i>Pneumocystis jirovecii</i> (PJP)
	<i>Pseudomonas aeruginosa</i>
	<i>Staphylococcus aureus</i>
	<i>Staphylococcus epidermidis</i>
	<i>Staphylococcus salivarius</i>
	<i>Streptococcus pneumoniae</i>
	<i>Streptococcus pyogenes</i>
	<i>Aspergillus</i> sp
	<i>Candida albicans</i>

Wet-Testing

In addition to the *in silico* analysis of the SARS-CoV-2 primers and probes for cross-reactivity, the analytical specificity of the Xpert Xpress CoV-2 plus test was evaluated by bench-testing a panel of 62 microorganisms comprising (4) human coronaviruses, (1) MERS-coronavirus, (1) SARS-coronavirus, (20) other respiratory viruses, (30) respiratory bacteria, (2) yeast strains, (3) fungal strain, and 1 human nasal wash fluid representing diverse microbial flora in the human respiratory tract.

The intact viruses were tested at concentrations of $\geq 10^5$ TCID₅₀/mL, $\geq 10^5$ CEID₅₀/mL, and $\geq 10^5$ copies/mL. Bacteria and yeast were tested at $\geq 10^6$ CFU/mL. The bacteria *Chlamydia pneumoniae* and *Mycoplasma pneumoniae* were tested at concentrations of $\geq 10^6$ IFU/mL and

$\geq 10^6$ CCU/mL, respectively. Live strains of *Coxiella burnetii*, *Lactobacillus reuteri*, and *Neisseria meningitidis* were not available, therefore genomic DNA at $\geq 10^6$ genome equivalent copies/mL was used. The *in vitro* transcribed (IVT) RNA sample was tested at $\geq 10^6$ genome equivalent/mL for human Coronavirus HKU1, and the genomic RNA sample was tested at $\geq 10^6$ genome equivalent/mL for SARS-coronavirus Urbani.

The microorganisms being evaluated for cross-reactivity were tested in groups or individually and spiked into negative simulated NPS/NS background matrix for testing. If a grouped microorganisms produced a SARS-CoV-2 positive result, then each member of the group was tested separately. If the individual non-target organism yielded a positive result, retesting was performed at lower concentrations until a concentration that no longer produce a false positive result was identified.

The results of the analytical specificity/exclusivity study demonstrate that the primer/probe sets included in the Xpert Xpress CoV-2 *plus* does not cross-react with the nucleic acids from non-intended respiratory microorganisms and phylogenetically related human coronavirus species. The one exception was the SARS-coronavirus, Urbani, which yielded the expected test result of SARS-CoV-2 POSITIVE. We expected cross-reactivity for the E gene target with the SARS-coronavirus Urbani strain which has the highly conserved E gene target shared among coronaviruses in the B lineage *Betacoronavirus*. Results are shown in **Table 12**.

Table 12. Analytical Specificity (Exclusivity) of the Xpert Xpress CoV-2 *plus* Test

Respiratory Microorganisms	Test Group	Tested Concentration	Final SARS-CoV-2 Result Call-Out	Number of Positive Results Obtained from Total Number of Replicates Tested		
				E	N2	RdRP
Human coronavirus, 229E	G1	1.1e5 TCID ₅₀ /mL	NEG	0/3	0/3	0/3
Human coronavirus, OC43		1.1e5 TCID ₅₀ /mL				
MERS-coronavirus		1.1e5 TCID ₅₀ /mL				
Human coronavirus, NL63	G2	1.1e5 TCID ₅₀ /mL	NEG	0/3	0/3	0/3
Human coronavirus, HKU1 ^a	G3	1.1e6 genome cp/mL	NEG	0/3	0/3	0/3
SARS-coronavirus, Urbani ^a	G4	1.1e6 genome cp/mL	POS	3/3	0/3	0/3
Influenza A H1N1 (pdm2009), Michigan/272/2017	G5	1.1e5 TCID ₅₀ /mL	NEG	0/3	0/3	0/3
Influenza B (Victoria Lineage), Hawaii/01/2018 (NA D197N)		1.1e5 TCID ₅₀ /mL				
RSV-A, Strain: 4/2015 Isolate #1		1.1e5 TCID ₅₀ /mL				
Adenovirus Type 1	G6	1.1e5 TCID ₅₀ /mL	NEG	0/3	0/3	0/3
Adenovirus Type 7A		1.1e5 TCID ₅₀ /mL				
Cytomegalovirus		1.1e5 TCID ₅₀ /mL				
Echovirus	G7	1.1e5 TCID ₅₀ /mL	NEG	0/3	0/3	0/3
Enterovirus, D68 strain US/KY/14-18953		1.1e5 TCID ₅₀ /mL				
Epstein Barr Virus (Human Herpes Virus 4 [Hhv-4])		1.1e5 cp/mL				
Herpes Simplex Virus (HSV) type 1		1.1e5 TCID ₅₀ /mL				
Human metapneumovirus (hMPV-5, type B1)		1.1e5 TCID ₅₀ /mL				
Measles		1.1e5 TCID ₅₀ /mL				
Mumps virus		1.1e5 TCID ₅₀ /mL				
Human parainfluenza Type 1	G8	1.1e5 TCID ₅₀ /mL	NEG	0/3	0/3	0/3
Human parainfluenza Type 2		1.1e5 TCID ₅₀ /mL				
Human parainfluenza Type 3		1.1e5 TCID ₅₀ /mL				
Human parainfluenza Type 4		1.1e5 TCID ₅₀ /mL				
Rhinovirus, Type 1A ^b		4.5e4 TCID ₅₀ /mL				
<i>Acinetobacter baumannii</i>	G9	1.1e6 CFU/mL	NEG	0/3	0/3	0/3
<i>Burkholderia cepacia</i>		1.1e6 CFU/mL				
<i>Candida albicans</i>		1.1e6 CFU/mL				

Respiratory Microorganisms	Test Group	Tested Concentration	Final SARS-CoV-2 Result Call-Out	Number of Positive Results Obtained from Total Number of Replicates Tested		
				E	N2	RdRP
<i>Candida parapsilosis</i>		1.1e6 CFU/mL				
<i>Bordetella pertussis</i>		1.1e6 CFU/mL				
<i>Chlamydia pneumoniae</i>		1.1e6 IFU/mL				
<i>Citrobacter freundii</i>		1.1e6 CFU/mL				
<i>Corynebacterium xerosis</i>	G10	1.1e6 CFU/mL	NEG	0/3	0/3	0/3
<i>Escherichia coli</i>		1.1e6 CFU/mL				
<i>Enterococcus faecalis</i>		1.1e6 CFU/mL				
<i>Hemophilus influenzae</i>		1.1e6 CFU/mL				
<i>Legionella spp.</i> ^d		1.1e6 CFU/mL				
<i>Moraxella catarrhalis</i>		1.1e6 CFU/mL				
<i>Mycobacterium tuberculosis (avirulent)</i>	G11	1.1e6 CFU/mL	NEG	0/3	0/3	0/3
<i>Mycoplasma pneumoniae</i>		1.1e6 CCU/mL				
<i>Neisseria mucosa</i>		1.1e6 CFU/mL				
<i>Propionibacterium acnes</i> (= <i>Cutibacterium acnes</i>) Z144		1.1e6 CFU/mL				
<i>Pseudomonas aeruginosa</i> , Z139		1.1e6 CFU/mL				
<i>Staphylococcus aureus</i>		1.1e6 CFU/mL				
<i>Staphylococcus epidermidis</i>	G12	1.1e6 CFU/mL	NEG	0/3	0/3	0/3
<i>Staphylococcus haemolyticus</i>		1.1e6 CFU/mL				
<i>Streptococcus agalactiae</i>		1.1e6 CFU/mL				
<i>Streptococcus pneumoniae</i>		1.1e6 CFU/mL				
<i>Streptococcus pyogenes</i>		1.1e6 CFU/mL				
<i>Streptococcus salivarius</i>		1.1e6 CFU/mL				
<i>Streptococcus sanguinis</i>		1.1e6 CFU/mL				
<i>Pneumocystis jirovecii</i> (PJP)		1.1e6 CFU/mL				
<i>Lactobacillus reuteri</i> , F275 ^e	G13	1.1e6 genome cp/mL	NEG	0/3	0/3	0/3
<i>Neisseria meningitides</i> ^e		1.1e6 genome cp/mL				
Pooled human nasal wash	G14	N/A	NEG	0/3	0/3	0/3
Influenza C (Taylor/1233/1947)	G15	1.1e5 CEID ₅₀ /mL	NEG	0/3	0/3	0/3
Rhinovirus, Type 1A ^b	G16	1.1e5 TCID ₅₀ /mL	NEG	0/3	0/3	0/3
<i>Aspergillus flavus</i>	G17	1.35e7 CFU/mL	NEG	0/3	0/3	0/3
<i>Aspergillus fumigatus</i>	G18	3.67e6 CFU/mL	NEG	0/3	0/3	0/3

Respiratory Microorganisms	Test Group	Tested Concentration	Final SARS-CoV-2 Result Call-Out	Number of Positive Results Obtained from Total Number of Replicates Tested		
				E	N2	RdRP
<i>Coxiella burnetii</i> °	G19	2.5e6 genome cp/mL	NEG	0/3	0/3	0/3
<i>Leptospira broomii</i>	G20	2.9e6 CFU/mL	NEG	0/3	0/3	0/3
<i>Parechovirus type 1</i>	G21	3.39e5 TCID ₅₀ /mL	NEG	0/3	0/3	0/3
<i>Fusobacterium necrophorum</i>	G22	6.52e6 CFU/mL	NEG	0/3	0/3	0/3
<i>Mycoplasma genitalium</i> °	G23	3.6e6 genome cp/mL	NEG	0/3	0/3	0/3

- Microorganism in the form of genomic RNA were tested in Tris-EDTA+ ([NH₄]₂SO₄) buffer using an assay definition file (ADF) without sample preparation.
- Rhinovirus, Type 1A was initially tested at 4.5 e4 TCID₅₀/mL in test group G8 using virus lot 325725. It was re-tested individually (G16) at a higher concentration of 1.1e5 TCID₅₀/mL using virus lot 326527.
- Microorganisms in the form of genomic DNA were tested in simulated NPS/NS background matrix using the ADF with full sample preparation.
- Legionella pneumophila* was tested in this study.

Microbial Interference

A microbial interference study was performed to assess the inhibitory effects of commensal microorganisms potentially encountered in upper respiratory tract specimens on the performance of the Xpert Xpress CoV-2 *plus* test. A panel of 18 commensal microorganisms, consisting of 15 viral strains and 3 bacterial strains was tested. Contrived samples consisted of SARS-CoV-2 virus seeded at 1.15x-3x LoD into simulated nasopharyngeal swab (NPS)/anterior nasal swab (NS) matrix in the presence of commensal virus and bacterial strains spiked at their respective concentrations listed in **Table 13**.

Replicates of 8 positive samples were tested with SARS-CoV-2 virus and each potential microbial interference strain combination. All 8 of 8 positive replicate samples were correctly identified as SARS-CoV-2 POSITIVE using the Xpert Xpress CoV-2 *plus* test. No interference by the listed above commensal viral or bacterial strains was reported at the concentrations tested.

Table 13: Microbial Interference Study Results

Commensal Strain	Concentration Tested	Number of Correct Test Results/ Number Tested
Adenovirus Type 1C ^a	1x10 ⁵ TCID ₅₀ /mL	8/8
Adenovirus Type 2C ^b	1x10 ⁵ TCID ₅₀ /mL	8/8 ^c
Adenovirus Type 3B ^b	1x10 ⁵ TCID ₅₀ /mL	8/8
Human Coronavirus-OC43 ^a	1x10 ⁵ copies/mL	8/8

Commensal Strain	Concentration Tested	Number of Correct Test Results/ Number Tested
Human Coronavirus-229E ^b	1x10 ⁵ TCID ₅₀ /mL	8/8
Human Coronavirus-NL63 ^b	1x10 ⁵ TCID ₅₀ /mL	8/8
Human Coronavirus-HKU1 ^b	1x10 ⁵ copies/mL	8/8
Metapneumovirus 5, Type B1 ^a	1x10 ⁵ TCID ₅₀ /mL	8/8
Parainfluenza Type 1 ^a	1x10 ⁵ TCID ₅₀ /mL	8/8
Parainfluenza Type 2 ^a	1x10 ⁵ TCID ₅₀ /mL	8/8
Parainfluenza Type 3 ^a	1x10 ⁵ TCID ₅₀ /mL	8/8
Rhinovirus Type 1A ^a	1x10 ⁵ PFU/mL	8/8
Influenza A H1N1 ^b	1x10 ⁵ CEID ₅₀ /mL	8/8
Influenza B ^b	1x10 ⁵ TCID ₅₀ /mL	8/8 ^d
Respiratory Syncytial Virus A ^b	1x10 ⁵ TCID ₅₀ /mL	8/8
<i>Haemophilus influenzae</i> ^a	1x10 ⁷ CFU/mL	8/8
<i>Staphylococcus aureus</i> ^a	1x10 ⁷ CFU/mL	8/8
<i>Staphylococcus epidermidis</i> ^a	1x10 ⁷ CFU/mL	8/8

a. These commensal strains were tested with SARS-CoV-2 at a concentration of 3x LoD.

b. These commensal strains were tested with SARS-CoV-2 at a concentration of 1.15x LoD.

c. Two of 8 replicates reported as ERROR. The runs were successfully repeated to obtain 8 valid replicates.

d. One of 8 replicates reported as ERROR. The run was successfully repeated to obtain 8 valid replicates.

Potentially Interfering Substances

Substances that could be present in the nasopharynx (or introduced during specimen collection and handling) and potentially interfere with accurate detection of SARS-CoV-2 were evaluated with direct testing on the Xpert Xpress CoV-2 *plus* test. Potentially interfering substances in the nasal passage and nasopharynx may include, but are not limited to: blood, nasal secretions or mucus, and nasal and throat medications (used to relieve congestion, nasal dryness, irritation, or asthma and allergy symptoms), as well as antibiotics and antivirals. Positive and negative samples were prepared in simulated nasopharyngeal swab (NPS)/ anterior nasal swab (NS) matrix. Negative samples (N = 8) were tested in the presence of each substance to determine the effect on the performance of the sample processing control (SPC). Positive samples (N = 8) were tested per substance with SARS-CoV-2 virus spiked at 3x LoD. The controls were samples with and without SARS-CoV-2 virus spiked at 3x LoD into simulated NPS/ NS matrix containing no potentially interfering substance.

The 23 potentially interfering substances, with active ingredients, that were evaluated are listed in **Table 14**. For substances that resulted in an INVALID test result, the concentration of the substance was reduced by half and re-tested. Interfering effects were observed for fluticasone propionate and mucin type I-S.

Table 14. Potentially Interfering Substances Tested

Substance ID	Substance/Class	Substance/Active Ingredient	Concentrations Tested
No substance	Control	Simulated NPS/NS Matrix	100% (v/v)
Afrin	Nasal Spray	Oxymetazoline (0.05%)	15% (v/v)
Albuterol Sulfate	Beta-adrenergic bronchodilator	Albuterol Sulfate (5mg/mL)	0.83 mg/mL (equivalent to 1 dose per day)
BD Universal Transport Medium	Transport Media	BD Universal Transport Medium	100% (v/v)
Blood	Blood	Blood (Human)	2% (v/v)
Copan Swab M	Transport Media	Copan Swab M	100% (v/v)
FluMist	FluMist®	Live intranasal influenza virus vaccine	6.7% (v/v)
Fluticasone Propionate Nasal Spray	Nasal corticosteroid	Fluticasone Propionate	5 µg/mL, 2.5 µg/mL
Ibuprofen	Analgesic (nonsteroidal anti-inflammatory drugs (NSAID))	Ibuprofen	21.9 mg/dL
Leukocytes	Leukocytes	Leukocytes (Human)	1 x 10 ⁶ cells/mL
Menthol	Throat lozenges, oral anesthetic, and analgesic	Benzocaine, Menthol	1.7 mg/mL
Mucin type II	Mucin	Purified Mucin protein (Porcine submaxillary gland, type II)	0.1% (w/v)
Mucin type I-S	Mucin	Purified Mucin protein (Bovine submaxillary gland, type I-S)	2.5 mg/mL, 1.25 mg/mL
Mupirocin	Antibiotic, nasal ointment	Mupirocin (2%)	10 mg/mL
Human peripheral blood mononuclear cells (PBMC)	Human peripheral blood mononuclear cells (PBMC)	Human peripheral blood mononuclear cells (PBMC)	1x10 ³ cells/µL
PHNY	Nasal Drops	Phenylephrine (1%)	15% (v/v)
Remel M4RT	Transport Media	Remel M4RT	100% (v/v)
Remel M5	Transport Media	Remel M5	100% (v/v)
Saline	Saline Nasal Spray	Sodium Chloride (0.65%)	15% (v/v)
Snuff	Tobacco	Nicotine	1% (w/v)
Tamiflu	Anti-viral drugs	Zanamivir	7.5 mg/mL
Tobramycin	Antibacterial, systemic	Tobramycin	4 µg/mL
Zicam	Nasal Gel	Luffa operculata, Galphimia glauca, Histaminum hydrochloricum Sulfur (0.05%)	15% (w/v)
Zinc	Zinc supplement	Zinc Gluconate	0.1 µg/mL

The results for the negative and positive samples are presented in **Table 15**.

The results from the SARS-CoV-2 negative samples showed that in the presence of fluticasone propionate nasal spray at 5 µg/mL and mucin type I-S at 2.5 mg/mL, 7 of 8 replicates correctly reported “**SARS-CoV-2 NEGATIVE**” test results, and 1 of 8 replicates reported “**INVALID**” test result for each substance. When the concentrations of fluticasone propionate nasal spray and mucin type I-S were reduced by half, 8/8 replicates provided valid test results and correctly reported “**SARS-CoV-2 NEGATIVE**” for each substance. No interference was observed at these lower concentrations.

The results from the SARS-CoV-2 positive samples showed that in the presence of fluticasone propionate nasal spray at 5 µg/mL, 7 of 8 replicates correctly reported “**SARS-CoV-2 POSITIVE**” test results and 1 of 8 replicate reported “**INVALID**” test result. When the concentration of the fluticasone propionate nasal spray was reduced to 2.5 µg/mL, 8/8 replicates provided valid test results and correctly reported “**SARS-CoV-2 POSITIVE**”. No interference was observed at this lower concentration.

Table 15. SARS-CoV-2 Negative Samples Tested in the Presence of Potentially Interfering Substances

Substance	Concentration Tested	Number SARS-CoV-2 Negative / Number of Valid Replicates Tested	Number SARS-CoV-2 Positive / Number of Valid Replicates Tested
SARS-CoV-2 Negative (No Substance) Control	100% (v/v)	8/8	Not Applicable
SARS-CoV-2 Positive (No Substance) Control	100% (v/v)	Not Applicable	8/8
Afrin	15% (v/v)	8/8	8/8
Albuterol Sulfate	0.83 mg/mL	8/8	8/8
BD Universal Transport Medium	100 (v/v)	8/8	8/8
Blood	2% (v/v)	8/8	8/8
Copan Swab M	100 (v/v)	8/8	8/8
FluMist®	6.7% (v/v)	8/8	8/8
Fluticasone Propionate Nasal Spray	5 µg/mL	7/8 ^a	7/8 ^a
	2.5 µg/mL	8/8	8/8
Ibuprofen	21.9 mg/dL	8/8	8/8
Leukocytes	1 x 10 ⁶ cells/mL	8/8	8/8
Menthol	1.7 mg/mL	8/8	8/8
Mucin (Type II)	0.1% (w/v)	8/8	8/8
Mucin (Type I-S)	2.5 mg/mL	7/8 ^a	8/8
	1.25 mg/mL	8/8	Not Applicable
Mupirocin	10 mg/mL	8/8	8/8
Human Peripheral Blood Mononuclear Cells (PBMC)	1 x 10 ³ cells/µL	8/8	8/8
PHNY	15% (v/v)	8/8	8/8

Substance	Concentration Tested	Number SARS-CoV-2 Negative / Number of Valid Replicates Tested	Number SARS-CoV-2 Positive / Number of Valid Replicates Tested
Remel M4RT	100% (v/v)	8/8	8/8
Remel M5	100% (v/v)	8/8	8/8
Saline	15% (v/v)	8/8	8/8
Snuff	1% (w/v)	8/8	8/8
Tamiflu	7.5 mg/mL	8/8	8/8
Tobramycin	4 µg/mL	8/8	8/8
Zicam	15% (w/v)	8/8	8/8
Zinc	0.1 µg/mL	8/8	8/8

- a. One of 8 replicates reported **INVALID** test result, indicating interference from the substance. The substance was subsequently tested with 8 replicates at half the initial concentration.

Carryover Contamination

A study was conducted to assess whether the single-use, self-contained Xpert Xpress CoV-2 *plus* cartridge prevents specimen and amplicon carryover by testing a negative sample immediately after testing of a very high positive sample in the same GeneXpert module. The negative sample used in this study consisted of simulated NPS/NS matrix and the positive sample consisted of high SARS-CoV-2 virus concentration (inactivated SARS-CoV-2 USA-WA1/2020 at 5e4 copies/mL) spiked into negative NPS/NS matrix. The negative sample was tested in a GeneXpert module at the start of the study. Following the initial testing of the negative sample, the high positive sample was processed in the same GeneXpert module immediately followed by another negative sample. This was repeated 20 times in the same module, resulting in 20 positives and 21 negatives for the module. The study was repeated using a second GeneXpert module for a total of 40 positive and 42 negative samples. All 40 positive samples were correctly reported as SARS-CoV-2 POSITIVE and all 42 negative samples were correctly reported as SARS-CoV-2 NEGATIVE with the Xpert Xpress CoV-2 *plus* test. No specimen or amplicon carry-over contamination was observed in this study.

Reproducibility

The reproducibility of the Xpert Xpress CoV-2 *plus* test on the GeneXpert Xpress System in a CLIA-Waived (CW) environment was established at three (3) sites using a 3-member panel including one negative sample, one low positive (~1.5x LoD) sample and one moderate positive (~3x LoD) sample. The negative sample consisted of simulated matrix without target microorganism or target RNA. The positive samples were contrived samples in a simulated matrix using inactivated NATtrol SARS-CoV-2 (ZeptoMetrix).

Testing was conducted over five (5) independent, non-consecutive days, using one (1) lot of Xpert Xpress CoV-2 *plus* cartridges at three (3) participating CW sites each with three (3) untrained/inexperienced operators to yield a total of 90 observations per panel member (3

Sites x 3 Operators x 1 Lot x 5 Days x 1 Run x 2 Replicates = 90 observations/panel member).

The percent agreement of the qualitative results for SARS-CoV-2 detection for each panel member analyzed by each of the 9 operators and each site is shown in **Table 16**. In addition, the overall percent agreement for each sample (% total agreement) and the two-sided Wilson Score confidence intervals (CI) are presented in the last column.

The results from the study are summarized in **Table 16**.

Table 14. Summary of the Reproducibility Results - % Agreement

Sample	Site 1				Site 2				Site 3				%Total Agreement and 95% CI by Panel Member
	Op1	Op2	Op3	Site	Op1	Op2	Op3	Site	Op1	Op2	Op3	Site	
Negative	100% (10/10)	100% (10/10)	100% (10/10)	100% (30/30)	100% (10/10)	100% (10/10)	100% (10/10)	100% (30/30)	100% (10/10)	100% (10/10)	100% (10/10)	100% (30/30)	100% (90/90) [95.9 % - 100%]
SARS-CoV-2 Low Pos	100% (10/10)	100% (10/10)	100% (10/10)	100% (30/30)	100% (10/10)	100% (10/10)	100% (10/10)	100% (30/30)	100% (10/10)	100% (10/10)	100% (10/10)	100% (30/30)	100% (90/90) [95.9% - 100%]
SARS-CoV-2 Mod Pos	100% (10/10)	100% (10/10)	100% (10/10)	100% (30/30)	100% (10/10)	100% (10/10)	100% (10/10)	100% (30/30)	100% (10/10)	100% (10/10)	100% (10/10)	100% (30/30)	100% (90/90) [95.9% - 100%]

Abbreviation: Op, operator; Pos, positive; Mod, moderate; CI, confidence interval

The evaluation of reproducibility the underlying analyte Ct values for E, N2 and RdRP obtained in the Xpert Xpress CoV-2 *plus* test was analyzed using nested Analysis of Variance (ANOVA). The mean, standard deviation (SD), and coefficient of variation (CV) between-sites, between-operators, between-lots, between-days, between-runs and within-run for each panel member are presented in **Table 17**.

Table 17. Summary of Nested ANOVA by Coefficient of Variation*

Sample	Analyte	Non-missing Value	Mean Ct	Site		Operator		Day		Error		Total	
				SD	CV	SD	CV	SD	CV	SD	CV	SD	CV
Negative	SPC	90	28.76	0.17	0.6	0.16	0.6	0.00	0.0	0.41	1.4	0.47	1.6
SARS-CoV-2 Low Pos	E	90	34.43	0.00	0.0	0.27	0.8	0.17	0.5	0.55	1.6	0.63	1.8
	N2	90	38.03	0.17	0.5	0.29	0.8	0.15	0.4	0.65	1.7	0.73	1.9
	RdRp	90	37.63	0.12	0.3	0.28	0.7	0.32	0.9	0.61	1.6	0.75	2.0
SARS-CoV-2 Mod Pos	E	90	32.19	0.04	0.1	0.21	0.7	0.24	0.7	0.36	1.1	0.49	1.5
	N2	90	35.84	0.12	0.3	0.23	0.7	0.15	0.4	0.45	1.3	0.55	1.5
	RdRp	90	35.48	0.12	0.4	0.26	0.7	0.18	0.5	0.42	1.2	0.54	1.5

*Coefficient of variation = SD/Mean Ct * 100%

1.4.2 Clinical Performance

The clinical performance of the Xpert Xpress CoV-2 plus test was evaluated in a multi-site, observational and method comparison study that included 22 geographically diverse CLIA-waived sites in the United States (US) using specimens collected from individuals showing signs and symptoms of respiratory infection. Of the 22 sites, 21 performed Xpert testing and specimen collection, and 1 site performed comparator and discrepant testing.

Specimens tested included prospective clinical NPS and NS specimens collected in UTM/VTM. Available demographic data from the individuals from whom the specimens were collected are presented in **Table 18**.

Table 18. Demographic Summary

Prospectively Collected Fresh Specimens from 2022	NPS (N=961)	NS (N=973)	Overall (N=1934)
Gender			
Female	569 (59.2%)	616 (63.3%)	1185 (61.3%)
Male	392 (40.8%)	357 (36.7%)	749 (38.7%)
Age Group (Years)			
≤5	8 (0.8%)	54 (5.5%)	62 (3.2%)
6-21	275 (28.6%)	285 (29.3%)	560 (29.0%)
22-59	537 (55.9%)	505 (51.9%)	1042 (53.9%)
≥60	141 (14.7%)	129 (13.3%)	270 (14.0%)
Race			

Prospectively Collected Fresh Specimens from 2022	NPS (N=961)	NS (N=973)	Overall (N=1934)
American Indian or Alaska Native	0 (0%)	1 (0.1%)	1 (0.1%)
Asian	27 (2.8%)	25 (2.6%)	52 (2.7%)
Black or African American	225 (23.4%)	217 (22.3%)	442 (22.9%)
White	656 (68.3%)	661 (67.9%)	1317 (68.1%)
Black or African American, White	4 (0.4%)	4 (0.4%)	8 (0.4%)
Other Mixed (each N≤5)	3 (0.3%)	2 (0.2%)	5 (0.3%)
Participant Declined to Answer or Unknown	46 (4.8%)	63 (6.5%)	109 (5.6%)
Ethnicity			
Hispanic	75 (7.8%)	81 (8.3%)	156 (8.1%)
Non-Hispanic	861 (89.6%)	857 (88.1%)	1718 (88.8%)
Participant Declined to Answer or Unknown	25 (2.6%)	35 (3.6%)	60 (3.1%)
COVID-19 Vaccination Status			
Vaccinated	721 (75.0%)	717 (73.7%)	1438 (74.4%)
Not Vaccinated	222 (23.1%)	239 (24.6%)	461 (23.8%)
Unknown	18 (1.9%)	17 (1.7%)	35 (1.8%)

Positive Percent Agreement (PPA) and Negative Percent Agreement (NPA) were determined when specimens were tested using Xpert Xpress CoV-2 *plus* side-by-side with a U.S. FDA-cleared molecular respiratory panel that included SARS-CoV-2 in a randomized and blinded fashion. Discrepant results between Xpert Xpress CoV-2 *plus* and the comparator were investigated using a U.S. FDA EUA SARS-CoV-2 molecular test. Additionally, non-determinate rate was determined for Xpert Xpress CoV-2 *plus*. Based on the specimens that yielded non-determinate results with Xpert Xpress CoV-2 *plus*, the non-determinate rate was 0.7% (7/961) with NPS and 0.4% (4/973) with NS specimens.

A total of 1783 specimens, including 883 NPS and 900 NS specimens that yielded valid results by both the Xpert Xpress CoV-2 *plus* and the U.S. FDA-cleared molecular respiratory panel were included in the clinical performance evaluation. The Xpert Xpress CoV-2 *plus* demonstrated a PPA and NPA of 98.2% and 99.1%, respectively in NPS specimens and 99.0% and 99.1%, respectively in NS specimens (**Table 19**).

Table 19. Xpert Xpress CoV-2 *plus* Performance Results

		Comparator								
		Overall			NPS			NS		
		Positive	Negative	Total	Positive	Negative	Total	Positive	Negative	Total
Xpert Xpress CoV-2 <i>plus</i>	Positive	215	14	229	112	7 ^b	119	103	7 ^d	110
	Negative	3	1551	1554	2 ^a	762	764	1 ^c	789	790
	Total	218	1565	1783	114	769	883	104	796	900
PPA (95% CI)		98.6% (96.0% - 99.5%)			98.2% (93.8% - 99.5%)			99.0% (94.8% - 99.8%)		
NPA (95% CI)		99.1% (98.5% - 99.5%)			99.1% (98.1% - 99.6%)			99.1% (98.2% - 99.6%)		
NPS: nasopharyngeal swab; NS: anterior nasal swab; PPA: Positive Percent Agreement; NPA: Negative Percent Agreement; CI: Confidence Interval										
a. Discrepant test results based on a U.S. FDA EUA SARS-CoV-2 molecular test: 1/2 SARS-CoV-2 positive; 1/2 SARS-CoV-2 negative										
b. Discrepant test results based on a U.S. FDA EUA SARS-CoV-2 molecular test: 7/7 SARS-CoV-2 negative										
c. Discrepant test results based on a U.S. FDA SARS-CoV-2 molecular test: 1/1 SARS-CoV-2 negative										
d. Discrepant test results based on a U.S. FDA SARS-CoV-2 molecular test: 1/7 SARS-CoV-2 positive; 5/7 SARS-CoV-2 negative; 1/7 invalid results										

1.5 Conclusions

The results of the analytical and clinical performance studies summarized above demonstrate that the Xpert Xpress CoV-2 *plus* test is substantially equivalent to the predicate device.