



June 5, 2025

Jiangsu Yuyue Medical Equipment & Supply Co., LTD.
% Yijie You
Manager
Qimmiq Medical Consulting Service Co., Ltd.
RM.406, Building C, Run Science Park, No.18 Shenzhou Road,
Huangpu
Guangzhou, Guangdong 510663
China

Re: K242736

Trade/Device Name: Portable Oxygen Concentrator (Spirit-3)
Regulation Number: 21 CFR 868.5440
Regulation Name: Portable Oxygen Generator
Regulatory Class: Class II
Product Code: CAW
Dated: April 1, 2025
Received: April 1, 2025

Dear Yijie You:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device" (<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

All medical devices, including Class I and unclassified devices and combination product device constituent parts are required to be in compliance with the final Unique Device Identification System rule ("UDI Rule"). The UDI Rule requires, among other things, that a device bear a unique device identifier (UDI) on its label and package (21 CFR 801.20(a)) unless an exception or alternative applies (21 CFR 801.20(b)) and that the dates on the device label be formatted in accordance with 21 CFR 801.18. The UDI Rule (21 CFR 830.300(a) and 830.320(b)) also requires that certain information be submitted to the Global Unique Device Identification Database (GUDID) (21 CFR Part 830 Subpart E). For additional information on these requirements, please see the UDI System webpage at <https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/unique-device-identification-system-udi-system>.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory->

[assistance/contact-us-division-industry-and-consumer-education-dice](#)) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Bradley Quinn
Assistant Director
DHT1C: Division of Anesthesia,
Respiratory, and Sleep Devices
OHT1: Office of Ophthalmic, Anesthesia,
Respiratory, ENT, and Dental Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

Indications for Use

Submission Number (if known)

K242736/S001

Device Name

Portable Oxygen Concentrator (Spirit-3)

Indications for Use (Describe)

The Portable Oxygen Concentrator, model: Spirit-3, is intended to provide supplemental oxygen in a home, institutional, or travel environment. And the device is used with adult only, not used with pediatrics, infant, or neonate patients, etc.

Type of Use (Select one or both, as applicable)



Prescription Use (Part 21 CFR 801 Subpart D)



Over-The-Counter Use (21 CFR 801 Subpart C)

CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services
Food and Drug Administration
Office of Chief Information Officer
Paperwork Reduction Act (PRA) Staff
PRASaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

510(K) Summary

1. 510(k) owner (Applicant)

Establishment Registration number:		3017825461
Name:		Jiangsu Yuyue Medical Equipment & Supply Co., LTD.
Address:		No.1 Baisheng Road Development Zone, Danyang, Jiangsu 212300 CHINA
Contact Person	Name:	Wei Yan
	Address:	No.1 Baisheng Road Development Zone, Danyang, Jiangsu 212300 CHINA
	TEL:	86-511-86900833
	Email:	weiyang@yuyue.com.cn

2. Contact Person (Authorized representative)

Name:	You Yijie
Address:	RM.406, Building C, Run Science Park, No.18 Shen Zhou Road, Huangpu, Guangzhou, Guangdong 510663 P.R. China
TEL:	(+86)020-82245821
FAX:	(+86)020-82245821
Email:	jet.you@qimmiq-med.com

Date prepared: Nov. 13, 2024

3. Device Information

Device Common Name:	Generator, Oxygen, Portable
Trade Name:	Portable Oxygen Concentrator
Model:	Spirit-3
Classification name:	Portable oxygen generator
Review panel:	Anesthesiology

Product code:	CAW
Regulation Class:	II
Regulation Number:	868.5540

4. Predicate Device and reference device Information

Predicate device:

510(k) submitter/holder:	Qingdao Kingon Medical Science and Technology Co., Ltd
510(K) Number:	K230702
Trade Name:	Portable Oxygen Concentrator
Model:	P2-K4, P2-S4
Classification name:	Generator, Oxygen, Portable
Review panel:	Anesthesiology
Product code:	CAW
Regulation Class:	II
Regulation Number:	868.5440

Reference device

510(k) submitter/holder:	Inogen, Inc.
510(K) Number:	K222086
Trade Name:	Inogen Rove 4 Portable Oxygen Concentrator
Model:	Inogen Rove 4
Classification name:	Generator, Oxygen, Portable
Review panel:	Anesthesiology
Product code:	CAW
Regulation Class:	II
Regulation Number:	868.5440

5. Device description

Portable Oxygen Concentrator, model: Spirit-3, is a portable oxygen generator that is intended to release oxygen for respiratory therapy by means of physical means (a molecular sieve). It supplies a pulsed high concentration of oxygen and is used with a nasal cannula to deliver oxygen from the concentrator to the patient. The Portable Oxygen Concentrator is small, portable and may be used in home, institutional, or travel environment.

Portable Oxygen Concentrator, model: Spirit-3, is capable of continuous use in a home, institutional, or travel environment. Power options include 100 – 240 V (50/ 60Hz) AC power supply, DC power supply and rechargeable lithium-ion battery.

The portable oxygen concentrator consists of two parts: an oxygen concentrator and accessories.

6. Indications for Use

The Portable Oxygen Concentrator, model: Spirit-3, is intended to provide supplemental oxygen in a home, institutional, or travel environment. And the device is used with adult only, not used with pediatrics, infant, or neonate patients, etc.

7. Comparison to predicate device

Comparison Analysis: The subject device is substantially equivalent to the predicate devices, Portable Oxygen Concentrator, Model: P2-K4, P2-S4 marketed by Qingdao Kingon Medical Science and Technology Co., Ltd, the substantial equivalence chart is provided as follows:

SE Comparisons	Subject device (Portable Oxygen Concentrator, model: Spirit-3)	Predicate device (Kingon Portable Oxygen Concentrator, Model: P2-K4, P2-S4)	Reference device (Inogen Rove™ 4 Portable Oxygen Concentrator)	Discussion of difference
510K Number	/	K230702	K222086	/
Classification	21CFR 868.5440	21CFR 868.5440	21CFR 868.5440	Same
Product Code	CAW	CAW	CAW	Same
FDA Class	II	II	II	Same
Indications for Use	The Portable Oxygen Concentrator (Model: Spirit-3) is intended to provide supplemental Oxygen in a home, institutional, or Travel environment. And the device is used with adult only, not used with pediatrics, infant, or	The Portable Oxygen Concentrator (Model: P2-K4, P2-S4) is intended to provide supplemental Oxygen in a home, institutional, or Travel environment.	The Inogen Rove 4 Portable Oxygen Concentrator provides a high concentration of supplemental oxygen to patients requiring respiratory therapy on a prescriptive basis. It may be used in home,	Same with the predicate device which is also used with adult only.

	neonate patients, etc.		institution and transport modalities.	
Model	Spirit-3	P2-K4, P2-S4	Inogen Rove 4	/
Environment of Use	Home, institutional, or travel environment.	home, institutional, or travel environment.	home, institution and transport modalities	Same with the predicate device
Patient Population	Adult patient only	Adult	Adult patient only	Same
Material of Patient contact components	Top cover: PC/PMMA Coextruded film Housing shell: ABS	Wiring cover: PC+ABS Intake hood: PC+ABS Nozzle fitting: Aluminum alloy Button panel: PET Main housing: PC+ABS	Not public	Different The material of patient contact components of subject device Spirit -3 is different from predicate device P2-K4 and P2-S4 (K230702), which affect safety of biocompatibility, while not affect effectiveness. Material of subject device has been implemented biocompatibility evaluation in agreement with recommendations according to FDA Guidance-Use of International Standard ISO 10993-1, "Biological evaluation

				of medical devices - Part 1: Evaluation and testing within a risk management process", and Cytotoxicity, Sensitization and Irritation meet the requirement of ISO 10993- 5, ISO 10993-10 and ISO 10993-23 respectively. Therefore, the difference of subject device with predicate device do not raise new questions of safety and effectiveness.
Duration and type of contact	Type of contact: surface device; Duration: permanent (> 30 d);	Type of contact: surface device; Duration: permanent (> 30 d);	Externally Communicating, Tissue, Permanent Duration (>30 days)	Same with the predicate device
Complete list of all the biocompatibility tests performed	ISO 10993- 5 tested for Cytotoxicity; ISO 10993-10 tested for Sensitization; ISO 10993-23 tested for Irritation; ISO 18562-2 tested for Particulate matter; ISO 18562-3 tested for Volatile	ISO 10993- 5 tested for Cytotoxicity; ISO 10993-10 tested for Sensitization; ISO 10993-23 tested for Irritation; ISO 18562-2 tested for Particulate matter; ISO 18562-3 tested for Volatile organic Compounds;	ISO 18562-2 tested for Particulate matter; ISO 18562-3 tested for Volatile organic Compounds;	Same with the predicate device

	organic Compounds;			
Single Patient, multi-use	Yes	Yes	Yes	Same
Patient Interface	Cannula Port	Cannula Port	Cannula Port	Same
Dimensions	22 × 8 × 22 (cm)	17×8.6×15.7(cm)	With 4-cell battery: 6 x 2.7 x 7.5 in (15.1 x 6.9 x 19.1 cm) With 8-cell battery: 6 x 2.7 x 7.8 in (15.1 x 6.9 x 19.7 cm)	Similar Dimension of subject device Spirit-3(22 × 8 × 22 (cm)) is a little larger than and similar with that of predicate device P2-K4 and P2-S4 (K230702)(17×8.6×15.7(cm)), does not raise new questions of safety and effectiveness.
Weight	2.4 kg (with battery)	1.5 kg (with battery)	With 4-cell battery: 3.0 pounds (1.4kg) With 8-cell battery: 3.3 pounds (1.5kg)	Similar The weight of subject device Spirit-3(2.4Kg) is similar with that of predicate device P2-K4 and P2-S4 (K230702)(1.5kg), does not raise new questions of safety and effectiveness.
Startup time	2 minutes	2 minutes	2 minutes	Same

Oxygen Concentration	90%-3%/+6% at all settings	90%-3%/+6% at all settings	90% - 3%/+6% at all settings	Same
Setting	adjustable in 1 increments from 1 to 4	adjustable in 1 increment from 1 to 4	adjustable in 1 increment from 1 to 4	Same
Pulse mode bolus size	37.5mL per breath at setting 4 with 20BPM	42 mL per breath at setting 4 with 20 BPM	42 mL per breath at setting 4 with 20 BPM	Different The Pulse mode bolus size of subject device Spirit-3(37.5ml at setting 4 with 20BPM) is a little smaller than 42.0 ml at setting 4 with 20BPM of predicate device P2-K4 and P2-S4 (K230702), but is covered by the subject device, therefore the difference of subject device with predicate device do not raise new questions of safety and effectiveness.
Principle of operation	by means of molecular sieve	by means of molecular sieve	by means of molecular sieve	Same
Type of Protection Against Electrical Shock	Class II	Class II	Not public	Same with the predicate device
Degree of Protection	Type BF	Type BF	Not public	Same with the predicate device

to Concentrator Components Against Electrical Shock:				
Degree of Protection to Concentrator Components Against Ingress of Water	IP22	IP22	Not public	Same with the predicate device
Rated breath rate	15 - 40 Breath per minute	10 - 40 Breath per minute	10 - 40 Breath per minute	Different Breath rate range of the predicate device P2-K4 and P2-S4 is 10 to 40 BPM, and subject device Spirit-3 is 15 to 40 BPM, the minimum Breath rate rang for ISO80601-2-67:2020 Clause 201.12.1.101 requirement to disclosure is 15BPM, and it is also covered by the predicate device range, therefore the difference of subject device with predicate device do not raise new questions of safety and effectiveness.

Work System	Work continuously	Work continuously	Work continuously	Same
User Interface	Buttons, LCD Display	Buttons, LCD Display	Buttons, LCD Display	Same
Power requirements	AC power supply: 100-240V AC, 50-60Hz, 2.0A; 12.0VDC 10.0A out DC Power Supply: DC 12-16V, 10A Battery: DC 14.4V	AC power supply: 100-240VAC ;50-60 Hz, 19VDC 5.26A out DC power supply: 12 - 16V DC input, 19V 6A DC output Battery: not available	AC power supply: 100 to 240 VAC, 50 to 60 Hz Autosensing 2.0 – 1.0A DC power supply: 13.5-15.0 VDC, 100W Max voltage: 12.0 to 16.8 VDC (+ 0.5) Battery: not available	Similar with the predicate device The subject device Spirit-3 has same AC power supply, DC power supply and battery with those of predicate device. And the subject device Spirit-3 is complied with IEC 60601-1, IEC 60601-1-2, ISO 80601-2-67 and ISO 80601-2-69. So the difference of subject device Spirit-3 with predicate device do not raise new questions of safety and effectiveness.
Maximum oxygen discharge pressure	150kPa(21.8psi)	110kPa(16.0psi)	151.68kPa(22psi)	Different from the predicate device, but similar with the reference device The Maximum oxygen discharge pressure of subject device Spirit-3 is 150kPa(21.8psi), is higher than 110Kpa(16.0psi) of predicate device P2-K4 and P2-S4

				(K230702), but close to 22 PSI of reference device Inogen Rove 4 (K222086). Therefore, the maximum oxygen discharge pressure of difference does not raise new questions of safety and effectiveness.
Inspiratory trigger sensitivity	-0.12 cmH2O	-0.12cmH20	-0.12cmH20	Same
Software	Embedded	Embedded	Embedded	Same
Acoustic Noise	55 dBA at 0.75 LPM	55.3dBA at 0.84LPM	Not public	Similar The acoustic Noise of subject device Spirit-3 is 55 dBA at 0.75 LPM, is very close to 55.3dBA at 0.84LPM of predicate device P2-K4 and P2-S4 (K230702). Therefore, the acoustic noise of subject device Spirit-3 does not raise new questions of safety and effectiveness.

Alarms	Empty battery Pressure failure including high pressure and low pressure No flow System over temperature Battery over temperature Compressor Failure Low oxygen concentration Low battery No breath detected Start-up failure	Battery empty Low pressure No pulse High Temp Compressor Failure Fan Failure Low Flow Low Battery No Breath Detected EEPROM Failure	Not public	Different The subject device Spirit-3 and predicate device P2-K4 and P2-S4 have difference in alarms, the predicate device has fan Failure and EEPROM Failure while the subject device has no such alarm for fan Failure or EEPROM Failure. If the fan failure of the subject causes the system temperature to rise, and then the system over temperature alarm is active, and they have the same alarms related to the basic function of safety and performance, such as pressure, flow, system temperature, oxygen concentration, compressor, battery, breath detection, which meet the requirement of clause 201.12.4.101 and 201.12.4.102 of ISO 80601-2-67: 2020, clause 201.11.8.101.2, 201.12.4.102 and 201.13.2.101 of ISO 80601-2-69: 2020.
--------	---	--	------------	--

				Hence, the difference between the subject device and the predicate device does not raise new issues of safety and effectiveness.
Status Indicator	Flow rates	Flow rates	Flow rates	Same
	Battery Condition	Battery Condition	Battery Condition	Same
	Alarms	Alarms	Alarms	Same
Battery Duration	5 hours at 0.21LPM	BA-K201: Up to 2.5 hours at 0.21 LPM; BA-K200/BA-S200: Up to 4.8hours at 0.21LPM; BA-K202/BA-S202: Up to 7.2 hours at 0.21LPM	Not public	Similar Battery Duration at 0.21 LPM of subject device Spirit-3 is about 5 hours, which is very close to battery BA-K200/BA-S200 equipped with predicate device(4.8 hours), and battery of subject device Spirit-3 meets IEC 62133-2:2017/AMD1:2021 requirement, the battery duration does not raise any issues in safety and effectiveness.

Operating Condition	Temperature: 10°C to 35°C Humidity: 15% to 75%, non-condensing Atmosphere pressure: 70kPa to 106 kPa Altitude: 0 to 3000 meters	Temperature: 5 to 40°C Humidity: 10% to 90%, non-condensing Atmosphere pressure: 70 to 106 kPa Altitude: 0 to 3048 meters	Temperature: 41 to 104°F (5 to 40°C) Humidity: 15% to 90%, non-condensing Atmosphere pressure: 0 to 10,000 ft (0 to 3048 meters)	Different The range of operating temperature, relative humidity, and altitude of subject device Spirit-3 is a little narrower than those of predicate device P2-K4 and P2-S4. The performance of subject device Spirit-3 under the operating ambient is in compliance with ISO 80601-2-69 requirement, and the operating ambient is described in labeling of subject device Spirit-3. So the difference doesn't raise new questions of safety and effectiveness.
Storage and Transportation Condition	Temperature: -20 to +60°C Humidity: ≤93%, non-condensing Atmosphere pressure: 70 to 106 kPa	Temperature: -20 to 70°C Humidity: 5% to 90%, non-condensing Atmosphere pressure: 70 to 106 kPa	Temperature: -13 to 158°F (-25 to 70°C) Humidity: 5% to 95%, non-condensing Store in a dry environment. Atmosphere pressure: Not public	Different The range of temperature (-20 to +60°C) of subject device Spirit-3 is a little narrower than predicate device P2-K4 and P2-S4(-20 to 70°C); The range of relative humidity (≤93%) of subject device Spirit-3 is a little wider than predicate

				device P2-K4 and P2-S4(5% to 90%); The range of atmospheric pressure (70 to 106 kPa) of subject device Spirit-3 is same with predicate device P2-K4 and P2-S4 (70 to 106 kPa); The Storage and transport of subject device Spirit-3 under the ambient is in compliance with ASTM D4169-22, and the Storage and transport ambient is described in labeling of subject device. So the difference doesn't raise new questions of safety and effectiveness.
Electrical Safety	IEC 60601-1	IEC 60601-1	IEC 60601-1	Same
Electromagnetic Compatibility	IEC 60601-1-2	IEC 60601-1-2	IEC 60601-1-2	Same
Standards Met	ANSI AAMI ES60601-1:2005/(R)2012 & A1:2012, C1:2009/(R)2012 & A2:2010/(R)2012 (Cons. Text) [Incl. AMD2:2021]	ANSI AAMI ES60601-1: 2005 /A1:2012 and A2:2020 IEC 60601-1-2: 2014 IEC 60601-1-11: 2015/A1:2020	IEC 60601-1:2012 IEC 60601-1-2: 2020 IEC 60601-1-6:2020	Different While both the subject and the predicate devices have met IEC 60601-1-2 and IEC 62133-2

	IEC 60601-1-2:2014 +A1:2020 IEC 60601-1-11:2015 +A1:2020 IEC 60601-1-8:2006 +A1:2012 +A2:2020 IEC60601-1-6:2010 +A1:2013 + A2:2020 ISO 80601-2-69: 2020 ISO 80601-2-67: 2020 ISO 18562-2: 2017 ISO 18562-3: 2017 IEC 62133-2: 2017/AMD1:2021 ISO 10993-5:2009 ISO 10993-10:2021 ISO 10993-23:2021	IEC 60601-1-8: 2006/ A1:2012 /A2:2020 IEC 60601-1-6: 2010/ A1:2013 /A2:2020 ISO 80601-2-69: 2020 ISO 80601-2-67: 2020 ISO 18562-2: 2017 ISO 18562-3: 2017 IEC 62133: 2017 ISO 10993-5:2009 ISO 10993-10:2021 ISO 10993-23:2021	IEC 60601-1-8:2012 IEC 60601-1-11:2015 ISO 80601-2-69:2020	requirements While both the subject and the predicate devices have met IEC 60601-1-2 and IEC 62133-2 requirements, the subject device Spirit-3 satisfied the requirements of the latest version of standards, therefore the difference between subject device Spirit-3 and predicate device does not raise new questions of safety and effectiveness.
--	---	---	--	--

8. Discussion of Non-Clinical Tests Performed for Safety and effectiveness are as follows

The device performed performance testing to demonstrate and support safety and effectiveness when compared to the predicate and the applicable standards.

Standards	Standards Name
ANSI AAMI ES60601-1:2005/(R)2012 & A1:2012, C1:2009/(R)2012 & A2:2010/(R)2012 (Cons. Text) [Incl. AMD2:2021]	Medical Electrical Equipment - Part 1: General Requirements For Basic Safety And Essential Performance
IEC 60601-1-2 Edition 4.1 2020-09 CONSOLIDATED VERSION	Medical Electrical Equipment -- Part 1-2: General Requirements For Basic Safety And Essential Performance -- Collateral Standard: Electromagnetic Disturbances -- Requirements And Tests
IEC TS 60601-4-2:2024	Medical electrical equipment -Part 4-2: Guidance and interpretation-Electromagnetic immunity: performance of medical electrical equipment and medical electrical systems
CISPR 25: 2021	Radio disturbance characteristics for the protection of receivers used on board vehicles, boats, and on devices - Limits and methods of measurement
IEC 60601-1-11 Edition 2.1 2020-07 CONSOLIDATED VERSION	Medical Electrical Equipment - Part 1-2: General Requirements For Basic Safety And Essential Performance - Collateral Standard: Electromagnetic Disturbances - Requirements And Tests
IEC 60601-1-8 Edition 2.2 2020-07 CONSOLIDATED VERSION	Medical electrical equipment - Part 1-8: General requirements for basic safety and essential performance - Collateral Standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems
IEC 60601-1-6:2005/AMD2:2020	Medical electrical equipment - Part 1-6:General requirements for basic safety and essential performance -Collateral standard: Usability
ISO 80601-2-69: 2020	Medical electrical equipment. Particular requirements for the basic safety and essential performance of oxygen

	concentrator equipment
ISO 80601-2-67: 2020	Medical electrical equipment. Particular requirements for basic safety and essential performance of oxygen-conserving equipment
ISO 18562-2: 2017	Biocompatibility evaluation of breathing gas pathways in healthcare applications - Part 3: Tests for emissions of volatile organic compounds
ISO 18562-3: 2017	Biocompatibility evaluation of breathing gas pathways in healthcare applications. Tests for emissions of volatile organic compounds (VOCs)
IEC 62133-2: 2017/AMD1:2021	Secondary cells and batteries containing alkaline or other non-acid electrolytes - Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications - Part 2: Lithium systems
ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
ISO 10993-10:2021	Biological evaluation of medical devices - Part 10: Tests for skin sensitization
ISO 10993-23:2021	Biological evaluation of medical devices - Part 23: Tests for irritation
IEC 62304:2006+A1:2015	Medical device software - Software life cycle processes

● **Electrical Safety and Electromagnetic Compatibility Testing**

The Device was tested and complied with the applicable requirements of the following standards: ANSI AAMI ES60601-1, IEC 60601-1-2, IEC 60601-1-11, IEC 60601-1-8, ISO 80601-2-69, ISO 80601-2-67, IEC TS 60601-4-2:2024, CISPR 25: 2021 demonstrated that the basic safety and performance of the device met the requirements.

● **Biocompatibility Testing**

According to ISO 10993-1:2018 and ISO 18562-1:2017, we performed the biocompatibility evaluation and test of skin surface direct-contacting components:

- (1) In vitro Cytotoxicity test according to ISO 10993-5:2009

- (2) Skin Sensitization test according to ISO 10993-10:2021 tested
- (3) Irritation test according to ISO 10993-23:2021 tested

And biocompatibility evaluation and test of the breathing gas pathway:

- (1) Performed test for emissions of particulate matter according to ISO 18562-2:2017.
- (2) Performed test for emissions of volatile organic compounds according to ISO 18562-3:2017.

- **Software Verification and Validation**

Software verification and validation was performed, and it was demonstrated that the software performs as intended according to the FDA Guidance for the Content of Premarket Submissions for Software Contained in Medical Devices.

The test results verified the device compliance to recognized consensus standards and therefore does not raise new questions of safety and effectiveness.

9. Concise summary for any performance testing

9.1 Summary of Performance testing – bench

The performance of device, such Sound Level, Oxygen Concentration, Flow Control Settings and Pulse Volumes, Maximum Outlet Pressure, Oxygen delivery and trigger sensitivity, was tested and complied with the applicable requirements of the following standards: ISO 80601-2-69, ISO 80601-2-67, demonstrated that the performance of the device met the requirements.

Battery safety was tested and verified via IEC 62133-2, and battery performance testing was performed which meet the requirement of performance.

9.2 Performance Testing – Animal

Not applicable to this submission for the Portable Oxygen Concentrator. The subject of this premarket submission, Portable Oxygen Concentrator, did not require animal studies and test results to support substantial equivalence.

9.3 Performance Testing – Clinical

Not applicable to this submission for the Portable Oxygen Concentrator. The subject of this premarket submission, Portable Oxygen Concentrator, did not require clinical studies to support substantial equivalence.

10. Substantial Equivalence Conclusion

The Portable Oxygen Concentrator, model: Spirit-3, has the same intended use and similar characteristics as the cleared predicate device Portable Oxygen Concentrator, Model: P2-K4, P2-S4.

Moreover, the portable oxygen concentrator, model Spirit-3, meets the safety and performance standards required for Portable Oxygen Concentrators, as verified by the bench tests included in this submission, which demonstrate that the differences that existed between Spirit-3 and P2-K4, P2-S4 do not raise any new questions of safety or effectiveness.

Thus, the Portable Oxygen Concentrator, model: Spirit-3 is Substantially Equivalent (SE) to the predicate device.