



January 23, 2026

Basic Medical Technology, Inc.  
Amy Gao  
RA Engineer  
5300 Concours St.  
Ontario, California 91764

Re: K253160

Trade/Device Name: Syntex Exam Gloves  
Regulation Number: 21 CFR 880.6250  
Regulation Name: Non-powdered patient examination glove  
Regulatory Class: Class I, reserved  
Product Code: LZA  
Dated: December 29, 2025  
Received: December 29, 2025

Dear Amy Gao:

We have reviewed your section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (the Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database available at <https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Additional information about changes that may require a new premarket notification are provided in the FDA guidance documents entitled "Deciding When to Submit a 510(k) for a Change to an Existing Device"

(<https://www.fda.gov/media/99812/download>) and "Deciding When to Submit a 510(k) for a Software Change to an Existing Device" (<https://www.fda.gov/media/99785/download>).

Your device is also subject to, among other requirements, the Quality System (QS) regulation (21 CFR Part 820), which includes, but is not limited to, 21 CFR 820.30, Design controls; 21 CFR 820.90, Nonconforming product; and 21 CFR 820.100, Corrective and preventive action. Please note that regardless of whether a change requires premarket review, the QS regulation requires device manufacturers to review and approve changes to device design and production (21 CFR 820.30 and 21 CFR 820.70) and document changes and approvals in the device master record (21 CFR 820.181).

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR Part 803) for devices or postmarketing safety reporting (21 CFR Part 4, Subpart B) for combination products (see <https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR Part 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR Parts 1000-1050.

All medical devices, including Class I and unclassified devices and combination product device constituent parts are required to be in compliance with the final Unique Device Identification System rule ("UDI Rule"). The UDI Rule requires, among other things, that a device bear a unique device identifier (UDI) on its label and package (21 CFR 801.20(a)) unless an exception or alternative applies (21 CFR 801.20(b)) and that the dates on the device label be formatted in accordance with 21 CFR 801.18. The UDI Rule (21 CFR 830.300(a) and 830.320(b)) also requires that certain information be submitted to the Global Unique Device Identification Database (GUDID) (21 CFR Part 830 Subpart E). For additional information on these requirements, please see the UDI System webpage at <https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/unique-device-identification-system-udi-system>.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance>) and CDRH Learn (<https://www.fda.gov/training-and-continuing-education/cdrh-learn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

**ALLAN GUAN -S**

For Bifeng Qian, M.D., Ph.D.  
Assistant Director  
DHT4C: Division of Infection  
Control Devices  
OHT4: Office of Surgical and  
Infection Control Devices  
Office of Product Evaluation and Quality  
Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K253160

Device Name

Syntex Exam Gloves

Indications for Use (Describe)

The glove is a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.

Type of Use (Select one or both, as applicable)

☐ Prescription Use (Part 21 CFR 801 Subpart D)

☒ Over-The-Counter Use (21 CFR 801 Subpart C)

**CONTINUE ON A SEPARATE PAGE IF NEEDED.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

# Basic Medical Technology Inc.

## 510(K) Summary

### 1. Submission Information

Date Prepared: Jan 22, 2026

Submission Format: Traditional 510(k)

510(K)#: K253160

### 2. Submitter Information

Applicant Name: Basic Medical Technology Inc.

Address: 5300 Concourse Ontario, CA 91764

Contact Person: Amy Gao

Tel : (909) 980-1678

Email: amygao@basicmedical.com

### 3. Device Information

Trade Name: Syntex Exam Gloves

Common Name: Synthetic Butadiene-acrylonitrile latex Powder-Free Exam Gloves

Classification Name: Non-Powdered Patient Examination Glove

Regulation: 21 CFR 880.6250

Product Code: LZA

Classification Panel: General Hospital

Device Class: Class I

### 4. Predicate Device Information

#### **Predicate Device:**

Applicant: Tianjin Aoshang Outdoor Equipment Co.,Ltd.

510(k) #: K220442

Trade Name: Synguard Nitrile Exam Glove

### 5. Device Description:

Syntex Exam Gloves is a Class I patient examination glove, that made from synthetic Butadiene-acrylonitrile latex. They are natural latex color, non-sterile, powder-free, ambidextrous with beaded cuff, and single use only, and come in different sizes- XS, S, M, L, XL and XXL.

The device meets all the specifications in ASTM D6319-19, Standard specification for Nitrile Examination Gloves. Additionally, the gloves have been tested for biocompatibility.

## Basic Medical Technology Inc.

### 6. Indications for Use:

The glove is a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner.

### 7. Comparison of Subject Device and Predicate Device:

**General Comparison Table:**

| Device                       | Proposed Device                                                                                                                                           | Predicate Device                                                                                                                                                                     | Result    |
|------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|
| <b>510K #</b>                | K253160                                                                                                                                                   | K220442                                                                                                                                                                              | -         |
| <b>Applicant Name</b>        | Basic Medical Technology Inc.                                                                                                                             | Tianjin Aoshang Outdoor Equipment Co.,Ltd.                                                                                                                                           | -         |
| <b>Product Name</b>          | Syntex Exam Gloves                                                                                                                                        | Synguard Nitrile Exam Glove                                                                                                                                                          | -         |
| <b>Product Code</b>          | LZA                                                                                                                                                       | LZA                                                                                                                                                                                  | Same      |
| <b>Regulation Number</b>     | 21 CFR 880.6250                                                                                                                                           | 21 CFR 880.6250                                                                                                                                                                      | Same      |
| <b>Device Classification</b> | Class I                                                                                                                                                   | Class I                                                                                                                                                                              | Same      |
| <b>Indications for use</b>   | The glove is a disposable device intended for medical purposes that is worn on the examiner's hand to prevent contamination between patient and examiner. | Synguard Nitrile Exam Glove is a disposable device intended for medical purposes that is worn on the examiner's hand or finger to prevent contamination between patient and examiner | Same      |
| <b>Powder free</b>           | Yes                                                                                                                                                       | Yes                                                                                                                                                                                  | Same      |
| <b>Material</b>              | synthetic Butadiene-acrylonitrile latex                                                                                                                   | Acrylonitrile, Butadiene                                                                                                                                                             | Same      |
| <b>Color</b>                 | Natural latex color                                                                                                                                       | Blue                                                                                                                                                                                 | Different |
| <b>Size</b>                  | XS, S, M, L, XL, XXL                                                                                                                                      | S, M, L, XL                                                                                                                                                                          | Different |
| <b>Sterile</b>               | Non-Sterile                                                                                                                                               | Non-Sterile                                                                                                                                                                          | Same      |
| <b>Single Use</b>            | Single use                                                                                                                                                | Single use                                                                                                                                                                           | Same      |
| <b>Design feature</b>        | Ambidextrous                                                                                                                                              | Ambidextrous                                                                                                                                                                         | Same      |

\*Different 1- Color

The color of the proposed device is different from the predicate device. However, the biocompatibility and specification performance tests have been performed on each of the proposed devices and the test results show no significant concerns. Therefore, the difference will not affect the safety and effectiveness of the proposed device.

\*Different 2- Size

The size of the proposed device is not exactly same as the predicate device. However, the size of the proposed device has been covered by ASTM D6319-19. The user can select appropriate model depended on size of user's hand. Therefore, the difference will not affect the safety and effectiveness of the proposed device.

**Technological Characteristic Comparison Table:**

| Characteristics | Proposed device | Predicate device | Result |
|-----------------|-----------------|------------------|--------|
|-----------------|-----------------|------------------|--------|

## Basic Medical Technology Inc.

|                                                                    |                         | K253160                                                                                       | K220442                                                   |           |
|--------------------------------------------------------------------|-------------------------|-----------------------------------------------------------------------------------------------|-----------------------------------------------------------|-----------|
| <b>Dimension</b><br><br><b>Reference standard</b><br>ASTM D6319-19 | Length                  | Min 220mm for XS, S<br>Min 230mm for M, L, XL, XXL                                            | Min 230mm for S, M<br>Min 240mm for L, XL                 | Different |
|                                                                    | Width                   | XS: 70±10 mm<br>S: 80±10 mm<br>M: 95±10 mm<br>L: 110±10 mm<br>XL: 120±10 mm<br>XXL: 130±10 mm | S: 84±5 mm<br>M: 94±5 mm<br>L: 110±10 mm<br>XL: 120±10 mm | Different |
|                                                                    | Thickness               | Palm: Min 0.05 mm<br>Finger: Min 0.05 mm                                                      | Palm: Min 0.05 mm<br>Finger: Min 0.05 mm                  | Same      |
| <b>Physical properties</b><br><br>ASTMD6319-19<br>ASTM D412-16     | <b>Before aging</b>     |                                                                                               |                                                           | Same      |
|                                                                    | Tensile strength        | Min 14MPa                                                                                     | Min 14MPa                                                 |           |
|                                                                    | Ultimate elongation     | Min 500%                                                                                      | Min 500%                                                  |           |
|                                                                    | <b>After aging</b>      |                                                                                               |                                                           |           |
|                                                                    | Tensile strength        | Min 14MPa                                                                                     | Min 14MPa                                                 |           |
|                                                                    | Ultimate elongation     | Min 400%                                                                                      | Min 400%                                                  |           |
| <b>Freedom from holes</b><br>ASTM D6319-19<br>ASTM D5151-19        |                         | No water leakage occurs<br>G-I, AQL2.5                                                        | No water leakage occurs<br>G-I, AQL2.5                    | Same      |
| <b>Residual Powder</b><br>ASTMD6319-19<br>ASTM D6124-06            |                         | ≤ 2 mg per glove                                                                              | ≤ 2 mg per glove                                          | Same      |
| <b>Biocompatibility</b>                                            | Skin Irritation         | No Irritation                                                                                 | No Irritation                                             | Same      |
|                                                                    | Sensitization           | No Sensitization                                                                              | No Sensitization                                          | Same      |
|                                                                    | Acute Systemic Toxicity | No Acute Systemic Toxicity                                                                    | No Acute Systemic Toxicity                                | Same      |

### 8. Summary of Non-Clinical Performance Data

Non-clinical tests were conducted to verify that the performance of the proposed device met all design specifications. The test results demonstrated that the proposed device met the performance criteria with the following standards:

- ASTM D6319-19, Standard Specification for Nitrile Examination Gloves for Medical Application.
- ASTM D5151-19, Standard Test Method for Detection of Holes in Medical Gloves.
- ASTM D6124-06, Standard Test Method for Residual Powder on Medical Gloves
- ASTM D412-16, Standard Test Methods for Vulcanized Rubber and Thermoplastic Elastomers—Tension
- ISO 10993-10: 2021 Biological evaluation of medical devices - Part 10: Tests for skin sensitization.

## Basic Medical Technology Inc.

- ISO10993-23:2021 Biological evaluation of medical devices - Part 23: Tests for irritation.
- ISO 10993-11:2017 Biological evaluation of medical devices — Part 11: Tests for systemic toxicity

| Test performed          | Methodology                    | Acceptance Criteria                                                                     | Result |
|-------------------------|--------------------------------|-----------------------------------------------------------------------------------------|--------|
| Freedom From Holes      | ASTM D6319-19<br>ASTM D5151-19 | Meet requirement inspection level G- 1, AQL 2.5 (ISO2859-1)                             | Pass   |
| Dimension-Length        | ASTM D6319-19                  | Minimum 220mm for size XS-S<br>Minimum 230mm for size M-XXL                             | Pass   |
| Dimension-Width         | ASTM D6319-19                  | XS: 70±10mm<br>S: 80±10mm<br>M: 95±10mm<br>L: 110±10mm<br>XL: 120±10mm<br>XXL: 130±10mm | Pass   |
| Dimension-Thickness     | ASTM D6319-19                  | Finger: 0.05mm (min)<br>Palm: 0.05mm (min)                                              | Pass   |
| Physical properties     | ASTM D6319-19<br>ASTM D412-16  | Tensile Strength (Min14 Mpa)<br>Elongation (Before Aging 500% and after aging 400%) Min | Pass   |
| Powder Residue          | ASTM D6319-19<br>ASTM D6124-06 | Not more than 2 mg per glove                                                            | Pass   |
| Skin Irritation         | ISO10993-23:2021               | Under the conditions of the study, not a primary skin irritant                          | Pass   |
| Skin Sensitization      | ISO10993-10:2021               | Under the conditions of the study, not a contact sensitizer                             | Pass   |
| Acute Systemic Toxicity | ISO 10993-11:2017              | Under the conditions of the study, no signs of acute systemic toxicity were observed    | Pass   |

### 9. Summary of Clinical Testing

Not applicable.

### 10. Conclusion

The conclusions drawn from the non-clinical tests demonstrate that the subject device is as safe, as effective, and performs as well as or better than the legally marketed predicate device.