

**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION
DECISION SUMMARY**

A. 510(k) Number: K163713

B. Purpose for Submission:

To obtain a substantial equivalence determination for the Sofia Lyme FIA for use with Sofia

C. Measurand:

IgM antibodies to *Borrelia burgdorferi* proteins

IgG antibodies to *Borrelia burgdorferi* proteins

D. Type of Test:

Fluorescent immunoassay (FIA), bi-directional lateral flow format

E. Applicant:

Quidel Corporation

F. Proprietary and Established Names:

Sofia Lyme FIA

Sofia Lyme Control Set

G. Regulatory Information:

1. Regulation section:

21 CFR 866.3830

2. Classification:

Class II

3. Product code:

LSR; Reagent, Borrelia Serological Reagent

JTT; Quality Control Material (assayed and unassayed)

4. Panel:

Immunology and Microbiology Devices

H. Intended Use:

1. Intended use(s):

Sofia Lyme FIA:

The Sofia Lyme FIA employs immunofluorescence for the rapid differential detection of human IgM and IgG antibodies to *Borrelia burgdorferi* from serum and plasma specimens from patients suspected of *B. burgdorferi* infection. This qualitative test is intended for use as an aid in the diagnosis of Lyme disease. A negative result does not preclude infection with *B. burgdorferi*. All positive results for IgM and/or IgG should be further tested by a corresponding second-tier western blot assay. Test results are to be used in conjunction with information obtained from the patient's clinical evaluation and other diagnostic procedures.

Sofia Lyme Control Set:

The Sofia Lyme Control Set is intended for use as assayed quality control materials to verify the performance of the Sofia Lyme FIA test system.

2. Indication(s) for use:

Same as Intended Use(s)

I. Device Description:

The Sofia Lyme FIA is an immunofluorescence-based, lateral flow assay for detection of IgM and/or IgG antibodies to *Borrelia burgdorferi* in patient specimens. Reagents for the assay are ready-to-use and provided in the kit.

The assay uses a bidirectional test strip format to detect both IgM and IgG antibodies to *B. burgdorferi*. One side of the test strip detects IgM antibodies to *B. burgdorferi* and the other side of the test strip detects IgG antibodies to *B. burgdorferi*.

To perform the test, the patient specimen (30 µL serum/plasma) is added to a pre-filled Reagent Vial. The diluted sample (100 µL) is then dispensed into the round sample well located near the center of the Test Cassette. The Test Cassette is loaded into the Sofia instrument for an automatically defined development time (WALK AWAY Mode) or pre-incubated on the bench top prior to loading into Sofia (READ NOW Mode). Sofia scans the test strip, analyzes the fluorescent signal, and then displays two (2) test results: IgM (positive, negative or invalid) and IgG (positive, negative or invalid). As an option, the results are automatically printed on the integrated printer.

Each Sofia Lyme FIA kit will contain one Positive and one Negative Control—each provided in separate dropper bottles. The external controls will be available separately as well in a Sofia Lyme Control Set. The Positive and Negative QC controls are formulated with patient Lyme IgM and IgG positive serum that are diluted into 1X PBS, and 0.3% Microcide is added to the solution as an antimicrobial. To use External Controls, two drops from the Negative or Positive Control are added to the round sample well located near the center of the cassette. The cassette is then treated in the same manner as a cassette containing a patient sample. The Positive and Negative Control each require their own cassette.

J. Substantial Equivalence Information:

1. Predicate device name(s): Vidas Lyme IgM
Vidas Lyme IgG

2. Predicate 510(k) number(s): K122979
K122986

3. Comparison with predicate:

Item Features	Proposed Device Sofia Lyme FIA with Sofia	Predicate IgM Device BioMérieux Vidas Lyme IgM (K122979)	Predicate IgG Device BioMérieux Vidas Lyme IgG (K122986)
Intended Use: Sofia Lyme FIA	Sofia Lyme FIA: The Sofia Lyme FIA employs immunofluorescence for the rapid differential detection of human IgM and IgG antibodies to <i>Borrelia burgdorferi</i> from serum and plasma specimens from patients suspected of <i>B. burgdorferi</i> infection. This qualitative test is intended for use as an aid in the diagnosis of Lyme disease. A negative result does not preclude infection with <i>B. burgdorferi</i> . All positive results for IgM and/or IgG should be further tested by a corresponding second- tier western blot assay. Test results are to be used in conjunction with information obtained from the patient's clinical evaluation and other diagnostic procedures.	The VIDAS Lyme IgM (LYM) assay is an automated qualitative enzyme immunoassay intended for use on the instruments of the VIDAS family in the presumptive detection of human IgM antibodies to <i>Borrelia burgdorferi</i> in human serum (plain or separation gel) or plasma (sodium heparin or lithium heparin). It should be used to test patients with a history and/or symptoms of infection with <i>B. burgdorferi</i> . All VIDAS Lyme IgM positive and equivocal specimens should be further tested with a Western Blot IgG assay to obtain supportive evidence of infection with <i>B. burgdorferi</i> .	The VIDAS Lyme IgG (LYG) assay is an automated qualitative enzyme immunoassay intended for use on the instruments of the VIDAS family in the presumptive detection of human IgG antibodies to <i>Borrelia burgdorferi</i> in human serum (plain or separation gel) or plasma (sodium heparin or lithium heparin). It should be used to test patients with a history and/or symptoms of infection with <i>B. burgdorferi</i> . All VIDAS Lyme IgG positive specimens should be further tested with a western blot IgG assay to obtain supportive evidence of infection with <i>B. burgdorferi</i> .

Item Features	Proposed Device Sofia Lyme FIA with Sofia	Predicate IgM Device BioMérieux Vidas Lyme IgM (K122979)	Predicate IgG Device BioMérieux Vidas Lyme IgG (K122986)
Intended Use: Sofia Lyme Control Set	Sofia Lyme Control Set: The Sofia Lyme Control Set is intended for use as assayed quality control materials to verify the performance of the Sofia Lyme FIA test system.		
Instrument	Sofia	VIDAS and miniVIDAS	VIDAS and miniVIDAS
Analyte	Human IgM and IgG antibodies against <i>B. burgdorferi</i> proteins	Human IgM antibodies against <i>B. burgdorferi</i> proteins	Human IgG antibodies against <i>B. burgdorferi</i> proteins
Automated Analysis	Yes	Yes	Yes
Read Results	Read results on instrument screen or print with optional printer	Result report is printed	Result report is printed
Read Result Time	10 Minutes	27 minutes	27 minutes
Specimen Types	Human serum and plasma	Human serum and plasma	Human serum and plasma
Qualitative	Yes	Yes	Yes
Test Principle	Immunofluorescence Device	Immunofluorescence Device	Immunofluorescence Device
Format	Lateral-flow Bi-directional Test Cassette	Enzyme-linked fluorescent assay (ELFA)	Enzyme-linked fluorescent assay (ELFA)
Antibodies Used	Monoclonal anti-human IgG and polyclonal anti- human IgM	Anti-human IgM antibodies	Anti-human IgG antibodies
Antigens Used	Recombinant Proteins and synthetic peptides of <i>B. burgdorferi</i>	Recombinant proteins of <i>B. burgdorferi</i>	Recombinant proteins of <i>B. burgdorferi</i>
Detection Method	Polystyrene microparticles dyed with Europium chelate	Alkaline phosphatase/4-MUP	Alkaline phosphatase/4-MUP
Storage	Room Temperature (15-30°C)	2-8°C	2-8°C
Running Buffer Solution	One pre-filled vial containing PBS	Sample diluent and wash buffer	Sample diluent and wash buffer

Item Features	Proposed Device Sofia Lyme FIA with Sofia	Predicate IgM Device BioMérieux Vidas Lyme IgM (K122979)	Predicate IgG Device BioMérieux Vidas Lyme IgG (K122986)
Quality Control Features	Built-in features include: - Built-in procedural control zone scanned by the analyzer to determine whether adequate flow occurred on the IgG side of the assay. - Built-in reference control line scanned by the analyzer to determine whether adequate flow occurred on the IgM side of the assay. - Analyzer prevents use of expired cartridge from being read by the reader - Cassette must be properly inserted	One positive and one negative control are included and must be tested after opening a new kit to monitor reagent failure.	One positive and one negative control are included and must be tested after opening a new kit to monitor reagent failure.

K. Standard/Guidance Document Referenced (if applicable):

Establishing the Performance Characteristics of *In Vitro* Diagnostic Devices for the Detection of Antibodies to *Borrelia burgdorferi*

L. Test Principle:

Immunofluorescence

M. Performance Characteristics (if/when applicable):

1. Analytical performance:

a. *Precision/Reproducibility:*

The precision study evaluated the within laboratory precision / repeatability of the Sofia Lyme FIA using serum samples.

Each negative or spiked positive serum sample was tested in duplicate with a total of 24 different runs (two runs per day over a total of 12 days) at one site (Quidel), utilizing two operators and two Sofia Analyzers during the study. Each sample was tested a total of 48 times over the course of the study. The testing days spanned across two calibration events. See summary tables below

Sofia Lyme FIA Precision – Within Run

Sample	IgM Positive		IgM % Positivity Total (n=96)	IgG Positive		IgG % Positivity Total (n=96)
	Run 1	Run 2		Run 1	Run 2	
Negative (C ₀)	0/48	0/48	0%	0/48	0/48	0%
High Negative (C ₅)	3/48	5/48	8.33%	2/48	5/48	7.29%
Low Positive (C ₉₅)	48/48	48/48	100%	48/48	46/48	97.92%
Moderate Positive (2-3X)	48/48	48/48	100%	48/48	48/48	100%

Sofia Lyme FIA Precision – Between Operators

Sample	% IgM Positivity		IgM % Positivity Total (n=96)	% IgG Positivity		IgG % Positivity Total (n=96)
	Op. 1	Op. 2		Op. 1	Op. 2	
Negative (C ₀)	0/48	0/48	0%	0/48	0/48	0%
High Negative (C ₅)	3/48	5/48	8.33%	3/48	4/48	7.29%
Low Positive (C ₉₅)	48/48	48/48	100%	47/48	47/48	97.92%
Moderate Positive (2-3X)	48/48	48/48	100%	48/48	48/48	100%

The reproducibility study demonstrated inter-laboratory reproducibility with a panel of test samples at various concentrations of IgM and IgG antibodies to *B. burgdorferi*. See summary of results below.

Sofia Lyme FIA IgM and IgG Reproducibility Study Inter-laboratory Agreement

Site	IgM				IgG			
	Negative	High Negative	Low Positive	Moderate Positive	Negative	High Negative	Low Positive	Moderate Positive
1	30/30	30/30	30/30	30/30	30/30	30/30	27/30	30/30
2	30/30	30/30	30/30	30/30	30/30	30/30	26/30	30/30
3	30/30	30/30	30/30	30/30	30/30	30/30	25/30	30/30
Total	90/90	90/90	90/90	90/90	90/90	90/90	78/90	90/90
% Overall Agreement (95% CI)	100% (95.1-100%)	100% (95.1-100%)	100% (95.1-100%)	100% (95.1-100%)	100% (95.1-100%)	100% (95.1-100%)	86.7% (78.0-92.4%)	100% (95.1-100%)

b. Linearity/assay reportable range: Not applicable

c. Traceability, Stability, Expected values (controls, calibrators, or methods):

Not applicable (N/A)

d. Detection limit: Not applicable

e. *Analytical specificity:*

This study evaluated the analytical specificity of the Sofia Lyme FIA using serum obtained from asymptomatic (healthy, normal) populations in both endemic and non-endemic regions. Individuals were split evenly between endemic and non-endemic populations. The specificity of the Sofia Lyme FIA and the predicate device is summarized below.

Summary of Analytical Specificity - IgM

Assay	N	Negative	Equivocal	Positive	% Specificity
Sofia Lyme FIA IgM	200	180	N/A	20	90.0%
Predicate IgM	200	181	16	3	90.5%

Summary of Analytical Specificity - IgG

Assay	N	Negative	Positive	% Specificity
Sofia Lyme FIA IgG	200	180	20	90.0%
Predicate Lyme IgG	200	197	3	98.5%

f. *Assay cut-off:* Not applicable

g. *Interfering Substances:*

A study was performed to assess whether common endogenous and exogenous substances would interfere with the Sofia Lyme FIA. The substances and concentrations tested are in the table below. There was no evidence of interference for the listed substances at the tested levels.

Substances Evaluated for Interference with Sofia Lyme FIA

Interfering Substance	Concentrations Tested
Bilirubin	15 mg/dL
Hemoglobin	20 g/dL
Lipids	750 mg/dL
Albumin	5.0 g/dL
Acetylsalicylic Acid	3.62 mmol/L
Amoxicillin	206 µmol/L
Azithromycin	15.3 µmol/L
Ceftriaxone	1460 µmol/L
Cefuroxime Axetil	1416 µmol/L
Doxycycline Hyclate	67.5 µmol/L
Erythromycin	81.6 µmol/L
Ibuprofen	2425 µmol/L
Minocycline	10.33 µmol/L
Penicillin G	33.67 µmol/L
Phenoxymethylpenicillin	14.27 µmol/L
Prednisolone	8.31 µmol/L
Tetracyclines	34 µmol/L

h. Cross-Reactivity:

The cross-reactivity of the Sofia Lyme FIA was evaluated with 17 disease state sample types for possible interference with the Sofia Lyme FIA assay. Samples were characterized and obtained from commercial vendors and did not include any first or second tier Lyme testing results. The total number of each disease state sample type, as well as the number of positives observed with the Sofia Lyme FIA is summarized in the table below.

Evaluation of Conditions that May Cross-React with Sofia Lyme FIA

Disease State Diagnosis	# of Samples	IgM Positive Results	IgG Positive Results
Anti-Nuclear Antibodies	10	4/10	1/10
Babesiosis	12	3/12	4/12
Chronic Fatigue Syndrome	12	2/12	2/12
Cytomegalovirus	11	2/11	1/11
Epstein Barr Virus	10	4/10	0/10
Fibromyalgia	10	1/10	0/10
<i>H. pylori</i>	10	4/10	0/10
HIV	11	2/11	0/11
Influenza	11	0/11	0/11
Leptospirosis	6	1/6	2/6
Lupus	20	8/20	4/20
Multiple Sclerosis	10	1/10	1/10
Parvovirus B19	15	4/15	1/15
Rheumatoid Factor	10	0/10	0/10
Rickettsia	3	0/3	0/3
Rocky Mountain Spotted Fever	10	0/10	0/10
Syphilis	28	8/28	4/28

i. Quality Controls/ Matrix Effects

The Sofia Lyme External positive and negative quality controls are individual bottles manufactured to generate positive IgG and IgM results, and negative IgG and IgM results. The Sofia Lyme QC materials simulate the composition of patient samples as closely as possible in order to minimize matrix effects and correctly reflect the expected performance with patient samples. The positive and negative QC controls are formulated with patient serum samples that are diluted into 1X PBS, filled into dropper bottles and added dropwise to the Sofia Lyme FIA test device. Unlike a patient sample, 0.3% Microcide is added to the external control solutions to act as an antimicrobial.

The purpose of this study was to evaluate matrix effects other than the analyte that may influence the measured value of the analyte in the positive and negative external quality controls. This protocol is based on the FDA Guidance for Industry and FDA Staff Assayed and Unassayed Quality Control Material: June 7, 2007.

This study addressed the following potential effects:

- Determine if there was any bias of the external control matrix relative to the natural sample. Demonstrate that the manufacturing formulation that includes 0.3% Microcide did not alter the properties of the sample and still reflects similar performance compared to natural human specimens.
- Demonstrate that the external quality control solution result output was similar to the sample result output regardless of the preparatory and operational steps that are required.
- Demonstrate the precision of the manufactured positive and negative controls over a span of 5 days.

Study results are summarized in the tables below.

External Controls – Evaluation of Matrix Effects

Matrix	Method	N	IgM % Positivity	IgG % Positivity
			Cutoff	Cutoff
Negative QC Control	Dropper	100	0%	0%
Negative Patient Sample	Pipette	100	0%	0%
Positive QC Control	Dropper	100	100%	100%
Positive Patient Sample	Pipette	100	100%	100%

External Controls – Precision

			# Positive Results		% Positivity	
QC Control	Date	n	IgM	IgG	IgM	IgG
Negative	8/22/2016	20	0/20	0/20	0%	0%
	8/23/2016	20	0/20	0/20	0%	0%
	8/24/2016	20	0/20	0/20	0%	0%
	8/25/2016	20	0/20	0/20	0%	0%
	8/26/2016	20	0/20	0/20	0%	0%
	Total	100	0/100	0/100	0%	0%
Positive	8/22/2016	20	20/20	20/20	100%	100%
	8/23/2016	20	20/20	20/20	100%	100%
	8/24/2016	20	20/20	20/20	100%	100%
	8/25/2016	20	20/20	20/20	100%	100%
	8/26/2016	20	20/20	20/20	100%	100%
	Total	100	100/100	100/100	100%	100%

In addition to the above studies, the external controls were used daily or as described in the package insert in the following studies, and performed as expected:

Blood Tube Matrix Comparison	Precision
Reproducibility	Cross-Reactivity
Interfering Substances	Read Now vs Walk Away
Specimen Volume control	Development Time
Diluted Sample Volume	

2. Comparison studies:

a. *Method comparison with predicate device: prospective*

A prospective study was performed on 678 samples to establish the clinical performance of the assay. Of these, 178 samples were prospectively-collected and 500 were frozen samples previously submitted for Lyme testing and sequentially retained without regard to results. At each of three laboratories, samples were tested in parallel using a commercially available Lyme IgM and IgG EIA method (predicate) and the Sofia Lyme FIA. Possible results for the predicate assay were negative, positive or equivocal. Possible results for the Sofia Lyme FIA were negative or positive. Positive Percent Agreement (PPA) was calculated between Sofia Lyme FIA positive results and predicate positive and equivocal results because the Centers for Disease Control (CDC) recommended 2-tier Lyme testing algorithm calls for first tier positive and equivocal results to be tested by a second tier western blot assay. Combined results from the three sites are shown below:

Lyme IgM Results (Tier 1) for Sofia Compared to the Predicate Test -Prospective Specimens

	Predicate Lyme IgM			% Agreement	95% CI
	Pos	Equivocal	Neg		
Sofia IgM Pos	37	20	67	PPA ² = 56.4% (57 / 101)	46.7% - 65.7%
Sofia IgM Neg	7	37 ¹	510	NPA = 88.4% (510 / 577)	85.5% - 90.8%
Total:	44	57	577		

¹Of the 37 specimens that were Sofia Lyme FIA negative and predicate equivocal, 34 were negative by western blot

² Additional Information: Four specimens that were negative by the predicate were positive by Sofia Lyme FIA and confirmed positive by western blot

Lyme IgG Results (Tier 1) for Sofia Compared to Predicate Test -Prospective Specimens

	Predicate Lyme IgG		% Agreement	95% CI
	Pos	Neg		
Sofia IgG Pos	47	54	PPA = 82.5% (47 / 57)	70.4% - 90.4%
Sofia IgG Neg	10	567	NPA = 91.3% (567 / 621)	88.8% - 93.3%
Total:	57	621		

Second-Tier Testing: In accordance with the CDC recommended 2-tier Lyme testing algorithm the Sofia Lyme FIA results and the predicate Lyme IgM and IgG positive and equivocal results were confirmed using a commercially available Lyme IgM or IgG western blot assay. A summary of second-tier testing results along with the PPA between Sofia Lyme FIA and the predicate for first and second tier testing are shown below.

Lyme IgM Results (Tier 2) for Sofia Compared to the Predicate Assay -Prospective Specimens

Assay	Tier 1 + or ±	IgM WB +	IgM WB -
Predicate IgM	101 ¹	37	63
Sofia IgM	124 ²	45	76
Predicate and Sofia IgM	57 ¹	33	23

¹ One sample had insufficient volume to be tested on WB insufficient

² Three samples had insufficient volume to be tested on WB insufficient

Lyme IgM PPA for Sofia Compared to the Predicate Assay -Prospective Specimens

	PPA	N	95% CI
1st Tier PPA	56.4%	57/101	46.7% - 65.7%
2nd Tier PPA	89.2%	33/37	74.7% - 96.3%%

Lyme IgG Results (Tier 2) for Sofia Compared to Predicate Test -Prospective Specimens

Assay	Tier 1 +	IgG WB +	IgG WB -
Predicate IgG	57	28	29
Sofia IgG	101 ¹	25	72
Predicate and Sofia IgG	47	25	22

¹ Four samples had insufficient volume to be tested on WB insufficient

Lyme IgG PPA for Sofia Compared to the Predicate Assay -Prospective Specimens

	PPA	N	95% CI
1st Tier PPA	82.5%	47/57	70.4% - 90.4%
2nd Tier PPA	89.3%	25/28	72.0% - 97.1%

b. Method comparison with predicate device: retrospective

Additionally, to supplement positive IgM and IgG results, 52 known positive samples were obtained from a commercial vendor, blinded and distributed to testing sites. Results of testing are summarized below:

Lyme IgM Results for Sofia Compared to Predicate Test - Retrospective Specimens

	Predicate Lyme IgM			% Agreement	95% CI
	Pos	Equivocal	Neg		
Sofia IgM Pos	17	11	4	PPA* = 90.3% (28 / 31)	74.3% - 97.4%
Sofia IgM Neg	0	3 ¹	17	NPA = 81.0% (17 / 21)	59.4% - 92.9%
Total:	17	14	21		

¹Of the three specimens that were Sofia Lyme FIA negative and predicate equivocal, 1 was negative by western blot

Lyme IgG Results for Sofia Compared to Predicate Test - Retrospective Specimens

	Predicate Lyme IgG		% Agreement	95% CI
	Pos	Neg		
Sofia IgG Pos	27	5	PPA = 100% (27 / 27)	89.2% - 100.0%
Sofia IgG Neg	0	20	NPA = 80.0% (20 / 25)	60.4% - 91.6%
Total:	27	25		

Second-Tier Testing: In accordance with the CDC recommended 2-tier Lyme testing algorithm the Sofia Lyme FIA results and the predicate Lyme IgM and IgG positive and equivocal results were confirmed using a commercially available Lyme IgM or IgG western blot method. A summary of second-tier testing results along with the PPA between Sofia Lyme FIA and the predicate for first and second tier testing are shown below.

IgM Second Tier Testing - Retrospective Specimens

Assay	Tier 1 + or ±	IgM WB +	IgM WB -
Predicate IgM	31	22	9
Sofia IgM	32 ¹	20	11
Predicate and Sofia IgM	28	20	8

¹One sample had insufficient volume to be tested on WB insufficient

Lyme IgM PPA for Sofia Compared to the Predicate Assay -Retrospective Specimens

	PPA	N	95% CI
1st Tier PPA	90.3%	28/31	74.3% - 97.4%
2nd Tier PPA	90.9%	20/22	71.0%-98.7%

IgG Second Tier Testing - Retrospective Specimens

Assay	Tier 1 +	IgG WB +	IgG WB -
Predicate IgG	27	8	19
Sofia IgG	32 ¹	8	22
Predicate and Sofia IgG	27	8	19

¹Two samples had insufficient volume to be tested on WB insufficient

Lyme IgG PPA for Sofia Compared to the Predicate Assay -Retrospective Specimens

	PPA	N	95% CI
1st Tier PPA	100.0%	27/27	85.2% - 100.0%
2nd Tier PPA	100.0%	8/8	62.8% - 100.0%

c. Matrix comparison:

A matrix comparison study was performed by testing replicates of donor blood. Ten donors were used for each concentration, with the exception of the negative IgG

sample, for which there were 21 donors. Five replicates from each donor were used.

The comparative matrix qualitative IgM and IgG results demonstrate that four (4) whole blood collection tube types, (serum (RED), serum (TIGER), sodium heparin, and lithium heparin tubes) are not significantly different when compared to the serum matrix IgM and IgG results.

The EDTA collection tube qualitative results were significantly different when compared to the serum matrix results. This difference was most pronounced in samples that were evaluated at the C5 level.

The IgM and IgG Matrix Equivalency summary tables for each collection tube including each level, the number and percentage (%) tested are shown below:

Matrix Comparison – Blood Collection Tubes - IgM

Whole Blood Collection Tube	Number and Percentage (%) of Samples			
	Negative	C5	C95	2-3X
Serum (RED)	0 of 50 (0%)	49 of 50 (98%)	50 of 50 (100%)	50 of 50 (100%)
Serum (Tiger)	0 of 50 (0%)	48 of 50 (96%)	50 of 50 (100%)	50 of 50 (100%)
EDTA	0 of 50 (0%)	40 of 50 (80%)	50 of 50 (100%)	50 of 50 (100%)
Sodium Heparin	0 of 50 (0%)	50 of 50 (100%)	50 of 50 (100%)	50 of 50 (100%)
Lithium Heparin	0 of 50 (0%)	46 of 50 (92%)	50 of 50 (100%)	50 of 50 (100%)

Matrix Comparison – Blood Collection Tubes - IgG

Whole Blood Collection Tube	Number and Percentage (%) of Samples			
	Negative	C5	C95	2-3X
Serum (RED)	0 of 105 (0%)	36 of 50 (72%)	44 of 50 (88%)	50 of 50 (100%)
Serum (Tiger)	0 of 105 (0%)	44 of 50 (88%)	45 of 50 (90%)	50 of 50 (100%)
EDTA	0 of 105 (0%)	50 of 50 (100%)	50 of 50 (100%)	50 of 50 (100%)
Sodium Heparin	0 of 105 (0%)	34 of 50 (68%)	46 of 50 (92%)	50 of 50 (100%)
Lithium Heparin	0 of 105 (0%)	43 of 50 (86%)	48 of 50 (96%)	50 of 50 (100%)

3. Clinical studies:

a. *Clinical Sensitivity:*

Well-characterized clinically or culture confirmed Lyme disease serum samples were tested with the Sofia Lyme FIA. The results were compared to those obtained with the predicate Lyme IgG and IgM assays. The samples include the following patient populations:

- Category 1a: Initial (acute) samples from patients with documented erythema migrans (EM) or culture positive <1 month after symptom onset.
- Category 1b: Initial (acute) samples from patients with documented erythema migrans (EM) or culture positive 1-2 months after symptom onset.
- Category 1c: Initial (acute) samples from patients with documented erythema migrans (EM) or culture positive 2-3 months after symptom onset.
- Category 2: Convalescent samples from patients with documented erythema migrans (EM) or culture positive 3-12 months after appearance of symptoms.
- Category 3: Late Lyme (>1 year) neurological or Arthritic.

Results are reported as Lyme IgG and/or Lyme IgM positive or negative by Sofia for each sample tested. Results are summarized for each test method as “Sofia % Sensitivity” and “Predicate Device % Sensitivity” including 95% Confidence Interval for each method. Percent Sensitivity and 95% Confidence Interval are calculated for each sample sub-group as well as for all groups combined.

IgM Results for each Clinical Characterization Category

Category	n	Sofia IgM				Predicate IgM				
		Pos	Neg	% Sens	95% CI	Pos	Equiv	Neg	%Sens ¹	95% CI
Acute, < 1 month, with EM	64	39	25	60.9%	48.7-72.0%	28	8	28	56.3%	44.1-67.7%
Acute, 1-2 months, with EM	4	3	1	75.0%	28.9-96.6%	2	0	2	50.0%	15.0-85.0%
Convalescent, 3-12 months, with EM	15	11	4	73.3%	47.6-89.5%	6	2	7	53.3%	30.1-75.2%
Late Lyme (>1 yr), Neuro or Arthritic	12	8	4	66.7%	38.8-86.5%	7	3	2	83.3%	54.0-96.5%
All Categories	95	61	34	64.2%	54.2-73.1%	43	13	39	58.9%	48.9-68.3%

¹Of the 13 samples that were “equivocal” by the predicate, 12 of 13 were negative by FDA approved western blot.

IgG Results for each Clinical Characterization Category

Category	n	Sofia IgG				Predicate IgG			
		Pos	Neg	% Sens	95%CI	Pos	Neg	% Sens	95%CI
Acute, < 1 month, with EM	64	50	14	78.1%	66.5-86.6%	28	36	43.8%	32.3-55.9%
Acute, 1-2 months, with EM	4	4	0	100%	54.3-100.0%	3	1	75.0%	28.9-96.6%
Convalescent, 3-12 months, with EM	15	10	5	66.7%	41.5-85.0%	6	9	40.0%	19.8-64.3%
Late Lyme (>1 yr), Neuro or Arthritic	12	12	0	100%	78.4-100.0%	10	2	83.3%	54.0-96.5%
All Categories	95	76	19	80.0%	70.8-86.9%	47	48	49.5%	39.6-59.4%

b. *Clinical specificity*: Not applicable

c. Other clinical supportive data (when a. and b. are not applicable):

CDC Lyme panel was provided and tested with the Sofia Lyme FIA. First tier testing on the panel provided by the CDC was completed with a total Lyme antibody assay. These assay results were reported as positive or negative only, the type of antibody (IgM or IgG) was not reported. As a result their ability to evaluate the performance of the Sofia Lyme FIA on the CDC provided panel is limited. Please note that provision of a panel for testing does not imply an endorsement of the assay by the CDC.

Sofia IgM % Agreement with expected CDC result

Group ID	N	Sofia IgM			Western Blot Lyme IgM		
		Pos	Neg	% Agreement	Pos	Neg	% Agreement
Negative Controls	190	33	157	82.6%	0	190	100.0%
Early Lyme EM Positive	60	49	11	81.7%	31	29	51.7%
Late Lyme	30	22	8	73.3%	17	13	56.7%

Sofia IgG % Agreement with expected CDC result

Group ID	N	Sofia IgG			Western Blot Lyme IgG		
		Pos	Neg	% Agreement	Pos	Neg	% Agreement
Negative Controls	190	25	165	86.8%	0	190	100.0%
Early Lyme EM Positive	60	49	11	81.7%	19	41	31.7%
Late Lyme	30	30	0	100.0%	26	4	86.7%

N. Proposed Labeling:

The labeling is sufficient and it satisfies the requirements of 21 CFR Part 809.10.

O. Conclusion:

The submitted information in this premarket notification supports a substantial equivalence determination.