



**510(k) SUBSTANTIAL EQUIVALENCE DETERMINATION
DECISION SUMMARY
ASSAY ONLY**

I Background Information:

A 510(k) Number

K220763

B Applicant

ALPCO

C Proprietary and Established Names

ALPCO Calprotectin Immunoturbidimetric Assay

D Regulatory Information

Product Code(s)	Classification	Regulation Section	Panel
NXO	Class II	21 CFR 866.5180 - Fecal Calprotectin Immunological Test System	IM - Immunology

II Submission/Device Overview:

A Purpose for Submission:

New device

B Measurand:

Fecal calprotectin

C Type of Test:

Quantitative, immunoturbidimetric assay

III Intended Use/Indications for Use:

A Intended Use(s):

See Indications for Use below.

B Indication(s) for Use:

The ALPCO Calprotectin Immunoturbidimetric Assay is an in-vitro diagnostic assay used for the quantitative measurement of human fecal calprotectin in human stool. The ALPCO Calprotectin Immunoturbidimetric Assay is intended for in-vitro diagnostic use as an aid in diagnosis of inflammatory bowel disease (IBD), specifically Crohn's disease (CD) and ulcerative colitis (UC), and as an aid in the differentiation of IBD from irritable bowel syndrome (IBS) in conjunction with other clinical and laboratory findings.

C Special Conditions for Use Statement(s):

Rx - For Prescription Use Only

D Special Instrument Requirements:

Beckman AU680 analyzer

IV Device/System Characteristics:

A Device Description:

The ALPCO Calprotectin Immunoturbidimetric Assay consists of the following materials:

Materials provided (100 tests):

- Reagent 1: Tris buffer solution with <0.1% sodium azide, ready to use, 1 x 22.5 mL
- Reagent 2: Suspension of latex particles (< 0.5%) coated with anti-calprotectin antibodies with <0.1% sodium azide, ready to use, 1 x 5.3 mL

Materials Required but not Provided:

The ALPCO Calprotectin Calibrator Set

5 levels: 50, 200, 500, 800, 1050 µg/g, 5 x 1 mL, sold separately

The Calprotectin Immunoturbidimetric Assay Control Set

2 levels: Control 1: 35 µg/g; Control 2: 135 µg/g), 2 x 1 mL, sold separately

Feces sample collection

- 1) Sample collection tube
- 2) Transport container

Feces preparation:

The ALPCO Calprotectin Immunoturbidimetric Assay can prepare fecal samples by ALPCO Easy Stool Extraction Device (sold separately) or manual weighing method.

Instrument:
Beckman AU680 analyzer

B Principle of Operation:

The ALPCO Calprotectin Immunoturbidimetric Assay is a particle-enhanced immunoturbidimetric assay. Calprotectin in the sample binds to anti-calprotectin antibodies immobilized on the latex particles, which causes agglutination. This agglutination is detected as an absorbance change and the amount of agglutination is directly proportional to the concentration of calprotectin in the sample. The concentration is determined by interpolation from a calibration curve prepared from calibrators of known concentration.

The fecal samples can be processed and extracted using ALPCO Easy Stool Extraction Device or manual weighing method.

The ALPCO Easy Stool Extraction Device should not be used for the extraction of liquid/watery stool samples. Refer to K191807 for the fecal specimen preparation using the ALPCO Easy Stool Extraction Device.

For manual weighing method, 50 to 100 mg of the stool sample is taken by the inoculation loop and placed it into the pre-weighed tube. The sample is weighed and Extraction Buffer is added (99 times the sample weight volume). Further, the sample is homogenized on a multi-tube vortex by vigorous shaking (at highest speed) for 30 minutes. Next, the sample is centrifuged for 5 minutes at 3000 x g. Finally, decant the supernatant into a fresh labeled tube and continue with the assay procedure or store the extracts at 2–8°C for up to 3 days or at -80°C for up to 14 days.

Results are reported in µg/g of calprotectin in stool. Interpretation of results is shown as follows:

Calprotectin Concentration	Interpretation	Follow-Up
< 50 µg/g	Normal	None
50–100 µg/g	Equivocal	Retest in 4–6 weeks
> 100 µg/g	Elevated	Repeat as clinically indicated

V Substantial Equivalence Information:

A Predicate Device Name(s):

ALPCO Calprotectin Chemiluminescence ELISA, ALPCO Easy Stool Extraction Device

B Predicate 510(k) Number(s):

K191807

C Comparison with Predicate(s):

Device & Predicate Device(s):	<u>K220763</u> (New Device)	<u>K191807</u> (Predicate)
Device Trade Name	ALPCO Calprotectin Immunoturbidimetric Assay	ALPCO Calprotectin Chemiluminescence ELISA
General Device Characteristic Similarities		
Intended Use / Indications For Use	The ALPCO Calprotectin Immunoturbidimetric Assay is an in-vitro diagnostic assay used for the quantitative measurement of human fecal calprotectin in human stool. The ALPCO Calprotectin Immunoturbidimetric Assay is intended for in-vitro diagnostic use as an aid in diagnosis of inflammatory bowel disease (IBD), specifically Crohn's disease (CD) and ulcerative colitis (UC), and as an aid in the differentiation of IBD from irritable bowel syndrome (IBS) in conjunction with other clinical and laboratory findings.	The ALPCO Calprotectin Chemiluminescence ELISA is an in vitro diagnostic chemiluminescent assay intended for the quantitative measurement of fecal calprotectin, a neutrophilic protein that is a marker of intestinal mucosal inflammation, in human stool. The ALPCO Calprotectin Chemiluminescence ELISA is intended for in vitro diagnostic use as an aid in the diagnosis of inflammatory bowel disease (IBD), specifically Crohn's disease (CD) and ulcerative colitis (UC), and as an aid in the differentiation of IBD from irritable bowel syndrome (IBS) in conjunction with other clinical and laboratory findings.
Analyte	Fecal Calprotectin	Same
Type of Test	Quantitative	Same
Sample Matrix	Human stool	Same
Extraction Procedure	ALPCO Easy Stool Extraction Device or manual weighting method	Same
Calibrator / Control Analyte	Native human calprotectin	Same
Results Interpretation	Normal: < 50 µg/g Equivocal: 50–100 µg/g Abnormal: >100 µg/g	Same
General Device Characteristic Differences		
Methodology	Immunoturbidimetric Assay	Chemiluminescence ELISA
Analytical Measuring Interval (AMI)	11.0–1000.0 µg/g	7.9–6000 µg/g
Assay Process	Automated	Manual ELISA
Sample Dilution	1:100	1:25000
Instrument	Beckman AU680 analyzer	Not applicable
Capture Antibody	Mouse monoclonal IgG antibodies immobilized on latex particles	Mouse monoclonal IgG antibodies immobilized on plate

Detection Antibody	Not applicable	Biotinylated mouse monoclonal anti-calprotectin antibody
Calibrators	Ready to use 5 levels: 50, 200, 500, 800, 1050 µg/g	Lyophilized 8 levels: 0, 5, 20, 40, 156, 625, 2500, 10000 µg/g
Controls	Ready to use 2 levels: Control 1: 35 µg/g Control 2: 135 µg/g	Lyophilized, 3 levels: Control 1: 20–40 µg/g Control 2: 80–125 µg/g Control 3: 1250–2500 µg/g
Shelf-Life (Kit)	2–8°C for 24 months	2–8°C for 18 months

VI Standards/Guidance Documents Referenced:

- CLSI EP05-A3, Evaluation of Precision of Quantitative Measurement Procedures, Approved Guideline – Third Edition
- CLSI EP06-Ed2, Evaluation of the Linearity of Quantitative Measurement Procedures – Second Edition
- CLSI EP07, Interference Testing in Clinical Chemistry – Third Edition
- CLSI EP17-A2, Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline – Second Edition
- CLSI EP28-A3c, Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory, Approved Guideline – Third Edition
- Class II Special Controls Guidance Document: Fecal Calprotectin Immunological Test Systems - Guidance for Industry and FDA Staff, July (2006)

VII Performance Characteristics (if/when applicable):

A Analytical Performance:

All results met the manufacturer's pre-determined acceptance criteria.

1. Precision/Reproducibility:

Within-Laboratory Precision

A study was conducted in accordance with CLSI EP05-A3 to evaluate repeatability and within-lab imprecision of the ALPCO Calprotectin Immunoturbidimetric Assay. The study was performed at one site/laboratory, using one reagent lot and calibrators, one Beckman Coulter AU680 analyzer, and one operator. Seven stool extracts were prepared using the ALPCO Easy Stool Extraction Device from seven human stool samples. The extracts were analyzed in duplicate per run, two runs per day, for 20 days to generate a total of 80 replicates per sample. The results (within-run, between-run, between-day, and within-laboratory) are shown in the following table:

Sample	N	Mean (µg/g)	Within-run		Between-run		Between-day		Within-Lab	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV
1	80	14.4	1.7	12.1	1.5	10.2	1.4	9.9	2.2	15.0
2	80	37.6	2.2	5.9	1.9	5.1	1.5	3.9	2.5	6.7
3	80	81.3	2.1	2.6	2.2	2.7	2.3	2.8	3.1	3.9
4	80	173.8	3.3	1.9	4.2	2.4	6.5	3.7	7.5	4.3
5	80	414.5	4.2	1.0	6.7	1.6	9.2	2.2	10.8	2.6
6	80	622.5	5.5	0.9	6.8	1.1	13.8	2.2	15.1	2.4
7	80	854.8	9.0	1.1	9.5	1.1	12.5	1.5	15.6	1.8

Lot-to-Lot Precision

The lot-to-lot imprecision of the ALPCO Calprotectin Immunoturbidimetric Assay was evaluated in accordance with CLSI EP05-A3. The study was performed using three reagent lots and one lot of calibrators, one AU680 analyzer, at one site, by one operator. Seven human stool samples containing various fecal calprotectin concentrations were extracted using the ALPCO Easy Stool Extraction Device and analyzed in five replicates, once per day, for 5 days to generate a total of 75 replicates per sample. The results are summarized in the following table:

Sample	N	Mean (µg/g)	Within-run		Between-Day		Between-Lot		Total	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV
1	75	14.3	1.8	12.3	1.5	10.3	1.6	11.3	2.8	19.6
2	75	38.4	1.9	4.9	1.4	3.6	1.8	4.6	2.9	7.6
3	75	86.0	2.6	3.1	2.5	2.9	4.0	4.6	5.4	6.2
4	75	173.4	3.6	2.1	3.9	2.3	5.0	2.9	7.3	4.2
5	75	416.5	5.0	1.2	7.8	1.9	10.1	2.4	13.8	3.3
6	75	626.3	7.1	1.1	9.5	1.5	11.4	1.8	16.5	2.6
7	75	869.8	7.5	0.9	8.9	1.0	11.3	1.3	16.3	1.9

Site-to-Site Reproducibility

The study was performed using one reagent lot and one lot of calibrators, three AU680 analyzers at three sites, and three operators (one operator per instrument/site). Seven stool samples containing various fecal calprotectin concentrations were extracted using the ALPCO Easy Stool Extraction Device and analyzed in five replicates per run, one run per day, for 5 days to generate a total of 75 replicates per sample. The results are reported in the following tables:

Sample	N	Mean (µg/g)	Within-run		Between-Day		Between-Site		Total	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV
1	75	11.2	1.8	15.7	1.7	14.9	1.8	16.2	3.0	27.1
2	75	35.0	2.1	6.0	1.7	4.8	3.2	9.2	4.2	11.9
3	75	77.3	3.2	4.1	2.6	3.3	3.3	4.3	5.3	6.8
4	75	168.5	2.9	1.7	6.1	3.6	6.7	4.0	9.5	5.6
5	75	445.7	4.7	1.0	5.8	1.3	9.6	2.2	12.2	2.7

Sample	N	Mean (µg/g)	Within-run		Between-Day		Between-Site		Total	
			SD	%CV	SD	%CV	SD	%CV	SD	%CV
6	75	643.8	6.4	1.0	10.9	1.7	8.9	1.4	15.5	2.4
7	75	782.2	14.3	1.8	15.2	1.9	21.2	2.7	29.7	3.8

Extraction Precision:

Refer to K191807 for extraction device precision using the ALPCO Easy Stool Extraction Device.

2. Linearity:

The linearity study was performed following CLSI EP06-Ed2 using two sample panels. Sample Panel 1 was prepared by dilution of a high native sample pool (576.3 µg/g) with a low native sample pool (9.2 µg/g) to obtain 11 dilutions. Sample Panel 2 was prepared by dilution of a high native sample pool (1285.4 µg/g) with a low native sample pool (10.8 µg/g) to obtain 11 dilutions. The samples were tested in triplicate, using one reagent lot. The expected value, predicted value, and deviation from predicted were calculated. The results are summarized in the following tables:

Sample Panel 1

Sample Dilution	Mean Measured Value (µg/g)	Expected Value (µg/g)	Predicted Value (µg/g)	Deviation (µg/g)	% Deviation
1	9.2	9.2	9.3	-0.1	-0.9%
2	60.5	65.9	66.3	-5.8	-8.8%
3	113.6	122.6	123.3	-9.7	-7.9%
4	166.1	179.3	180.4	-14.3	-7.9%
5	221.6	236.0	237.4	-15.9	-6.7%
6	275.1	292.8	294.5	-19.4	-6.6%
7	341.5	349.5	351.5	-10.0	-2.9%
8	417.0	406.2	408.6	8.4	2.1%
9	471.9	462.9	465.6	6.3	1.3%
10	543.0	519.6	522.7	20.3	3.9%
11	576.3	576.3	579.7	-3.5	-0.6%

Sample Panel 2

Sample Dilution	Mean Measured Value (µg/g)	Expected Value (µg/g)	Predicted Value (µg/g)	Deviation (µg/g)	% Deviation
1	10.4	10.8	10.5	0.0	-0.4%
2	95.4	100.0	97.0	-1.6	-1.7%
3	185.3	199.4	193.5	-8.2	-4.2%
4	290.3	304.0	294.8	-4.6	-1.5%
5	398.8	399.6	387.6	11.2	2.9%
6	487.3	501.5	486.5	0.8	0.2%
7	604.5	603.5	585.4	19.2	3.3%
8	716.3	702.3	681.2	35.1	5.1%
9	798.8	801.1	777.0	21.8	2.8%

Sample Dilution	Mean Measured Value (µg/g)	Expected Value (µg/g)	Predicted Value (µg/g)	Deviation (µg/g)	% Deviation
10	887.9	899.8	872.8	15.1	1.7%
11	995.9	998.6	968.6	27.3	2.8%

Linear Regression Summary

Dilution Panel	Range (µg/g)	Slope (95%CI)	Intercept (95%CI)	R ²	Recovery Range (%)
Sample Panel 1	9.2–576.3	1.03 (0.99–1.07)	-12.74 (-27.18–1.71)	0.997	91.2–103.9
Sample Panel 2	10.4–995.9	1.01 (0.99–1.03)	-6.95 (-18.40–4.49)	0.999	93.1–102.2
Combined	9.2–995.9	1.04 (1.02–1.06)	-12.19 (-20.28– -4.10)	0.999	91.2–103.9

Linearity was supported for the analytical measuring interval from 11.0 µg/g to 1000.0 µg/g.

Spiking Recovery:

The recovery of the ALPCO Calprotectin Immunoturbidimetric Assay was evaluated using seven extracted stool samples containing various concentrations of calprotectin across the analytical measuring interval of the assay. Purified native calprotectin diluted in standard diluent was used as the spiking material in a proportion of nine parts of sample to one part spiking material. Recovery was calculated compared to the baseline result. Percent recoveries ranged from 95.5% to 106.5%.

Hook Effect

Hook effect was evaluated by testing samples above the analytical measuring interval: from 2,000 to 24,000 µg/g. One lot of reagents was used; five samples were run in duplicate. The results showed no hook effect up to 18,000 µg/g.

3. Analytical Specificity/Interference:

The interference testing of the ALPCO Calprotectin Immunoturbidimetric Assay was evaluated following CLSI guideline EP07-Ed3. The study was performed using three samples: one high positive stool sample (477–488 µg/g, one moderately high positive stool sample (261–287 µg/g), and one low positive stool sample (41–46 µg/g). Interfering substances were spiked into each stool sample at 10% of the total specimen volume; a series of endogenous and exogenous interfering substances of different concentrations were tested. The solvent used to create the interferent stock was spiked into each stool sample at 10% of the total specimen volume to generate a control sample. Samples were tested in duplicate. Percent recovery was calculated to evaluate the differences between blank and control samples. No interference (% recovery range within 90–110%) of the ALPCO Calprotectin Immunoturbidimetric Assay was found to the concentration level tested for the following substances:

Interference Substance	Spiking Concentration
Oral Pharmaceuticals	
Iron (II) sulfate	0.11 mg/50 mg stool
Prednisone	0.31 mg/50 mg stool
Azathioprine	0.19 mg/50 mg stool
Mesalamine; 5-ASA	5.21 mg/50 mg stool
Lansoprazol	0.18 mg/50 mg stool
Vancomycin	2.0 mg/50 mg stool
Sulfamethoxazole	1.6 mg/50 mg stool
Trimethoprim	0.35 mg/50 mg stool
Ciprofloxacin	1.25 mg/50 mg stool
Nutritional Supplements	
Vitamin E (DL- α Tocopherol Acetate)	0.3 mg/50 mg stool
Multi-Vitamins (A, B1, B2, B3, B5, B6, B8, B9, B12, C, D, and E) and minerals	1.06 mg/50 mg stool
Vitamin B7 (biotin)	1750 ng/50 mg stool
Microorganisms	
<i>Escherichia coli</i>	1.5×10^7 cfu/mL
<i>Salmonella enterica subsp. enterica</i>	1.5×10^7 cfu/mL
<i>Klebsiella pneumoniae subsp. pneumonia</i>	1.5×10^7 cfu/mL
<i>Citrobacter freundii</i>	1.5×10^7 cfu/mL
<i>Shigella flexneri</i>	1.5×10^7 cfu/mL
<i>Yersinia enterocolitica subsp. enterocolitica</i>	1.5×10^7 cfu/mL
Other	
Human anti-mouse antibody (HAMA)	0.241 μ g/50 mg stool
Hemoglobin	7.5 mg/50 mg stool
Rheumatoid Factor (RF)	500 IU/50 mg stool

4. Assay Reportable Range:

The analytical measuring interval (AMI) of the ALPCO Calprotectin Immunoturbidimetric Assay is 11.0–1000.0 μ g/g.

5. Traceability, Stability, Expected Values (Controls, Calibrators, or Methods):

Traceability:

No international reference materials or reference measurement procedures are currently available for calprotectin. The ALPCO Calprotectin Immunoturbidimetric Assay is traceable to an internal standard that was value assigned using the ALPCO Calprotectin Chemiluminescent ELISA (K191807) predicate device. Calprotectin calibrators are prepared by spiking human recombinant calprotectin in a buffered matrix with stabilizers and preservatives to the target calprotectin concentrations.

Stability:

a. Shelf-Life Stability

The real-time and accelerated reagent stability studies were performed using two native stool sample pools and two controls. For real-time stability, three reagent lots of the ALPCO Calprotectin Immunoturbidimetric Assay (reagents, calibrators, and controls) were stored at 2–8°C, and the samples were tested every 4 months. For accelerated stability, three lots of the reagents were stored at 25°C, and 37°C, and the samples were tested every 7 days up to 35 days. Percent (%) recovery was calculated for each sample at each testing timepoint compared to the value at Day 0 stored under the same condition. The results support that the ALPCO Calprotectin Immunoturbidimetric Assay is stable at 2–8°C for 24 months.

b. Reagent On-Board Stability

Reagent on-board stability was tested using one lot of the ALPCO Calprotectin Immunoturbidimetric Assay reagents on one instrument (AU680 chemistry analyzer) using one calibration event. Two native sample pools and two controls were tested every 7 days using this lot of reagents stored on-board. The percent (%) recovery was evaluated for each sample compared to its value at Day 0. The results support the ALPCO Calprotectin Immunoturbidimetric Assay is stable when stored on-board for 30 days.

c. Transport Stability

The Calprotectin Immunoturbidimetric Assay (reagents, calibrators, and controls) transport stability was performed using one lot of reagents, one lot of calibrators, one instrument, and one lot of controls. Two sample pools and two controls were tested in two replicates. The reagents, calibrators, and controls were shipped using standards shipping packaging with two ice packs via commercial shipping for 2 days during peak summer conditions, and upon arrival the assay was stored at 30°C for an additional 14 days. Percent (%) recovery was calculated to evaluate the difference between transported and not transported reagents. Reagent percent (%) recovery ranged 91.0%–103.4%; for the calibrators: 94.9%–104.8%; and for the controls: 98.3%–101.7%.

d. Sample Stability

Refer to K191807

6. Detection Limit:

Limit of Blank (LoB):

The LoB of the ALPCO Calprotectin Immunoturbidimetric Assay y was determined in accordance with CLSI EP17-A2. Four blank stool samples were extracted and tested in replicates of five on two reagent lots across 3 days. The LoB was estimated as the 95th percentile for the distribution of blank sample results using nonparametric method. The LoB was determined to be 2.2 µg/g.

Limit of Detection (LoD):

The LoD of the ALPCO Calprotectin Immunoturbidimetric Assay was determined in accordance with CLSI EP17-A2. Four stool samples containing a low level of calprotectin

were extracted and tested in replicates of five with two reagent lots across 4 days. The reported LoD measurement is the maximum value of the individual LoDs obtained for each reagent lot. The LoD was determined to be 3.9 µg/g.

Limit of Quantitation (LoQ):

The LoQ of the ALPCO Calprotectin Immunoturbidimetric Assay was determined in accordance with CLSI EP17-A2. Five stool samples containing a low level of calprotectin were tested in replicates of five on two reagent lots across 4 days. The LoQ estimate for each reagent lot was determined as the measurand concentration interpolated from the intersection of the precision profile curve with the total within-laboratory imprecision goal of 20 %CV. The maximum value obtained for each reagent lot was reported as the LoQ. The LoQ was determined to be 11.0 µg/g.

7. Assay Cut-Off:

The cut-offs of the ALPCO Calprotectin Immunoturbidimetric Assay are presented in the following table:

Calprotectin Concentration	Interpretation	Follow-Up
< 50 µg/g	Normal	None
50–100 µg/g	Equivocal	Retest in 4–6 weeks
> 100 µg/g	Elevated	Repeat as clinically indicated

B Comparison Studies:

1. Method Comparison with Predicate Device:

A method comparison study was performed by testing 306 stool sample extracts in singlicate using the ALPCO Calprotectin Immunoturbidimetric Assay and the predicate device. Of the 306 samples, 168 tested within the analytical measuring intervals of both assays and further evaluated by qualitative agreement and regression analysis. Among the 168 samples, 12 samples were native samples spiked with purified native human calprotectin to supplement high values. The comparison results were summarized in the following tables:

ALPCO Calprotectin Immunoturbidimetric Assay	Predicate Device			Total
	Elevated (>100 µg/g)	Equivocal (50–100 µg/g)	Normal (<50 µg/g)	
Elevated (>100 µg/g)	80	6	3	89
Equivocal (50–100 µg/g)	6	18	7	31
Normal (<50 µg/g)	0	8	40	48
Total	86	32	50	168

Equivocal results considered as elevated:

Negative Percent Agreement (NPA)	80.00% (95%CI: 68.28–89.97%)
Positive Percent Agreement (PPA)	93.22% (95% CI: 87.08–97.03%)
Total Percent Agreement (TPA)	89.29% (95%CI: 83.60–93.53%)

Equivocal results considered as normal:

Negative Percent Agreement (NPA)	89.02% (95% CI: 80.18–94.86%)
Positive Percent Agreement (PPA)	93.02% (95% CI: 84.43–97.40%)
Total Percent Agreement (TPA)	91.07% (95% CI: 85.70–94.92%)

The 168 samples within the analytical measuring intervals of both assays were evaluated by Passing-Bablok regression. The results are summarized in the table below:

N	Range (µg/g)	Slope (95%CI)	Intercept (95%CI)	R ²
168	11.0–975.3	0.98 (0.95–1.02)	-0.52 (-4.74–3.23)	0.99

2. Sample Extraction Method Comparison:

Extraction method comparison was performed to demonstrate equivalency of the ALPCO Easy Stool Extraction Device and manual weighing method. A total of 60 native stool samples were extracted using the ALPCO Easy Stool Extraction Device and manual weighing method, respectively. The paired samples were run using one lot of the ALPCO Calprotectin Immunoturbidimetric Assay reagent, calibrators, and controls. Test results of 53 paired samples across the analytical measuring interval were evaluated by Passing-Bablok regression analysis. The results are summarized in the following table:

N	Range (µg/g)	Slope (95%CI)	Intercept (95%CI)	R ²
53	11.2–780.8	1.01 (0.99–1.03)	0.71 (-0.41–2.30)	0.99

C Clinical Studies:

a. *Clinical Sensitivity and specificity*

The clinical performance of the Calprotectin Immunoturbidimetric Assay was evaluated. A total of 349 stool specimens from the intended use population were prospectively collected at 15 sites. All enrolled patients were scheduled to undergo colonoscopy and were further screened to determine eligibility based upon the below inclusion/exclusion criteria:

- *Inclusion criteria:*
 ≥ 22 years of age; possible inflammatory bowel disease or irritable bowel syndrome with symptoms such as abdominal pain, diarrhea, altered appetite, weight loss, or anemia; able to provide a stool sample within 24–48 hours prior to colonoscopy preparation or 3–14 days following colonoscopy; able to understand the study and the tasks required and sign the Informed Consent Form.

- *Exclusion criteria:*
Unable or unwilling to provide a stool specimen; known intestinal cancer; receiving chemotherapy or systemic immunosuppressive drugs; previously diagnosed IBD managed with immunomodulators, 5-ASA or biologic therapies or who have undergone a surgical resection or diversion procedure; known intestinal infection; known upper GI disease (i.e., esophagitis or gastritis); have taken proton pump inhibitors (PPIs) or H2-receptor antagonists for the previous 2 weeks; have taken NSAIDs (including aspirin) on seven or more days during the 2 weeks before providing the sample.

Clinical diagnoses were established by the clinical investigator/gastroenterologist:

- IBD diagnosis was based on endoscopy results and/or histology of biopsies taken during the endoscopy.
- IBS diagnosis was based on the Rome IV criteria and confirmed by negative endoscopy including the colon and terminal ileum.
- Subjects were diagnosed with “Other GI conditions” when they did not meet the diagnostic criteria for IBD or IBS.

The clinical diagnosis of the samples in the study were summarized in the following table:

Clinical Diagnosis	Number of Subjects
Inflammatory Bowel Disease (IBD)	63
Ulcerative Colitis (UC)	29
Crohn’s Disease (CD)	24
Indeterminant/Undefined	10
Irritable Bowel Syndrome (IBS)	105
Other GI conditions	181
Total	349

The stool samples were de-identified and were tested using the ALPCO Calprotectin Immunoturbidimetric Assay at ALPCO site and at two external sites. The results are summarized in the following table:

ALPCO Calprotectin Immunoturbidimetric Assay	Clinical Diagnosis			Total
	IBD	IBS	Other GI	
Elevated (>100 µg/g)	48	4	3	55
Equivocal (50–100 µg/g)	9	6	6	21
Normal (<50 µg/g)	6	95	172	273
Total	63	105	181	349

IBD vs. non-IBD

The performance of the Calprotectin Immunoturbidimetric Assay as an aid in the diagnosis of IBD was evaluated by calculate the estimates of the clinical sensitivity, clinical specificity, positive predictive value (PPV), and negative predictive value (NPV). The results are summarized in the following tables:

Equivocal results as Elevated

ALPCO Calprotectin Immunoturbidimetric Assay	Clinical Diagnosis		Total
	IBD	Non-IBD	
Elevated (≥ 50 $\mu\text{g/g}$)	57	19	76
Normal (< 50 $\mu\text{g/g}$)	6	267	273
Total	63	286	349
Sensitivity: 90.5% (95% CI: 80.7–95.6%) Specificity: 93.4% (95% CI: 89.9–95.7%) PPV: 75.0% (95% CI: 65.8–82.4%) NPV: 97.8% (95% CI: 95.4–99.0%)			

Equivocal results as Normal:

ALPCO Calprotectin Immunoturbidimetric Assay	Clinical Diagnosis		Total
	IBD	Non-IBD	
Elevated (> 100 $\mu\text{g/g}$)	48	7	55
Normal (≤ 100 $\mu\text{g/g}$)	15	279	294
Total	63	286	349
Sensitivity: 76.2% (95% CI: 64.4–85.0%) Specificity: 97.6% (95% CI: 94.8–98.9%) PPV: 87.3% (95% CI: 76.5–93.5%) NPV: 94.9% (95% CI: 91.8–96.9%)			

IBD vs. IBS

The performance of the Calprotectin Immunoturbidimetric Assay as an aid to differentiate the IBD and IBS was evaluated by calculate the estimates of the clinical sensitivity, clinical specificity, positive predictive value (PPV), and negative predictive value (NPV). The results are summarized in the following tables:

Equivocal results as Elevated:

ALPCO Calprotectin Immunoturbidimetric Assay	Clinical Diagnosis		Total
	IBD	IBS	
Elevated (≥ 50 $\mu\text{g/g}$)	57	10	67
Normal (< 50 $\mu\text{g/g}$)	6	95	101
Total	63	105	168
Sensitivity: 90.5% (95% CI: 80.7–95.6%) Specificity: 90.5% (95% CI: 83.4–94.7%) PPV: 85.1% (95% CI: 75.9–91.2%) NPV: 94.1% (95% CI: 88.1–97.1%)			

Equivocal results as Normal:

ALPCO Calprotectin Immunoturbidimetric Assay	Clinical Diagnosis		Total
	IBD	IBS	
Elevated (> 100 $\mu\text{g/g}$)	48	4	52
Normal (≤ 100 $\mu\text{g/g}$)	15	101	116
Total	63	105	168

Sensitivity: 76.2% (95% CI: 64.4–85.0%) Specificity: 96.2% (95% CI: 90.6–98.5%) PPV: 92.3% (95% CI: 82.0–96.9%) NPV: 87.1% (95% CI: 81.2–91.3%)
--

2. Other Clinical Supportive Data (When 1 Is Not Applicable):

Not applicable

D Clinical Cut-Off:

Same as assay cut-off

E Expected Values/Reference Range:

A study was performed using 120 stool samples to establish the reference range for the ALPCO Calprotectin Immunoturbidimetric Assay following CLSI EP28-A3c. The stool samples were obtained from 120 asymptomatic individuals with no abdominal complaints and no history of IBS, IBD or other chronic intestinal disorders. The cohort contained 65 females and 55 males, 22 to 82 years old with the mean age of 43.8 years.

Among 120 samples, the value of stool calprotectin level for 104 samples were below 11 µg/g, the lower end of the AMI. Two samples were tested with elevated value of stool calprotectin, i.e., 111.2 and 130.3 µg/g, and one with equivocal, value 68.5 µg/g. The results are summarized in the following table:

	Calprotectin Concentration (µg/g)
Range	<3.9–130.3
Mean	8.6
Median	3.9
Lower Limit (90% CI)	3.9 (3.9–3.9)
Upper Limit (90% CI)	67.5 (20.0–130.3)

VIII Proposed Labeling:

The labeling supports the finding of substantial equivalence for this device.

IX Conclusion:

The submitted information in this premarket notification is complete and supports a substantial equivalence decision.