

K122711

**510(k) Summary**

**Submission Date:** Prepared September 4, 2012, updated September 17, 2013

**Submitter Information:**

*Company Name:* Nu Skin Enterprises, Inc.

*Company Address:* 75 West Center Street  
Provo, UT 84601

*Contact Person:* Tyler Whitehead

SEP 17 2013

**Device Information:**

*Trade Name:* Facial Spa

*Common Name:* Facial Toning Device

*Classification Name:* Transcutaneous electrical nerve stimulator 21 CFR 882.5890

*Regulation Class:* Class II

*Product Code:* NFO

**Predicate Devices:** NuFace® Plus (K103472)  
Carol Cole Company  
Class II

**Device Description:** The Facial Spa is a hand-held, battery-powered device used with conductive treatment gel to stimulate the face with microcurrents.

**Intended Use:** The Facial Spa is intended for the stimulation of healthy facial skin.

**Indications for Use:** The Facial Spa is indicated for facial stimulation for over-the-counter cosmetic use.

**Compliance Data:** The Facial Spa conforms to the applicable requirements of the following standards:

- IEC 60601-1 Issued: 1998/01/01 Ed: 2 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance; Amendment 1-1991, Amendment 2-1995, Corrigendum – 1995

- IEC 60601-1-2 Issued 2012/06/01 Ed: 3 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic compatibility – Requirements and tests
- IEC 60601-2-10 Issued:1987/01/01 Ed:1 Medical Electrical Equipment Part 2: Particular Requirements for the Safety of Nerve and Muscle Stimulators
- IEC 60601-1-4 Issued: 2000/04/01 Ed:1.1 Medical Electrical Equipment - Part 1-4: General Requirements for Safety – Collateral Standard: Programmable Electrical Medical Systems Consolidated with Amd. 1:1999
- ISO ISO10993-10 / Ed.3, Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization
- ISO 10993-12 Biological evaluation of medical devices - Part 12: Sample Preparation and Reference Material

**Comparison to Predicate Device:**

The Facial Spa is identical to its predicate device in intended use and indications for use. It uses the same technological principle as the predicate device to accomplish its intended use, namely transcutaneous electrical stimulation. Both are hand-held, battery-powered devices. The present device raises no new questions of safety or effectiveness. A technical comparison of the two devices has been completed and supports a substantial equivalence determination.

**Table: Comparison of Nu Skin Facial Spa and Carol Cole NuFace® Plus**

<b>Features</b>	<b>Nu Skin Facial Spa (Subject Device)</b>	<b>Carol Cole NuFace® Plus (Predicate Device)</b>
<b>Similarities</b>		
Regulation Number	21 C.F.R. § 882.5890	21 C.F.R. § 882.5890
Regulation Name	Transcutaneous Electrical Nerve Stimulator	Transcutaneous Electrical Nerve Stimulator
Regulatory Class	Class II	Class II
Product Code	NFO	NFO
Intended Use	Stimulate the face	Stimulate the face

Indications for Use	Over-the-Counter Cosmetic Use	Over-the-Counter Cosmetic Use
Power Source	Battery operated	Battery operated
Average DC current through electrodes when device is on but no pulse being applied	0	0
Number of output channels	1	1
Microprocessor controlled	Yes (PLC)	Yes
Automatic no-load trip	Yes	Yes
Automatic shut off	Yes	Yes
User over-ride control	Yes	Yes
Indicator Display		
<i>On/off</i>	Yes	Yes
<i>Low-battery</i>	Yes	Yes
On Time (seconds)	Constant	Constant
Off Time (seconds)	None	None
Compliance with Voluntary Standards	IEC 60601-1 IEC 60601-1-2	IEC 60601-1 IEC 60601-1-2
Compliance with 21 C.F.R. § 898	Yes	Yes
<b>Differences</b>		
Modes	One to Five	One
Maximum Output Voltage	@ 500 Ω Mode 1: 188 mV (214 mV) Mode 2: 189 mV (225 mV) Mode 3: 185 mV (220 mV) Mode 4: 184 mV (212 mV) Mode 5: 184 mV (213 mV)	@ 500 Ω 137 mV
	@ 2k Ω Mode 1: 722 mV (1527 mV) Mode 2: 746 mV (812 mV) Mode 3: 735 mV (815 mV) Mode 4: 740 mV (802 mV) Mode 5: 741 mV (801 mV)	@ 2k Ω 769 mV
	@ 10 k Ω Mode 1: 2.626 V (2.978 V) Mode 2: 3.787 V (4.011 V) Mode 3: 3.747 V (4.014 V) Mode 4: 3.503 V (3.980 V) Mode 5: 3.729 V (3.980 V)	@ 10 k Ω 3.82 V

Maximum Output Current	@ 500 $\Omega$ Mode 1: 376 $\mu$ A (427 $\mu$ A) Mode 2: 379 $\mu$ A (449 $\mu$ A) Mode 3: 371 $\mu$ A (440 $\mu$ A) Mode 4: 368 $\mu$ A (424 $\mu$ A) Mode 5: 369 $\mu$ A (427 $\mu$ A)	@ 500 $\Omega$ 274 $\mu$ A
	@ 2 k $\Omega$ Mode 1: 361 $\mu$ A (763 $\mu$ A) Mode 2: 373 $\mu$ A (406 $\mu$ A) Mode 3: 367 $\mu$ A (407 $\mu$ A) Mode 4: 370 $\mu$ A (401 $\mu$ A) Mode 5: 370 $\mu$ A (400 $\mu$ A)	@ 2 k $\Omega$ 387 $\mu$ A
	@ 10 k $\Omega$ Mode 1: 362 $\mu$ A (398 $\mu$ A) Mode 2: 379 $\mu$ A (401 $\mu$ A) Mode 3: 375 $\mu$ A (401 $\mu$ A) Mode 4: 350 $\mu$ A (398 $\mu$ A) Mode 5: 373 $\mu$ A (398 $\mu$ A)	@ 10 k $\Omega$ 383 $\mu$ A
Waveform	Direct Current	Pulsed Monophasic
Weight	120 g (approx. 4 oz)	9 oz
Dimensions	<i>Facial Spa with Large Conductor:</i> 143 x 31.2 x 67.5 mm (5.6" x 1.2" x 2.7")  <i>Facial Spa with Small Conductor:</i> 136 x 31.2 x 67.5 mm (5.4" x 1.2" x 2.7")	3" x 5.25" x 1.25"
Housing material and construction	ABS, Chrome plated ABS, rubber (elastomer)	Thermo Plastic

**Conclusion:**

The Facial Spa is substantially equivalent to the predicate device, as the devices share a common intended use, and technological differences between the Facial Spa and the predicate do not raise new questions of safety or effectiveness.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

September 17, 2013

NuSkin Enterprises, Inc.  
C/O Paul D. Rubin  
Ropes & Gray LLP  
One Metro Center  
700 12<sup>th</sup> Street, Suite 900  
Washington, DC 20005

Re: K122711

Trade/Device Name: NuSkin Facial Spa  
Regulation Number: 21 CFR 882.5890  
Regulation Name: Transcutaneous electrical nerve stimulator  
Regulatory Class: Class II  
Product Code: NFO  
Dated: September 4, 2013  
Received: September 4, 2013

Dear Mr. Rubin,

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA).

You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set

forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Joyce M. Whang -S**

for Victor Krauthamer, Ph.D.  
Acting Director  
Division of Neurological and Physical  
Medicine Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known): K122711

Device Name: NuSkin Facial Spa

Indications For Use:

The Facial Spa is indicated for facial stimulation for over-the-counter cosmetic use.

Prescription Use \_\_\_\_\_  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use  X   
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF  
NEEDED)

---

Concurrence of Center for Devices and Radiological Health (CDRH)

Joyce M. Whang -S