

510(K) Summary

Submitter: Shaser, Inc.
10 Maguire Road
Building 1
Lexington, MA 02421

JUN 13 2013

Contact: Anthony Burns
Sr. Director of Regulatory Affairs

Date Summary Prepared: December 31, 2012

Device Trade Name: Shaser V-MINI Hair Removal System

Common Name: Light based hair removal device

Classification Name: Powered Light Based Non-Laser Surgical Instrument with Thermal Effect. Product Code ONF

Equivalent Device: Shaser HRS2 Hair Removal System

Device Description: Over-The-Counter, Cordless, Rechargeable, Personal Light-Based Hair Removal System For Permanent Hair Reduction

Emission activation is by fingerswitch. Includes limited life treatment head and battery charger. Overall weight of the device is 0.45 Kg, and the size is 2.1 x 0.6 x 0.8 cm (HxWxD).

Charger Electrical Requirement is 115 VAC, 15A, 50-60 Hz, single phase.

Intended Use: The SHASER V-MINI hair removal system is an over the counter device intended for removal of unwanted hair. It is also intended for permanent reduction in unwanted hair. Permanent hair reduction is defined as the long-term stable reduction in the number of hairs regrowing when measured at 6, 9, and 12 months after the completion of a treatment regimen.

Comparison: The Shaser V-MINI Hair Removal System has the same indication for uses, the same principle of operation, the same pulse energy range and wavelength range as the predicate device.

Nonclinical Performance Data: none

Clinical Performance Data: Label comprehension and usability test of consumers' ability to understand the instructions for use and to evaluate their ability to use the device safely in a simulated OTC home-use environment.

- 150 study subject were tested for label comprehension and 123 study subjects tested for usability. Both test populations included low literacy subjects.

The results of the two tests confirms sufficient label comprehension and safe and appropriate use of the device.

Conclusion: The Shaser V-MINI Hair Removal System is a safe and effective device for the intended uses.

Additional Information: None



Food and Drug Administration
10903 New Hampshire Avenue
Document Control Center - WQ66-G609
Silver Spring, MD 20993-002

June 13, 2013

Shaser, Inc.
% Ms. Sarah Hackett
10 Maguire Road, Suite 120
Building One
Lexington, Massachusetts 02421

Re: K130015

Trade/Device Name: Shaser V-Mini Hair Removal System
Regulation Number: 21 CFR 878.4810
Regulation Name: Laser surgical instrument for use in general and
Plastic surgery and in dermatology
Regulatory Class: Class II
Product Code: ONF
Dated: May 06, 2013
Received: May 08, 2013

Dear Ms. Hackett:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA).

You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please go to <http://www.fda.gov/AboutFDA/CentersOffices/CDRH/CDRHOices/ucm115809.htm> for the Center for Devices and Radiological Health's (CDRH's) Office of Compliance. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

David Krause -S

for Mark N. Melkerson
Acting Director
Division of Surgical Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications for Use

510(k) Number (if known): K130015

Device Name: Shaser V-MINI Hair Removal System

Indications For Use:

The SHASER HRS3 hair removal system is an over the counter device intended for removal of unwanted hair. It is also intended for permanent reduction in unwanted hair. Permanent hair reduction is defined as the long-term stable reduction in the number of hairs regrowing when measured at 6, 9, and 12 months after the completion of a treatment regimen.

Prescription Use
(Part 21 CFR 801 Subpart D)

OR

Over-The-Counter Use √
(Part 21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Neil R Ogden
2013.06.11 16:24:42 -04'00'

(Division Sign-Off) for MXM
Division of Surgical Devices
510(k) Number K130015