



AcceleDent Aura™  
Special 510(k) NOTIFICATION

K13 0643

APR 23 2013

## 510(k) Summary

AcceleDent® Aura

April 11, 2013

This summary of Special 510(k) substantial equivalence determination is being submitted in accordance with the requirements of 21 CFR part 807.92.

510(k) Submitter: OrthoAccel® Technologies, Inc.  
6575 West Loop South, Suite 200  
Bellaire, TX 77401  
Phone: 832-803-0339  
Fax: 713-583-9972

Contact: Zaffer Syed

Device trade name: **AcceleDent® Aura**  
Common name: Orthodontic plastic bracket  
Classification name: Orthodontic Vibratory Accessory  
Regulation number: 21 CFR 872.5470  
Classification: Class II  
Panel: Dental  
Product code: OYH  
Predicate Device: AcceleDent®, 510(k) Number K110661

### Device Description:

AcceleDent® Aura is an orthodontic accessory for the treatment of tooth malocclusion. It is used as an adjunctive therapy for patients with orthodontic appliances such as braces to help facilitate tooth movement. AcceleDent® Aura should be used by patients for twenty minutes per day in conjunction with standard orthodontic treatment.

AcceleDent® Aura includes the Activator, Mouthpiece and Travel Case. The Activator and connected Mouthpiece are used by patients to provide a light vibration to the teeth – the Activator vibrates at a 0.25 N (25 grams) force level and 30 Hz frequency for 20 minutes when turned-on; the vibration is transmitted from the Activator through the Mouthpiece to the patient's teeth as they lightly bite down on the Mouthpiece.

The Activator and Mouthpiece assembly is light, comfortable, hands-free, and can be used while multi-tasking or while engaged in a variety of other daily activities. The Travel Case is an enclosure that may be used to



AcceleDent Aura™  
Special 510(k) NOTIFICATION

---

conceal, protect, and keep the AcceleDent® Aura Activator and Mouthpiece clean while not in use. The device includes a USB port, which can connect directly into a computer or power supply to recharge the battery. The USB port can also be connected to a computer to display usage data. A USB Cable and Power Adaptor are included to complete the system.

AcceleDent® Aura is a modified version of the predicate FDA-cleared AcceleDent® [(510(k) Number K110661)]. The fundamental scientific technology of delivering therapeutic vibrations to teeth and the intended use have not changed with the subject AcceleDent® Aura device. A comparison table of the AcceleDent® Aura to the predicate AcceleDent® is shown below.

**Comparison of AcceleDent® and AcceleDent® Aura**

	<b>AcceleDent® (Predicate Device)</b>	<b>AcceleDent® Aura (Subject Device)</b>
<b>Indication for Use</b>	AcceleDent® is an orthodontic accessory intended for use during orthodontic treatment. It is used in conjunction with orthodontic appliances such as braces and helps facilitate minor anterior tooth movement.	AcceleDent® Aura is an orthodontic accessory intended for use during orthodontic treatment. It is used in conjunction with orthodontic appliances such as braces and helps facilitate minor anterior tooth movement.
<b>Regulation Number</b>	21 CFR 872.5470	21 CFR 872.5470
<b>Product Code</b>	OYH	OYH
<b>Device Class</b>	II	II
<b>510(k) Number</b>	K110661	K130643
<b>For use with orthodontics</b>	Yes	Yes
<b>Material</b>	Elastomer	Elastomer
<b>Duration of Use</b>	20 minutes per day during Orthodontic treatment	20 minutes per day during Orthodontic treatment
<b>Power Source</b>	Lithium Polymer Battery	Lithium Polymer Battery
<b>Output Force</b>	0.25 N (25 g)	0.25 N (25 g)
<b>Frequency</b>	30 Hz	30 Hz
<b>Weight – Activator (grams)</b>	65	33
<b>Dimensions – Activator (HxLxW –mm-)</b>	76x41x28	79x36x30
<b>Rechargeable</b>	Yes	Yes
<b>Automatic Timer</b>	Yes	Yes
<b>Usage Data</b>	Yes – Displayed on charging port	Yes – Displayed on PC (FastTrac Usage Report)
<b>Activator</b>	Yes	Yes
<b>Mouthpiece Attachments</b>	Same	Same
<b>Storage</b>	Travel Shell	Travel Case
<b>USB Connector</b>	No	Yes
<b>Charging Port</b>	Yes	No
<b>Travel Shell</b>	Yes	No
<b>Travel Case</b>	No	Yes
<b>LED Battery Charge Indicator</b>	Yes	Yes
<b>Audible On/Off Switch</b>	No	Yes
<b>Shelf Life</b>	2.0 years	2.0 years
<b>Battery Life</b>	2.5 years	2.0 years*



AcceleDent Aura™

Special 510(k) NOTIFICATION

---

*\*Device deactivates after equivalent of 18 months of 20 minute daily sessions.*

**Intended Use:**

AcceleDent® Aura is intended for use during orthodontic treatment. It is used in conjunction with orthodontic appliances such as braces and helps facilitate minor anterior tooth movement.

**Substantial Equivalence Discussion:**

A thorough design control review of the modifications made in the AcceleDent® Aura has demonstrated that the device will perform safely and as intended when compared to the AcceleDent®. As a result, the AcceleDent® Aura is substantially equivalent to the predicate AcceleDent®.

**Performance Testing:**

Risk Analysis was performed according to ISO 14971 and the results supported implementation of the modifications. Verification testing to IEC 60601 and device specifications supplemented the analysis performed. Results of the analysis and testing showed that performance of the modifications met device specifications and further demonstrated substantial equivalence.

**Clinical Testing:**

No further clinical data were collected for AcceleDent® Aura.

**Conclusions:**

Based on the information provided in this premarket notification, it is concluded that the AcceleDent® Aura is safe and effective and is substantially equivalent to the predicate device which was the original AcceleDent®.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

April 23, 2013

Mr. Zaffer Syed  
Director, Clinical & Regulatory Affairs  
OrthoAccel® Technologies, Incorporated  
6575 West Loop South, Suite 200  
BELLAIRE TX 77401

Re: K130643  
Trade/Device Name: AcceleDent® Aura  
Regulation Number: 21 CFR 872.5470  
Regulation Name: Orthodontic Plastic Bracket  
Regulatory Class: II  
Product Code: OYH  
Dated: March 21, 2013  
Received: March 26, 2013

Dear Mr. Syed:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Kwame O. Ulmer for  
-S 

Anthony D. Watson, B.S., M.S., M.B.A.  
Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

OrthoAccel

AcceleDent Aura™  
Special 510(k) NOTIFICATION

---

**Indication for Use**

510(k) Number (if known): ~~Not Assigned~~ K130643

Device Name: AcceleDent® Aura

Indication For Use:

AcceleDent® Aura is intended for use during orthodontic treatment. It is used in conjunction with orthodontic appliances such as braces and helps facilitate minor anterior tooth movement.

Prescription Use   X    
(21 CFR Part 801 Subpart D)

AND/OR

Over the Counter Use \_\_\_\_\_  
(21 CFR Part 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE; CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office Device Evaluation (ODE)

 Mary S. Runner -S  
2013.04.23 12:11:09  
-04'00'

Division Sign-Off

510(k) K130643