



June 19, 2018

Johnson & Johnson Consumer, Inc.  
Cindy Abraham  
Associate Director, Regulatory Affairs  
7050 Camp Hill Road  
Fort Washington, Pennsylvania 19034

Re: K180856

Trade/Device Name: Neutrogena Light Therapy Aging Mask+  
Regulation Number: 21 CFR 878.4810  
Regulation Name: Laser Surgical Instrument For Use In General And Plastic Surgery And In Dermatology  
Regulatory Class: Class II  
Product Code: OHS  
Dated: March 21, 2018  
Received: April 2, 2018

Dear Cindy Abraham:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good

manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<https://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/>) and CDRH Learn (<http://www.fda.gov/Training/CDRHLearn>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<http://www.fda.gov/DICE>) for more information or contact DICE by email ([DICE@fda.hhs.gov](mailto:DICE@fda.hhs.gov)) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

Jennifer R. Stevenson -

S3

For Binita S. Ashar, M.D., M.B.A., F.A.C.S.

Director

Division of Surgical Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

## Indications for Use

510(k) Number (if known)

K180856

Device Name

Neutrogena Light Therapy Aging Mask+

Indications for Use (Describe)

The Neutrogena Light Therapy Aging Mask+ is an over the counter device that is indicated for the treatment of full face wrinkles.

Type of Use (Select one or both, as applicable)

☐ Prescription Use (Part 21 CFR 801 Subpart D)

☒ Over-The-Counter Use (21 CFR 801 Subpart C)

### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
[PRASStaff@fda.hhs.gov](mailto:PRASStaff@fda.hhs.gov)

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

## 5. 510(k) SUMMARY

510(k) NUMBER: K180856

510(k) SUBMITTER: Johnson & Johnson Consumer, Inc.  
7050 Camp Hill Road  
Fort Washington, PA 19034

CONTACT: Cindy R. Abraham  
Associate Director, Regulatory Affairs  
cabraha3@its.jnj.com

DATE PREPARED: March 21, 2018

PROPRIETARY NAME: Neutrogena® Light Therapy Aging Mask+

PANEL: General & Plastic Surgery

REGULATION NUMBER: CFR Title 21, 878.4810

CLASSIFICATION: Class II

PRODUCT CODES: OHS

COMMON NAME: Light-based Over The Counter Wrinkle Reduction

### 5.1 Predicate Devices [807.92(a)(3)]

- LG Beauty LED MASK, K170984
- Rejuvalite MD, K133896
- Pro X OTC 5 Light Therapy, K140471

### 5.2 Device Description [807.92(a)(4)]

The Neutrogena Light Therapy Aging Mask+ uses known LED light therapy technology for the treatment of facial wrinkles. The device emits a combination of red light and infrared light. Users place the lightweight mask over the face and press the “On” button on the controller module to start treatment. The device will automatically turn off after each 10-minute treatment.

**5.3 Intended Use and Indications for Use [807.92(a)(5)]**

The Neutrogena Light Therapy Aging Mask+ is an over the counter device that is indicated for the treatment of full face wrinkles.

**5.4 Comparison of Technological Characteristics with the Predicate Devices [807.92(a)(6)]**

The Neutrogena Light Therapy Aging Mask+ is substantially equivalent to the predicate devices in that it has the same intended use and technological characteristics, as demonstrated in **Table 5.1**. Performance testing performed on the Neutrogena Light Therapy Aging Mask+ is sufficient to demonstrate that the subject device is as safe and effective as the legally marketed predicate devices. The technological and labeling differences do not raise new or different questions about safety or effectiveness. The Neutrogena Light Therapy Aging Mask+ is substantially equivalent to the predicate devices.

**Table 5.1 Technological Characteristics**

| <b>Device Feature</b>                  | <b>Johnson &amp; Johnson<br/>Neutrogena Light<br/>Therapy Aging<br/>Mask+<br/>Subject Device</b>                                  | <b>LG Electronics<br/>LG Beauty LED<br/>Mask<br/>Predicate 1</b>                                                          | <b>Trophy<br/>Rejuvalite MD<br/>Predicate 2</b>                                                                      | <b>La Lumiere, LLC<br/>Pro X OTC 5 Light<br/>Therapy<br/>Predicate 3</b>                                                                                                                               |
|----------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 510(k) Number                          | Not yet assigned                                                                                                                  | K170984                                                                                                                   | K133896                                                                                                              | K140471                                                                                                                                                                                                |
| Classification                         | 21 CFR 878.4810                                                                                                                   |                                                                                                                           |                                                                                                                      |                                                                                                                                                                                                        |
| Common Name                            | Light-based over-the-counter wrinkle reduction                                                                                    |                                                                                                                           |                                                                                                                      |                                                                                                                                                                                                        |
| Intended Use                           | Light-based treatment to reduce wrinkles on the body in general or specific anatomical area                                       |                                                                                                                           |                                                                                                                      |                                                                                                                                                                                                        |
| Product Code                           | OHS                                                                                                                               |                                                                                                                           |                                                                                                                      |                                                                                                                                                                                                        |
| Rx and/or OTC                          | OTC                                                                                                                               |                                                                                                                           |                                                                                                                      |                                                                                                                                                                                                        |
| Use Environment                        | Home-Use                                                                                                                          |                                                                                                                           |                                                                                                                      |                                                                                                                                                                                                        |
| Anatomical Location                    | Full face                                                                                                                         |                                                                                                                           |                                                                                                                      |                                                                                                                                                                                                        |
| Indications for Use                    | The Neutrogena Light Therapy Aging Mask+ is an over the counter device that is indicated for the treatment of full face wrinkles. | The LG BEAUTY LED MASK is an over the counter device that is intended for the use in the treatment of full face wrinkles. | The Rejuvalite MD is an Over the Counter device that is intended for the use in the treatment of full face wrinkles. | The Pro X OTC 5 is an Over-The-Counter device intended for use in the treatment of facial wrinkles. It is for people with wrinkles on their face and who have Fitzpatrick skin types I, II and/or III. |
| Power Source                           | Nickel-Metal Hydride rechargeable batteries                                                                                       | Li-Ion rechargeable batteries                                                                                             | AC to DC                                                                                                             | AA alkaline batteries                                                                                                                                                                                  |
| Energy Source                          | Light-emitting diodes                                                                                                             |                                                                                                                           |                                                                                                                      |                                                                                                                                                                                                        |
| Wavelength (nm)                        | Red 620-640nm                                                                                                                     | Red 637nm                                                                                                                 | Red 600, 622, 660nm                                                                                                  | Red 620-640nm                                                                                                                                                                                          |
|                                        | IR 820-880nm                                                                                                                      | IR 854nm                                                                                                                  | IR 860nm                                                                                                             | IR 820-880nm                                                                                                                                                                                           |
| Intensity (mW/cm <sup>2</sup> )        | 1.32 mW/cm <sup>2</sup>                                                                                                           | 4.08 mW/cm <sup>2</sup>                                                                                                   | 2.35 mW/cm <sup>2</sup>                                                                                              | 0.76 mW/cm <sup>2</sup>                                                                                                                                                                                |
| Treatment Regimen                      | 10 minutes/day for 60 sessions                                                                                                    | 9 minutes/day, 5 days a week for 8 weeks                                                                                  | 3 minutes/day, 5 days a week for 8 weeks                                                                             | 15 minutes/day for 60 sessions                                                                                                                                                                         |
| Total Energy Dose (J/cm <sup>2</sup> ) | 47.58 J/cm <sup>2</sup>                                                                                                           | 88.21 J/cm <sup>2</sup>                                                                                                   | 16.94 J/cm <sup>2</sup>                                                                                              | 41.23 J/cm <sup>2</sup>                                                                                                                                                                                |

## **5.5 Performance Data [807.92(b)(1), (b)(2)]**

Performance bench and safety testing include

- Performance verification
- Human Factors Evaluation
- Biocompatibility testing in accordance with requirements put forth in ISO 10993-1
- Electrical Safety and Electromagnetic Compatibility testing

The verification, electrical safety testing, electromagnetic compatibility testing, and human factors data presented in this 510(k) submission demonstrate the Neutrogena Light Therapy Aging Mask+ device met the established specifications for its intended use and environment of use. In addition, the testing demonstrated that the subject device does not raise new or different questions of safety or effectiveness when compared to the legally-marketed predicate devices.

## **5.6 Conclusion [807.92(b)(3)]**

The basis for substantial equivalence for the Neutrogena Light Therapy Aging Mask+ device and the predicate devices is performance data and conformity with recognized standards. Performance bench verification and validation demonstrate that the subject device performs comparably to the predicate devices that are marketed for the same intended use. Based on the performance testing and the similarities of the indications for use and the technological characteristics, it can be concluded that the Neutrogena Light Therapy Aging Mask+ device is substantially equivalent to the predicate devices.