

April 25, 2023

Shenzhen Xiazhifeng Electronic Co., Ltd.
% Tracy Che
Registration engineer
Feiying Drug & Medical Consulting Technical Service Group
Rm 2401 Zhenye International Business Center, No. 3101-90,
Qianhai Road
Shenzhen, Guangdong 518052
China

Re: K230549

Trade/Device Name: IPL Hair Removal Device, Model(s): FY-B505AL, FY-B505AG, FY-B505BL,

FY-B505BG, FY-B507AL, FY-B507AG, FY-B507BL, FY-B507BG, FY-B508AL, FY-B508AG, FY-B508BL, FY-B508BG, FY-B509AL, FY-B509AG,

FY-B509BL, FY-B509BG

Regulation Number: 21 CFR 878.4810

Regulation Name: Laser Surgical Instrument For Use In General And Plastic Surgery And In

Dermatology

Regulatory Class: Class II

Product Code: OHT

Dated: February 28, 2023 Received: February 28, 2023

#### Dear Tracy Che:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. Although this letter refers to your product as a device, please be aware that some cleared products may instead be combination products. The 510(k) Premarket Notification Database located at <a href="https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm">https://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfpmn/pmn.cfm</a> identifies combination product submissions. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of

K230549 - Tracy Che Page 2

Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the <u>Federal Register</u>.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803) for devices or postmarketing safety reporting (21 CFR 4, Subpart B) for combination products (see <a href="https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products">https://www.fda.gov/combination-products/guidance-regulatory-information/postmarketing-safety-reporting-combination-products</a>); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820) for devices or current good manufacturing practices (21 CFR 4, Subpart A) for combination products; and, if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <a href="https://www.fda.gov/medical-devices/medical-device-safety/medical-device-reporting-mdr-how-report-medical-device-problems">https://www.fda.gov/medical-device-problems</a>.

For comprehensive regulatory information about medical devices and radiation-emitting products, including information about labeling regulations, please see Device Advice (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance</a>) and CDRH Learn (<a href="https://www.fda.gov/training-and-continuing-education/cdrh-learn">https://www.fda.gov/training-and-continuing-education/cdrh-learn</a>). Additionally, you may contact the Division of Industry and Consumer Education (DICE) to ask a question about a specific regulatory topic. See the DICE website (<a href="https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice">https://www.fda.gov/medical-devices/device-advice-comprehensive-regulatory-assistance/contact-us-division-industry-and-consumer-education-dice</a>) for more information or contact DICE by email (DICE@fda.hhs.gov) or phone (1-800-638-2041 or 301-796-7100).

Sincerely,

# Jianting Wang -S

Jianting Wang
Acting Assistant Director
DHT4A: Division of General Surgery Devices
OHT4: Office of Surgical
and Infection Control Devices
Office of Product Evaluation and Quality
Center for Devices and Radiological Health

Enclosure

# DEPARTMENT OF HEALTH AND HUMAN SERVICES Food and Drug Administration

## **Indications for Use**

Form Approved: OMB No. 0910-0120

Expiration Date: 06/30/2023 See PRA Statement below.

510(k) Number (if known) K230549

Device Name

IPL Hair Removal Device

Model(s): FY-B505AL, FY-B505AG, FY-B505BL, FY-B505BG, FY-B507AL, FY-B507AG, FY-B507BL, FY-B507BG, FY-B508AL, FY-B508AG, FY-B508BL, FY-B508BG, FY-B509AL, FY-B509BL, FY-B509BG

Indications for Use (Describe)

The IPL Hair Removal Device is an over-the-counter device, intended for removal of unwanted body and/or facial hair.

Type of Use (Select one or both, as applicable)	
Prescription Use (Part 21 CFR 801 Subpart D)	Over-The-Counter Use (21 CFR 801 Subpart C)

#### CONTINUE ON A SEPARATE PAGE IF NEEDED.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

#### \*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\*

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services Food and Drug Administration Office of Chief Information Officer Paperwork Reduction Act (PRA) Staff PRAStaff@fda.hhs.gov

"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."

# 510(k) Summary

# "510(k) Summary" as required by 21 CFR Part 807.92.

#### I. Submitter

Shenzhen Xiazhifeng Electronic Co., Ltd.

201, No. 39, Longfeng Road, Xinsheng Community, Longgang Street, Longgang District,

Shenzhen, Guangdong, China

Post code: 518116

Tel.: +86 0755-28992929

James Li

General Manager

Tel: +86 13760395988

Email: 2850977367@qq.com Date of preparation: 2023-4-6

#### II. Device

Name of Device: IPL Hair Removal Device

Model(s): FY-B505AL, FY-B505AG, FY-B505BL, FY-B505BG, FY-B507AL, FY-B507AG, FY-B507BL, FY-B507BG, FY-B508AL, FY-B508AG, FY-B508BL, FY-B508BG, FY-B509AL,

FY-B509AG, FY-B509BL, FY-B509BG

Common or Usual Name: Light Based Over-The-Counter Hair Removal

Regulation Name: Laser surgical instrument for use in general and plastic surgery and in

dermatology

Regulatory Class: II Product Code: OHT

Regulation Number: 21 CFR 878.4810

#### **III. Predicate Devices**

#### Primary predicate device:

Manufacturer Predicate Device		510(k) Number	Approval Date
Shenzhen Junbobeauty		K220669	May 16, 2022
Technology Co., Ltd	HANDSET (IPL-666)	1422000)	141dy 10, 2022

### Predicate device:

<u>Manufacturer</u>	turer Predicate Device 510(k) Number App		Approval Date
Shenzhen Mareal Tech Co., Ltd	Home use hair removal device (T4, T5, T8, T4-01, T5-01, T8-01, T7)	K220248	July 1, 2022

#### **IV. Device Description**

The IPL Hair Removal Device is a personal, light-based, hair reduction device. The device provides hair reduction using IPL technology. Of which, the Device includes FY-B505AL, FY-B505AG, FY-B505BL, FY-B505BG, FY-B507AL, FY-B507AG, FY-B507BL, FY-B507BG, FY-B508AL, FY-B508AG, FY-B508BL, FY-B508BG, FY-B509AL, FY-B509AG, FY-B509BL, FY-B509BG sixteen models. All have adopted the same structure design, consisting of IPL Hair Removal Device main unit and power adapter two parts, and one non-removable lamp head (light-emitting treatment window) located in the main body which is the source of optical radiation, namely a Xenon flashlamp. Meanwhile, the device is powered from power adapter via an external power. The difference of all models is mainly appearance, cooling function and enclosure color, letter A in the model name means with skin cooling function, letter B means without skin cooling function, letter L means the enclosure color is light green, letter G is deep green, which do not affect the intended use.

#### V. Indications for Use

The IPL Hair Removal Device is an over-the-counter device, intended for removal of unwanted body and/or facial hair.

#### VI. Comparison of Technological Characteristics With the Predicate Device

The IPL Hair Removal Device has the same intended use, mode of action and similar operational characteristics as the predicate device. Any minor differences between the subject device and the listed predicate devices do no raise any issues of safety or efficacy. Performance data supports that the device is safe and as effective as the predicate devices for its intended use. Therefore, the IPL Hair Removal Device may be found substantially equivalent to its predicate devices.

IPL Hair Removal Device is compared with the following Predicate Devices in terms of intended use, design, material, specifications, and performance:

Comparison Elements	Subject Device	Primary Predicate <u>Device</u>	Predicate Device	Remark
510(k) Number	Pending	K220669	K220248	/
Trade name	IPL Hair Removal Device (FY-B505AL, FY- B505AG, FY- B505BL, FY- B505BG, FY- B507AL, FY- B507AG, FY-	REMOVAL HANDSET	Home use hair removal device (T4, T5, T8, T4-01, T5-01, T8-01, T7)	/

Manufacturer	B507BL, FY- B507BG, FY- B508AL, FY- B508AG, FY- B508BL, FY- B508BG, FY- B509AL, FY- B509AG, FY- B509BL, FY- B509BG) Shenzhen Xiazhifeng	Shenzhen Junbobeauty	Shenzhen Mareal Tech	
Regulation number	Electronic Co., Ltd. 21 CFR 878.4810	Technology Co., Ltd 21 CFR 878.4810	Co., Ltd 21 CFR 878.4810	Same
Product code	OHT	OHT	OHT	Same
Device classification	Class II	Class II	Class II	Same
Indication for use/ Intended use	The IPL Hair Removal Device is an over-the-counter device, intended for removal of unwanted body and/or facial hair.	IPL HAIR REMOVAL HANDSET is an over-the-counter device intended for removal of unwanted body and/or facial hair.	The Home use hair removal device is an over-the-counter device intended for removal of unwanted hair such as but not limited to small areas such as underarm and facial hair below the chin line and large areas such as legs, in patients with Fitzpatrick Skin Phototypes I-V.	Same
Prescription or OTC	ОТС	OTC	ОТС	Same
Applicable skin	Fitzpatrick Skin Types I-V	Fitzpatrick Skin Phototypes I-V	Fitzpatrick skin types I - V	Same
Treatment area	Small areas such as underarm, bikini line. Large areas such as legs, arms.	Used on facial hair below the chin line, arms, legs, underarms, bikini line.	Small areas such as underarm and facial hair below the chin line. Large areas such as legs.	Similar
Device design				
Source energy	Supplied by external adapter	Supplied by external adapter	Supplied by external adapter	Same

			T	1
Power supply	100-240V AC Input 12V3A DC Output	100~240V AC Input 12V3A DC Output	Unknown	Same
Dimension	36*61.5*183 mm	124*83*48.5mm	Unknown	Different
Sterilization	Not required	Not required	Not required	Same
Output specific	cation			
Light source	Intense Pulsed Light	Intense Pulsed Light	Intense Pulsed Light	Same
Energy medium	Xenon Flashlamp	Xenon Quartz Tube	Xenon Arc Flashlamp	Similar
Wavelength range	510-1100nm	470-1100nm	530-1100nm	Similar
Energy density	$1.1 \sim 3.0 \text{ J/cm}^2$	1.3-2.49 J/cm <sup>2</sup>	2.0~4.3 J/cm <sup>2</sup>	Similar
Output energy	4.0 J∼10.8 J	Level 1: 3.92J Level 2: 4.72J Level 3: 5.62J Level 4: 6.49J Level 5: 7.48J	7.8-15.39 J	Similar
Spot size	3.6 cm <sup>2</sup>	3 cm <sup>2</sup>	3.5cm², 3.9cm²	Similar
Pulse duration	5.0±2.0 ms	11.5-15 ms	5.0~9.6 ms	Similar
Pulsing control	Finger switch	Finger switch	Finger switch	Same
Delivery device	Direct illumination to tissue	Direct illumination to tissue	Direct illumination to tissue	Same
Output intensity level	5 Levels	5 levels	5 levels	Same
Software/ Firmware/ Microprocess or Control?	Yes	Yes	Yes	Same
Additional featu	ures			1
Electrical safety	IEC 60601-1 IEC 60601-1-2 IEC 60601-1-11 IEC 60601-2-83	IEC 60601-1 IEC 60601-1-2 IEC 60601-1-11 IEC 60601-2-83	IEC 60601-1 IEC 60601-1-2 IEC 60601-1-6 IEC 60601-1-11 IEC 60601-2-57 IEC 60601-2-83	Same
Eye safety	IEC 62471	IEC 62471	IEC 62471	Same
Biocompatibil ity	ISO 10993-5 ISO 10993-10 ISO 10993-23	ISO 10993-5 ISO 10993-10	ISO 10993-5 ISO 10993-10	Same

#### VII.Performance Data

The following performance data were provided in support of the substantial equivalence determination.

#### 1) Biocompatibility Testing

The biocompatibility evaluation for the body-contacting components of the IPL Hair Removal Device was conducted in accordance with the "Use of International Standard ISO 10993-1, 'Biological Evaluation of Medical Devices –Part 1: Evaluation and Testing Within a Risk Management Process, Document Issued on September 4, 2020", as recognized by FDA. The following testing was performed to, and passed, including:

- ➤ ISO 10993-5:2009, Biological evaluation of medical devices –Part 5: Tests for in vitro cytotoxicity
- ➤ ISO 10993-10:2021, Biological evaluation of medical devices —Part 10: Tests for skin sensitization
- ➤ ISO 10993-23:2021, Biological evaluation of medical devices –Part 23: Tests for skin irritation

#### 2) Electrical Safety and EMC

Electrical safety and EMC testing was performed to, and passed, as per the following standards:

- ➤ IEC 60601-1-2 Medical electrical equipment –Part 1-2: General requirements for basic safety and essential performance –Collateral standard: electromagnetic compatibility
- ➤ ANSI AAMI ES60601-1 Medical electrical equipment Part 1: General requirements for basic safety and essential performance
- ➤ IEC 60601-1-11 Medical Electrical Equipment –Part 1: General Requirements for Basic Safety and Essential Performance –Collateral Standard: Requirements for Medical Electrical Equipment and Medical Electrical Systems Used in the Home Healthcare Environment
- ➤ IEC 60601-2-83 Medical electrical equipment Part 2-83: Particular requirements for the basic safety and essential performance of home light therapy equipment

#### 3) Eye Safety

➤ IEC 62471 Photobiological safety of lamps and lamp systems

#### 4) Software Verification and Validation

Software documentation consistent with *moderate level* of concern was submitted in this 510(k). System validation testing presented in this 510(k) demonstrated that all software requirement specifications are met and all software hazards have been mitigated to acceptable risk levels.

#### 5) Usability

The product usability has been evaluated and verified according to the following FDA guidance

➤ Applying Human Factors and Usability Engineering to Medical Devices, issued on FEBRUARY 2016

# 6) Performance Testing

- ➤ Product output energy density has been tested and verified according to product specification and the test results have confirmed that this parameter complied with specification requirement.
- ➤ Product pulse duration has been tested and verified according to product specification and the test results have confirmed that this parameter complied with specification requirement.

### VIII. Conclusions

Based on the above analysis and non-clinical tests performed, it can be concluded that the subject device IPL Hair Removal Device is as safe, as effective, and performs as well as the legally marketed predicate devices.