

K081669

**510(k) SUMMARY**

**DenTek Oral Care Inc.'s  
Improved Comfort-Fit NightGuard**

**JUL 10 2008**

**Submitter's Name, Address, Telephone Number, Contact Person and Date Prepared**

Howard M. Holstein, Esq.  
Hogan & Hartson LLP  
555 13th Street NW  
Washington, DC 20004

Phone: (202) 637-5813  
Facsimile: (202) 637-5910

Date Prepared: June 13, 2008

**Name of Device and Name/Address of Sponsor**

Improved Comfort-Fit NightGuard

DenTek Oral Care, Inc.  
307 Excellence Way  
Maryville, TN 37801  
Phone: (865) 983-1300  
Facsimile: (865) 983-2444

**Common or Usual Name**

NightGuard

**Classification Name**

Mouthguard, Over-the-Counter

**Classification Product Code**

OBR

**Predicate Devices**

DenTek Oral Care Inc.'s Comfort Fit NightGuard (K072147)

## **Purpose of the Special 510(k) notice.**

The Improved Comfort-Fit NightGuard is a modification to DenTek Oral Care Inc.'s Comfort Fit NightGuard (K072147).

## **Intended Use**

DenTek's Improved Comfort-Fit NightGuard is indicated for use for protection against bruxism or nighttime teeth grinding. It is intended to reduce damage to the teeth and to prevent the noise associated with bruxing or grinding.

## **Technological Characteristics**

The Improved Comfort-Fit NightGuard is a posterior-only occlusion nightguard, consisting of two bite pads connected by a buccal retention strap. The Improved Comfort-Fit consists entirely of ELVAX, a thermoplastic material. The bite pads move along the buccal strap in order to adjust to the individual user needs, with the strap always contained within the wings of the bite pads. There are 5 positions of adjustability for each molar pad.

## **Substantial Equivalence**

DenTek's Improved Comfort-Fit has the same intended use and similar indications, principles of operation, and technological characteristics as DenTek's Comfort Fit. The minor differences made in the device for patient comfort do not raise any new questions of safety or effectiveness. Thus, the Improved Comfort-Fit is substantially equivalent to its predicate devices.



JUL 10 2008

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

DenTek Oral Care, Incorporated  
C/O Mr. Howard M. Holstein, Esq  
Regulatory Counsel  
Hogan & Hartson LLP  
555 Thirteenth Street, NW  
Washington, DC 20004

Re: K081669  
Trade/Device Name: Improved Comfort-Fit NightGuard  
Regulation Number: None  
Regulation Name: None  
Regulatory Class: Unclassified  
Product Code: OBR  
Dated: June 13, 2008  
Received: June 13, 2008

Dear Mr. Holstein:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (240) 276-0115. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/industry/support/index.html>.

Sincerely yours,



Chiu Lin, Ph.D.

Director

Division of Anesthesiology, General Hospital,

Infection Control and Dental Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

K081669

### Indications for Use Statement

510(k) Number (if known): \_\_\_\_\_

Device Name: Improved Comfort-Fit NightGuard

Indications for Use:

The DenTek Improved Comfort-Fit NightGuard is indicated for use for protection against bruxism or nighttime teeth grinding. It is intended to reduce damage to the teeth and to prevent the noise associated with bruxing or grinding.

Prescription Use \_\_\_\_\_  
(Per 21 C.F.R. 801.109)

AND/OR

Over-The-Counter Use X  
(Per 21 C.F.R. 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE -- CONTINUE ON ANOTHER  
PAGE IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Robert S. Betz DDS for Dr. Susan Runner  
(Division Sign-Off)  
Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K081669