

NOV 5

K 99 6/70

10(k) Summary of Safety and Effectiveness
in Accordance with SDMA '90

December 16, 1994

Braun of America, Inc

Twelfth
Philadelphia, PA
(215) 691-5400

CONTACT: Ms. S. Alsberge, Regulatory Affairs Associate

PRODUCT NAME: Subcutaneous Injection Sets

MODEL NAME: Comfort-Q Subcutaneous Injection Sets

CLASSIFICATION: I.D. NAME:

General Hospital
Class 1: 800 FPA, Intravascular Administration Set
21 CFR 800.5440

SUBSTANTIAL EQUIVALENCE TO

510 (k) number	Name	Applicant
K8801941	MICROCATH INTRAVENOUS CATHETER PLACEMENT UNIT AND CATHETER	BURRON MEDICAL INC.
K8873885	SUBCUTANEOUS INFUSION SET	PHARMACIA DELTEC

DEVICE DESCRIPTION:

B. Braun intends to introduce into commerce the Comfort-Q Subcutaneous Injection Set. The Comfort-Q Subcutaneous Injection Set is a 23 gauge (0.9mm) catheter over 26 gauge introducer with 42 in. (106.7cm.) microbore tubing. Adhesive is enclosed. These all have similar designs and performance characteristics. The intended use for infusion of fluids or other therapeutic agents is as follows.

ATTACHMENT 2

ATT

MATERIAL :

B. Braun of America certifies that the biocompatibility tests recommended in the Tripartite Guidance for this category of contact duration will be completed for all the materials listed above for use in the manufacture of the device. All the materials listed above have also past requirements for USP Class 6 cytotoxicity, Hemolysis, and USP Physio-chemical testing.

SUBSTANTIAL EQUIVALENCE:

From a regulatory standpoint, the Comfort-Q Subcutaneous Injection Sets are equivalent in materials, form, and intended use to the Pharmacia Deltec (Subcutaneous Infusion Sets) and in materials to the Burron Medical, Inc. (Microcath Intravenous Catheter and Placement Unit and Catheter). There are no new issues of safety or effectiveness raised by the Comfort-Q Subcutaneous Injection Sets.

SAFETY AND EFFECTIVENESS:

All finished products are tested and must meet all required release specifications before distribution. The array of testing required for release include, but are not limited to; sterility, pyrogenicity (endotoxin/ LAL Method), physical testing, visual examination (in process and finished product).

The physical testing is defined by Quality Control Test Procedure documents. These tests are established testing procedures and parameters which conform to the product design specifications. The physical testing for the The Comfort-Q Subcutaneous Injection Sets is defined in detail in the "Device Master Records".

The testing instruction records for each of the individually required procedures are approved, released, distributed and revised in accordance with document control GMP's.