

K954788

Oct. 28, 1996

5- 510(k) SUMMARY as required by 807.92 (c)

Submitter's name EUREDUC
Submitter's address 35, rue des Broderies 78310 COIGNIERES FRANCE
Phone # (33) 1 34 61 01 61 Fax # (33) 1 34 61 40 61
Name of contact person: Mr. Didier CAUZOT
Device Trade name: EUREDUC TP05
Device common name: Pneumatic Drainage System
Classification name: Sleeve Limb Compressible 74JOW CFR Section 870 58 00
Predicate device name: JOBST EXTREMITY PUMP(R) SYSTEM 7500
Predicate device 510(k): K882683

EUREDUC TP05 is comprised of:

- a pneumatic unit in a metal casing with
 - a microprocessor controlled pump;
 - a control panel
- a range of garments connectable to the unit:
 - leg boots;
 - arm sleeves;

The EUREDUC TP05 is intended for:

- Lymphedema of lower extremities; post surgical or congenital;
- Postmastectomy lymphedema ;
- Pre, intra, post phases of surgery.

... / ...