

JUN 19 1996

K961870

9. **SUMMARY OF SAFETY AND EFFECTIVENESS**

TSST-1, Toxin testing was conducted; results are in the range of values of commercially sold tampons and there are no significant differences between Playtex Tampons #12496 (Slimfit Tampon) and the predicate device (Exhibit 4).

Syngyna test results are reported in Section 6 - Regular and Super absorbency ranges are within 6-9 gms and 9-12 gms respectively.

Dioxin Analysis - Not detectable at limits of detection of HRGG/HRMS Method 8290 (Triangle Labs, North Carolina). Summary provided in Exhibit 5.