## **SECTION 3** quantex C4 - 510(k) SUMMARY (Summary of Safety and Effectiveness)

Submitted by:

Betty Lane

Director, Regulatory Affairs

Instrumentation Laboratory Company

113 Hartwell Avenue

Lexington, MA 02173

phone (617) 861-4182

fax (617) 861-4464

1C963198

**Contact Person:** 

Betty Lane

phone (617) 861-4182

**Summary Prepared:** 

May 29, 1996

Name of the device:

quantex C4:

8 x 6 mL anti-human C4

P/N 3000-22137

2 x 100 mL Buffer

P/N 3000-22130

Classification name(s):

866.5240 Complement components immunological test systems

Class II

**Identification of predicate device(s):** 

IL Test™ C4

#K920010/B

## **Description of the device/intended use(s):**

This product permits the quantitative in vitro diagnostic determination of Complement 4 in serum and plasma on the ILab Clinical Chemistry System by turbidimetric immunoassay method. C4 is one of a group of serum proteins which destroy infectious agents. Measurement of this protein aids in the diagnosis of immunologic disorders, especially those associated with deficiencies of C4.

## Statement of how the Technological Characteristics of the Device compare to the Predicate device:

Biokit's quantex C4 is substantially equivalent in performance, intended use, safety and effectiveness to the predicate, IL Test™ C4 -- 510(k) #K920010/B.

## **Summary of Performance Data:**

A comparative performance study was performed with IL Test<sup>TM</sup> C4 on a Monarch Clinical Chemistry System and Biokit's quantex C4 on an ILab Clinical Chemistry System. Method comparison studies evaluated C4 levels of thirty-eight serum samples on the respective instruments. The correlation (r) in these studies was 0.9791.