

K962394

Implant Integration Systems, Inc.

6161 Clark Road, Suite 8 ♦ Paradise, CA 95969 USA ♦ (916) 872-1020 ♦ Fax: (916) 877-7555

OCT 11 1996

Summary of Safety and Effectiveness Information

The Implant Integration Systems Surgical Drills are intended to be used in the placement of the "Dual Integrator" Dental Implant (K944225).

This product is similar in design, composition (titanium nitride coated stainless steel) and function to the Implant Innovations, Inc. Surgical Drills (K891615), the Impl-Med Corp. Titanium/Stainless Steel Dental Implant Drills (K902158) and the Steri-oss Drills and Taps (K915733).



Cory D. Hanosh, Sales Manager

6-19-96

Date