

K962474

**510(k) SUMMARY
FOR
OLYMPUS FG GRASPING FORCEPS**

DEC 20 1996

Device Name: Olympus FG-4/14/20/21/22/26/32/38/40/##-1 Grasping Forceps (Urology)

Common/Usual Name: Grasping Forceps

Classification Name: Endoscope and Accessories
21 CFR 876.1500, Class II

Predicate Device: Olympus Pre-Amendment and 510(k) cleared grasping forceps
Olympus FG-Series Rat tooth Grasping Forceps (K955051)

Contact Person: Barry Sands
Olympus America, Inc.
Endoscope Division
2 Corporate Center Drive
Melville, NY 11747-3157
(516) 844-5474

Summary Preparation Date: June 17, 1996

Statement of Intended Use: The Olympus FG-6/7/8/9/42/47-1 grasping forceps are intended to be used to grasp tissue and/or retrieve foreign body, and excised tissue endoscopically.

Device Description: The Olympus FG-6/7/8/9/42/47-1 grasping forceps are designed for grasping tissue and/or retrieving foreign body, and excised tissue under endoscopic visualization. These forceps consist of a flexible shaft and a proximal control handle. Operation of the proximal control handle actuates the distal tip grasping jaws.