

NOV 27 1996

K964502

**A Summary of Safety and Effectiveness Information**

**Name of Device:** Unitron Industries, Sound F/XP <sup>+</sup>4

**Type of Device:** Behind-The-Ear Hearing Aid. Substantially equivalent to other behind-the-ear hearing aids.

**Intended Use:** To amplify and transmit sound via air conduction to the ear.

**Features:** Three fitter controls - Cross over frequency, low channel gain, and threshold knee point.

Two channel compression circuitry.

Class D circuitry for excellent battery life.

Telecoil with pre-amplifier.

**Assembly:** Assembled from standard components that are used by Unitron and other hearing aid manufacturers.

**Technical Characteristics:** Technical specifications comply with ANSI Standard S 3.22-1987.

**Fit:** Hearing aid adjustment dictated by individual audiogram.

**Controls:** Volume control, similar to other devices. Switch for selecting microphone, telecoil and off.

**Power:** Standard hearing aid battery - size 13.

A user's guide and other information is supplied with each hearing aid.