ECS Division 3726 E. Miraloma Avenue Anaheim, CA 92806 Telephone (714) 630-6400 Toll Free (800) 854-6737



510(k) Summary

Submitter

Bio-Rad Laboratories, ECS Division 3726 E. Miraloma Avenue Anaheim, CA 92806 (714)630-6400 Fax (714)666-1383

Contact Person

Elizabeth Platt

Date of Summary Preparation

November 27, 1996

Device (Trade & Common Name)

Lyphochek Coagulation Control

Classification Name

Class II, 81GGN CFR 864.5425: Plasma Coagulation Control

Devices to Which Substantial Equivalence is Claimed

ThromboScreen Coagulation Control Plasma Pacific Hemostasis K820450

Statement of Intended Use

Lyphochek Coagulation Control is intended for use as a quality control plasma to monitor the precision of citrated coagulation systems.



ECS Division 3726 E. Miraloma Avenue Anaheim, CA 92806 Telephone (714) 630-6400 Toll Free (800) 854-6737

Description of the Device

Lyphochek Coagulation Control is prepared from human plasma with added constituents of non-human origin and purified biochemicals. The control is provided in lyophilized form for increased stability.

Statement of How Technological Characteristics Compare to Substantial Equivalence Device

A table is provided below comparing the similarities between the Bio-Rad Lyphochek Coagulation Control and the devices to which substantial equivalence is claimed.

	Bio-Rad Lyphochek Coagulation Control	Pacific Hemostasis ThromboScreen Coagulation Control Plasma
Intended Use	A quality control plasma to monitor the precision of citrated coagulation systems.	A control in prothrombin time and partial thromboplastin time determinations
Form	Lyophilized	Lyophilized
Matrix	Human plasma	Human Plasma
Storage	2 - 8°C	2 - 8°C
Open Vial Claim	16 hours at 2-25°C	8 hours at 2-8°C
Levels	Level 1, 2 and 3	Level 1, 2 and 3
Assays	prothrombin time (PT),	prothrombin time (PT),
	activated partial	activated partial
	thromboplastin time (APTT),	thromboplastin time (APTT),
	fibrinogen (Level 1)	fibrinogen (Level 1)