

**PREMARKET NOTIFICATION 510(k)**  
**Cordis Corporation**  
**Avanti + Catheter Sheath Introducer**

K970392

APR 24 1997

**SUMMARY OF SAFETY AND EFFECTIVENESS**

**I. General Provisions**

Common or Usual Name: Catheter Sheath Introducer or Sheath Introducer System

Proprietary Name: Cordis Avanti + Catheter Sheath Introducer System

**II. Name of Predicate Devices**

Brite Tip Catheter Sheath Introducer, K 9954595, October 02, 1995

Avanti™ Catheter Sheath Introducer, K 945616, November 14, 1994

**III. Classification**

Class II

**IV. Performance Standards**

Performance standards have not been established by the FDA under section 514 of the Food, Drug and Cosmetic Act.

**V. Intended Use and Device Description**

The Avanti + Catheter Sheath Introducer is indicated for use in arterial and venous procedures requiring percutaneous introduction of catheters and other intravascular devices.

As with other currently marketed Cordis CSIs mentioned in this submission, these devices provide vascular access for various intravascular devices through the valve while simultaneously maintaining hemostasis. Infusion of fluids into the vasculature and withdrawal of blood from the vasculature are possible using the sheath sideport.

**VI. Biocompatibility**

All appropriate biocompatibility tests were previously performed on the materials used for the Cordis Avanti + Catheter Sheath Introducer. No new tests were performed, since all materials had been successfully tested on previously concurred devices.

**PREMARKET NOTIFICATION 510(k)**  
**Cordis Corporation**  
**Avanti + Catheter Sheath Introducer**

**VII. Summary of Substantial Equivalence**

The Cordis Avanti + Catheter Sheath Introducer is similar in design, construction, indication for use, and performance characteristics to other commercially available sheath introducers.