

K972537

MAR 28 1997

A SUMMARY OF SAFETY AND EFFECTIVENESS INFORMATION

Name of Device: Rexton's Model 48 and 49PP

Type of Device: Body Aid – Substantially equivalent to other body aids.

Intended Use: To amplify and transmit sound to the ear.

MAR 28 1997

Features: Both the 48 and the 49PP are small slim-lined body aids. Their unique shape fits snugly and inconspicuously into any pocket or other places, e.g. brassiere etc. The 48 and the 49PP both have audio input that allows the connection to external sources. Also both have 3 position switch for off-telecoil-microphone.

Assembly: Assembled from standard components which are widely utilized by other hearing aid manufacturers.

Technical Characteristics: Technical specifications comply with S3.22-1987 ANSI Specifications. Preliminary data sheets for the 48 and 49PP are enclosed.

Fit: The frequency response of this product is dictated by the individual audiogram from each client.

Power: Standard battery size AA.

A user's manual and other information is supplied with each hearing aid (enclosed).