



510(k) Summary

SEP | 9 1997

Contact/Submitter:

Janell Colley, Regulatory Affairs Associate

Date Prepared:

September 8, 1997

Trade Name:

The Amplatz Thrombectomy Device

Common Name:

Thrombectomy Catheter

Classification Name:

Percutaneous Catheter (per 21 CFR 870.1250)

Predicate Devices:

MICROVENA Amplatz Thrombectomy Device

Device Description:

The Amplatz Thrombectomy Device (ATD) is a percutaneous, rotational thrombectomy device consisting of a polyurethane catheter. Inside the catheter, a distal rotor housing contains a small diameter, recessed impeller, attached to a drive shaft. The shaft is connected to a disposable, high speed, air driven motor. An infusion line with luer connector, attached to the proximal motor assembly, allows for infusion of saline through the catheter. A foot pedal/regulator assembly is required to operate the air motor.

Intended Use:

The modified Amplatz Thrombectomy Device (ATD) is intended for use in the mechanical dissolution of acute and sub-acute thrombus within dialysis fistulae.

Functional/Safety

Testing:

The modified Amplatz Thrombectomy Device (ATD) device has successfully undergone functional and safety testing of new design features, including joint strengths, flow rate tests, burst tests, flexibility tests, device life tests,

maceration capability tests, and material biocompatibility tests.

Conclusion:

The modified Amplatz Thrombectomy Device (ATD) is substantially

equivalent to the predicate device.





Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

SEP 1 9 ""

Ms. Janell Colley Microvena Corporation 1861 Buerkle Road White Bear Lake, Minnesota 55110-5246

Re: K971984

Amplatz Thrombectomy Device Regulatory Class: II (two) Product Code: 74 MCW Dated: September 5, 1997

Received: September 9, 1997

Dear Ms. Colley:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4639. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "http://www.fda.gov/cdrh/dsmamain.html".

Thomas J. Cellelon

Thomas J. Callahan, Ph.D.

Director

Division of Cardiovascular,

Respiratory, and Neurological Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure

Indications for Use

510(k) Number (if known):	K971984
Device Name:	Amplatz Thrombectomy Device
Indications for Use:	The intended use for the device is the mechanical dissolution of acute and sub-acute thrombus within dialysis fistulae.
Division Sign-Off) Control of Cardiovascular, Respiratory, Courological Devices 510(k) Number Control of Cardiovascular, Respiratory,	
(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)	
Concurrence of CDRH, Office of Device Evaluation (ODE)	
Prescription Use	OR Over-the-Counter Use
(Per 21 CFR 801.109)	(Optional Format 1-2-96)