



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

JAN 22 1998

Cynthia C. Knapp
Director of Lab Services
AccuMed International, Inc.
29299 Clemens Road, Suite 1-K
Westlake, OH 44145

Re: K972576
Trade Name: Gram Positive Autoidentification Plates, Sensititre AP90
Regulatory Class: I
Product Code: LQL
Dated: January 7, 1998
Received: January 12, 1998

Dear Ms. Knapp:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.


Page 2

Under the Clinical Laboratory Improvement Amendments of 1988 (CLIA-88), this device may require a CLIA complexity categorization. To determine if it does, you should contact the Centers for Disease Control and Prevention (CDC) at (770)488-7655.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll free number (800) 638-2041 or at (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>"

Sincerely yours,

A handwritten signature in black ink that reads "Steven Gutman". The signature is written in a cursive, slightly slanted style.

Steven I. Gutman, M.D., M.B.A.
Director
Division of Clinical Laboratory Devices
Office of Device Evaluation
Center for Devices and Radiological Health

Enclosure

510(k) Number (if known): K972576

Device Name: Gram Positive AutoIdentification Panel, Sensititre® AP90

Indications For Use:

The Sensititre AP90 panel is an *in vitro* diagnostic product intended to be used in conjunction with the Sensititre system for the automated identification of clinically significant gram positive organisms.

This 510(k) is for gram positive identification. Reagents are listed as follows:

REAGENTS & PLATE TEST LAYOUT

| | 1 | 2 | 3 | 4 |
|---|----------|-----------------------|----------|----------|
| A | UREA | ESCULIN | ARGININE | FR13 |
| B | FR16 | RHAMNOSE | FR15 | MANNITOL |
| C | FR14 | TREHALOSE | FR17 | MALTOSE |
| D | FR18 | FR19 | FR20 | FR21 |
| E | FR22 | FR23 | GLYCEROL | FR24 |
| F | FR25 | GLUCOSE | SUCROSE | FR26 |
| G | FR27 | B-METHYL GLUCOSIDE | FR28 | FR29 |
| H | SORBITOL | FR30 | FR31 | FR32 |

The above 4 by 8 layout is repeated 3 times across the plate.

FR-Fluorescent reagent.

The Urease test requires an oil overlay.

Organisms with indications for testing include:

Enterococcus avium
Enterococcus durans
Enterococcus faecalis
Enterococcus faecium
Erysipelothrix rhusiopathiae
Listeria species
Micrococcus species
Rhodococcus equi
Staphylococcus aureus
Staphylococcus capitis
Staphylococcus epidermidis
Staphylococcus haemolyticus

Staphylococcus hominis
Staphylococcus saprophyticus
Staphylococcus warneri
Streptococcus agalactiae
Streptococcus milleri group
Streptococcus Group G
Streptococcus mitis
Streptococcus pneumoniae
Streptococcus pyogenes
Streptococcus uberis

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED).

Concurrence of CDRH, Office of Device Evaluation (ODE)

Sally J. Seligson for J. Tschand
(Division Sign-Off)

Division of Clinical Laboratory Devices

510(k) Number K972576

Prescription Use ☒
(Per 21 CFR 801.109)

OR

Over-The-Counter Use ☐

(Optional Format 1-2-96)