



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Yasuharu Ishii
Manager
Shimadzu Medical Systems
20101 South Vermont Ave.
Torrance, California 90502-1328

MAR - 4 1998

Re: K973802
Shimasonic SDU-450 Linear/Convex Diagnostic
Ultrasound System (Echo View)
Dated: January 20, 1998
Received: January 21, 1998
Regulatory class: II
21 CFR 892.1560/21 CFR 892.1570
Procode: 90 IYO/Procode: 90 ITX

Dear Mr. Ishii:

We have reviewed your section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug and Cosmetic Act. You may, therefore, market the device, subject to the general controls provisions Act (Act). The general controls provisions of the Act include requirements for registration, listing of devices, good manufacturing practices, labeling, and prohibitions against misbranding and adulteration.

This determination of substantial equivalence applies to the following transducers intended for use with the Shimasonic SDU-450 System, as described in your premarket notification:

Transducer Model Number

L070-075U
L040-100U
EL033-075UA
L110-035U
VA13R-035U
VA13R-050U
VA15R-035U
VA15R-050U
VA40R-035U
VA57R-0375U
TV11R-055U
TV13R-050U
L072-050U
UB10R-065U

If your device is classified (see above) into either class II (Special Controls) or class III, (Premarket Approval) it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic QS inspections, the FDA will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, the Food and Drug Administration (FDA) may publish further announcements concerning your device in the Federal Register.

Please note: this response to your premarket notification does not affect any obligation you may have under sections 531 and 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Please be advised that the determination above is based on the fact that no medical devices have been demonstrated to be safe and effective for in vitro fertilization or percutaneous umbilical blood sampling, nor have any devices been marketed for these uses in interstate commerce prior to May 28, 1976, or reclassified into class I (General Controls) or class II (Special Controls). FDA considers devices specifically intended for in vitro fertilization and percutaneous umbilical blood sampling to be investigational, and subject to the provision of the investigational device exemptions (IDE) regulations, 21 CFR, Part 812. Therefore, your product labeling must be consistent with FDA's position on this use.

This determination of substantial equivalence is granted on the condition that prior to shipping the first device, you submit a postclearance special report. This report should contain complete information, including acoustic output measurements based on production line devices, requested in Appendix G, (enclosed) of the Center's September 30, 1997 "Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers." If the special report is incomplete or contains unacceptable values (e.g., acoustic output greater than approved levels), then the 510(k) clearance may not apply to the production units which as a result may be considered adulterated or misbranded.

The special report should reference the manufacturer's 510(k) number. It should be clearly and prominently marked "ADD-TO-FILE" and should be submitted in duplicate to:

Food and Drug Administration
Center for Devices and Radiological Health
Document Mail Center (HFZ-401)
9200 Corporate Boulevard
Rockville, Maryland 20850

This letter will allow you to begin marketing your device as described in your premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus permits your device to proceed to market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4591. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Page 3 - Yasuharu Ishii

If you have any questions regarding the content of this letter, please contact **Rodrigo C. Perez** at (301) 594-1212.

Sincerely yours,

for David L. Lippman

Lillian Yin, Ph.D.
Director, Division of Reproductive,
Abdominal, Ear, Nose and Throat,
and Radiological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosures

510(k) Number (if known) : K973802

Device Name : Diagnostic Ultrasound System SDU-450, L070-075U

Fill out one form for each ultrasound system or transducer.

Indications for use: Diagnostic ultrasound imaging or Doppler analysis
of the human body as follows:

Modes of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify)**	Other (Specify)
Ophthalmic										
Fetal										
Abdominal										
Intra-operative (Specify)										
Intra-operative Neurological										
Pediatric										
Small Organ (Specify) *		X	X						X	
Neonatal Cephalic										
Adult Cephalic										
Cardiac Adult										
Cardiac Pediatric										
Trans-esophageal										
Trans-rectal										
Trans-vaginal										
Intra-luminal										
Trans-utethral										
Peripheral vessel		X	X						X	
Laparoscopic										
Musculo-skeletal Conventional		X	X						X	
Musculo-skeletal Superficial										
Other (Specify)										

Additional Comments:

* Breast, Thyroid, Testes

** B/M mode

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, ENT,
and Radiological Devices

510(k) Number

K973802

510(k) Number (if known) : K973802

Device Name : Diagnostic Ultrasound System SDU-450, L040-100U

Fill out one form for each ultrasound system or transducer.

Indications for use: Diagnostic ultrasound imaging or Doppler analysis
of the human body as follows:

Modes of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify)**	Other (Specify)
Ophthalmic										
Fetal										
Abdominal										
Intra-operative (Specify)										
Intra-operative Neurological										
Pediatric										
Small Organ (Specify) *		X	X						X	
Neonatal Cephalic										
Adult Cephalic										
Cardiac Adult										
Cardiac Pediatric										
Trans-esophageal										
Trans-rectal										
Trans-vaginal										
Intra-luminal										
Trans-utethral										
Peripheral vessel		X	X						X	
Laparoscopic										
Musculo-skeletal Conventional		X	X						X	
Musculo-skeletal Superficial										
Other (Specify)										


Additional Comments:

* Breast, Thyroid, Testes

** B/M mode

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)

Division of Reproductive, Abdominal, ENT,
and Radiological Devices

510(k) Number K973802

Prescription Use (Per 21 CFR 801.109)

Ultrasound Device Indications Statement Page 3 of 14

510(k) Number (if known) : K973802

Device Name : Diagnostic Ultrasound System SDU-450, EL033-075UA

Fill out one form for each ultrasound system or transducer.

Indications for use: Diagnostic ultrasound imaging or Doppler analysis
of the human body as follows:

Modes of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify)**	Other (Specify)
Ophthalmic										
Fetal										
Abdominal										
Intra-operative (Specify)										
Intra-operative Neurological										
Pediatric										
Small Organ (Specify) *		X	X						X	
Neonatal Cephalic										
Adult Cephalic										
Cardiac Adult										
Cardiac Pediatric										
Trans-esophageal										
Trans-rectal										
Trans-vaginal										
Intra-luminal										
Trans-utethral										
Peripheral vessel		X	X						X	
Laparoscopic										
Musculo-skeletal Conventional		X	X						X	
Musculo-skeletal Superficial										
Other (Specify)										


Additional Comments:

* Breast, Thyroid, Testes

** B/M mode

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)

Division of Reproductive, Abdominal, ENT,
and Radiological Devices

510(k) Number K973802

Prescription Use (Per 21 CFR 801.109)

510(k) Number (if known) : K973802

Device Name : Diagnostic Ultrasound System SDU-450, L110-035U

Fill out one form for each ultrasound system or transducer.

Indications for use: Diagnostic ultrasound imaging or Doppler analysis
of the human body as follows:

Modes of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify)**	Other (Specify)
Ophthalmic										
Fetal		X	X						X	
Abdominal		X	X						X	
Intra-operative (Specify)										
Intra-operative Neurological										
Pediatric										
Small Organ (Specify) *										
Neonatal										
Cephalic										
Adult Cephalic										
Cardiac Adult										
Cardiac Pediatric										
Trans-esophageal										
Trans-rectal										
Trans-vaginal										
Intra-luminal										
Trans-utethral										
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Specify)										


Additional Comments:

** B/M mode

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)


(Division Sign-Off)
Division of Reproductive, Abdominal, ENT,
and Radiological Devices
510(k) Number K973802

510(k) Number (if known) : K973802

Device Name : Diagnostic Ultrasound System SDU-450, VA13R-035U

Fill out one form for each ultrasound system or transducer.

Indications for use: Diagnostic ultrasound imaging or Doppler analysis
of the human body as follows:

Modes of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify)**	Other (Specify)
Ophthalmic										
Fetal										
Abdominal		X	X						X	
Intra-operative (Specify)										
Intra-operative Neurological										
Pediatric										
Small Organ (Specify) *										
Neonatal Cephalic										
Adult Cephalic										
Cardiac Adult		X	X						X	
Cardiac Pediatric										
Trans-esophageal										
Trans-rectal										
Trans-vaginal										
Intra-luminal										
Trans-utethral										
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Specify)										

Additional Comments:

** B/M mode

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, ENT, and Radiological Devices

510(k) Number K973802

510(k) Number (if known) : K973802

Device Name : Diagnostic Ultrasound System SDU-450, VA13R-050U

Fill out one form for each ultrasound system or transducer.

Indications for use: Diagnostic ultrasound imaging or Doppler analysis
of the human body as follows:

Modes of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify)**	Other (Specify)
Ophthalmic										
Fetal										
Abdominal		X	X						X	
Intra-operative (Specify)										
Intra-operative Neurological										
Pediatric										
Small Organ (Specify) *										
Neonatal Cephalic										
Adult Cephalic										
Cardiac Adult		X	X						X	
Cardiac Pediatric		X	X						X	
Trans-esophageal										
Trans-rectal										
Trans-vaginal										
Intra-luminal										
Trans-utethral										
Peripheral vessel		X	X						X	
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Specify)										

Additional Comments:

** B/M mode

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division of Reproductive, Abdominal, ENT,
and Radiological Devices

510(k) Number K973802

Ultrasound Device Indications Statement Page 7 of 14

510(k) Number (if known) : K973802

Device Name : Diagnostic Ultrasound System SDU-450, VA15R-035U

Fill out one form for each ultrasound system or transducer.

Indications for use: Diagnostic ultrasound imaging or Doppler analysis
of the human body as follows:

Modes of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify)**	Other (Specify)
<i>Ophthalmic</i>										
<i>Fetal</i>		X	X						X	
<i>Abdominal</i>		X	X						X	
<i>Intra-operative (Specify)</i>										
<i>Intra-operative Neurological</i>										
<i>Pediatric</i>										
<i>Small Organ (Specify) *</i>										
<i>Neonatal Cephalic</i>										
<i>Adult Cephalic</i>										
<i>Cardiac Adult</i>		X	X						X	
<i>Cardiac Pediatric</i>										
<i>Trans-esophageal</i>										
<i>Trans-rectal</i>										
<i>Trans-vaginal</i>										
<i>Intra-luminal</i>										
<i>Trans-utethral</i>										
<i>Peripheral vessel</i>										
<i>Laparoscopic</i>										
<i>Musculo-skeletal Conventional</i>										
<i>Musculo-skeletal Superficial</i>										
<i>Other (Specify)</i>										


Additional Comments:

** B/M mode

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)


 (Division Sign-Off)
 Division of Reproductive, Abdominal, ENT,
 and Radiological Devices
 510(k) Number K973802

510(k) Number (if known) : K973802

Device Name : Diagnostic Ultrasound System SDU-450, VA15R-050U

Fill out one form for each ultrasound system or transducer.

Indications for use: Diagnostic ultrasound imaging or Doppler analysis
of the human body as follows:

Modes of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify)**	Other (Specify)
Ophthalmic										
Fetal		X	X						X	
Abdominal		X	X						X	
Intra-operative (Specify)										
Intra-operative Neurological										
Pediatric										
Small Organ (Specify) *										
Neonatal Cephalic										
Adult Cephalic										
Cardiac Adult		X	X						X	
Cardiac Pediatric		X	X						X	
Trans-esophageal										
Trans-rectal										
Trans-vaginal										
Intra-luminal										
Trans-utethral										
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Specify)										


Additional Comments:

** B/M mode

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)


(Division Sign-Off)
Division of Reproductive, Abdominal, ENT,
and Radiological Devices
510(k) Number K973802

510(k) Number (if known) : K973802

Device Name : Diagnostic Ultrasound System SDU-450, VA40R-035U

Fill out one form for each ultrasound system or transducer.

Indications for use: Diagnostic ultrasound imaging or Doppler analysis
of the human body as follows:

Modes of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify) **	Other (Specify)
Ophthalmic										
Fetal		X	X						X	
Abdominal		X	X						X	
Intra-operative (Specify)										
Intra-operative Neurological										
Pediatric										
Small Organ (Specify) *										
Neonatal Cephalic										
Adult Cephalic										
Cardiac Adult		X	X						X	
Cardiac Pediatric		X	X						X	
Trans-esophageal										
Trans-rectal										
Trans-vaginal										
Intra-luminal										
Trans-utethral										
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Specify)										

Additional Comments:

** B/M mode

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)

(Division Sign-Off)

Division of Reproductive, Abdominal, ENT,
and Radiological Devices

510(k) Number

K973802

510(k) Number (if known) : K973802

Device Name : Diagnostic Ultrasound System SDU-450, VA57R-0375U

Fill out one form for each ultrasound system or transducer.

Indications for use: Diagnostic ultrasound imaging or Doppler analysis
of the human body as follows:

Modes of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify) **	Other (Specify)
Ophthalmic										
Fetal		X	X						X	
Abdominal		X	X						X	
Intra-operative (Specify)										
Intra-operative Neurological										
Pediatric										
Small Organ (Specify) *										
Neonatal Cephalic										
Adult Cephalic										
Cardiac Adult		X	X						X	
Cardiac Pediatric		X	X						X	
Trans-esophageal										
Trans-rectal										
Trans-vaginal										
Intra-luminal										
Trans-utethral										
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Specify)										


Additional Comments:

** B/M mode

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)


(Division Sign-Off)
Division of Reproductive, Abdominal, ENT,
and Radiological Devices
510(k) Number K973802

510(k) Number (if known) : K973802

Device Name : Diagnostic Ultrasound System SDU-450, TV11R-055U

Fill out one form for each ultrasound system or transducer.

Indications for use: Diagnostic ultrasound imaging or Doppler analysis
of the human body as follows:

Modes of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify) **	Other (Specify)
Ophthalmic										
Fetal										
Abdominal										
Intra-operative (Specify)										
Intra-operative Neurological										
Pediatric										
Small Organ (Specify) *										
Neonatal Cephalic										
Adult Cephalic										
Cardiac Adult										
Cardiac Pediatric										
Trans-esophageal										
Trans-rectal		X	X						X	
Trans-vaginal		X	X						X	
Intra-luminal										
Trans-utethral										
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Specify)										

Additional Comments:

** B/M mode

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

(Division Sign-Off)

Division of Reproductive, Abdominal, ENT,
and Radiological Devices

510(k) Number K973802

510(k) Number (if known) : K973802

Device Name : Diagnostic Ultrasound System SDU-450, TV13R-050U

Fill out one form for each ultrasound system or transducer.

Indications for use: Diagnostic ultrasound imaging or Doppler analysis
of the human body as follows:

Modes of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify)**	Other (Specify)
Ophthalmic										
Fetal										
Abdominal										
Intra-operative (Specify)										
Intra-operative Neurological										
Pediatric										
Small Organ (Specify) *										
Neonatal Cephalic										
Adult Cephalic										
Cardiac Adult										
Cardiac Pediatric										
Trans-esophageal										
Trans-rectal										
Trans-vaginal		X	X						X	
Intra-luminal										
Trans-utethral										
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Specify)										


Additional Comments:

** B/M mode

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)


(Division Sign-Off)
Division of Reproductive, Abdominal, ENT,
and Radiological Devices
510(k) Number K973802

510(k) Number (if known) : K973802

Device Name : Diagnostic Ultrasound System SDU-450, L072-050U

Fill out one form for each ultrasound system or transducer.

Indications for use: Diagnostic ultrasound imaging or Doppler analysis
of the human body as follows:

Modes of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify) **	Other (Specify)
Ophthalmic										
Fetal		X	X						X	
Abdominal										
Intra-operative (Specify)										
Intra-operative Neurological										
Pediatric										
Small Organ (Specify) *										
Neonatal Cephalic										
Adult Cephalic										
Cardiac Adult										
Cardiac Pediatric										
Trans-esophageal										
Trans-rectal										
Trans-vaginal										
Intra-luminal										
Trans-utethral										
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional		X	X						X	
Musculo-skeletal Superficial										
Other (Specify)										

Additional Comments:

** B/M mode

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

David A. Keyman
(Division Sign-Off)

Division of Reproductive, Abdominal, ENT,
and Radiological Devices

510(k) Number K973802

Prescription Use (Per 21 CFR 801.109)

Ultrasound Device Indications Statement Page 14 of 14

510(k) Number (if known) : K973802

Device Name : Diagnostic Ultrasound System SDU-450, UB10R-065U

Fill out one form for each ultrasound system or transducer.

Indications for use: Diagnostic ultrasound imaging or Doppler analysis
of the human body as follows:

Modes of Operation

Clinical Application	A	B	M	PWD	CWD	Color Doppler	Power (Amplitude) Doppler	Color Velocity Imaging	Combined (Specify)**	Other (Specify)
Ophthalmic										
Fetal										
Abdominal										
Intra-operative (Specify)										
Intra-operative Neurological										
Pediatric										
Small Organ (Specify) *										
Neonatal Cephalic										
Adult Cephalic										
Cardiac Adult										
Cardiac Pediatric										
Trans-esophageal										
Trans-rectal		X	X						X	
Trans-vaginal										
Intra-luminal										
Trans-utethral										
Peripheral vessel										
Laparoscopic										
Musculo-skeletal Conventional										
Musculo-skeletal Superficial										
Other (Specify)										


Additional Comments:

** B/M mode

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use (Per 21 CFR 801.109)


 (Division Sign-Off)
 Division of Reproductive, Abdominal, ENT,
 and Radiological Devices
 510(k) Number K973802