

OCT 28 1998

**Class II 510(K) Summary**  
**Shelhigh *No-React*® VASCUPATCH™**

This summary of the 510(k) information is being submitted as required by section 807.92(a).

**I. Proprietary and Common Name:**

Proprietary name: Shelhigh *No-React*® *VascuPatch*™

Common name: Cardiovascular patch

**II. Regulatory Class:**

Class II device

**III. Intended Use**

For repair of carotid endarterectomy and vascular patch grafting.

**IV. Product Description**

The Shelhigh *No-React*® *VascuPatch*™ is made of glutaraldehyde fixed bovine pericardium, rinsed with the detoxification process *No-React*®. It is stored in 2% Benzyl alcohol like the Shelhigh *No-React*® pericardial patch.

The material exhibits good tensile strength, shrink temperature, excellent biocompatibility and suture retention. It is soft and pliable making it convenient to implant.

Bovine pericardial material has been used successfully as a tissue patch for pericardial closure. Glutaraldehyde processed bovine pericardium has a long history of success as a permanently implanted material.

**V. Substantial Equivalence**

The Shelhigh *No-React*® *VascuPatch*™ is equivalent to the Knitted MiniCrimp Vascupatch manufactured by Meadox Medical, Inc. #K905496 and identical to the Shelhigh *No-React*® pericardial patch K974914 currently manufactured by Shelhigh Inc.

**VI. Comparison with Predicate Device**

The Shelhigh *No-React*® *VascuPatch*™, like the Shelhigh *No-React*® pericardial patch, is a glutaraldehyde cross-linked bovine pericardium membrane which

exhibits substantially equivalent physical/mechanical properties as measured by suture retention, tensile strength, and shrink temperature. They both have identical flexibility, wall thickness and both are stored in bezyle alcohol.

## **VII. Nonclinical / Animal Tests**

An animal study was conducted to evaluate Shelhigh *No-React*® vs. the conventional glutaraldehyde treated patch, the Shelhigh *No-React*® patch shows higher level of biocompatibility.

## **VIII. Conclusions**

The non clinical /Animal testing data showed that the Shelhigh *No-React*® *VascuPatch*<sup>TM</sup>. has higher level of cytocompatibility when compared with the conventional glutaraldehyde treated and storage. The *VascuPatch*<sup>TM</sup> is identical to the Shelhigh *No-React*® *pericardial patch*, thus is substantially equivalent to the predicate device.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

OCT 28 1998

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Shlomo Gabbay, M.D.  
Chief Scientific Advisor  
Shelhigh, Inc.  
P.O. Box 884  
Millburn, NJ 07041

Re: K982810  
Shelhigh No-React® VASCUPATCH™  
Regulatory Class: II  
Product Code: DXZ  
Dated: July 15, 1998  
Received: August 11, 1998

Dear Dr. Gabbay:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. All references to the No-React® treatment or to animal studies regarding said treatment in the labeling or promotional materials/advertisements must be immediately followed by a statement that no clinical data are available which evaluate the long term impact of the detoxification treatment in humans.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in

Page 2 - Shlomo Gabbay, M.D.

the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4648. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597, or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, reading "Thomas J. Callahan". The signature is fluid and cursive, with the first name "Thomas" being the most prominent part.

Thomas J. Callahan, Ph.D.

Director

Division of Cardiovascular, Respiratory,  
and Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Page 1 of 1

510(k) Number (if known): K982810

Device Name: Shelhigh No-React® VASCUPATCH™

Indications For Use:

For repair of carotid endarterectomy and vascular patch grafting.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

[Signature]  
(Division Sign-Off)

Division of Cardiovascular, Respiratory,  
and Neurological Devices

510(k) Number K 982810

Prescription Use ✓  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_

(Optional Format 1-2-96)