

OCT 21 1998

**I. 510(k) SUMMARY OF SAFETY AND EFFECTIVENESS**

Greiner Mediatech, Inc. („Greiner“) is submitting a 510(k) pre-market notification for its Greiner MiniCollect® Serum Gel tube with clot activator and separation gel for capillary blood collection. The Greiner MiniCollect® is a non-sterile, non-evacuated blood collection device containing a clot activator and an inert acrylic barrier material, intended to collect, transport, separate and process capillary blood for testing serum.

Greiner is claiming substantial equivalence to Becton Dickinson's Microtainer® Brand tubes with gel separator. Both blood collection tubes have the same intended use and are made out of the same material, polypropylene plastic. The Greiner MiniCollect® cap is made from rubber and has "cross-cuts" to allow for direct collection into and sampling from the tube without having to remove the cap. Becton Dickinson's Microtainer® Brand tube caps are polypropylene stoppers. The equivalency of assay results for both tubes was evaluated by testing paired samples collected in Greiner MiniCollect® Serum Gel tubes and Becton Dickinson Microtainer® Brand tubes with gel separator. Test results from paired samples for 23 analytes and 13 hormones were evaluated demonstrating good correlation.

Greiner's 510(k) has been submitted on August 27, 1998 by Doug Harris, Managing Director, Greiner Mediatech, Inc., 260 Gateway Drive, Suite 17A, Bel Air, MD 21014 (410/836-8228).



Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

OCT 21 1998

Douglas L. Harris  
Managing Director  
Greiner Meditech, Inc.  
260 Gateway Drive, Suite 17A  
P.O. Box 943  
Bel Air, Maryland 21014

Re: K982998  
MiniCollect® Serum Gel Blood Collection Tube  
Regulatory Class: II  
Product Code: JKA  
Dated: August 27, 1998  
Received: August 27, 1998

Dear Mr. Harris:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.


Page 2

Under the Clinical Laboratory Improvement Amendments of 1988 (CLIA-88), this device may require a CLIA complexity categorization. To determine if it does, you should contact the Centers for Disease Control and Prevention (CDC) at (770) 488-7655.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Steven I. Gutman, M.D., M.B.A.  
Director  
Division of Clinical  
Laboratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health


Enclosure

**STATEMENT OF INDICATIONS FOR USE**

510(k) Number (if known): K982998

Device Name: MiniCollect® Serum Gel Blood Collection Tube

Indications for Use: To collect, transport, separate and process capillary blood for testing serum in the clinical laboratory.

  
(Division Sign-Off)  
Division of Clinical Laboratory Devices  
510(k) Number K982998

Prescription Use X


Over-The-Counter Use \_\_\_\_\_

**STATEMENT OF INDICATIONS FOR USE**

510(k) Number (if known): K982998

Device Name: MiniCollect® Serum Gel Blood Collection Tube

Indications for Use: To collect, transport, separate and process capillary blood for testing serum in the clinical laboratory.

  
(Division Sign-Off)  
Division of Clinical Laboratory Devices  
510(k) Number K982998

Prescription Use X

Over-The-Counter Use \_\_\_\_\_