

DEC 21 1998

K983640



Summary of Safety & Effectiveness
SYNCHRON® Systems Lipid Calibrator

1.0 **Submitted By:**

Lucinda Stockert
Staff Regulatory Specialist, Product Submissions
Beckman Coulter, Inc.
200 S. Kraemer Blvd., W-104
Brea, California 92822-8000
Telephone: (714) 961-3777
FAX: (714) 961-4123

2.0 **Date Submitted:**

October 15, 1998

3.0 **Device Name(s):**

3.1 **Proprietary Names**

SYNCHRON® Systems Lipid Calibrator

3.2 **Classification Name**

Calibrator, (21 CFR §862.1150)

4.0 **Predicate Device(s):**

SYNCHRON Systems Reagent	Predicate	Manufacturer	Docket Number
SYNCHRON® Systems Lipid Calibrator	Beckman™ VIGIL Lipid Control	Beckman Coulter, Inc.	K974452

5.0 **Description:**

The SYNCHRON Systems Lipid Calibrator Set is a two level ready-to-use human serum-based liquid calibrator set manufactured by Beckman Coulter, Inc. Each kit contains 3 X 2 mL bottles of a specific level of calibrator (identified as Level 1 and Level 2). Once opened, the calibrators are stable for 65 days when stored at +2°C to +8°C.

5.0 Intended Use:

The SYNCHRON® Systems Lipid Calibrator is intended for use with the SYNCHRON Systems for the calibration of direct HDL Cholesterol reagent.

7.0 Comparison to Predicate(s):

Similarities to the Predicates

Reagent	Aspect/Characteristic	Comments
SYNCHRON Systems Lipid Calibrator	Source Material: Defibrinated human plasma spiked with human lipids and stabilized by freezing.	Same as Beckman™ VIGIL Lipid Control
	Storage Temperature (-15°C to -20°C)	
	Liquid, ready-to-use form	
	Value Assignment Methodology	

Differences from the Predicate

Reagent	Aspect/Characteristic	Comments
SYNCHRON Systems Lipid Calibrator	Intended Use:	SYNCHRON Systems Lipid Calibrator is intended for use in calibration of SYNCHRON Systems HDL Cholesterol Reagent. Beckman™ VIGIL Lipid Control is intended for use in monitoring the reliability of automated <i>in vitro</i> diagnostic assays of HDL Cholesterol.
	Analytes	SYNCHRON Systems Lipid Calibrator contains HDL Cholesterol. Beckman™ VIGIL Lipid Control contains cholesterol, HDL cholesterol, triglycerides, apolipoprotein A-1, apolipoprotein B.
	Levels of Analyte	SYNCHRON Systems Lipid Calibrator: 2 levels Beckman™ VIGIL Lipid Control: 4 levels

8.0 Summary of Performance Data:

The data in the Premarket Notification on safety and effectiveness supports a finding of substantial equivalence to products already in commercial distribution. Stress stability and open bottle stability studies of the lipid calibrator support the Beckman stability claim of 24 months and 65 days for open bottles.

This summary of safety and effectiveness is being submitted in accordance with the requirements of the Safe Medical Device Act of 1990 and the implementing regulation 21 CFR 807.92.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

DEC 21 1998

Food and Drug Administration
2098 Gaither Road
Rockville MD 20850

Ms. Lucinda Stockert
Staff Regulatory Specialist, Product Submissions
Beckman Coulter, Inc.
200 S. Kraemer Blvd., W-104
Brea, CA 92822-8000

Re: K983640

Trade Name: SYNCHRON® Systems Lipid Calibrator
Regulatory Class: II
Product Code: 75 JIS
Dated: October 15, 1998
Received: October 16, 1998

Dear Ms. Stockert:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597, or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in dark ink, reading "Steven Gutman". The signature is written in a cursive, slightly slanted style.

Steven I. Gutman, M.D., M.B.A.
Director
Division of Clinical
Laboratory Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K983640

Device Name: **SYNCHRON® Systems**
Lipid Calibrator

Indications for Use:

The SYNCHRON® Systems Lipid Calibrator, used in conjunction with SYNCHRON® HDL Cholesterol reagent, is intended for use on Beckman's SYNCHRON Systems for the calibration of HDL Cholesterol test systems.

Clinical Significance:

A calibrator is a device intended for medical purposes for use in a test system to establish points of reference that are used in the determination of values in the measurement of substances in human specimens.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓
(per 21 CFR 801.109)

OR

Over-the-Counter Use _____
Optional Format 1-2-96

Veronica J. Caluin for GWM
(Division Sign-Off)

Division of Clinical Laboratory Devices

510(k) Number K983640