

4/16/99

K990928



**Bio-Rad  
Laboratories**

*Diagnostics Group  
9500 Jeronimo Road  
Irvine, California 92618-2017  
Telephone: (949) 598-1200*

---

## **510(k) Summary**

### Submitter

Bio-Rad Laboratories  
9500 Jeronimo Road  
Irvine, CA  
(949)598-1285  
Fax (949)598-1555

### Contact Person

Elizabeth Platt

### Date of Summary Preparation

March 18, 1999

### Device (Trade & Common Name)

Lyphochek Urine Metals Control

### Classification Name

Class I, 91DIE  
CFR 862.3280: Materials, Heavy Metals Control

### Devices to Which Substantial Equivalence is Claimed

Lyphochek Urine Metals Control  
Bio-Rad Laboratories  
Irvine, California  
K891691/A

### Statement of Intended Use

Lyphochek Urine Metals Control is intended for use as an assayed quality control urine to monitor the precision of laboratory testing procedures for the analytes listed in the package insert.



Description of the Device

Lyphocek Urine Metals Control is prepared from human urine with added trace elements, heavy metals and other chemical compounds. The control is provided in lyophilized form for increased stability.

Statement of How Technological Characteristics Compare to Substantial Equivalence Device

A table is provided below comparing the similarities between the Bio-Rad Lyphocek Urine Metals Control and the devices to which substantial equivalence is claimed.

	Bio-Rad Lyphocek Urine Metals Control <b>New Submission</b>	Bio-Rad Lyphocek Urine Metals Control <b>Prior Submission</b>
Intended Use	An assayed quality control urine to monitor the precision of laboratory testing procedures for the analytes listed in the package insert.	An assayed quality control urine to monitor the precision of laboratory testing procedures for the analytes listed in the package insert.
Form	Lyophilized	Lyophilized
Matrix	Human urine.	Human urine.
Storage	2-8°C	2-8°C
Open Vial Claim	5 days when stored tightly capped at 2-8°C.	5 days when stored tightly capped at 2-8°C.
Differences	<ul style="list-style-type: none"><li>• Domestic and International products have been combined.</li><li>• Same claims as the substantially equivalent device plus claims for: Hippuric Acid, Mandelic Acid, Pentachlorophenol, Phenol, Phenylglyoxylic Acid, and Trichloroacetic Acid.</li></ul>	<ul style="list-style-type: none"><li>• Domestic and International products are separate.</li><li>• Refer to substantially equivalent package insert.</li></ul>



APR 16 1999

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Elizabeth Platt  
Regulatory Affairs Supervisor  
Bio-Rad Laboratories  
Diagnostics Group  
9500 Jeronimo Road  
Irvine, California 92618-2017

Re: K990928  
Trade Name: Lyphocek Urine Metals Control  
Regulatory Class: I  
Product Code: DIE  
Dated: March 18, 1999  
Received: March 19, 1999

Dear Ms. Platt:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.


Page 2

Under the Clinical Laboratory Improvement Amendments of 1988 (CLIA-88), this device may require a CLIA complexity categorization. To determine if it does, you should contact the Centers for Disease Control and Prevention (CDC) at (770) 488-7655.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification"(21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597, or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Steven I. Gutman, M.D, M.B.A.  
Director  
Division of Clinical  
Laboratory Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number: K990928

Device Name: Lyphochek Urine Metals Control


Indications for Use:

Lyphochek Urine Metals Control is intended for use as an assayed quality control urine to monitor the precision of laboratory testing procedures for the analytes listed in the package insert.

(PLEASE DO NOT WRITE BELOW THE LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

(Concurrence of CDRH, Office of Device Evaluation)



(Division Sign-Off)  
Division of Clinical Laboratory Devices  
510(k) Number \_\_\_\_\_

Prescription Use

OR Over-The Counter Use \_\_\_\_\_