

JUL 29 1999

MiniMed Inc.
Premarket [510(k)] Notification
MiniMed® Quick-serter™ Model 395



K992300
MiniMed®

Part C. 510(k) Summary

Submitter: MiniMed® Inc. 12744 San Fernando Rd., Sylmar, California 91342

Contact: Don Selvey, Regulatory Affairs (818) 362-5958, ext. 3011; (480) 704-8070 (v/f)

Name of Device: MiniMed Quick-serter™ infusion set insertion system


Predicate Device: MiniMed Sof-serter™ infusion set insertion system

Description of the New Device: The MiniMed Quick-serter is a manually operated, spring-loaded infusion set insertion device. It is similar in features and intended use to the MiniMed Sof-serter infusion set insertion system. The device has been designed for use exclusively with the Maersk Medical A/S Contour infusion sets. The device is contraindicated for use with other infusion sets.

Intended Use of the New Device: The MiniMed Quick-serter infusion set insertion system is intended to help make insertion of the Contour infusion sets simpler and with minimal discomfort. Use of the device may improve the user's consistency of infusion set insertion.

Comparison of the Technological Features of the New Device and Predicate Device: Technologically, both insertion devices are spring-driven devices that require the user to load the infusion set into a carrier, compress a spring, then activate the device by depressing the release or trigger buttons. Both devices are made of plastic with metal springs. The differences between the new device and the predicate device are limited to differences required for compatibility with the specific infusion set to be inserted. These modifications do not affect the safety or effectiveness of the device.

Signed,


Don Selvey
Senior Regulatory Affairs Specialist
Department of Clinical Research and Regulatory Affairs
MiniMed Inc.

7-6-99
date

®MiniMed is a Registered Trademarks of MiniMed Inc.
™Sof-serter and Quick-serter are Trademarks of MiniMed Inc.

000004



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUL 29 1999

Mr. Don Selvey
Senior Regulatory Affairs Specialist
Department of Clinical and Regulatory Affairs
MiniMed® Technologies, Incorporated
12744 San Fernando Road
Sylmar, California 91342

Re: K992300
Trade Name: MiniMed® Quick-serter™ Infusion Set
Insertion System, Model 395
Regulatory Class: II
Product Code: KZH
Dated: July 7, 1999
Received: July 8, 1999

Dear Mr. Selvey:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Page 2 - Mr. Selvey

Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control,
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure



INDICATIONS FOR USE

510(k) Number:

Device Name: MiniMed Quick-serter infusion set insertion system

Indications for Use:

The Quick-serter infusion set insertion system (Quick-serter), model 395, is indicated for use as an aid for insertion of the Maersk Medical A/S Contour infusion sets.

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use _____
(Per 21 CFR 801.109)

or

Over-the-Counter Use _____ ✓

Deborah Hubbard for Pat Cicerone

(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number K992300

000007