MAR - 2 2000 K993117

September 1, 1999

[1] 510(k) Summary of Safety and Effectiveness Information

[2] Safeskin Corporation12671 High Bluff DriveSan Diego, CA 92130

Telephone: 619-794-8111 Fax: 619-350-2382

Contact: Eugene V. Goorchenko

Telephone: 619-509-7010 Fax: 619-350-2382

[3] Trade Names: "Safeskin PFE-XTRA Latex Examination Glove"
Common Name: Patient Examination Glove
Classification Name: Patient Examination Glove

- [4] The predicate is a Safeskin powder free examination glove which meets all of the requirements of ASTM D 3578-00, Standard Specification for Rubber Examination Gloves for Medical Application.
- [5] The Safeskin PFE-XTRA Examination Glove will meet all the current specifications for ASTM D 3578-99.
- [6] The Safeskin PFE-XTRA Examination Glove is a disposable device intended to be worn by healthcare and similar personnel to prevent contamination between such personnel and the patient and for use with chemotherapeutic agents.
- [7] The Safeskin PFE-XTRA Examination Glove possesses the following technological characteristics (as compared to ASTM or equivalent standards):

<u>Characteristics</u> <u>Standards</u>

Dimensions Meets ASTM D 3578-99

Physical Properties Meets ASTM D 3578-99

Freedom from pinholes Meets ASTM D 3578-99

Meets ASTM D 5151

Powder-Free Meets ASTM D 6124

4 mg/glove maximum

Biocompatability

Primary Skin Irritation in Rabbits Passes
Guinea Pig Sensitization Passes

- [8] The performance test data that support a determination of substantial equivalence are described above.
- [9] Clinical data are not needed for examination gloves.
- [10] It can be concluded that the Safeskin PFE-XTRA Examination Glove is safe and effective for use with chemotherapeutic agents and will perform according to the glove performance standards referenced in Section 7 above, thereby meeting ASTM standards, FDA requirements, and the labeling claims for the product. Consequently, this device is substantially equivalent to currently marketed devices.



MAR - 2 2000

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Mr. Eugene V. Goorchenko Director of Regulatory Affairs Safeskin Corporation 12671 High Bluff Drive San Diego, California 92130

Re: K993117

Trade Name: Safeskin Pfe-xtra® Latex Examination Glove

with Protein Labeling Claim (50 microgram

or less)

Regulatory Class: I Product Code: LYY

Dated: December 29, 1999 Received: December 30, 1999

Dear Mr. Goorchenko:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any

the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4690. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "http://www.fda.gov/cdrh/dsmamain.html".

Sincerely yours

Timothy A. Ulatowski

Director

Division of Dental, Infection Control, and General Hospital Devices Office of Device Evaluation Center for Devices and

Radiological Health

Enclosure

l'agcoi
INDICATIONS FOR USE
Applicant: Safeskin Corporation
510(k) Number (if known): K993117
Device Name:
Device Name: Safeskin PFE-XTRA Powder-Free Latex Examination Glove
Protein Claim 50 mcg or less
Indications For Use:
A medical examination glove intended to be worn on the hands of
A medical examination grove intended to be worn on the hands of
healthcare personnel and similar personnel to prevent contamination
between such personnel and the patient and for use with chemotherapeutic agents.
(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)
Comment of CDDH OF CD in File (CDD)
Concurrence of CDRH, Office of Device Evaluation (ODE)
(Jun Silvan
(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices
510(k) Number 49311
Confidential Use
THE REPORT OF THE PROPERTY OF

(Ose as. I Formet 1-2-96)