

**Summary of 510(k) Submission**

---

<b>Name and address of submitter</b>	Vistakon, Johnson & Johnson Vision Products, Inc. 4500 Salisbury Road, Suite 300 Jacksonville, Florida 32216 Contact: Martine D. Martino Phone: (904) 443-1808 Date of Original: April 22, 1999 Date Modified: May 26, 1999
<hr/>	
<b>Identification of device</b>	<ul style="list-style-type: none"><li>• The trade name is VISTAKON (lenefilcon A) soft (hydrophilic) contact lenses, clear and visibility tint, with UV blocker, for daily wear</li><li>• The common or usual name is Soft (hydrophilic) Contact Lens (daily wear).</li><li>• The FDA Classification is Class II.</li></ul>
<hr/>	
<b>Predicate devices</b>	The predicate device is VISTAKON (lenefilcon A) soft (hydrophilic) contact lenses, clear and visibility tint, with UV blocker, for daily wear covered under K983912.
<hr/>	
<b>Description of device</b>	The device description does not change from that cleared under K983912.

---

*Continued on next page*

## Summary of 510(k) Submission, Continued

---

**Intended use** The reason for the 510(k) is to add the following Indications Statement (see **bold text**) for VISTAKON (lenefilcon A) soft (hydrophilic) contact lenses, clear and visibility tint, with UV blocker, for daily wear. The revised Indications Statement is included with this 510(k):

The VISTAKON Contact Lens (spherical) is indicated for daily wear for the correction of refractive ametropia (myopia and hyperopia) in aphakic or non-aphakic persons with non-diseased eyes who may have 1.00 D or less of astigmatism.

The VISTAKON BIFOCAL Contact Lens is indicated for daily wear for the correction of distance and near vision in presbyopic, aphakic or non-aphakic persons with non-diseased eyes who may have 0.75 D or less of astigmatism.

The VISTAKON TORIC Contact Lens is indicated for daily wear for the correction of visual acuity in aphakic or non-aphakic persons with non-diseased eyes that are hyperopic or myopic and may have 10.00 D or less of astigmatism.

The VISTAKON TORIC BIFOCAL Contact Lens is indicated for daily wear for the correction of distance and near vision in presbyopic aphakic or non-phakic persons with non-diseased eyes who may have 10.00D of astigmatism or less.

**VISTAKON UV Blocking Contact Lenses help protect against transmission of harmful UV radiation to the cornea and into the eye.**

Eye care practitioners may prescribe the lens for either single-use disposable wear or for frequent/ planned replacement wear, with cleaning, disinfection and scheduled replacement. When prescribed for frequent/planned replacement wear, the lens may be disinfected using a chemical disinfection system only.

---

*Continued on next page*

## Summary of 510(k) Submission, Continued

---

<b>Characteristics</b>	The characteristics do not change. They are the same as previously submitted in K983912.
<b>Non-clinical studies</b>	<ul style="list-style-type: none"><li>• Non-clinical studies (microbiology, toxicology, chemistry, shelf-life, and leachability) on the lens material were not conducted because the lens material, lenefilcon A, does not change.</li><li>• Non-clinical studies were conducted to support the following modified Actions in the labeling.<ul style="list-style-type: none"><li>• The average UV blocking for VISTAKON (lenefilcon A) Contact Lenses with UV Blocker for UVA is 95% and for UVB is 99.9%.</li></ul></li></ul>
<b>Clinical studies</b>	This 510(k) describes a labeling modification: an additional Indications Statement. There is no change in lens material, the manufacturing process, nor the parameters and properties, therefore, the clinical data previously submitted in K983912 supports the clinical safety of the subject device.
<b>Conclusions drawn from studies</b>	Additional studies were not conducted, therefore, the conclusions drawn from studies previously submitted in K983912 support the non-clinical and clinical safety of the subject device.

---



JUN 8 1999

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Martin D. Martino  
Regulatory Affairs Specialist  
Vistakon  
Johnson & Johnson Vision Products, Inc.  
4500 Salisbury Road  
Suite 300  
Jacksonville, FL 32216

Re: K991163  
Trade Name: VISTAKON (lenefilcon A) Soft (hydrophilic) Contact Lens, Clear and  
Visibility Tint, with UV Blocker, for daily Wear  
Regulatory Class: II  
Product Code: 86 LPL  
Dated: April 6, 1999  
Received: April 7, 1999

Dear Ms. Martino:


We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4613. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

A handwritten signature in cursive script, reading "A. Ralph Rosenthal".

A. Ralph Rosenthal, M.D.

Director

Division of Ophthalmic Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

K991163

## Indications Statement

510(k) Number: K991163

Device Name: VISTAKON (lenefilcon A) Soft (hydrophilic) Contact Lens Clear and Visibility Tint with UV Blocker

The VISTAKON Contact Lens (spherical) is indicated for daily wear for the correction of refractive ametropia (myopia and hyperopia) in aphakic or non-aphakic persons with non-diseased eyes who may have 1.00 D or less of astigmatism.

The VISTAKON BIFOCAL Contact Lens is indicated for daily wear for the correction of distance and near vision in presbyopic, aphakic or non-aphakic persons with non-diseased eyes who may have 0.75 D or less of astigmatism.

The VISTAKON TORIC Contact Lens is indicated for daily wear for the correction of visual acuity in aphakic or non-aphakic persons with non-diseased eyes that are hyperopic or myopic and may have 10.00 D or less of astigmatism.

The VISTAKON TORIC BIFOCAL Contact Lens is indicated for daily wear for the correction of distance and near vision in presbyopic aphakic or non-phakic persons with non-diseased eyes who may have 10.00D of astigmatism or less.

VISTAKON UV Blocking Contact Lenses help protect against transmission of harmful UV radiation to the cornea and into the eye.

Eye care practitioners may prescribe the lens for either single-use disposable wear or for frequent/ planned replacement wear, with cleaning, disinfection and scheduled replacement. When prescribed for frequent/planned replacement wear, the lens may be disinfected using a chemical disinfection system only.

PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)



(Division Sign-Off)

Division of Ophthalmic Devices



Prescription Use ✓

510(k) Number K991163

OR

Over the Counter