JUL 26 1999 June 9, 1999

K991951

[1] 510(k) Summary of Safety and Effectiveness Information

Safeskin Corporation [2] 12671 High Bluff Drive San Diego, CA 92130

Telephone:

619-509-7010

Fax:

619-350-2382

Contact:

Eugene V. Goorchenko

Telephone:

619-509-7010

Fax:

619-350-2382

[3] Trade Name: Safeskin Synthetic Powder-Free Examination Gloves

Common Name:

Patient Examination Gloves

Classification Name: Patient Examination Gloves

[4] The predicate is a powder-free synthetic (vinyl) examination glove which meets all of the requirements of ASTM D 5250-99.

- [5] The Safeskin Synthetic Powder-Free Examination Glove will meet all the current specifications for ASTM D 5250-99, vinyl examination gloves.
- [6] Safeskin Synthetic Powder-Free Examination Gloves are disposable medical gloves intended to be worn by healthcare and similar personnel to prevent contamination between such personnel and the patient.
- [7] Safeskin® Synthetic Powder-Free Examination Gloves possess the following technological characteristics (as compared to ASTM or equivalent standards):

Characteristics Standards

**Dimensions** Meets ASTM D 5250

**Physical Properties** Meets ASTM D 5250

Freedom from pinholes Meets ASTM D 5250

Meets ASTM D 5151

Powder-Free Meets ASTM D 6124

2 mg/glove maximum

## Biocompatability

Primary Skin Irritation in Rabbits Passes
Guinea Pig Sensitization Passes

- [8] The performance test data that support a determination of substantial equivalence are described above.
- [9] Clinical data are not needed for surgical gloves.
- [10] It can be concluded that the Safeskin® Synthetic Powder-Free Examination Glove will perform according to the glove performance standards referenced in Section 7 above and therefore will meet ASTM standards, FDA requirements, and the labeling claims for the product. Consequently, this device is substantially equivalent to currently marketed devices.



JUL 26 1999

Food and Drug Administration 9200 Corporate Boulevard Rockville MD 20850

Mr. Eugene V. Goorchenko Director of Regulatory Affairs Safeskin Corporation 12671 High Bluff Drive San Diego, California 92130

Re: K991951

Trade Name: Safeskin Synthetic Vinyl Powder-Free Polymer

Coated Examination Gloves (green, white)

Regulatory Class: I Product Code: LYZ Dated: June 9, 1999 Received: June 10, 1999

Dear Mr. Goorchenko:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "http://www.fda/gov/cdrh/dsmamain.html".

Sincerely yours

Timothy A. Ulatowski

Director

Division of Dental, Infection Control, and General Hospital Devices Office of Device Evaluation Center for Devices and Radiological Health

Enclosure

## INDICATIONS FOR USE

Applicant:

Safeskin Corporation

510(k) Number:

K991951

Device Name: Viny Synthetic powder-free examination glove (Ving), polymer Coaled)

Indications for Use:

A medical glove intended to be worn on the hands of healthcare and similar personnel to prevent contamination between such personnel and the patient.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH Office of Device Evaluation (ODE)

Over-The-Counter X

(Division Sign-Off)

Division of Dental, Infection Control,

9