

SEP 10 1999

K992210

August 25, 1999

- [1] 510(k) Summary of Safety and Effectiveness Information
- [2] Safeskin Corporation  
12671 High Bluff Drive  
San Diego, CA 92130
- Telephone: 619-794-8111  
Fax: 619-350-2382
- Contact: Eugene V. Goorchenko  
Telephone: 619-509-7010  
Fax: 619-350-2382
- [3] Trade Name: ShieldMaster Powdered Purple Nitrile Gloves  
Common Name: Medical Gloves  
Classification Name: Patient Examination Gloves
- [4] The predicate is a Safeskin blue Nitrile Glove which meets all of the requirements of ASTM D 3578-99, Standard Specification for Rubber Examination Gloves (with the exception of elongation).
- [5] The ShieldMaster Powdered Purple Nitrile Glove will meet all the current specifications for ASTM D 3578-99 (with the exception of elongation).
- [6] ShieldMaster Powdered Purple Nitrile Gloves are disposable devices intended to be worn by healthcare and similar personnel to prevent contamination between such personnel and the patient.
- [7] ShieldMaster Powdered Purple Nitrile Gloves possess the following technological characteristics (as compared to ASTM or equivalent standards):

Characteristics

Dimensions

Physical Properties

Freedom from pinholes

Standards

Meets ASTM D 3578-99

Meets ASTM D 3578-99 and  
ASTM D 6319-99

Meets ASTM D 3578-99  
Meets ASTM D 5151

## Biocompatibility

Primary Skin Irritation in Rabbits

Passes

Guinea Pig Sensitization

Passes

- [8] The performance test data that support a determination of substantial equivalence are described above.
- [9] Clinical data are not needed for medical gloves.
- [10] It can be concluded that the ShieldMaster Powdered Purple Nitrile Glove will perform according to the glove performance standards referenced in Section 7 above and therefore will meet ASTM standards, FDA requirements, and the labeling claims for the product. Consequently, this device is substantially equivalent to currently marketed devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

SEP 10 1999

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Eugene V. Goorchenko  
Director of Regulatory Affairs  
Safeskin Corporation  
12671 High Bluff Drive  
San Diego, California 92130

Re: K992210  
Trade Name: Shieldmaster Powdered Purple Nitrile  
Examination Gloves  
Regulatory Class: I  
Product Code: LZA  
Dated: June 30, 1999  
Received: July 1, 1999

Dear Mr. Goorchenko:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of


Page 2 -Mr. Goorchenko

the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4690. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski  
Director  
Division of Dental, Infection Control,  
and General Hospital Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

K992210

FDA 510(k), Premarket Notification

Safeskin Corporation

### INDICATIONS FOR USE

Applicant: Safeskin Corporation

510(k) Number:

Device Name: Powdered Purple Nitrile Medical Glove

Examination  
^

Indications for Use:

A medical glove intended to be worn on the hands of healthcare and similar personnel to prevent contamination between such personnel and the patient.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE  
IF NEEDED)

Concurrence of CDRH Office of Device Evaluation (ODE)

Over-The-Counter X

Chin S. Lin  
(Division Sign-Off)  
Division of Dental, Infection Control,  
and General Hospital Devices  
510(k) Number K992210