

DEC 29 1999

510(k) Premarket Notification  
October 20, 1999

K993575  
**510(k) Summary**

**Trade Name:** ReSolve™ QuickAnchor

**Sponsor:** Mitek Products  
60 Glacier Drive  
Westwood, MA 02090  
Registration #1221934

**Contact:** Paula E. Bulger  
Regulatory Affairs Project Manager  
Mitek Products  
60 Glacier Drive  
Westwood, MA 02090  
Phone: (781) 251-2700  
Fax: (781) 461-9166

**Device Generic Name:** Staple, Fixation, Bone

**Classification:** According to Section 513 of the Federal Food, Drug, and Cosmetic Act, the device classification is Class II.

**Product Code:** JDR (21 CFR 888.3030)

**Predicate Devices:** K970896 - Mitek Panalok Anchor  
K962511, K982420 - Mitek Micro Anchor  
K992611 - Mitek Rotator Cuff QuickAnchor Plus

**Product Description:** The device described in this 510(k) is a sterile implant used to anchor or lock suture within pre-drilled bone sites and firmly secure soft tissue to bone.

**Indications for Use:**

The ReSolve™ QuickAnchor is used for the fixation of absorbable monofilament surgical suture to bone.

This product is intended for the following indications:

Cranio/Maxillofacial:

Repair, reconstruction or reattachment of tendons, ligaments, muscles and soft tissue flaps to the parietal, temporal ridge, frontal, mandible, maxilla, zygoma and periorbital bones of the skull.

**Safety and Performance:**

The following safety and performance data has been provided to support substantial equivalence of ReSolve™ QuickAnchor:

**Performance testing:** Pullout force (preserved human cadaver skull)  
Strength comparison (ReSolve vs. Mitek Micro Anchor)  
Skull thickness measurements

**Conclusion:**

Based on safety and performance data, similarities in design, operating principle, materials, biocompatibility and sterilization method, the ReSolve™ QuickAnchor has been shown to be substantially equivalent to predicate devices under the Federal Food, Drug and Cosmetic Act.

000249



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

DEC 29 1999

Ms. Paula E. Bulger  
Regulatory Affairs Project Manager  
Mitek Products  
Ethicon, Inc.  
60 Glacier Drive  
Westwood, MA 02090

Re: K993575

Trade Name: Resolve™ QuickAnchor (Cranio/Maxillofacial)  
Regulatory Class: II  
Product Code: DZL  
Dated: October 20, 1999  
Received: October 21, 1999

Dear Ms. Bulger:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any

Page 2 - Ms. Bulger

obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski *for*  
Director  
Division of Dental, Infection Control  
and General Hospital Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

Page 1 of 1

510(k) Number (if known): K993575

Device Name: ReSolve™ QuickAnchor

Indications for Use:

The ReSolve™ QuickAnchor is used for the fixation of absorbable monofilament polydioxanone surgical suture to bone. This product is intended for the following indications:

Cranio/Maxillofacial:

Repair, reconstruction or reattachment of tendons, ligaments, muscles and soft tissue flaps to the parietal, temporal ridge, frontal, mandible, maxilla, zygoma and periorbital bones of the skull.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓  
(Per 21 CFR 801.109)

OR

Over-the -Counter Use \_\_\_\_\_

(Division Sign-Off) Pamela Scott for Susan Runner  
Division of Dental, Infection Control,  
and General Hospital Devices  
510(k) Number K993575

000001